


# LEYBOLD®


## PHYSICS EXPERIMENTS



PHYSICS

- MECHANICS
- HEAT
- ELECTRICITY
- ELECTRONICS
- OPTICS
- ATOMIC AND NUCLEAR PHYSICS
- SOLID-STATE PHYSICS

# PHYSICS EXPERIMENTS IN TRUSTED LEYBOLD QUALITY



Experiments have become an indispensable part of education. Indispensable because the combination of theoretical knowledge with experimental learning sessions ensures sustainable and successful learning.

We provide a wide range of high-quality experiments from all areas of physics. Our proven LEYBOLD quality guarantees durability and safety during the complete experiment setup.

MORE THAN  
500 EXPERIMENTS  
IN VARIOUS  
PHYSICS RANGES



FURTHER EXPERIMENTS  
ARE AVAILABLE AT  
OUR WEBSITE UNDER

[WWW.LEYBOLD-SHOP.COM](http://WWW.LEYBOLD-SHOP.COM)


These are available in several versions (e.g. with or without PC support) which can be adapted to the time requirements and student's level of knowledge. The experiments are supported by manuals, which includes clear step by step instructions, sample measurements and safety advices.

We will advise you personally and demonstrate our experiments for you.

Have fun experimenting!

# CONTENT

TRADITION AND INNOVATION	II-III
MOBILE-CASSY 2	IV-V
X-RAY FUNDAMENTAL & PROFESSIONAL	VI-VII
PHOTONICS	VIII
LOGIC BOARDS	IX
DOCUMENT CENTER	X-XI
OVERVIEW OF EXPERIMENTS	XII-XIII
HOW TO USE THIS CATALOGUE	XIV



P1	MECHANICS	1
P2	HEAT	63
P3	ELECTRICITY	87
P4	ELECTRONICS	149
P5	OPTICS	165
P6	ATOMIC AND NUCLEAR PHYSICS	229
P7	SOLID-STATE PHYSICS	275
	REGISTER	293

# TRADITION AND INNOVATION

EXPERIMENTS FOR STUDENT PRACTICALS AND  
DEMONSTRATIONS FOR MORE THAN 160 YEARS



1898



1973



1991



2004

The LD DIDACTIC Group is a world-leading manufacturer of high-quality scientific and technical training systems.

Our teaching systems have been making a decisive contribution to the transfer of knowledge in schools, further education colleges and universities, and also in the industrial sector, for generations.

Founded in 1850 in Cologne, LD DIDACTIC can look back on 160 years of company history. This former subsidiary of LEYBOLD, a Cologne company of long-standing tradition, sees itself as an innovative supplier of quality products and sells these products and complete solutions under the brand names LEYBOLD, ELWE Technik and FEEDBACK.

## LEYBOLD®



2014

Our many years of experience and our innovative technical potential combined with the close cooperation with teachers and trainers from the relevant fields enable us to offer targeted solution to our customers and at the same time make complex topics in biology, chemistry, physics and technology transparent for the student.

DEMO-MULTIMETER



## EXPERIMENTATION WITH PASSION

Our staff develop with great passion tailor-made solutions for our customers. Tailor-made because they are matched to school type and teaching content, but also because they have been developed in close cooperation with teachers and trainers.

Our aim is to make the everyday planning and execution of lessons easier for teachers. We also aspire to teach complex topics to students clearly and transparently and to awaken their enthusiasm for scientific subjects.

## PREMIUM PRODUCTS „made in Germany“

We manufacture high-quality products and complete solutions for our customers daily at two sites in Germany.

Our tried and trusted LEYBOLD quality is recognised not only in Germany, but also worldwide, and is our guarantee for technically demanding products of the highest quality standard.



## VISIT US

In our showrooms in Huerth and Mannheim you are able to experience our products and complete solutions first hand.

We cordially invite you to visit us.

# MOBILE-CASSY 2

THE ONLY ONE MEASURING DEVICE YOU WILL NEED FOR MEASURED VALUES IN PHYSICS

# LEYBOLD®

## MOBILE-CASSY 2

(524 005)

### THE MEASURING DEVICE FOR STUDENT EXPERIMENTS IN NATURAL SCIENCES

The new generation of our Mobile-CASSY impresses through its


- universal connection possibilities,
- intuitive operation,
- fast recording of measured values and
- its graphics capability

This results in a wide range of applications from student experiments to demonstrations.



Set-up of the experiment: Compound pendulum with rotation sensor S

## THE CORE TECHNICAL DATA





HIGH-CONTRAST DISPLAY SHOWN AT ORIGINAL SIZE AND RESOLUTION

From the display icons, above the measured value, the measurement channels are shown and can be selected. Whilst the icons below on the right are accessing the set-up menus ▲ and other displays through the touch wheel ⌚.

Graphics display:	9 cm (3.5"), colour QVGA
Operation:	large capacitive touch wheel (42 mm)
Resolution:	12 bit
Integrated measurement ranges Voltage:	$\pm 0.1/\pm 0.3/\pm 1/\pm 3/\pm 10/\pm 30$ V
Integrated measurement ranges Current:	$\pm 0.03/\pm 0.1/\pm 0.3/\pm 1/\pm 3$ A
Integrated measurement ranges Temperature:	-200 ... +200 °C / -200 ... +1200 °C
Sensor connections:	2 (CASSY sensors and sensors M)
Sampling rate:	max. 100,000 values/second
Time resolution of the timer inputs:	20 nanoseconds
Loudspeaker:	integrated for key tones and GM counter tube (can be disabled as required)
Data storage device:	integrated micro SD card for more than a thousand measurement files and screen shots, optionally also via a USB stick
Battery capacity:	14 watt-hours (AA size, replaceable)

# THE UNIVERSAL CONNECTION POSSIBILITIES



LEYBOLD®

- 1 2 CASSY SENSORS**  
all CASSY sensors and sensor boxes are supported
- 2 VOLTAGE AND CURRENT**  
directly via 4 mm safety sockets
- 3 2 NEW SENSORS M**  
e.g. for the light barriers in the Advanced Science Kit Set MEC 6
- 4 NiCr-Ni TEMPERATURE PROBE**  
via a type K socket
- 5 USB STICK**  
for the simple transfer of measured data and screen shots
- 6 PC**  
via USB cable with full CASSY Lab 2 support
- 7 CHARGER (SUPPLIED)**  
with status display
- 8 KENSINGTON LOCK**  
as anti-theft protection
- 9 WIRELESS LAN (OPTIONAL)**  
for wireless data transfer e.g. to a PC with CASSY Lab 2 or a tablet

# FOR EVERY REQUIREMENT AND EVERY BUDGET

The LEYBOLD X-ray system has a modular structure and enables the individual configuration of the separate appliances, so that you only buy what you actually need.

In addition to the basic equipment, you can choose your accessories for basic experiments (FUNDAMENTAL Experiments) or advanced applications (PROFESSIONAL Experiments) depending on the experiment requirements.



## BASIC EQUIPMENT

### X-RAY APPARATUS

The X-ray apparatus is available in two variants – as a basic apparatus or as a complete apparatus with a Mo tube, goniometer and NaCl monocrystal. If you wish to use other tubes, the X-ray basic apparatus is the most flexible solution.

You can extend the X-ray apparatus with a drawer for your accessories irrespective of this.

### Goniometer

No matter whether you are interested in Bragg spectra, X-ray energy spectra or computed tomography, you will be happy with the precision and high resolution of the goniometer.



### TUBES

In addition to the Mo tube, there are other tubes, which are more suitable for special areas of application, e.g. Cu tube for Debye-Scherrer diagrams, Ag tube for X-ray fluorescence due to its high energy K-lines, W or Au tubes for radiation and computed tomography due to their high intensity.






## X-RAY ENERGY DETECTOR

The X-ray energy detector exposes energy-dispersive X-ray spectra with the CASSY system. Using the X-ray energy spectra, various chemical elements can be easily distinguished by means of their characteristic X-ray radiation and their mass fraction can also be determined. This also confirms the Compton effect.

# FUNDAMENTAL EXPERIMENTS

## COMPUTED TOMOGRAPHY MODULE

The X-ray apparatus has a laterally installed fluorescent screen, on which X-rays can be seen directly. The computed tomography module captures this visible X-ray image and the software provided there controls the rotation of the object in the X-ray apparatus through 360° and carries out the back projection of the X-ray image for the real-time 3D reconstruction.



## BRAGG SPECTRA

With the complete apparatus with a Mo tube, you have everything you need for capturing your first Bragg spectrum. Other available monocrystals and/or X-ray tubes offer several possible variations.

# PROFESSIONAL EXPERIMENTS

## COMPUTED TOMOGRAPHY PRO

If the resolution of the computed tomography module is insufficient, the X-ray image sensor provides the solution with its megapixel resolution and its 12-bit grey scale. With this, voxel resolutions up to an edge length of below 50  $\mu\text{m}$  are possible. The X-ray image sensor also provides a comfortable and fast solution for Laue diagrams.



## HD accessory

With the HD accessory consisting of a high-resolution collimator and counter tube holder with narrow gaps and new software, the angular resolution of the goniometer increases to 0.01°. Bragg spectra with a 4-fold resolution are possible with this.


## GOLD TUBE

The gold tube is the Leybold tube with the highest intensity. It is particularly suitable for capturing X-ray images, Laue diagrams or computed tomograms. In addition to the Tungsten tube, it is also the only tube whose Bragg spectra consist of L-lines.

# PHOTONICS

## EDUCATIONAL KITS

- Laser and advanced optics
- For higher education and university level
- Wide range of experiments with educational manuals
- Sophisticated topics, easy to setup





## BASICS


- Gas, Fibre and Solid-state Lasers
- open frame setup
- hands-on experiments

### SHORT LIST OF TOPICS

- He-Ne Laser and Nd:YAG Laser from setup to physics
- Semi conductor lasers
- Interferometers (Michelson, Mach-Zehnder, Fabry-Perot)
- Heterodyne Interferometer
- Non-linear optics
- Special relativity (Laser Gyroscope)



FOR MORE INFORMATION ON OUR PHOTONICS RANGE,  
PLEASE REFER TO CHAPTER P5.8 (STARTING PAGE 206)  
IN THIS CATALOG.



# LOGIC BOARDS

**LEYBOLD®**

THE ENTRY INTO THE DIGITAL ELECTRONICS!

## LOGIC BOARD 1 (571 401)

Introducing the basic logic gates (AND, OR, NOT, NAND, XOR) used in digital electronics. These are used to investigate the laws of logical operations (de Morgan's law, associative law and distributive law) and non-feedback logic circuits (switch networks). Finally simple flip-flop circuits with feedback are assembled to study storage of information.

Switch states are indicated by means of an LED at each output.

### GATES:

- AND
- OR
- NOT
- NAND
- XOR

### FLIP-FLOPS:

- RS-flip-flop
- D-flip-flop
- RC module for construction of a multivibrator

## LOGIC BOARD 2 (571 402)

The second board is used for advanced topics. The adders are investigated as practical examples of combinatorial logic (logic circuits without feedback). Various flip-flop circuits add students' knowledge on circuits with feedback like shift registers or latches. Applications of digital technology will be investigated, e.g. multiplexing, demultiplexing and the topics of digital-to analogue and analogue-to-digital conversion are covered.

### GATES:

- NOT
- NAND

### FLIP-FLOPS:


- RS-flip-flop
- D-flip-flop

### ADDITIONAL:

- RC module for construction of a multivibrator
- Adders
- AD converter/DA converter
- 7 segment display

FOR MORE INFORMATION ON OUR LOGIC BOARDS, PLEASE REFER TO CHAPTER P4.5.1 TO 4.5.3 IN THIS CATALOG.

# EXPERIMENT INSTRUCTIONS



## THE DOCUMENT CENTER IS THE ELECTRONIC LEYBOLD LIBRARY

- For information on experiments for students and demonstration experiments
- For operation manuals



DOWNLOAD THE  
DOCUMENT CENTER  
FREE OF CHARGE:  
[WWW.LD-DIDACTIC.COM](http://WWW.LD-DIDACTIC.COM)

# IN THE DOCUMENT CENTER


## THE DOCUMENT CENTER OFFERS

- Leaflets (demonstration instructions) with link to CASSY LAB 2 and Spektralab
- Experiments for students as interactive pdf files:
  - Easy-to-understand worksheets for students
  - Complete information with experiment results for teachers
  - Ability to go from student to teacher version and back with one mouse click
  - Student documents can be filled out on the computer and stored or printed out as a protocol
- Sorted into literature packages – facilitates and encourages the compilation of own test series
- Free-of-charge online update of literature packages following initial acquisition

## HOW DOES IT WORK?

The literature packages are clearly displayed in a table of contents, which is structured to guide you to the document you need. The more literature packages you have installed, the more entries are listed in the table of contents. Once the system has been installed, the documents can be set to update automatically if desired.

The convenient, fault-tolerant search function helps you find the right document rapidly.



Student version with typewriter tool for filling out protocols on the computer

Teacher version with an example of the solution and notes about the experiments

## P1 MECHANICS

page 1

### P1.1 Measuring methods

Measuring lengths, measuring volume and density, determining the gravitational constant

page 3

### P1.2 Forces

Static effects of forces, force as vector, lever, block and tackle, inclined plane, friction

page 6

### P1.3 Translational motions of a mass point

One-dimensional motions on Fletcher's trolley and on the linear air track, conservation of linear momentum, free fall, angled projection, two-dimensional motions on the air table

page 12

### P1.4 Rotational motions of a rigid body

Rotational motions, conservation of angular momentum, centrifugal force, Motions of a gyroscope, moment of inertia, conservation of Energy

page 26

## P2 HEAT

page 63

### P2.1 Thermal expansion

Thermal expansion of solids and liquids, thermal anomaly of water

page 65

### P2.2 Heat transfer

Thermal conductivity, solar collector

page 68

### P2.3 Heat as a form of energy

Mixing temperatures, heat capacities, converting mechanical and electrical energy into heat

page 70

### P2.4 Phase transitions

Latent heat and vaporization heat, measuring vapor pressure, critical temperature

page 74

## P3 ELECTRICITY

page 87

### P3.1 Electrostatics

Basic experiments on electrostatics, Coulomb's law, field lines and equipotential lines, effects of force in an electric field, charge distributions on electrical conductors, definition of capacitance, plate capacitor

page 89

### P3.2 Fundamentals of electricity

Charge transfer with drops of water, Ohm's law, Kirchhoff's laws, circuits with electrical measuring instruments, conducting electricity by means of electrolysis, experiments on electrochemistry

page 104

### P3.3 Magnetostatics

Basic experiments on magnetostatics, magnetic dipole moment, effects of force in a magnetic field, Biot-Savart's law

page 111

### P3.4 Electromagnetic induction

Voltage impulse, induction in a moving conductor loop, induction by means of a variable magnetic field, eddy currents, transformer, measuring the earth's magnetic field

page 115

## P4 ELECTRONICS

page 149

### P4.1 Components and basic circuits

Current and voltage sources, special resistors, diodes, diode circuits, transistors, transistor circuits, optoelectronics

page 151

### P4.2 Operational amplifier

Internal design of an operational amplifier, operational amplifier circuits

page 159

### P4.3 Open- and closed-loop control

Closed-loop control

page 161

### P4.5 Digital electronics

Simple combinations, logic circuits, analog inputs and outputs

page 162

## P5 OPTICS

page 165

### P5.1 Geometrical optics

Reflection and refraction, laws of imaging, image distortion, optical instruments

page 167

### P5.2 Dispersion and chromatics

Refractive index and dispersion, color mixing, absorption spectra, reflection spectra

page 171

### P5.3 Wave optics

Diffraction, two-beam interference, Newton's rings, Michelson interferometer, other types of interferometers, white-light reflection holography, transmission holography

page 176

### P5.4 Polarization

Basic experiments, birefringence, optical activity and polarimetry, Kerr effect, Pockels effect, Faraday effect

page 189

## P6 ATOMIC AND NUCLEAR PHYSICS

page 229

### P6.1 Introductory experiments

Oil-spot experiment, Millikan experiment, specific electron charge, Planck's constant, dual nature of wave and particle, Paul trap

page 231

### P6.2 Atomic shell

Balmer series of hydrogen, emission and absorption spectra, inelastic collisions of electrons, Franck-Hertz experiment, electron spin resonance, normal Zeeman effect, optical pumping (anomalous Zeeman effect)

page 239

### P6.3 X-rays physics

Detection of X-rays, attenuation of X-rays, Physics of the atomic shell, X-ray energy spectroscopy, structure of X-ray spectrums, Compton effect at X-rays, X-ray tomography

page 250

### P6.4 Radioactivity

Detecting radioactivity, Poisson distribution, radioactive decay and half-life, attenuation of  $\alpha$ -,  $\beta$ - and  $\gamma$  radiation

page 261

## P7 SOLID-STATE PHYSICS

page 275

### P7.1 Properties of crystals

Crystal structure, x-ray scattering, elastic and plastic deformation

page 277

### P7.2 Conduction phenomena

Hall effect, electrical conductivity in solids, photoconductivity, luminescence, thermoelectricity, superconductivity

page 281

### P7.3 Magnetism

Dia-, para- and ferromagnetism, ferromagnetic hysteresis

page 288

### P7.4 Scanning probe microscopy

Scanning tunneling microscope

page 290

**P1.5 Oscillations**

Simple and compound pendulum, harmonic oscillations, torsion pendulum, coupling of oscillations

page 32

**P1.6 Wave mechanics**

Transversal and longitudinal waves, wave machine, circularly polarized waves, propagation of water waves, interference of water waves, resonances

page 39

**P1.7 Acoustics**

Sound waves, oscillations of a string, wavelength and velocity of sound, reflection of ultrasonic waves, interference of ultrasonic waves, Acoustic Doppler effect, fourier analysis, ultrasound in media

page 45

**P1.8 Aero- and hydrodynamics**

Barometric measurements, buoyancy, viscosity, surface tension, introductory experiments on aerodynamics, measuring air resistance, measurements in a wind tunnel

page 55

**P2.5 Kinetic theory of gases**

Brownian motion of molecules, gas laws, specific heat of gases, real gases

page 77

**P2.6 Thermodynamic cycle**

Hot-air engine: qualitative and quantitative experiments, heat pump

page 81

**P3.5 Electrical machines**

Basic experiments on electrical machines, electric generators, electric motors, three-phase machines

page 122

**P3.6 DC and AC circuits**

Circuit with capacitor, circuit with coil, impedances, measuring-bridge circuits, measuring AC voltages and AC currents, electrical work and power, electromechanical devices

page 126

**P3.7 Electromagnetic oscillations and waves**

Electromagnetic oscillator circuit, decimeter-range waves, propagation of decimeter-range waves along lines, microwaves, propagation of microwaves along lines, directional characteristic of dipole radiation

page 134

**P3.8 Free charge carriers in a vacuum**

Tube diode, tube triode, Maltese-cross tube, Perrin tube, Thomson tube

page 140

**P3.9 Electrical conduction in gases**

Spontaneous and non-spontaneous discharge, gas discharge at reduced pressure, cathode rays and canal rays

page 145

**P5.5 Light intensity**

Quantities and measuring methods of lighting engineering, laws of radiation

page 195

**P5.6 Velocity of light**

Measurement according to Foucault/Michelson, measuring with short light pulses, measuring with an electronically modulated signal

page 198

**P5.7 Spectrometer**

Prism spectrometer, grating spectrometer

page 202

**P5.8 Photonics**

Basic Optics, optical applications, optical imaging and colour, laser basics, solid state laser, optical fibres, technical applications

page 206

**P6.5 Nuclear physics**

Demonstrating paths of particles, Rutherford scattering, Nuclear magnetic resonance,  $\alpha$  spectroscopy,  $\gamma$  spectroscopy, Compton effect, properties of radiation particles

page 265

**P6.6 Quantum physics**

Quantum optics, particles

page 272

**P7.5 Applied solid-state physics**

X-ray fluorescence analysis

page 292

# HOW TO USE THIS CATALOGUE

## PAGE STRUCTURE

Section  
Subsection

MECHANICS  
TRANSLATIONAL MOTIONS OF A MASS POINT

Experiment set-up  
Topic

P1.3.4  
CONSERVATION OF LINEAR  
MOMENTUM

Experiment

Each experiment is identified by a P and a four-digit number.

P1.3.4.5  
Newton's third law and laws of collision  
- Recording and evaluating with two ultrasonic motion sensors and CASSY



Newton's third law and laws of collision - Recording and evaluating with two ultrasonic motion sensors and CASSY (P1.3.4.5)

Short description of the experiment

Equipment list

Display of available  
CASSY® device

Experiment results

CASSY® marking


Cat. No.	Description	P1.3.4.5
337 130	Track, 1.5 m	1
337 110	Trolley	2
337 114	Additional weights, pair	1
337 473	Impact spring, soft	1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
524 0701	Ultrasonic motion sensor S	2
501 11	Extension cable, 15 pin	2
301 25	Support block	2
340 89ET5	Coupling plugs, 4 mm, set of 5	1
additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

In the experiment P1.3.4.5 the way of the two trolleys are measured during an elastic or inelastic collision. This is done with two ultrasonic sensors.  $v(t)$  and  $a(t)$ -diagrams are determined by the Sensor-CASSY. The forces  $F_1(t)$  and  $F_2(t)$  can be calculated from the accelerations  $a_1(t)$  and  $a_2(t)$  and the masses  $m_1$  and  $m_2$  of the trolleys. It is confirmed that  $F_1(t) = -F_2(t)$  applies during the collision. It follows the conservation of momentum during the whole collision process:

$$\Delta p_1 + \Delta p_2 = \int F_1(t) dt + \int F_2(t) dt = \int (F_1(t) + F_2(t)) dt = 0$$

In addition, the center of mass motion  $s_3(t)$  is observed and it can be shown that the speed  $v_3(t)$  of the center of gravity during the whole collision process is constant. This shows alternative the conservation of momentum during the whole collision process:

$$v_3 = \frac{m_1 \cdot v_1 + m_2 \cdot v_2}{m_1 + m_2} = \frac{p}{m_1 + m_2} = \text{const.}$$



Path, velocity and momentum transfer during the collision

CASSY®

WWW.LEYBOLD-SHOP.COM



# P1 MECHANICS



P1.1	MEASURING METHODS	3
P1.2	FORCES	6
P1.3	TRANSLATIONAL MOTIONS OF A MASS POINT	12
P1.4	ROTATIONAL MOTIONS OF A RIGID BODY	26
P1.5	OSCILLATIONS	32
P1.6	WAVE MECHANICS	39
P1.7	ACOUSTICS	45
P1.8	AERO- AND HYDRODYNAMICS	55

# P1 MECHANICS



## P1.1 MEASURING METHODS

P1.1.1 Measuring lengths	3
P1.1.2 Measuring volume and density	4
P1.1.3 Determining the gravitational constant	5

## P1.2 FORCES

P1.2.1 Static effects of forces	6
P1.2.2 Force as vector	7
P1.2.3 Lever	8
P1.2.4 Block and tackle	9
P1.2.5 Inclined plane	10
P1.2.6 Friction	11

## P1.3 TRANSLATIONAL MOTIONS OF A MASS POINT

P1.3.2 One-dimensional motions on Fletcher's trolley	12-14
P1.3.3 One-dimensional motions on the linear air track	15-17
P1.3.4 Conservation of linear momentum	18-20
P1.3.5 Free fall	21-22
P1.3.6 Angled projection	23
P1.3.7 Two-dimensional motions on the air table	24-25

## P1.4 ROTATIONAL MOTIONS OF A RIGID BODY

P1.4.1 Rotational motions	26
P1.4.2 Conservation of angular momentum	27
P1.4.3 Centrifugal force	28
P1.4.4 Motions of a gyroscope	29
P1.4.5 Moment of inertia	30
P1.4.6 Conservation of Energy	31

## P1.5 OSCILLATIONS

P1.5.1 Simple and compound pendulum	32-33
P1.5.2 Harmonic oscillations	34
P1.5.3 Torsion pendulum	35-36
P1.5.4 Coupling of oscillations	37-38

## P1.6 WAVE MECHANICS

P1.6.1 Transversal and longitudinal waves	39
P1.6.2 Wave machine	40
P1.6.3 Circularly polarized waves	41
P1.6.4 Propagation of water waves	42
P1.6.5 Interference of water waves	43
P1.6.6 Resonances	44

## P1.7 ACOUSTICS

P1.7.1 Sound waves	45
P1.7.2 Oscillations of a string	46
P1.7.3 Wavelength and velocity of sound	47-49
P1.7.4 Reflection of ultrasonic waves	50
P1.7.5 Interference of ultrasonic waves	51
P1.7.6 Acoustic Doppler effect	52
P1.7.7 Fourier analysis	53
P1.7.8 Ultrasound in media	54

## P1.8 AERO- AND HYDRODYNAMICS

P1.8.1 Barometric measurements	55
P1.8.2 Bouyancy	56
P1.8.3 Viscosity	57-58
P1.8.4 Surface tension	59
P1.8.5 Introductory experiments on aerodynamics	60
P1.8.6 Measuring air resistance	61
P1.8.7 Measurements in a wind tunnel	62

P1.1.1

MEASURING LENGTHS

P1.1.1.1

Using a caliper gauge with vernier

P1.1.1.2

Using a micrometer screw

P1.1.1.3

Using a spherometer to determine bending radii



Left to right: Using a caliper gauge with vernier (P1.1.1.1), Using a micrometer screw (P1.1.1.2), Using a spherometer to determine bending radii (P1.1.1.3)

Cat. No.	Description	P1.1.1	P1.1.1.2	P1.1.1.3
311 54	Precision vernier callipers	1		
311 83	Precision micrometer		1	
550 35	Copper resistance wire, 0.2 mm diam., 100 m		1	
550 39	Brass resistance wire, 0.5 mm diameter, 50 m		1	
311 86	Spherometer			1
460 291	Plane mirror, 11.5 cm x 10 cm			1
662 092	Cover slips			1
664 154	Watch glass dish, 80 mm diam.			1
664 157	Watch glass dish, 125 mm diam.			1


The caliper gauge, micrometer screw and spherometer are precision measuring instruments; their use is practiced in practical measuring exercises.

In the experiment P1.1.1.1, the caliper gauge is used to determine the outer and inner dimensions of a test body. The vernier scale of the caliper gauge increases the reading accuracy to 1/20 mm.

Different wire gauges are measured in the experiment P1.1.1.2. In this exercise a fundamental difficulty of measuring becomes apparent, namely that the measuring process changes the measurement object. Particularly with soft wire, the measured results are too low because the wire is deformed by the measurement.

The experiment P1.1.1.3 determines the bending radii  $R$  of watch-glasses using a spherometer. These are derived on the basis of the convexity height  $h$  at a given distance  $r$  between the feet of the spherometer, using the formula

$$R = \frac{r^2}{2h} + \frac{h}{2}$$



Vertical section through the measuring configuration with spherometer  
Left: object with convex surface, Right: Object with concaves surface

### P1.1.2

#### MEASURING VOLUME AND DENSITY

##### P1.1.2.1

Determining the volume and density of solids

##### P1.1.2.2

Determining the density of liquids using the plumb bob

##### P1.1.2.3

Determining the density of liquids using the pycnometer after Gay-Lussac

##### P1.1.2.4

Determining the density of air



Determining the density of air (P1.1.2.4)

Cat. No.	Description	P1.1.2.1	P1.1.2.2	P1.1.2.3	P1.1.2.4
362 04	Overflow vessel	1			
590 08ET2	Measuring cylinders, 100 ml, set of 2	1			
590 06	Plastic beaker	1			
309 48ET2	Fishing line, set of 2	1			
311 54	Precision vernier callipers	1			
315 05	Single-pan suspension balance 311	1		1	1
352 52	Steel balls, 30 mm, set of 6	1			
361 63	Cubes (2x) and ball (1x)	1			
590 33	Gauge blocks, set of 2	1			
309 42	Colouring, red, 10 g	1			
362 025	Plumb bob		1		
315 011	Hydrostatic balance		1		
315 31	Set of weights, 10 mg to 200 g		1		
382 21	Stirring thermometer, -30...+110 °C		1	1	
665 754	Measuring cylinder, 100 ml, with plastic base		2	2	
671 9720	Ethanol, denaturated, 1 l		1	1	
666 145	Gay-Lussac pycnometer, 50 ml			1	
379 07	Sphere with 2 stopcocks, glass, 1 l				1
667 072	Support ring for round flask, 250 ml, cork				1
375 58	Hand vacuum pump				1

Depending on the respective aggregate state of a homogeneous substance, various methods are used to determine its density

$$\rho = \frac{m}{V}$$

$m$ : mass,  $V$ : volume

The mass and volume of the substance are usually measured separately. To determine the density of solid bodies, a weighing is combined with a volume measurement. The volumes of the bodies are determined from the volumes of liquid which they displace from an overflow vessel. In the experiment P1.1.2.1, this principle is tested using regular bodies for which the volumes can be easily calculated from their linear dimensions.

To determine the density of liquids, the plumb bob is used in the experiment P1.1.2.2. The measuring task is to determine the densities of water-ethanol mixtures. The Plumb bob determines the density from the buoyancy of a body of known volume in the test liquid.

To determine the density of liquids, the pycnometer after Gay-Lussac is used in the experiment P1.1.2.3. The measuring task is to determine the densities of water-ethanol mixtures. The pycnometer is a pear-shaped bottle in which the liquid to be investigated is filled for weighing. The volume capacity of the pycnometer is determined by weighing with a liquid of known density (e.g. water).

In the experiment P1.1.2.4, the density of air is determined using a sphere of known volume with two stop-cocks. The weight of the enclosed air is determined by finding the difference between the overall weight of the air-filled sphere and the empty weight of the evacuated sphere.



### P1.1.3

#### DETERMINING THE GRAVITATIONAL CONSTANT

##### P1.1.3.1

Determining the gravitational constant with the gravitation torsion balance after Cavendish - Measuring the excursion with a light pointer

Determining the gravitational constant with the gravitation torsion balance after Cavendish - Measuring the excursion with a light pointer (P1.1.3.1)

Cat. No.	Description	P1.1.3.1
332 101	Gravitation torsion balance	1
471 830	He-Ne Laser, linearly polarised	1
313 05	Table stop-clock	1
311 77	Steel tape measure, 2 m	1
300 02	Stand base, V-shaped, small	1
301 03	Rotatable clamp	1
301 01	Leybold multiclamp	1
300 42	Stand rod, 47 cm, 12 mm diam.	1

The heart of the gravitation torsion balance after Cavendish is a light-weight beam horizontally suspended from a thin torsion band and having a lead ball with the mass  $m_2 = 15 \text{ g}$  at each end. These balls are attracted by the two large lead spheres with the mass  $m_1 = 1.5 \text{ kg}$ . Although the attractive force


$$F = G \cdot \frac{m_1 \cdot m_2}{r^2}$$

$r$ : distance between sphere midpoints

is less than  $10^{-9} \text{ N}$ , it can be detected using the extremely sensitive torsion balance. The motion of the small lead balls is observed and measured using a light pointer. Using the curve over time of the motion, the mass  $m_1$  and the geometry of the arrangement, it is possible to determine the gravitational constant  $G$  using either the end-deflection method or the acceleration method.

In the end-deflection method, a measurement error of less than 5 % can be achieved through careful experimenting. The gravitational force is calculated from the resting position of the elastically suspended small lead balls in the gravitational field of the large spheres and the righting moment of the torsion band. The righting moment is determined dynamically using the oscillation period of the torsion pendulum. The acceleration method requires only about 1 min. observation time. The acceleration of the small balls by the gravitational force of the large spheres is measured, and the position of the balls as a function of time is registered.

In the experiment P1.1.3.1, the light pointer is a laser beam which is reflected in the concave reflector of the torsion balance onto a scale. Its position on the scale is measured manually point by point as a function of time.



Oscillations of the gravitation torsion balance around the final equilibrium position  $S_{II}$

### P1.2.1

#### STATIC EFFECTS OF FORCES

P1.2.1.1  
Expansion of a helical spring

P1.2.1.2  
Bending of a leaf spring



Left: Expansion of a helical spring (P1.2.1.1), Right: Bending of a leaf spring (P1.2.1.2)

Cat. No.	Description	P1.2.1.1	P1.2.1.2
352 07ET2	Helical spring, 10 N/m, set of 2	1	
352 08ET2	Helical spring, 25 N/m, set of 2	1	
340 85	Weights, each 50 g, set of 6	1	1
301 21	Stand base MF	2	2
301 27	Stand rod, 50 cm, 10 mm diam.	2	2
301 26	Stand rod, 25 cm, 10 mm diam.	1	1
301 25	Support block	1	
311 77	Steel tape measure, 2 m	1	1
301 29	Pointer, pair	1	1
340 811ET2	Plug-in axles, set of 2	1	
352 051ET2	Leaf springs, set of 2		1
666 615	Universal bosshead		1
686 50ET5	Metal plates, set of 5		1
309 48ET2	Fishing line, set of 2		1


Forces can be recognized by their effects. Thus, static forces can e.g. deform a body. It becomes apparent that the deformation is proportional to the force acting on the body when this force is not too great.

The experiment P1.2.1.1 shows that the extension  $s$  of a helical spring is directly proportional to the force  $F_s$ . Hooke's law applies:

$$F_s = -D \cdot s$$

$D$ : spring constant

The experiment P1.2.1.2 examines the bending of a leaf spring arrested at one end in response to a known force generated by hanging weights from the free end. Here too, the deflection is proportional to the force acting on the leaf spring.



Schematic diagram of bending a leaf spring



P1.2.2  
FORCE AS VECTOR

P1.2.2.1  
Composition and resolution of forces

Composition and resolution of forces (P1.2.2.1\_b)

Cat. No.	Description	P1.2.2.1 (b)
301 301	Adhesive magnetic board	1
314 215	Circular dynamometer, 5 N	2
301 331	Magnetic base with hook	1
352 08ET2	Helical spring, 25 N/m, set of 2	1
311 77	Steel tape measure, 2 m	1
342 61	Weights, 50 g, set of 12	1
301 300	Demonstration panel frame	1


The nature of force as a vectorial quantity can be easily and clearly verified in experiments on the adhesive magnetic board. The point of application of all forces is positioned at the midpoint of the angular scale on the adhesive magnetic board, and all individual forces and the angles between them are measured. The underlying parallelogram of forces can be graphically displayed on the adhesive magnetic board to facilitate understanding.

In experiment P1.2.2.1, a force  $F$  is compensated by the spring force of two dynamometers arranged at angles  $\alpha_1$  and  $\alpha_2$  with respect to  $F$ . The component forces  $F_1$  and  $F_2$  are determined as a function of  $\alpha_1$  and  $\alpha_2$ . This experiment verifies the relationships

$$F = F_1 \cdot \cos \alpha_1 + F_2 \cdot \cos \alpha_2$$

and

$$0 = F_1 \cdot \sin \alpha_1 + F_2 \cdot \sin \alpha_2$$



Parallelogram of forces

### P1.2.3

#### LEVER

##### P1.2.3.1

One-sided and two-sided lever

##### P1.2.3.2

Wheel and axle as a lever with unequal sides



Wheel and axle as a lever with unequal sides (P1.2.3.2)

Cat. No.	Description	P1.2.3.1	P1.2.3.2
342 60	Lever, 1 m	1	
342 61	Weights, 50 g, set of 12	1	1
314 45	Spring balance, 2 N	1	1
314 46	Spring balance, 5 N	1	1
300 02	Stand base, V-shaped, small	1	1
301 01	Leybold multiclamp	1	1
300 42	Stand rod, 47 cm, 12 mm diam.	1	1
342 75	Multiple pulley and moment disc		1

In physics, the law of levers forms the basis for all forms of mechanical transmission of force. This law can be explained using the higher-level concept of equilibrium of angular momentum.

The experiment P1.2.3.1 examines the law of levers:

$$F_1 \cdot x_1 = F_2 \cdot x_2$$

for one-sided and two-sided levers. The object is to determine the force  $F_1$  which maintains a lever in equilibrium as a function of the load  $F_2$ , the load arm  $x_2$  and the power arm  $x_1$ .

The experiment P1.2.3.2 explores the equilibrium of angular momentum using a wheel and axle. This experiment broadens the understanding of the concepts force, power arm and line of action, and explicitly proves that the absolute value of the angular momentum depends only on the force and the distance between the axis of rotation and the line of action.



One-sided and two-sided lever (P1.2.3.1)





P1.2.4  
BLOCK AND TACKLE

P1.2.4.1  
Fixed pulley, loose pulley and block and tackle as simple machines


Fixed pulley, loose pulley and block and tackle as simple machines (P1.2.4.1)

Cat. No.	Description	P1.2.4.1
342 28	Block and tackle D	1
315 36	Set of weights, 0.1 kg to 2 kg	1
300 01	Stand base, V-shaped, large	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
300 44	Stand rod, 100 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
314 181	Precision dynamometer, 20 N	1
341 65	Pulley	2*

\* additionally recommended

The fixed pulley, loose pulley and block and tackle are classic examples of simple machines. Experiments with these machines represent the most accessible introduction to the concept of work in mechanics.

In the experiment P1.2.4.1, the block and tackle is set up on the lab bench using a stand base. The block and tackle can be expanded to three pairs of pulleys and can support loads of up to 20 N. The pulleys are mounted virtually friction-free in ball bearings.



Setup with block and tackle

### P1.2.5

#### INCLINED PLANE

##### P1.2.5.1

Inclined plane: force along the plane and force normal to the plane

##### P1.2.5.2

Determining the coefficient of static friction using the inclined plane



*Inclined plane: force along the plane and force normal to the plane (P1.2.5.1)*

Cat. No.	Description	P1.2.5.1	P1.2.5.2
341 21	Inclined plane	1	1
314 141	Precision dynamometer, 1 N	1	
342 10	Wooden blocks for friction experiments, pair		1
311 77	Steel tape measure, 2 m		1


The motion of a body on an inclined plane can be described most easily when the force exerted by the weight  $G$  on the body is vectorially decomposed into a force  $F_1$  along the plane and a force  $F_2$  normal to the plane. The force along the plane acts parallel to a plane inclined at an angle  $\alpha$ , and the force normal to the plane acts perpendicular to the plane. For the absolute values of the forces, we can say:

$$F_1 = G \cdot \sin \alpha \quad \text{and} \quad F_2 = G \cdot \cos \alpha$$

This decomposition is verified in the experiment P1.2.5.1. Here, the two forces  $F_1$  and  $F_2$  are measured for various angles of inclination  $\alpha$  using precision dynamometers.

The experiment P1.2.5.2 uses the dependency of the force normal to the plane on the angle of inclination for quantitative determination of the coefficient of static friction  $\mu$  of a body. The inclination of a plane is increased until the body no longer adheres to the surface and begins to slide. From the equilibrium of the force along the plane and the coefficient of static friction

$$F_1 = \mu \cdot F_2 \quad \text{we can derive} \quad \mu = \tan \alpha$$



*Calculating the coefficient of static friction (P1.2.5.2)*

P1.2.6  
FRICTION

P1.2.6.1  
Static friction, sliding friction  
and rolling friction



Static friction, sliding friction and rolling friction (P1.2.6.1)

Cat. No.	Description	P1.2.6.1
315 36	Set of weights, 0.1 kg to 2 kg	1
300 40	Stand rod, 10 cm, 12 mm diam.	6
314 47	Spring balance, 10 N	1
342 10	Wooden blocks for friction experiments, pair	1

In discussing friction between solid bodies, we distinguish between static friction, sliding friction and rolling friction. Static friction force is the minimum force required to set a body at rest on a solid base in motion. Analogously, sliding friction force is the force required to maintain a uniform motion of the body. Rolling friction force is the force which maintains the uniform motion of a body which rolls on another body.

To begin, the experiment P1.2.6.1 verifies that the static friction force  $F_H$  and the sliding friction force  $F_G$  are independent of the size of the contact surface and proportional to the resting force  $G$  on the base surface of the friction block. Thus, the following applies:


$$F_H = \mu_H \cdot G \text{ and } F_G = \mu_G \cdot G$$

The coefficients  $\mu_H$  and  $\mu_G$  depend on the material of the friction surfaces. The following relationship always applies:

$$\mu_H > \mu_G$$

To distinguish between sliding and rolling friction, the friction block is placed on top of multiple stand rods laid parallel to each other. The rolling friction force  $F_R$  is measured as the force which maintains the friction block in a uniform motion on the rolling rods. The sliding friction force  $F_G$  is measured once more for comparison, whereby this time the friction block is pulled over the stand rods as a fixed base (direction of pull = direction of rod axes). This experiment confirms the relationship:

$$F_G > F_R$$




Comparison of sliding (point) and rolling friction (triangle)

### P1.3.2

#### ONE-DIMENSIONAL MOTIONS ON FLETCHER'S TROLLEY

##### P1.3.2.1


Path-time diagram of straight motion - Recording the time with an electronic stopclock



Path-time diagram of straight motion - Recording the time with an electronic stopclock (P1.3.2.1\_c)

Cat. No.	Description	P1.3.2.1 (G)
337 130	Track, 1.5 m	1
337 110	Trolley	1
337 114	Additional weights, pair	1*
315 410	Slotted mass hanger, 10 g, small	1
315 418	Slotted weight, 10 g, grey	4
309 48ET2	Fishing line, set of 2	1
337 462	Combination light barrier	2
337 463	Holder for combination spoked wheel	1
337 464	Combination spoked wheel	1
575 471	Counter S	1
501 16	Multi-core cable, 6-pole, 1.5 m	2


\* additionally recommended



Path-time diagram of a linear motion (P1.3.2.1)

Fletcher's trolley is the classical experiment apparatus for investigating linear translational motions. The trolley has a ball bearing, his axles are spring-mounted and completely immersed in order to prevent being overloaded. The wheels are designed in such a way that the trolley centers itself on the track and friction at the wheel flanks is avoided.

Using extremely simple means, the experiment P1.3.2.1 makes the definition of the velocity  $v$  as the quotient of the path difference  $\Delta s$  and the corresponding time difference  $\Delta t$  directly accessible to the students. The path difference  $\Delta s$  is read off directly from a scale on the track. The electronic measurement of the time difference is started and stopped using two light barriers. To enable investigation of uniformly accelerated motions, the trolley is connected to a thread which is laid over a pulley, allowing various weights to be suspended.



**P1.3.2**  
**ONE-DIMENSIONAL MOTIONS ON FLETCHER'S TROLLEY**

P1.3.2.2  
Path-time diagram of straight motion - Recording and evaluating with CASSY

P1.3.2.3  
Definition of the Newton as a unit of force - Recording and evaluating with CASSY

P1.3.2.4  
Path-time diagram of straight motion - Recording and evaluating with VideoCom

Path-time diagram of straight motion - Recording and evaluating with CASSY (P1.3.2.2\_a)

Cat. No.	Description	P1.3.2.2 (a)	P1.3.2.3 (a)	P1.3.2.4
337 130	Track, 1.5 m	1	1	1
337 110	Trolley	1	1	1
337 114	Additional weights, pair	1*		1*
315 410	Slotted mass hanger, 10 g, small	1		1
315 418	Slotted weight, 10 g, grey	4		4
309 48ET2	Fishing line, set of 2	1	1	1
337 462	Combination light barrier	1	1	
337 464	Combination spoked wheel	1	1	1
683 41	Holding magnet for track	1	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	<b>1</b>	
524 220	CASSY Lab 2	1	1	
524 074	Timer S	1	1	
501 16	Multi-core cable, 6-pole, 1.5 m	1	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1	
337 115	Newton weights		1	
337 463	Holder for combination spoked wheel			1
337 47USB	VideoCom			1
300 59	Camera tripod			1
501 38	Connecting lead, 32 A, 200 cm, black			4
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

\* additionally recommended

The experiment P1.3.2.2 looks at motion events which can be transmitted to the combination spoked wheel by means of a thin thread on Fletcher's trolley. The combination spoked wheel serves as an easy-running deflection pulley. The signals of the combination light barrier are recorded by the computer-assisted measuring system CASSY and converted to a path-time diagram. As this diagram is generated in real

time while the experiment is running, the relationship between the motion and the diagram is extremely clear.

In the experiment P1.3.2.3, a calibrated weight exercises an accelerating force of 1 N on a trolley with the mass 1 kg. As one might expect, CASSY shows the value


$$a = 1 \frac{m}{s^2}$$

for the acceleration. At the same time, this experiment verifies that the trolley is accelerated to the velocity

$$v = 1 \frac{m}{s}$$

in the time 1 s.

The single-line CCD video camera VideoCom represents in the experiment P1.3.2.4 a new, easy-to-use method for recording one-dimensional motions. It illuminates one or more bodies in motion with LED flashes and images them on the CCD line with 2048-pixel resolution via a camera lens (CCD: charge-coupled device). A piece of retro-reflecting foil is attached to each of the bodies to reflect the light rays back to the lens. The current positions of the bodies are transmitted to the computer up to 160 times per second via the USB interface. The automatically controlled flash operates at speeds of up to 1/800 s, so that even a rapid motion on the linear air track or any other track can be sharply imaged. The software supplied with VideoCom represents the entire motion of the bodies in the form of a path-time diagram, and also enables further evaluation of the measurement data.



Comparison of sliding (point) and rolling friction (triangle)

### P1.3.2

#### ONE-DIMENSIONAL MOTIONS ON FLETCHER'S TROLLEY

##### P1.3.2.5

Accelerated motions with an acceleration sensor




Accelerated motions with an acceleration sensor (P1.3.2.5)

Cat. No.	Description	P1.3.2.5
337 130	Track, 1.5 m	1
337 110	Trolley	1
337 114	Additional weights, pair	1
524 018	Pocket-CASSY 2 Bluetooth	1
524 019	Rechargeable battery for Pocket-CASSY 2 Bluetooth	1
524 0031	Bluetooth dongle	1
524 220	CASSY Lab 2	1
524 0424	3D acceleration sensor S	1
300 761	Support blocks, set of 6	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

Fletcher's trolley is the classical experiment apparatus for investigating linear translational motions. The trolley has a ball bearing, his axles are spring-mounted and completely immersed in order to prevent being overloaded. The wheels are designed in such a way that the trolley centers itself on the track and friction at the wheel flanks is avoided.

In the experiment P1.3.2.5 the acceleration of a cart is recorded from within the system using a wireless, bluetooth connected and battery powered Pocket-CASSY equipped with a 3D acceleration sensor. The gravitational forces and accelerations are measured in different situations.




Acceleration, velocity and path of the trolley

P1.3.3

ONE-DIMENSIONAL MOTIONS ON THE LINEAR AIR TRACK

P1.3.3.1

Path-time diagram of straight motion - Recording the time with forked light barrier




Path-time diagram of straight motion - Recording the time with forked light barrier (P1.3.3.1\_a)

Cat. No.	Description	P1.3.3.1 (a)	P1.3.3.1 (d)
337 501	Air track	1	1
337 53	Air supply	1	1
667 823	Power controller	1	1
311 02	Metal rule, 1 m	1	1
337 46	Fork-type light barrier	1	2
501 16	Multi-core cable, 6-pole, 1.5 m	1	2
524 013	Sensor-CASSY 2	1	
524 220	CASSY Lab 2	1	
524 074	Timer S	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	
575 471	Counter S		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	

The advantage of studying linear translational motions on the linear air track is that interference factors such as frictional forces and moments of inertia of wheels do not occur. The sliders on the linear air track are fitted with an interrupter flag which obscures a light barrier. By adding additional weights, it is possible to double and triple the masses of the sliders.

Using extremely simple means, the experiment P1.3.3.1 (a) makes the definition of the velocity  $v$  as the quotient of the path difference  $\Delta s$  and the corresponding time difference  $\Delta t$  directly accessible to the students. The path difference  $\Delta s$  is read off directly from a scale on the track. The electronic measurement of the time difference is started by switching off the holding magnet. The instantaneous velocity of the slider can also be calculated from the obscuration time of a forked light barrier and the width of the interrupter flag. To enable investigation of uniformly accelerated motions, the slider is connected to a thread which is laid over a pulley, allowing weights to be suspended.

Using extremely simple means, the experiment P1.3.3.1 (d) makes the definition of the velocity  $v$  as the quotient of the path difference  $\Delta s$  and the corresponding time difference  $\Delta t$  directly accessible to the students. The path difference  $\Delta s$  is read off directly from a scale on the track. The electronic measurement of the time difference is started by using a light barrier. The instantaneous velocity of the slider can also be calculated from the obscuration time of a forked light barrier and the width of the interrupter flag. To enable investigation of uniformly accelerated motions, the slider is connected to a thread which is laid over a pulley, allowing weights to be suspended.



Path-time diagram for uniform motion (P1.3.3.1)

### P1.3.3

#### ONE-DIMENSIONAL MOTIONS ON THE LINEAR AIR TRACK

P1.3.3.4  
Path-time and velocity-time diagrams of straight motion - Recording and evaluating with CASSY

P1.3.3.5  
Uniformly accelerated motion with reversal of direction - Recording and evaluating with CASSY

P1.3.3.6  
Kinetic energy of a uniformly accelerated mass - Recording and evaluating with CASSY



Path-time and velocity-time diagrams of straight motion - Recording and evaluating with CASSY (P1.3.3.4)

Cat. No.	Description	P1.3.3.4-6
337 501	Air track	1
337 53	Air supply	1
667 823	Power controller	1
337 462	Combination light barrier	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>
524 220	CASSY Lab 2	1
524 074	Timer S	1
501 16	Multi-core cable, 6-pole, 1.5 m	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The computer-assisted measurement system CASSY is particularly suitable for simultaneously measuring transit time  $t$ , path  $s$ , velocity  $v$  and acceleration  $a$  of a slider on the linear air track. The linear motion of the slider is transmitted to the motion sensing element by means of a lightly tensioned thread; the signals of the motion sensing element are matched to the CASSY measuring inputs by the Timer S.

The PC supports significantly the evaluation of the measured data. Alternatively the measured data can be exported in table format to external data evaluation.

The object of the experiment P1.3.3.4 is to study uniform and uniformly accelerated motions on the horizontally aligned linear air track.

In the experiment P1.3.3.5 the patch, velocity and acceleration of a slider is record, which moves uphill on an inclined plane, then stops, moves downhill, reflected elastically at the lower end and oscillated several times back and forth.

The experiment P1.3.3.6 records the kinetic energy


$$E = \frac{m}{2} \cdot v^2$$

of a uniformly accelerated slider of the mass  $m$  as a function of the time and compares it with the work

$$W = F \cdot s$$

which the accelerating force  $F$  has performed. This verifies the relationship

$$E(t) = W(t)$$



Path-time, velocity-time and acceleration-time diagram





**P1.3.3**  
**ONE-DIMENSIONAL MOTIONS ON THE LINEAR AIR TRACK**

P1.3.3.7  
Confirming Newton's first and second laws for linear motions - Recording and evaluating with VideoCom

P1.3.3.8  
Uniformly accelerated motion with reversal of direction - Recording and evaluating with VideoCom

P1.3.3.9  
Kinetic energy of a uniformly accelerated mass - Recording and evaluating with VideoCom

Confirming Newton's first and second laws for linear motions - Recording and evaluating with VideoCom (P1.3.3.7)

Cat. No.	Description	P1.3.3.7-9
337 501	Air track	1
337 53	Air supply	1
667 823	Power controller	1
337 47USB	VideoCom	1
300 59	Camera tripod	1
311 02	Metal rule, 1 m	1
501 38	Connecting lead, 32 A, 200 cm, black	4
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The object of the experiment P1.3.3.7 is to study uniform and uniformly accelerated motions of a slider on the linear air track and to represent these in a path-time diagram. The software also displays the velocity  $v$  and the acceleration  $a$  of the body as a function of the transit time  $t$ , and the further evaluation verifies Newton's equation of motion

$$F = m \cdot a$$

$F$ : accelerating force

$m$ : mass of accelerated body

In the experiment P1.3.3.8 the patch, velocity and acceleration of a slider is record, which moves uphill on an inclined plane, then stops, moves downhill, reflected elastically at the lower end and oscillated several times back and forth.

The experiment P1.3.3.9 records the kinetic energy


$$E = \frac{m}{2} \cdot v^2$$

of a uniformly accelerated slider of the mass  $m$  as a function of the time and compares it with the work

$$W = F \cdot s$$

which the accelerating force  $F$  has performed. This verifies the relationship

$$E(t) = W(t)$$



Investigating uniformly accelerated motions with VideoCom

### P1.3.4

#### CONSERVATION OF LINEAR MOMENTUM

##### P1.3.4.1

Energy and linear momentum in elastic collision - Measuring with two forked light barriers

##### P1.3.4.2

Energy and linear momentum in inelastic collision - Measuring with two forked light barriers

##### P1.3.4.3

Rocket principle: conservation of momentum and reaction



Energy and linear momentum in elastic collision - Measuring with two forked light barriers (P1.3.4.1\_c)

Cat. No.	Description	P1.3.4.1 (c)	P1.3.4.2 (c)	P1.3.4.3
337 130	Track, 1.5 m	1	1	
337 110	Trolley	2	2	
337 114	Additional weights, pair	1	1	
337 112	Impact spring for track	1		
337 462	Combination light barrier	2	2	
575 451	Counter P	1	1	
501 16	Multi-core cable, 6-pole, 1.5 m	2	2	
337 501	Air track			1
337 53	Air supply			1
667 823	Power controller			1
337 561	Jet slider with dynamometric device			1
314 081	Precision dynamometer, 0.01 N			1

The use of a track makes possible superior quantitative results when verifying the conservation of linear momentum in an experiment.

In the experiments P1.3.4.1 and P1.3.4.2, the obscuration times  $\Delta t_i$  of two light barriers are measured, e.g. for two bodies on a linear track before and after elastic and inelastic collision. These experiments investigate collisions between a moving body and a body at rest, as well as collisions between two bodies in motion. The velocities

$$v_i = \frac{d}{\Delta t_i}$$

$d$ : width of interrupter flags

the momentum values

$$p_i = m_i \cdot v_i$$

$m_i$ : masses of bodies

and the energies

$$E_i = \frac{1}{2} \cdot m_i \cdot v_i^2$$

of the bodies before and after collision can be calculated and compared.

In the experiment P1.3.4.3, the recoil force on a jet slider is measured for different nozzle cross-sections using a sensitive dynamometer in order to investigate the relationship between repulsion and conservation of linear momentum.



**P1.3.4**  
**CONSERVATION OF LINEAR MOMENTUM**

P1.3.4.4  
Newton's third law and laws of collision - Recording and evaluating with VideoCom

Newton's third law and laws of collision - Recording and evaluating with VideoCom (P1.3.4.4\_a)

Cat. No.	Description	P1.3.4.4 (a)
337 130	Track, 1.5 m	1
337 110	Trolley	2
337 114	Additional weights, pair	1
337 112	Impact spring for track	2
337 47USB	VideoCom	1
300 59	Camera tripod	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The single-line CCD video camera is capable of recording pictures at a rate of up to 160 pictures per second. This time resolution is high enough to reveal the actual process of a collision (compression and extension of springs) between two bodies on the track. In other words, VideoCom registers the positions  $s_1(t)$  and  $s_2(t)$  of the two bodies, their velocities  $v_1(t)$  and  $v_2(t)$  as well as their accelerations  $a_1(t)$  and  $a_2(t)$  even during the actual collision. The energy and momentum balance can be verified not only before and after the collision, but also during the collision itself.

The experiment P1.3.4.4 records the elastic collision of two bodies with the masses  $m_1$  and  $m_2$ . The evaluation shows that the linear momentum

$$p(t) = m_1 \cdot v_1(t) + m_2 \cdot v_2(t)$$

remains constant during the entire process, including the actual collision. On the other hand, the kinetic energy

$$E(t) = \frac{m_1}{2} \cdot v_1^2(t) + \frac{m_2}{2} \cdot v_2^2(t)$$


reaches a minimum during the collision, which can be explained by the elastic strain energy stored in the springs. This experiment also verifies Newton's third law in the form

$$m_1 \cdot a_1(t) = -m_2 \cdot a_2(t)$$

From the path-time diagram, it is possible to recognize the time  $t_0$  at which the two bodies have the same velocity

$$v_1(t_0) = v_2(t_0)$$

and the distance  $s_2 - s_1$  between the bodies is at its lowest. At time  $t_0$ , the acceleration values (in terms of their absolute values) are greatest, as the springs have reached their maximum tension.



Confirmation of Newton's third law

### P1.3.4

#### CONSERVATION OF LINEAR MOMENTUM

##### P1.3.4.5

Newton's third law and laws of collision  
- Recording and evaluating with two ultrasonic motion sensors and CASSY



Newton's third law and laws of collision – Recording and evaluating with two ultrasonic motion sensors and CASSY (P1.3.4.5)


Cat. No.	Description	P1.3.4.5
337 130	Track, 1.5 m	1
337 110	Trolley	2
337 114	Additional weights, pair	1
337 473	Impact spring, soft	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>
524 220	CASSY Lab 2	1
524 0701	Ultrasonic motion sensor S	2
501 11	Extension cable, 15 pin	2
301 25	Support block	2
340 89ET5	Coupling plugs, 4 mm, set of 5	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

In the experiment P1.3.4.5 the way of the two trolleys are measured during an elastic or inelastic collision. This is done with two ultrasonic sensors.  $v(t)$  and  $a(t)$ -diagrams are determined by the Sensor-CASSY. The forces  $F_1(t)$  and  $F_2(t)$  can be calculated from the accelerations  $a_1(t)$  and  $a_2(t)$  and the masses  $m_1$  and  $m_2$  of the trolleys. It is confirmed that  $F_1(t) = -F_2(t)$  applies during the collision. It follows the conservation of momentum during the whole collision process:


$$\Delta p_1 + \Delta p_2 = \int F_1(t) dt + \int F_2(t) dt = \int (F_1(t) + F_2(t)) dt = 0$$

In addition, the center of mass motion  $s_3(t)$  is observed and it can be showed that the speed  $v_3(t)$  of the center of gravity during the whole collision process is constant. This shows alternative the conservation of momentum during the whole collision process:

$$v_3 = \frac{m_1 \cdot v_1 + m_2 \cdot v_2}{m_1 + m_2} = \frac{p}{m_1 + m_2} = \text{const.}$$



Path, velocity and momentum transfer during the collision



## P1.3.5

### FREE FALL

#### P1.3.5.1

Free fall: time measurement with the contact plate and the counter S

Free fall: time measurement with the contact plate and the counter S (P1.3.5.1)

Cat. No.	Description	P1.3.5.1
336 23	Contact plate, large	1
336 21	Holding magnet	1
336 25	Holding magnet adapter with a release mechanism	1
575 471	Counter S	1
301 21	Stand base MF	2
301 26	Stand rod, 25 cm, 10 mm diam.	3
300 46	Stand rod, 150 cm, 12 mm diam.	1
301 01	Leybold multiclamp	2
311 23	Scale with pointers	1
501 25	Connecting lead, 32 A, 50 cm, red	1
501 26	Connecting lead, 32 A, 50 cm, blue	1
501 35	Connecting lead, 32 A, 200 cm, red	1
501 36	Connecting lead, 32 A, 200 cm, blue	1

To investigate free fall, a steel ball is suspended from an electromagnet. It falls downward with a uniform acceleration due to the force of gravity


$$F = m \cdot g$$

$m$ : mass of ball,  $g$ : gravitational acceleration

as soon as the electromagnet is switched off. The friction of air can be regarded as negligible as long as the falling distance, and thus the terminal velocity, are not too great; in other words, the ball falls freely.

In the experiment P1.3.5.1, electronic time measurement is started as soon as the ball is released through interruption of the magnet current. After traveling a falling distance  $h$ , the ball falls on a contact plate, stopping the measurement of time  $t$ . The measurements for various falling heights are plotted as value pairs in a path-time diagram. As the ball is at rest at the beginning of timing,  $g$  can be determined using the relationship

$$h = \frac{1}{2} g \cdot t^2$$



Path-time diagram of the free fall of the ball

### P1.3.5

#### FREE FALL

##### P1.3.5.3

Free fall: multiple measurements with the g ladder


##### P1.3.5.4

Free fall: Recording and evaluating with VideoCom



Free fall: multiple measurements with the g ladder (P1.3.5.3)

Cat. No.	Description	P1.3.5.3	P1.3.5.4
529 034	g ladder	1	
337 46	Fork-type light barrier	1	
501 16	Multi-core cable, 6-pole, 1.5 m	1	
524 013	Sensor-CASSY 2	1	
524 220	CASSY Lab 2	1	
524 074	Timer S	1	
337 47USB	VideoCom		1
300 59	Camera tripod		1
337 472	Falling body		1
336 21	Holding magnet		1
300 01	Stand base, V-shaped, large		1
300 46	Stand rod, 150 cm, 12 mm diam.		1
300 41	Stand rod, 25 cm, 12 mm Ø		1
301 01	Leybold multiclamp		1
501 38	Connecting lead, 32 A, 200 cm, black		4
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1



Path-time and acceleration-time diagrams of the falling body

The disadvantage of preparing a path-time diagram by recording the measured values point by point is that it takes a long time before the dependency of the result on experiment parameters such as the initial velocity or the falling height becomes apparent. Such investigations become much simpler when the entire measurement series of a path-time diagram is recorded in one measuring run using the computer.

In the experiment P1.3.5.3, a ladder with several rungs falls through a forked light barrier, which is connected to the CASSY computer interface device to measure the obscuration times. This measurement is equivalent to a measurement in which a body falls through multiple equidistant light barriers. The height of the falling body corresponds to the rung width. The measurement data are recorded and evaluated using CASSY Lab. The instantaneous velocities are calculated from the obscuration times and the rung width and displayed in a velocity-time diagram  $v(t)$ . The measurement points can be described by a straight line


$$v(t) = v_0 + g \cdot t$$

$g$ : gravitational acceleration

whereby  $v_0$  is the initial velocity of the ladder when the first rung passes the light barrier.

In the experiment P1.3.5.4, the motion of a falling body is tracked as a function of time using the single-line CCD camera VideoCom and evaluated using the corresponding software. The measurement series is displayed directly as the path-time diagram  $h(t)$ . This curve can be described by the general relationship

$$s = v_0 \cdot t + \frac{1}{2} g \cdot t^2$$



### P1.3.6

#### ANGLED PROJECTION

##### P1.3.6.1

Point-by-point recording of the projection parabola as a function of the speed and angle of projection

##### P1.3.6.2

Principle of superposing: comparison of inclined projection and free fall


Principle of superposing: comparison of inclined projection and free fall (P1.3.6.2)

Cat. No.	Description	P1.3.6.1	P1.3.6.2
336 56	Projection apparatus, large	1	1
301 06	Bench clamp	2	2
311 77	Steel tape measure, 2 m	1	
300 76	Laboratory stand II	1	
311 22	Vertical rule	1	
300 11	Saddle base	1	
649 42	Tray, 552 mm x 197 mm x 48 mm	1	1
688 108	Quartz sand, 1 kg	1	1
336 21	Holding magnet		1
521 231	Low-voltage power supply, 3/6/9/12 V		1
311 02	Metal rule, 1 m		1
300 44	Stand rod, 100 cm, 12 mm diam.		1
301 07	Simple bench clamp		1
501 26	Connecting lead, 32 A, 50 cm, blue		1
501 35	Connecting lead, 32 A, 200 cm, red		1
501 36	Connecting lead, 32 A, 200 cm, blue		1

The trajectory of a ball launched at a projection angle  $\alpha$  with a projection velocity  $v_0$  can be reconstructed on the basis of the principle of superposing. The overall motion is composed of a motion with constant velocity in the direction of projection and a vertical falling motion. The superposition of these motions results in a parabola, whose height and width depend on the angle and velocity of projection.

The experiment P1.3.6.1 measures the trajectory of the steel ball point by point using a vertical scale. Starting from the point of projection, the vertical scale is moved at predefined intervals; the two pointers of the scale are set so that the projected steel ball passes between them. The trajectory is a close approximation of a parabola. The observed deviations from the parabolic form may be explained through friction with the air.

In the experiment P1.3.6.2, a second ball is suspended from a holding magnet in such a way that the first ball would strike it if propelled in the direction of projection with a constant velocity. Then, the second ball is released at the same time as the first ball is projected. We can observe that, regardless of the launch velocity  $v_0$  of the first ball, the two balls collide; this provides experimental confirmation of the principle of superposing.



Schematic diagram comparing angled projection and free fall (P1.3.6.2)

### P1.3.7

#### TWO-DIMENSIONAL MOTIONS ON THE AIR TABLE

P1.3.7.1  
Uniform linear motion and uniform circular motion

P1.3.7.2  
Uniformly accelerated motion

P1.3.7.3  
Two-dimensional motion on an inclined plane

P1.3.7.4  
Two-dimensional motion in response to a central force

P1.3.7.5  
Superposing translational and rotational motion on a rigid body



Uniform linear motion and uniform circular motion (P1.3.7.1)

Cat. No.	Description	P1.3.7.1-3	P1.3.7.4	P1.3.7.5
337 801	Air table	1	1	1
352 10	Helical spring, 3 N/m		1	

The air table makes possible recording of any two-dimensional motions of a slider for evaluation following the experiment. To achieve this, the slider is equipped with a recording device which registers the position of the slider on metallized recording paper every 20 ms.

The aim of the experiment P1.3.7.1 is to examine the instantaneous velocity of straight and circular motions. In both cases, their absolute values can be expressed as

$$v = \frac{\Delta s}{\Delta t}$$

where  $\Delta s$  is the straight path traveled during time  $\Delta t$  for linear motions and the equivalent arc for circular motions.

In the experiment P1.3.7.2, the slider without an initial velocity moves on the air table inclined by the angle  $\alpha$ . Its motion can be described as a one-dimensional, uniformly accelerated motion. The marked positions permit plotting of a path-time diagram from which we can derive the relationship

$$s = \frac{1}{2} \cdot a \cdot t^2 \quad \text{where} \quad a = g \cdot \sin \alpha$$


In the experiment P1.3.7.3, a motion „diagonally upward“ is imparted on the slider on the inclined air table, so that the slider describes a parabola. Its motion is uniformly accelerated in the direction of inclination and virtually uniform perpendicular to this direction.

The aim of the experiment P1.3.7.4 is to verify Kepler's law of areas. Here, the slider moves under the influence of a central force exerted by a centrally mounted helical screw. In the evaluation, the area

$$\Delta A = |r \times \Delta s|$$


“swept“ due to the motion of the slider in the time  $\Delta t$  is determined from the radius vector  $r$  and the path section  $\Delta s$  as well as from the angle between the two vectors.

The experiment P1.3.7.5 investigates simultaneous rotational and translational motions of one slider and of two sliders joined together in a fixed manner. One recorder is placed at the center of gravity, while a second is at the perimeter of the “rigid body“ under investigation. The motion is described as the motion of the center of gravity plus rotation around that center of gravity.



Relationship between circular motion and (tangential) uniform rectilinear motion





**P1.3.7**  
TWO-DIMENSIONAL MOTIONS  
ON THE AIR TABLE

- P1.3.7.6  
Two-dimensional motion of two elastically coupled bodies
- P1.3.7.7  
Experimentally verifying the equality of a force and its opposing force
- P1.3.7.8  
Elastic collision in two dimensions
- P1.3.7.9  
Inelastic collision in two dimensions

*Elastic collision in two dimensions (P1.3.7.8)*

Cat. No.	Description	P1.3.7.6-9
337 801	Air table	1

The air table is supplied complete with two sliders. This means that this apparatus can also be used to investigate e.g. two-dimensional collisions. In the experiment P1.3.7.6, the motions of two sliders which are elastically coupled by a rubber band are recorded. The evaluation shows that the common center of gravity moves in a straight line and a uniform manner, while the relative motions of the two sliders show a harmonic oscillation. In the experiment P1.3.7.7, elastically deformable metal rings are attached to the edges of the sliders before the start of the experiment. When the two rebound, the same force acts on each slider, but in the opposite direction. Therefore, regardless of the masses  $m_1$  and  $m_2$  of the two sliders, the following relationship applies for the total two-dimensional momentum

$$m_1 \cdot v_1 + m_2 \cdot v_2 = 0$$


The experiments P1.3.7.8 and P1.3.7.9 investigate elastic and inelastic collisions between two sliders. The evaluation consists of calculating the total two-dimensional momentum

$$p = m_1 \cdot v_1 + m_2 \cdot v_2$$

and the total energy

$$E = \frac{m_1}{2} \cdot v_1^2 + \frac{m_2}{2} \cdot v_2^2$$

both before and after collision.



*Elastic impact. Addition of the impulse vectors.*

### P1.4.1

#### ROTATIONAL MOTIONS

##### P1.4.1.1

Path-time diagrams of rotational motions - Time measurement with the counter

##### P1.4.1.2

Path-time diagrams of rotational motions - Recording and evaluating with CASSY



Path-time diagrams of rotational motions - Time measurement with the counter (P1.4.1.1\_a)

Cat. No.	Description	P1.4.1.1 (a)	P1.4.1.2
347 23	Rotation model	1	1
337 46	Fork-type light barrier	1	
575 471	Counter S	1	
501 16	Multi-core cable, 6-pole, 1.5 m	1	1
300 76	Laboratory stand II	1	1
301 07	Simple bench clamp	1	1
337 462	Combination light barrier		1
524 013	Sensor-CASSY 2		1
524 220	CASSY Lab 2		1
524 074	Timer S		1
336 21	Holding magnet		1
300 41	Stand rod, 25 cm, 12 mm Ø		1
300 11	Saddle base		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

The low-friction Plexiglas disk of the rotation model is set in uniform or uniformly accelerated motion for quantitative investigations of rotational motions. Forked light barriers are used to determine the angular velocity; their light beams are interrupted by a 10° flag mounted on the rotating disk. When two forked light barriers are used, measurement of time  $t$  can be started and stopped for any angle  $\varphi$  (optional possible). This experiment determines the mean velocity

$$\omega = \frac{\varphi}{t}$$

If only one forked light barrier is available, the obscuration time  $\Delta t$  is measured, which enables calculation of the instantaneous angular velocity

$$\omega = \frac{10^\circ}{\Delta t}$$

The use of the computer-assisted measured-value recording system CASSY facilitates the study of uniform and uniformly accelerated rotational motions. A thread stretched over the surface of the rotation model transmits the rotational motion to the motion sensing element whose signals are adapted to the measuring inputs of CASSY by a box. In the experiment P1.4.1.1, the angular velocity  $\omega$  and the angular acceleration  $\alpha$  are recorded analogously to acceleration in translational motions. Both uniform and uniformly accelerated rotational motions are investigated. The results are graphed in a velocity-time diagram  $\omega(t)$ . In the case of a uniformly accelerated motion of a rotating disk initially at rest, the angular acceleration can be determined from the linear function

$$\omega = \alpha \cdot t$$

The topic of the experiment P1.4.1.2 are homogeneous and constantly accelerated rotational motions, which are studied on the analogy of homogeneous and constantly accelerated translational motions.



## P1.4.2

### CONSERVATION OF ANGULAR MOMENTUM

#### P1.4.2.1

Conservation of angular momentum in elastic rotational collision

#### P1.4.2.2

Conservation of angular momentum in inelastic rotational collision

Conservation of angular momentum in elastic rotational collision (P1.4.2.1\_b)

Cat. No.	Description	P1.4.2.1-2 (b)
347 23	Rotation model	1
524 431	Light barrier M	2
524 005	Mobile-CASSY 2	1

Torsion impacts between rotating bodies can be described analogously to one-dimensional translational collisions when the axes of rotation of the bodies are parallel to each other and remain unchanged during the collision. This condition is reliably met when carrying out measurements using the rotation model. The angular momentum is specified in the form

$$L = I \cdot \omega$$

$I$ : moment of inertia,  $\omega$ : angular velocity

The principle of conservation of angular momentum states that for any torsion impact of two rotating bodies, the quantity

$$L = I_1 \cdot \omega_1 + I_2 \cdot \omega_2$$

before and after impact remains the same.

The experiments P1.4.2.1 and P1.4.2.2 investigate the nature of elastic and inelastic torsion impact. Using two light barriers and the measuring system CASSY, the obscuration times of two interrupter flags are registered as a measure of the angular velocities before and after torsion impact.

### P1.4.3

#### CENTRIFUGAL FORCE

##### P1.4.3.3

Centrifugal force of an orbiting body  
- Measuring with the central force apparatus and CASSY



Centrifugal force of an orbiting body - Measuring with the central force apparatus and CASSY (P1.4.3.3)

Cat. No.	Description	P1.4.3.3
524 068	Centrifugal force apparatus S	1
521 49	AC/DC power supply, 0...12 V	1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
524 074	Timer S	1
337 46	Fork-type light barrier	1
501 16	Multi-core cable, 6-pole, 1.5 m	1
301 06	Bench clamp	1
300 02	Stand base, V-shaped, small	1
300 40	Stand rod, 10 cm, 12 mm diam.	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The centrifugal force apparatus S enables experimental investigation of the centrifugal force  $F$  as a function of the rotating mass  $m$ , the distance  $r$  of the mass from the centre of rotation and the angular velocity  $\omega$ , thus making it possible to confirm the relation

$$F = m \cdot \omega^2 \cdot r$$

$r$ : radius of orbit,  $\omega$ : angular velocity

for the centrifugal force.

In the centrifugal force apparatus S, the centrifugal force  $F$  acting on a rotating mass  $m$  is transmitted via a lever with ball-and-socket joint and a push pin in the axis of rotation to a leaf spring, whose deflection is measured electrically by means of a bridge-connected strain gauge. In the measuring range relevant for the experiment, the deformation of the leaf spring is elastic and thus proportional to the force  $F$ .


In the experiment P1.4.3.3, the relationship

$$F \propto \omega^2$$

is derived directly from the parabolic shape of the recorded curve  $F(\omega)$ . To verify the proportionalities

$$F \propto r, \quad F \propto m$$

curves are recorded and evaluated for different orbit radii  $r$  and various masses  $m$ .



Centrifugal force  $F$  as a function of the angular velocity  $\omega$

P1.4.4  
MOTIONS OF A GYROSCOPE

P1.4.4.3  
Precession of a gyroscope

P1.4.4.4  
Nutation of a gyroscope



Precession of a gyroscope (P1.4.4.3)

Cat. No.	Description	P1.4.4.3	P1.4.4.4
348 20	Gyroscope	1	1
342 61	Weights, 50 g, set of 12	1	
524 082	Rotary motion sensor S	1	1
337 468	Reflection light barrier	1	1
590 021	Double spring clip	1	1
524 074	Timer S	1	1
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

Gyroscopes generally execute extremely complex motions, as the axis of rotation is supported at only one point and changes directions constantly. We distinguish between the precession and the nutation of a gyroscope.


The aim of the experiment P1.4.4.3 is to investigate the precession of a gyroscope. The precession frequency  $f_p$  is measured by means of the rotary motion sensor S, the rotary frequency  $f$  of the gyroscope disk by means of the reflection light barrier, each in combination with CASSY. The dependence of the precession frequency  $f_p$  on the applied force, i.e. the torque  $M$  and the rotary frequency  $f$  is determined quantitatively. The relationship

$$\omega_p = \frac{M}{I} \cdot \frac{1}{\omega}$$

applies for the corresponding angular frequencies  $\omega_p$  and  $\omega$  and for a known moment of inertia  $I$  of the gyroscope around its axis of symmetry. In the experiment P1.4.4.4, the nutation of a force-free gyroscope is investigated. The nutation frequency  $f_N$  is measured by means of the rotary motion sensor S, the rotary frequency  $f$  of the gyroscope disk by means of the reflection light barrier, each in combination with CASSY. The dependence of the nutation frequency  $f_N$  on the rotary frequency  $f$  is determined quantitatively. The relationship

$$\omega_N = \frac{I \cdot \omega}{I_{\perp}}$$

applies for the corresponding angular frequencies  $\omega_N$  and  $\omega$  and for known moments of inertia  $I$  of the gyroscope around its axis of symmetry (rotational axis of the gyroscope disk) and  $I_{\perp}$  around the pivotal point (point of support) of the axis.



Precession of the gyroscope

### P1.4.5

#### MOMENT OF INERTIA

##### P1.4.5.1

Definition of moment of inertia

##### P1.4.5.2

Moment of inertia and body shape

##### P1.4.5.3

Confirming Steiner's theorem



Definition of moment of inertia (P1.4.5.1), Moment of inertia and body shape (P1.4.5.2), Confirming Steiner's theorem (P1.4.5.3)

Cat. No.	Description	P1.4.5.1	P1.4.5.2	P1.4.5.3
347 801	Torsion axle	1	1	1
300 02	Stand base, V-shaped, small	1	1	1
313 07	Hand-held stop watch I, mechanical	1	1	1
347 81	Cylinders for torsion axle, set		1	
347 82	Sphere for the torsion axle		1	
347 83	Circular disc for the torsion axle			1

For any rigid body whose mass elements  $m_i$  are at a distance of  $r_i$  from the axis of rotation, the moment of inertia is

$$I = \sum_i m_i \cdot r_i^2$$

For a particle of mass  $m$  in an orbit with the radius  $r$ , we can say

$$I = m \cdot r^2$$

The moment of inertia is determined from the oscillation period of the torsion axle on which the test body is mounted and which is elastically joined to the stand via a helical spring. The system is excited to harmonic oscillations. For a known directed angular quantity  $D$ , the oscillation period  $T$  can be used to calculate the moment of inertia of the test body using the equation

$$I = D \cdot \left( \frac{T}{2\pi} \right)^2$$

In the experiment P1.4.5.1, the moment of inertia of a "mass point" is determined as a function of the distance  $r$  from the axis of rotation. In this experiment, a rod with two weights of equal mass is mounted transversely on the torsion axle. The centers of gravity of the two weights have the same distance  $r$  from the axis of rotation, so that the system oscillates with no unbalanced weight.

The experiment P1.4.5.2 compares the moments of inertia of a hollow cylinder, a solid cylinder and a solid sphere. This measurement uses two solid cylinders with equal mass but different radii. Additionally, this experiment examines a hollow cylinder which is equal to one of the solid cylinders in mass and radius, as well as a solid sphere with the same moment of inertia as one of the solid cylinders.

The experiment P1.4.5.3 verifies Steiner's law using a flat circular disk. Here, the moments of inertia  $I_A$  of the circular disk are measured for various distances  $a$  from the axis of rotation, and compared with the moment of inertia  $I_S$  around the axis of the center of gravity. This experiment confirms the relationship

$$I_A - I_S = M \cdot a^2$$



Steiner's theorem (P1.4.5.3)



P1.4.6  
CONSERVATION OF ENERGY

P1.4.6.1  
Maxwell's wheel

P1.4.6.2  
Maxwell's wheel - Recording  
and evaluating with ultrasonic  
motion sensor

Maxwell's wheel (P1.4.6.1)

Cat. No.	Description	P1.4.6.1	P1.4.6.2
331 22	Maxwell's wheel	1	1
337 46	Fork-type light barrier	1	
501 16	Multi-core cable, 6-pole, 1.5 m	1	
575 471	Counter S	1	
336 25	Holding magnet adapter with a release mechanism	1	
311 23	Scale with pointers	1	1
300 11	Saddle base	1	1
301 25	Support block	1	
301 21	Stand base MF	2	2
301 27	Stand rod, 50 cm, 10 mm diam.	2	3
300 44	Stand rod, 100 cm, 12 mm diam.	2	2
301 01	Leybold multiclamp	4	7
524 0701	Ultrasonic motion sensor S		1
501 11	Extension cable, 15 pin		1
524 005	Mobile-CASSY 2		1

The law of conservation of energy states that the total amount of energy in an isolated system remains constant over time. Within this system the energy can change form, for instance potential in kinetic energy. In the daily experience (also during experiments) energy apparently is lost. The reason for this is a change to an energy form which is not considered like the friction.

Experiment P1.4.6.1 is used to examine the conservation of energy at the Maxwell's wheel. During the experiment potential energy  $E_{pot}$  is transformed to kinetic energy  $E_{kin}$  due a translational motion ( $E_{trans}$ ) and a rotational motion ( $E_{rot}$ ). For different heights times and velocities are measured. From the data one can determine the inertia of the Maxwell's wheel. With a known inertia, one can calculate the gravitational acceleration.

In the experiment P1.4.6.2 the conservation of energy at the Maxwell's wheel is studied. The position and velocity, required for calculation of kinetic energy, is measured with an ultrasonic motion sensor.

## P1.5.1

### SIMPLE AND COMPOUND PENDULUM

P1.5.1.1  
Determining the gravitational  
acceleration with a simple pendulum

P1.5.1.2  
Determining the acceleration of  
gravity with a reversible pendulum



Determining the gravitational acceleration with a simple pendulum (P1.5.1.1\_a)

Cat. No.	Description	P1.5.1.1 (a)	P1.5.1.2 (a)
346 39	Ball with pendulum suspension	1	
313 07	Hand-held stop watch I, mechanical	1	1
311 77	Steel tape measure, 2 m	1	1
346 111	Reversible pendulum		1

A simple, or "mathematic" pendulum is understood to be a point-shaped mass  $m$  suspended on a massless thread with the length  $s$ . For small deflections, it oscillates under the influence of gravity with the period

$$T = 2\pi \cdot \sqrt{\frac{s}{g}}$$


Thus, a mathematic pendulum could theoretically be used to determine the gravitational acceleration  $g$  precisely through measurement of the oscillation period and the pendulum length.

In the experiment P1.5.1.1, the ball with pendulum suspension is used to determine the gravitational acceleration. As the mass of the ball is much greater than that of the steel wire on which it is suspended, this pendulum can be considered to be a close approximation of a mathematic pendulum. Multiple oscillations are recorded to improve measuring accuracy. For gravitational acceleration, the error then depends essentially on the accuracy with which the length of the pendulum is determined.

The reversible pendulum used in the experiment P1.5.1.2 has two edges for suspending the pendulum and two sliding weights for "tuning" the oscillation period. When the pendulum is properly adjusted, it oscillates on both edges with the same period

$$T_0 = 2\pi \cdot \sqrt{\frac{s_{\text{red}}}{g}}$$

and the reduced pendulum length  $s_{\text{red}}$  corresponds to the very precisely known distance  $d$  between the two edges. For gravitational acceleration, the error then depends essentially on the accuracy with which the oscillation period  $T_0$  is determined.



Measurement diagram for reversible pendulum (P1.5.1.2)





**P1.5.1**  
**SIMPLE AND COMPOUND PENDULUM**

P1.5.1.3  
Oscillations of a rod pendulum

P1.5.1.4  
Dependency of period of the oscillation of a rod pendulum on the amplitude

P1.5.1.5  
Determination of the acceleration due to gravity on earth by means of a rod pendulum

P1.5.1.6  
Pendulum with changeable acceleration due to gravity (variable g-pendulum)

Oscillations of a rod pendulum (P1.5.1.3)


Cat. No.	Description	P1.5.1.3-5	P1.5.1.6
346 20	Physical pendulum	1	1
524 082	Rotary motion sensor S	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	<b>1</b>
524 220	CASSY Lab 2	1	1
301 21	Stand base MF	2	2
301 26	Stand rod, 25 cm, 10 mm diam.	1	2
301 27	Stand rod, 50 cm, 10 mm diam.		1
301 01	Leybold multiclamp		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

In the experiment P1.5.1.3, the oscillation of a rod pendulum, i.e. an simple physical pendulum, is investigated. Using the rotary motion sensor S the oscillation of the pendulum is recorded as a function of time. Angle  $\alpha(t)$ , velocity  $\omega(t)$  and acceleration  $a(t)$  are compared. In addition, the effective length of the pendulum is determined from the measured oscillation period  $T$ .

In the experiment P1.5.1.4, the dependence of the period  $T$  on the amplitude  $A$  of an oscillation is investigated. For small deflections the oscillation of a pendulum is approximately harmonic and the period is independent from the amplitude. For high deflections this approximation is no longer satisfied: the higher the amplitude is the larger the period.

In experiment P1.5.1.5, the rod pendulum is applied as reversible pendulum. The value of the acceleration due to gravity is determined. The pendulum is set up at two pivot points at opposite sides of the rod. The position of two sliding weights influences the period. When the pendulum is properly adjusted, it oscillates on both edges with the same period  $T$ . The effective pendulum length  $l_e$  corresponds to the distance  $d$  between the two pivot points. The acceleration due to gravity is calculated from the effective pendulum length  $l_e$  and the period  $T$ .

In the experiment P1.5.1.6, a pendulum with variable acceleration due to gravity (variable g pendulum) is assembled and investigated. The oscillation plane is tilted. Therefore, the acceleration due to gravity is reduced. This leads to different oscillation periods depending on the tilt. In the experiment the dependence of the period on the tilt angle is determined. Additionally, the acceleration due to gravity on different celestial bodies is simulated.



Oscillations of the pendulum

## P1.5.2

### HARMONIC OSCILLATIONS

#### P1.5.2.1

Oscillations of a spring pendulum –  
Recording the path, velocity and  
acceleration with CASSY

#### P1.5.2.2

Determining the oscillation period  
of a spring pendulum as a function  
of the oscillating mass



Oscillations of a spring pendulum – Recording the path, velocity and acceleration with CASSY (P1.5.2.1)

Cat. No.	Description	P1.5.2.1-2
352 10	Helical spring, 3 N/m	1
342 61	Weights, 50 g, set of 12	1
336 21	Holding magnet	1
337 462	Combination light barrier	1
337 464	Combination spoked wheel	1
524 074	Timer S	1
501 16	Multi-core cable, 6-pole, 1.5 m	1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
300 01	Stand base, V-shaped, large	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
300 46	Stand rod, 150 cm, 12 mm diam.	1
301 01	Leybold multiclamp	2
301 08	Clamp with hook	1
309 48ET2	Fishing line, set of 2	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

When a system is deflected from a stable equilibrium position, oscillations can occur. An oscillation is considered harmonic when the restoring force  $F$  is proportional to the deflection  $x$  from the equilibrium position.

$$F = D \cdot x$$

$D$ : directional constant


The oscillations of a spring pendulum are often used as a classic example of this.

In the experiment P1.5.2.1, the harmonic oscillations of a spring pendulum are recorded as a function of time using the motion transducer and the computer-assisted measured value recording system CASSY. In the evaluation, the oscillation quantities path  $x$ , velocity  $v$  and acceleration  $a$  are compared on the screen. These can be displayed either as functions of the time  $t$  or as a phase diagram.


The experiment P1.5.2.2 records and evaluates the oscillations of a spring pendulum for various suspended masses  $m$ . The relationship

$$T = 2\pi \cdot \sqrt{\frac{D}{m}}$$

for the oscillation period is verified.



Path, velocity and acceleration of the pendulum



**P1.5.3**  
**TORSION PENDULUM**

P1.5.3.1  
Free rotational oscillations -  
Measuring with a hand-held stopclock

P1.5.3.2  
Forced rotational oscillations -  
Measuring with a hand-held stopclock

Free rotational oscillations - Measuring with a hand-held stopclock (P1.5.3.1)

Cat. No.	Description	P1.5.3.1	P1.5.3.2
346 00	Torsion pendulum	1	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1	1
313 07	Hand-held stop watch I, mechanical	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	2
562 793	Plug-in power supply for torsion pendulum		1
531 120	Multimeter LDanalog 20		1

The torsion pendulum after Pohl can be used to investigate free or forced harmonic rotational oscillations. An electromagnetic eddy current brake damps these oscillations to a greater or lesser extent, depending on the set current. The torsion pendulum is excited to forced oscillations by means of a motor-driven eccentric rod.

The aim of the experiment P1.5.3.1 is to investigate free harmonic rotational oscillations of the type

$$\varphi(t) = \varphi_0 \cdot \cos \omega t \cdot e^{-\delta \cdot t} \text{ where } \omega = \sqrt{\omega_0^2 - \delta^2}$$

$\omega_0$ : characteristic frequency of torsion pendulum

To distinguish between oscillation and creepage, the damping constant  $\delta$  is varied to find the current  $I_0$  which corresponds to the aperiodic limiting case. In the oscillation case, the angular frequency  $\omega$  is determined for various damping levels from the oscillation period  $T$  and the damping constant  $\delta$  by means of the ratio

$$\left| \frac{\varphi_{n+1}}{\varphi_n} \right| = e^{-\delta \cdot \frac{T}{2}}$$

of two sequential oscillation amplitudes. Using the relationship


$$\omega^2 = \omega_0^2 - \delta^2$$

we can determine the characteristic frequency  $\omega_0$ .

In the experiment P1.5.3.2, the torsion pendulum is excited to oscillations with the frequency  $\omega$  by means of a harmonically variable angular momentum. To illustrate the resonance behavior, the oscillation amplitudes determined for various damping levels are plotted as a function of  $\omega^2$  and compared with the theoretical curve

$$\varphi_0 = \frac{M_0}{I} \cdot \frac{1}{\sqrt{(\omega^2 - \omega_0^2)^2 + \delta^2 \cdot \omega^2}}$$

$I$ : moment of inertia of torsion pendulum



Resonance curves for two different damping constants (P1.5.3.2)

## P1.5.3


### TORSION PENDULUM

#### P1.5.3.3

Free rotational oscillations -  
Recording with CASSY

#### P1.5.3.4

Forced harmonic and chaotic  
rotational oscillations -  
Recording with CASSY



Free rotational oscillations - Recording with CASSY (P1.5.3.3)

Cat. No.	Description	P1.5.3.3	P1.5.3.4
346 00	Torsion pendulum	1	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	<b>1</b>
524 220	CASSY Lab 2	1	1
524 082	Rotary motion sensor S	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	2
562 793	Plug-in power supply for torsion pendulum		1
531 120	Multimeter LDanalog 20		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

The computer-assisted CASSY measured-value recording system is ideal for recording and evaluating the oscillations of the torsion pendulum. The numerous evaluation options enable a comprehensive comparison between theory and experiment. Thus, for example, the recorded data can be displayed as path-time, velocity-time and acceleration-time diagrams or as a phase diagram (path-velocity diagram).

The aim of the experiment P1.5.3.3 is to investigate free harmonic rotational oscillations of the general type


$$\varphi(t) = (\varphi(0) \cdot \cos \omega t + \dot{\varphi}(0) \cdot \sin \omega t) \cdot e^{-\delta t}$$

$$\text{where } \omega = \sqrt{\omega_0^2 - \delta^2}$$

where  $\omega_0$ : characteristic frequency of torsion pendulum

This experiment investigates the relationship between the initial deflection  $\varphi(0)$  and the initial velocity  $\dot{\varphi}(0)$ . In addition, the damping constant  $\delta$  is varied in order to find the current  $I_0$  which corresponds to the aperiodic limiting case.

To investigate the transition between forced harmonic and chaotic oscillations, the linear restoring moment acting on the torsion pendulum is deliberately altered in the experiment P1.5.3.4 by attaching an additional weight to the pendulum. The restoring moment now corresponds to a potential with two minima, i.e. two equilibrium positions. When the pendulum is excited at a constant frequency, it can oscillate around the left minimum, the right minimum or back and forth between the two minima. At certain frequencies, it is not possible to predict when the pendulum will change from one minimum to another. The torsion pendulum is then oscillating in a chaotic manner.



Potential energy of double pendulum with and without additional mass



**P1.5.4**  
**COUPLING OF OSCILLATIONS**

P1.5.4.1  
Coupled pendulum - Measuring with a hand-held stopclock

P1.5.4.2  
Coupled pendulum - Recording and evaluating with VideoCom

Coupled pendulum - Recording and evaluating with VideoCom (P1.5.4.2)

Cat. No.	Description	P1.5.4.1	P1.5.4.2
346 45	Double pendulum	1	1
300 02	Stand base, V-shaped, small	2	2
300 44	Stand rod, 100 cm, 12 mm diam.	2	2
300 42	Stand rod, 47 cm, 12 mm diam.	1	1
301 01	Leybold multiclamp	4	4
460 97	Metal rule, 0.5 m	1	1
309 48ET2	Fishing line, set of 2	1	1
313 07	Hand-held stop watch I, mechanical	1	
337 47USB	VideoCom		1
300 59	Camera tripod		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

Two coupled pendulums oscillate in phase with the angular frequency  $\omega_+$  when they are deflected from the equilibrium position by the same amount. When the second pendulum is deflected in the opposite direction, the two pendulums oscillate in phase opposition with the angular frequency  $\omega_-$ . Deflecting only one pendulum generates a coupled oscillation with the angular frequency

$$\omega = \frac{\omega_+ + \omega_-}{2}$$

in which the oscillation energy is transferred back and forth between the two pendulums. The first pendulum comes to rest after a certain time, while the second pendulum simultaneously reaches its greatest amplitude. Then the same process runs in reverse. The time from one pendulum stand still to the next is called the beat period  $T_S$ . For the corresponding beat frequency, we can say

$$\omega_s = \omega_+ - \omega_-$$

The aim of the experiment P1.5.4.1 is to observe in-phase, phase-opposed and coupled oscillations. The angular frequencies  $\omega_+$ ,  $\omega_-$ ,  $\omega_s$  and  $\omega$  are calculated from the oscillation periods  $T_+$ ,  $T_-$ ,  $T_S$  and  $T$  measured using a stopclock and compared with each other.

In the experiment P1.5.4.2, the coupled motion of the two pendulums is investigated using the single-line CCD camera VideoCom. The results include the path-time diagrams  $s_1(t)$  and  $s_2(t)$  of pendulums 1 and 2, from which the path-time diagrams  $s_+(t) = s_1(t) + s_2(t)$  of the purely in-phase motion and  $s_-(t) = s_1(t) - s_2(t)$  of the purely opposed-phase motion are calculated. The corresponding characteristic frequencies are determined using - Fourier transforms. Comparison identifies the two characteristic frequencies of the coupled oscillations  $s_1(t)$  and  $s_2(t)$  as the characteristic frequencies  $\omega_+$  of the function  $s_+(t)$  and  $\omega_-$  of the function  $s_-(t)$ .



Phase shift of coupled oscillation - recorded with VideoCom (P1.5.4.2)

## P1.5.4


### COUPLING OF OSCILLATIONS

#### P1.5.4.3

Coupling of longitudinal and rotational oscillations with the helical spring after Wilberforce

#### P1.5.4.4

Coupled pendulum - Recording and evaluating with CASSY



Coupled pendulum - Recording and evaluating with CASSY (P1.5.4.4)

Cat. No.	Description	P1.5.4.3	P1.5.4.4
346 51	Wilberforce's pendulum	1	
311 22	Vertical rule	1	
300 11	Saddle base	1	
313 17	Hand-held stop-watch II, mechanical	1	
346 03	Pendulums with axle, pair		1
340 85	Weights, each 50 g, set of 6		1
314 04ET5	Support clip, for plugging in, set of 5		1
352 10	Helical spring, 3 N/m		1
579 43	Motor and tachogenerator, STE 4/19/50		2
524 013	Sensor-CASSY 2		1
524 220	CASSY Lab 2		1
301 25	Support block		2
301 26	Stand rod, 25 cm, 10 mm diam.		1
301 27	Stand rod, 50 cm, 10 mm diam.		2
301 21	Stand base MF		2
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		2
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1



Coupling of longitudinal and rotational oscillations with the helical spring after Wilberforce (P1.5.4.3)

Wilberforce's pendulum is an arrangement for demonstrating coupled longitudinal and rotational oscillations. When a helical spring is elongated, it is always twisted somewhat as well. Therefore, longitudinal oscillations of the helical screw always excite rotational oscillations also. By the same token, the rotational oscillations generate longitudinal oscillations, as torsion always alters the spring length somewhat. The characteristic frequency  $f_l$  of the longitudinal oscillation is determined by the mass  $m$  of the suspended metal cylinder, while the characteristic frequency  $f_R$  of the rotational oscillation is established by the moment of inertia  $I$  of the metal cylinder. By mounting screwable metal disks on radially arranged threaded pins, it becomes possible to change the moment of inertia  $I$  without altering the mass  $m$ .


The first step in the experiment P1.5.4.3 is to match the two frequencies  $f_l$  and  $f_R$  by varying the moment of inertia  $I$ . To test this condition, the metal cylinder is turned one full turn around its own axis and raised 10 cm at the same time. When the frequencies have been properly matched, this body executes both longitudinal and rotational oscillations which do not affect each other. Once this has been done, it is possible to observe for any deflection how the longitudinal and rotational oscillations alternately come to a standstill. In other words, the system behaves like two classical coupled pendulums.

Two coupled pendulums swing in experiment P1.5.4.4 in phase with a frequency  $f_1$  when they are deflected from the rest position by the same distance. When the second pendulum is deflected in the opposite direction, the two pendulums oscillate in opposing phase with the frequency  $f_2$ . Deflecting only one pendulum generates a coupled oscillation with the frequency

$$f_n = \frac{f_1 + f_2}{2}$$

in which oscillation energy is transferred back and forth between the two pendulums. The first pendulum comes to rest after a certain time, while the second pendulum simultaneously reaches its greatest amplitude. The time from one standstill of a pendulum to the next is called  $T_s$ . For the corresponding beat frequency, we can say

$$f_s = |f_1 - f_2|$$



## P1.6.1

### TRANSVERSAL AND LONGITUDINAL WAVES

#### P1.6.1.1

Standing transversal waves on a thread

#### P1.6.1.2

Standing longitudinal waves on a helical spring

Transversal and longitudinal waves (P1.6.1)

Cat. No.	Description	P1.6.1.1	P1.6.1.2
686 57ET5	Rubber cords, 3 m, set of 5	1	1
301 21	Stand base MF	2	2
301 26	Stand rod, 25 cm, 10 mm diam.	1	1
301 27	Stand rod, 50 cm, 10 mm diam.	2	1
666 615	Universal bosshead	1	
301 25	Support block	1	1
314 04ET5	Support clip, for plugging in, set of 5	1	1
579 42	Motor with rocker, STE 2/19	1	1
522 621	Function generator S 12	1	1
301 29	Pointer, pair	1	1
311 77	Steel tape measure, 2 m	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1
352 07ET2	Helical spring, 10 N/m, set of 2		1
352 08ET2	Helical spring, 25 N/m, set of 2		1

A wave is formed when two coupled, oscillating systems sequentially execute oscillations of the same type. The wave can be excited e.g. as a transversal wave on an elastic string or as a longitudinal wave along a helical spring. The propagation velocity of an oscillation state - the phase velocity  $v$  - is related to the oscillation frequency  $f$  and the wavelength  $\lambda$  through the formula

$$v = \lambda \cdot f$$

When the string or the helical spring is fixed at both ends, reflections occur at the ends. This causes superposing of the "outgoing" and reflected waves. Depending on the string length  $s$ , there are certain frequencies at which this superposing of the waves forms stationary oscillation patterns - standing waves. The distance between two oscillation nodes or two antinodes of a standing wave corresponds to one half the wavelength. The fixed ends correspond to oscillation nodes. For a standing wave with  $n$  oscillation antinodes, we can say

$$s = n \cdot \frac{\lambda_n}{2}$$

This standing wave is excited with the frequency

$$f_n = n \cdot \frac{v}{2s}$$

The experiment P1.6.1.1 examines standing string waves. The relationship

$$f_n \sim n$$

is verified.

The experiment P1.6.1.2 looks at standing waves on a helical spring. The relationship

$$f_n \sim n$$

is verified. Two helical springs with different phase velocities  $v$  are provided for use.

P1.6.2

WAVE MACHINE

P1.6.2.1  
Wavelength, frequency and phase velocity of travelling waves



Wavelength, frequency and phase velocity of travelling waves (P1.6.2.1)

Cat. No.	Description	P1.6.2.1
401 20	Wave machine, basic module 1	2
401 22	Drive module for wave machine	1
401 23	Damping module for wave machine	1
401 24	Brake unit for wave machine	2
521 231	Low-voltage power supply, 3/6/9/12 V	1
521 25	Transformer, 2...12 V, 120 W	1
313 07	Hand-held stop watch I, mechanical	1
311 77	Steel tape measure, 2 m	1
501 451	Connecting leads, 19 A, 50 cm, black, pair	1
501 461	Connecting leads, 19 A, 100 cm, black, pair	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1

The "modular wave machine" equipment set enables us to set up a horizontal torsion wave machine, while allowing the size and complexity of the setup within the system to be configured as desired. The module consists of 21 pendulum bodies mounted on edge bearings in a rotating manner around a common axis. They are elastically coupled on both sides of the axis of rotation, so that the deflection of one pendulum propagates through the entire system in the form of a wave.


The aim of the experiment P1.6.2.1 is to explicitly confirm the relationship

$$v = \lambda \cdot f$$

between the wavelength  $\lambda$ , the frequency  $f$  and the phase velocity  $v$ . A stopclock is used to measure the time  $t$  required for any wave phase to travel a given distance  $s$  for different wavelengths; these values are then used to calculate the phase velocity

$$v = \frac{s}{t}$$

The wavelength is then "frozen" using the built-in brake, to permit measurement of the wavelength  $\lambda$ . The frequency is determined from the oscillation period measured using the stopclock. When the recommended experiment configuration is used, it is possible to demonstrate all significant phenomena pertaining to the propagation of linear transversal waves. In particular, these include the excitation of standing waves by means of reflection at a fixed or loose end.



Relationship between the frequency and the wavelength of a propagating wave





### P1.6.3

#### CIRCULARLY POLARIZED WAVES

##### P1.6.3.1

Investigating circularly polarized thread waves in the experiment setup after Melde

##### P1.6.3.2

Determining the phase velocity of circularly polarized thread waves in the experiment setup after Melde

Investigating circularly polarized thread waves in the experiment setup after Melde (P1.6.3.1)

Cat. No.	Description	P1.6.3.1	P1.6.3.2
401 03	Vibrating thread apparatus	1	1
311 77	Steel tape measure, 2 m	1	1
451 281	Stroboscope		1
315 05	Single-pan suspension balance 311		1

The experiment setup after Melde generates circularly polarized string waves on a string with a known length  $s$  using a motordriven eccentric. The tensioning force  $F$  of the string is varied until standing waves with the wavelength

$$\lambda_n = \frac{2s}{n}$$

$n$ : number of oscillation nodes

appear.

In the experiment P1.6.3.1, the wavelengths  $\lambda_n$  of the standing string waves are determined for different string lengths  $s$  and string masses  $m$  at a fixed excitation frequency and plotted as a function of the respective tensioning force  $F_m$ . The evaluation confirms the relationship

$$\lambda \propto \sqrt{\frac{F}{m^*}}$$

with the mass assignment

$$m^* = \frac{m}{s}$$


$m$ : string mass,  $s$ : string length

In the experiment P1.6.3.2, the same measuring procedure is carried out, but with the addition of a stroboscope. This is used to determine the excitation frequency  $f$  of the motor. It also makes the circular polarization of the waves visible in an impressive manner when the standing string wave is illuminated with light flashes which have a frequency approximating that of the standard wave. The additional determination of the frequency  $f$  enables calculation of the phase velocity  $c$  of the string waves using the formula

$$c = \lambda \cdot f$$

as well as quantitative verification of the relationship

$$c = \sqrt{\frac{F}{m^*}}$$



Wavelength  $\lambda$  of thread waves as a function of the tension force  $F$ , the thread length  $s$  and thread density  $m^*$  (P1.6.3.1)

## P1.6.4

### PROPAGATION OF WATER WAVES

#### P1.6.4.1

Exciting circular and straight water waves

#### P1.6.4.2

Huygens' principle in water waves

#### P1.6.4.3

Propagation of water waves in two different depths

#### P1.6.4.4

Refraction of water waves

#### P1.6.4.5

Doppler effect in water waves

#### P1.6.4.6

Reflection of water waves at a straight obstacle

#### P1.6.4.7

Reflection of water waves at curved obstacles



Exciting circular and straight water waves (P1.6.4.1)

Cat. No.	Description	P1.6.4.1	P1.6.4.2	P1.6.4.3	P1.6.4.4-7
401 501	Ripple tank D	1	1	1	1
313 033	Electronic time clock P	1			
311 77	Steel tape measure, 2 m	1		1	

Fundamental concepts of wave propagation can be explained particularly clearly using water waves, as their propagation can be observed with the naked eye.

The experiment P1.6.4.1 investigates the properties of circular and straight waves. The wavelength  $\lambda$  is measured as a function of each excitation frequency  $f$  and these two values are used to calculate the wave velocity

$$v = f \cdot \lambda$$

The aim of the experiment P1.6.4.2 is to verify Huygens' principle. In this experiment, straight waves strike an edge, a narrow slit and a grating. We can observe a change in the direction of propagation, the creation of circular waves and the superposing of circular waves to form one straight wave.

The experiments P1.6.4.3 and P1.6.4.4 aim to study the propagation of water waves in different water depths. A greater water depth corresponds to a medium with a lower refractive index  $n$ . At the transition from one "medium" to another, the law of refraction applies:

$$\frac{\sin \alpha_1}{\sin \alpha_2} = \frac{\lambda_1}{\lambda_2}$$


$\alpha_1, \alpha_2$ : angles with respect to axis of incidence in zone 1 and 2

$\lambda_1, \lambda_2$ : wavelength in zone 1 and 2

A prism, a biconvex lens and a biconcave lens are investigated as practical applications for water waves.

The experiment P1.6.4.5 observes the Doppler effect in circular water waves for various speeds  $u$  of the wave exciter.

The experiments P1.6.4.6 and P1.6.4.7 examine the reflection of water waves. When straight and circular waves are reflected at a straight wall, the "wave beams" obey the law of reflection. When straight waves are reflected by curved obstacles, the originally parallel wave rays travel in either convergent or divergent directions, depending on the curvature of the obstacle. We can observe a focusing to a focal point, respectively a divergence from an apparent focal point, just as in optics.



Convergent beam path behind a biconvex lens (P1.6.4.4)



## P1.6.5

### INTERFERENCE OF WATER WAVES

#### P1.6.5.1

Two-beam interference of water waves

#### P1.6.5.2

Lloyd's experiment on water waves

#### P1.6.5.3

Diffraction of water waves at a slit and at an obstacle

#### P1.6.5.4

Diffraction of water waves at a multiple slit

#### P1.6.5.5

Standing water waves in front of a reflecting barrier

Two-beam interference of water waves (P1.6.5.1)

Cat. No.	Description	P1.6.5.1-4	P1.6.5.5
401 501	Ripple tank D	1	1
311 77	Steel tape measure, 2 m		1

Experiments on the interference of waves can be carried out in an easily understandable manner, as the diffraction objects can be seen and the propagation of the diffracted waves observed with the naked eye.


In the experiment P1.6.5.1, the interference of two coherent circular waves is compared with the diffraction of straight waves at a double slit. The two arrangements generate identical interference patterns.

The experiment P1.6.5.2 reproduces Lloyd's experiment on generating two-beam interference. A second wave coherent to the first is generated by reflection at a straight obstacle. The result is an interference pattern which is equivalent to that obtained for two-beam interference with two discrete coherent exciters.

In the experiment P1.6.5.3, a straight wave front strikes slits and obstacles of various widths. A slit which has a width of less than the wavelength acts like a point-shaped exciter for circular waves. If the slit width is significantly greater than the wavelength, the straight waves pass the slit essentially unaltered. Weaker, circular waves only propagate in the shadow zones behind the edges. When the slit widths are close to the wavelength, a clear diffraction pattern is formed with a broad main maximum flanked by lateral secondary maxima. When the waves strike an obstacle, the two edges of the obstacle act like excitation centers for circular waves. The resulting diffraction pattern depends greatly on the width of the obstacle.

The object of the experiment P1.6.5.4 is to investigate the diffraction of straight water waves at double, triple and multiple slits which have a fixed slit spacing  $d$ . This experiment shows that the diffraction maxima become more clearly defined for an increasing number  $n$  of slits. The angles at which the diffraction maxima are located remain the same.

The experiment P1.6.5.5 demonstrates the generation of standing waves by means of reflection of water waves at a wall arranged parallel to the wave exciter. The standing wave demonstrates points at regular intervals at which the crests and troughs of the individual traveling and reflected waves cancel each other out. The oscillation is always greatest at the midpoint between two such nodes.



Diffraction of water waves at a narrow obstacle (P1.6.5.3)

### P1.6.6

#### RESONANCES

##### P1.6.6.1

Demonstration of resonance vibrations

##### P1.6.6.2

Demonstration of standing waves on a thread

##### P1.6.6.3

Demonstration of standing waves on a helical spring



Demonstration of resonance vibrations (P1.6.6.1)

Cat. No.	Description	P1.6.6.1	P1.6.6.2	P1.6.6.3
346 54	Leaf spring resonance model	1		
587 09	Vibration generator	1	1	1
522 561	Function generator P	1	1	1
501 33	Connecting lead, 32 A, 100 cm, black	2	2	2
309 50	Demonstration cord		1	
340 921ET2	Pulley 100 mm Ø, plug-in, pair, set of 2		1	
683 10	Weight, 0.1 kg		1	
300 01	Stand base, V-shaped, large		2	1
300 44	Stand rod, 100 cm, 12 mm diam.		2	1
301 01	Leybold multiclamp		2	1
301 25	Support block		1	
352 11	Helical spring, 2.7 N/m			1
301 08	Clamp with hook			1

The study of resonance phenomena is an essential topic not just physics but also in applied sciences. A resonance occurs when the excitation frequency is equal to the natural frequency of the oscillation system.

In the experiment P1.6.6.1 the leaf spring model is stimulated with different frequencies  $f$  and thus examined for resonances. Hence the wavelength  $\lambda$  of the vibrations can be identified.

$$\lambda = \frac{c}{f}$$

$c$ : propagation speed

In the experiment P1.6.6.2 a transversal wave is generated. By moving a cord up and down, the wave will be reflected at the fixed end and runs back on the cord. For defined frequencies, certain points will be fixed (nodal) while others (antinodes) oscillate with high amplitude.

In the experiment P1.6.6.3, a longitudinal wave at a helical spring is observed. By moving one end up and down and the other end is fixed, reflects the wave there and runs back. Again for defined frequencies, certain points of the spring remain at rest while others oscillate.



P1.7.1  
SOUND WAVES

P1.7.1.2  
Acoustic beats - Displaying on the oscilloscope

P1.7.1.3  
Acoustic beats - Recording with CASSY

Acoustic beats - Recording with CASSY (P1.7.1.3\_a)

Cat. No.	Description	P1.7.1.2	P1.7.1.3 (a)
414 72	Resonance tuning forks, pair	1	1
586 26	Multi-purpose microphone	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	
575 35	Adapter, BNC/4 mm, 2-pole	1	
300 11	Saddle base	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>		<b>1</b>
524 220	CASSY Lab 2		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

Acoustics is the study of sound and all its phenomena. This discipline deals with both the generation and the propagation of sound waves.

The experiment P1.7.1.2 demonstrates the wave nature of sound. Here, acoustic beats are investigated as the superposing of two sound waves generated using tuning forks with slightly different frequencies  $f_1$  and  $f_2$ . The beat signal is received via a microphone and displayed on the oscilloscope. By means of further (mis-) tuning of one tuning fork by moving a clamping screw, the beat frequency

$$f_s = f_2 - f_1$$

is increased, and the beat period (i. e. the interval between two nodes of the beat signal)


$$T_s = \frac{1}{f_s}$$

is reduced.

In the experiment P1.7.1.3, the acoustic beats are recorded and evaluated via the CASSY computer interface device. The individual frequencies  $f_1$  and  $f_2$ , the oscillation frequency  $f$  and the beat frequency  $f_s$  are determined automatically and compared with the calculated values

$$f = \frac{f_1 + f_2}{2}$$

$$f_s = f_2 - f_1$$



Acoustic beats

## P1.7.2

### OSCILLATIONS OF A STRING

#### P1.7.2.1

Determining the oscillation frequency of a string as a function of the string length and tension



Determining the oscillation frequency of a string as a function of the string length and tension (P1.7.2.1)

Cat. No.	Description	P1.7.2.1
414 01	Monochord	1
314 201	Precision dynamometer, 100.0 N	1
524 013	Sensor-CASSY 2	1
524 074	Timer S	1
524 220	CASSY Lab 2	1
337 46	Fork-type light barrier	1
501 16	Multi-core cable, 6-pole, 1.5 m	1
300 02	Stand base, V-shaped, small	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

In the fundamental oscillation, the string length  $s$  of an oscillating string corresponds to half the wavelength. Therefore, the following applies for the frequency of the fundamental oscillation:

$$f = \frac{c}{2s}$$

where the phase velocity  $c$  of the string is given by

$$c = \sqrt{\frac{F}{A \cdot \rho}}$$


$F$ : tensioning force,  $A$ : area of cross-section,  $\rho$ : density

In the experiment P1.7.2.1, the oscillation frequency of a string is determined as a function of the string length and tensioning force. The measurement is carried out using a forked light barrier and the computer-assisted measuring system CASSY, which is used here as a high-resolution stop-clock. The aim of the evaluation is to verify the relationships

$$f \propto \sqrt{F}$$

and

$$f \propto \frac{1}{s}$$



Frequency  $f$  as a function of the string length  $s$



### P1.7.3 WAVELENGTH AND VELOCITY OF SOUND

P1.7.3.1  
Kundt's tube: determining the wavelength of sound with the cork-powder method

P1.7.3.2  
Determining the wavelength of standing sound waves

Kundt's tube: determining the wavelength of sound with the cork-powder method (P1.7.3.1)

Cat. No.	Description	P1.7.3.1	P1.7.3.2
413 01	Kundt's tube	1	
460 97	Metal rule, 0.5 m	1	
586 26	Multi-purpose microphone		1
587 08	Broad-band speaker		1
522 621	Function generator S 12		1
587 66	Reflection plate		1
300 11	Saddle base		3
311 77	Steel tape measure, 2 m		1
531 120	Multimeter LDanalog 20		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1

Just like other waves, reflection of sound waves can produce standing waves in which the oscillation nodes are spaced at

$$d = \frac{\lambda}{2}$$

Thus, the wavelength  $\lambda$  of sound waves can be easily measured at standing waves.

The experiment P1.7.3.1 investigates standing waves in Kundt's tube. These standing waves are revealed in the tube using cork powder which is stirred up in the oscillation nodes. The distance between the oscillation nodes is used to determine the wavelength  $\lambda$ .

In the experiment P1.7.3.2, standing sound waves are generated by reflection at a barrier. This setup uses a function generator and a loudspeaker to generate sound waves in the entire audible range. A microphone is used to detect the intensity minima, and the wavelength  $\lambda$  is determined from their spacings.




Determining the wavelength of standing sound waves (P1.7.3.2)

P1.7.3

WAVELENGTH AND  
VELOCITY OF SOUND

P1.7.3.3  
Determining the velocity of sound  
in air as a function of the temperature

P1.7.3.4  
Determining the velocity of sound  
in gases



Determining the velocity of sound in air as a function of the temperature (P1.7.3.3)

Cat. No.	Description	P1.7.3.3	P1.7.3.4
413 60	Apparatus for sound velocity	1	1
516 249	Stand for tubes and coils	1	1
587 07	Tweeter	1	1
586 26	Multi-purpose microphone	1	1
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
524 034	Timer box	1	1
524 0673	NiCr-Ni adapter S, type K	1	
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K	1	
521 25	Transformer, 2...12 V, 120 W	1	
300 11	Saddle base	2	2
460 97	Metal rule, 0.5 m	1	1
501 44	Connecting leads, 19 A, 25 cm, red/blue, pair	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	1
660 999	Minican pressurised gas canister, carbon dioxide		1
660 984	Minican pressurised gas canister, helium		1
660 985	Minican pressurised gas canister, neon		1
660 980	Fine regulating valve for minican gas canisters		1
667 194	Silicone tubing, 7 mm diam., 1 m		1
604 481	Rubber tubing, 1 m x 4 mm diam., DIN 12865		1
604 510	Tubing connector, 4...15 mm		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

Sound waves demonstrate only slight dispersion, i.e. group and phase velocities demonstrate close agreement for propagation in gases. Therefore, we can determine the velocity of sound  $c$  as simply the propagation speed of a sonic pulse. In ideal gases, we can say

$$c = \sqrt{\frac{p \cdot \kappa}{\rho}} \quad \text{where } \kappa = \frac{C_p}{C_v}$$

$p$ : pressure,  $\rho$ : density,  $\kappa$ : adiabatic coefficient

$C_p$ ,  $C_v$ : specific heat capacities

The experiment P1.7.3.3 measures the velocity of sound in the air as a function of the temperature  $\vartheta$  and compares it with the linear function resulting from the temperature-dependency of pressure and density

$$c(\vartheta) = c(0) + 0.6 \cdot \frac{\vartheta}{^\circ\text{C}} \frac{m}{s}$$

The value  $c(0)$  determined using a best-fit straight line and the literature values  $p(0)$  and  $\rho(0)$  are used to determine the adiabatic coefficient  $\kappa$  of air according to the formula

$$\kappa = \frac{c(0)^2 \cdot \rho(0)}{p(0)}$$

The experiment P1.7.3.4 determines the velocity of sound  $c$  in carbon dioxide and in the inert gases helium and neon. The evaluation demonstrates that the great differences in the velocities of sound of gases are essentially due to the different densities of the gases. The differences in the adiabatic coefficients of the gases are comparatively small.





### P1.7.3 WAVELENGTH AND VELOCITY OF SOUND

P1.7.3.5  
Determining the velocity of  
sound in solids

Determining the velocity of sound in solids (P1.7.3.5)

Cat. No.	Description	P1.7.3.5
413 651	Metal rods, 1.5 m, set of 3	1
300 46	Stand rod, 150 cm, 12 mm diam.	1
587 251	Rochelle salt crystal (piezo-electric element)	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>
524 220	CASSY Lab 2	1
301 07	Simple bench clamp	1
501 38	Connecting lead, 32 A, 200 cm, black	2
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

In solid bodies, the velocity of sound is determined by the modulus of elasticity  $E$  and the density  $\rho$ . For the velocity of sound in a long rod, we can say

$$c = \sqrt{\frac{E}{\rho}}$$


In the case of solids, measurement of the velocity of sound thus yields a simple method for determining the modulus of elasticity.

The object of the experiment P1.7.3.5 is to determine the velocity of sound in aluminum, copper, brass and steel rods. This measurement exploits the multiple reflections of a brief sound pulse at the rod ends. The pulse is generated by striking the top end of the rod with a hammer, and initially travels to the bottom. The pulse is reflected several times in succession at the two ends of the rod, whereby the pulses arriving at one end are delayed with respect to each other by the time  $\Delta t$  required to travel out and back. The velocity of sound is thus

$$c = \frac{2s}{\Delta t}$$

$s$ : length of rod

To record the pulses, the bottom end of the rod rests on a piezoelectric element which converts the compressive oscillations of the sound pulse into electrical oscillations. These values are recorded using the CASSY system for computer-assisted measured-value recording.



Velocity of sound in different materials

## P1.7.4

### REFLECTION OF ULTRASONIC WAVES

#### P1.7.4.1

Reflection of planar ultrasonic waves at a plane surface

#### P1.7.4.2

Principle of an echo sounder



Reflection of planar ultrasonic waves at a plane surface (P1.7.4.1)

Cat. No.	Description	P1.7.4.1	P1.7.4.2
416 000	Ultrasonic transducer, 40 kHz	2	2
416 014	Generator, 40 kHz	1	1
416 015	AC amplifier	1	1
389 241	Concave mirror	1	
416 020	Sensor holder for concave mirror	1	
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	1
575 24	Screened cable, BNC/4 mm	1	2
460 43	Small optical bench	2	
460 40	Swivel joint with protractor scale	1	
587 66	Reflection plate	1	1
300 01	Stand base, V-shaped, large	1	
300 02	Stand base, V-shaped, small	2	
300 40	Stand rod, 10 cm, 12 mm diam.	1	
301 27	Stand rod, 50 cm, 10 mm diam.	1	
300 41	Stand rod, 25 cm, 12 mm Ø	1	
301 01	Leybold multiclamp	2	
666 615	Universal bosshead	1	
361 03	Spirit level	1	
311 77	Steel tape measure, 2 m	1	
300 42	Stand rod, 47 cm, 12 mm diam.		1
300 11	Saddle base		3
311 02	Metal rule, 1 m		1

When investigating ultrasonic waves, identical, and thus interchangeable transducers are used as transmitters and receivers. The ultrasonic waves are generated by the mechanical oscillations of a piezoelectric body in


the transducer. By the same token, ultrasonic waves excite mechanical oscillations in the piezoelectric body.

The aim of the experiment P1.7.4.1 is to confirm the law of reflection "angle of incidence = angle of reflection" for ultrasonic waves. In this setup, an ultrasonic transducer as a point-type source is set up in the focal point of a concave reflector, so that flat ultrasonic waves are generated. The flat wave strikes a plane surface at an angle of incidence  $\alpha$  and is reflected there. The reflected intensity is measured at different angles using a second transducer. The direction of the maximum reflected intensity is defined as the angle of reflection  $\beta$ .

The experiment P1.7.4.2 utilizes the principle of an echo sounder to determine the velocity of sound in the air, as well as to determine distances. An echo sounder emits pulsed ultrasonic signals and measures the time at which the signal reflected at the boundary surface is received. For the sake of simplicity, the transmitter and receiver are set up as nearly as possible in the same place. When the velocity of sound  $c$  is known, the time difference  $t$  between transmission and reception can be used in the relationship

$$c = \frac{2s}{t}$$

to determine the distance  $s$  to the reflector or, when the distance is known, the velocity of sound.



Echo sounder signal



## P1.7.5 INTERFERENCE OF ULTRASONIC WAVES

- P1.7.5.1  
Beating of ultrasonic waves
- P1.7.5.2  
Interference of two ultrasonic beams
- P1.7.5.3  
Diffraction of ultrasonic waves at a single slit
- P1.7.5.4  
Diffraction of ultrasonic waves at a double slit, a multiple slit and a grating

Beating of ultrasonic waves (P1.7.5.1)

Cat. No.	Description	P1.7.5.1	P1.7.5.2	P1.7.5.3	P1.7.5.4
416 000	Ultrasonic transducer, 40 kHz	3	3	2	2
416 015	AC amplifier	1	1	1	1
416 014	Generator, 40 kHz	2	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analog	1			
575 24	Screened cable, BNC/4 mm	1			
300 11	Saddle base	3	2		
311 902	Rotating platform with motor drive		1	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>		<b>1</b>	<b>1</b>	<b>1</b>
524 220	CASSY Lab 2		1	1	1
524 031	Current source box		1	1	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A		1	1	1
501 031	Connecting lead, protected, 8 m, screened		1	1	1
311 77	Steel tape measure, 2 m		1	1	
300 01	Stand base, V-shaped, large		1	1	1
300 02	Stand base, V-shaped, small		1	1	1
300 41	Stand rod, 25 cm, 12 mm Ø		1	1	1
300 42	Stand rod, 47 cm, 12 mm diam.		1	1	1
301 01	Leybold multiclamp		1	1	1
500 424	Connecting lead 19 A, 50 cm, black		1	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		2	2	2
416 020	Sensor holder for concave mirror			1	1
416 021	Frame with holder			1	1
416 030	Grating and slit for ultrasonics experiments			1	1
389 241	Concave mirror			1	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1	1	1

Experiments on the interference of waves can be carried out in a comprehensible manner using ultrasonic waves, as the diffraction objects are visible with the naked eye. In addition, it is not difficult to generate coherent sound beams.

In the experiment P1.7.5.1, beating of ultrasonic waves is investigated using two transducers which are operated using slightly different frequencies  $f_1$  and  $f_2$ . The signal resulting from the superposing of the two individual signals is interpreted as an oscillation with the periodically varying amplitude

$$A(t) \sim \cos(\pi \cdot (f_2 - f_1) \cdot t)$$

The beat frequency  $f_3$  determined from the period  $T_3$  between two beat nodes and compared with the difference  $f_2 - f_1$ .

In the experiment P1.7.5.2, two identical ultrasonic transducers are operated by a single generator. These transducers generate two coherent ultrasonic beams which interfere with each other. The interference pattern corresponds to the diffraction of flat waves at a double slit when the two transducers are operated in phase. The measured intensity is thus greatest at the diffraction angles  $\alpha$  where

$$\sin \alpha = n \cdot \frac{\lambda}{d} \quad \text{where } n = 0, \pm 1, \pm 2, \dots$$

$\lambda$ : wavelength,  $d$ : spacing of ultrasonic transducers


The experiments P1.7.5.3 and P1.7.5.4 use an ultrasonic transducer as a point-shaped source in the focal point of a concave reflector. The flat ultrasonic waves generated in this manner are diffracted at a single slit, a double slit and a multiple slit. An ultrasonic transducer and the slit are mounted together on the turntable for computer-assisted recording of the diffraction figures. This configuration measures the diffraction at a single slit for various slit widths  $b$  and the diffraction at multiple slits and gratings for different numbers of slits  $N$ .

P1.7.6

ACOUSTIC DOPPLER EFFECT

P1.7.6.1  
Investigating the Doppler effect  
with ultrasonic waves

P1.7.6.2  
Investigating the Doppler effect  
with ultrasonic waves - Recording  
and evaluating with CASSY



Investigating the Doppler effect with ultrasonic waves (P1.7.6.1)

Cat. No.	Description	P1.7.6.1	P1.7.6.2
416 000	Ultrasonic transducer, 40 kHz	2	2
416 015	AC amplifier	1	1
416 014	Generator, 40 kHz	1	1
501 031	Connecting lead, protected, 8 m, screened	1	1
501 644	Two-way adapters, black, set of 6	1	1
685 44ET4	Batteries, 1.5 V (AA), set of 4	1	1
337 07	Trolley with electric drive	1	1
460 81	Precision metal rail, 100 cm	2	2
460 85	Rail connector	1	1
460 88	Feet for metal rails, pair	1	1
460 95ET5	Clamp riders, set of 5	1	1
416 031	Accessories for acoustic Doppler effect	1	1
575 471	Counter S	1	
575 214	Oscilloscope 30 MHz, two-channel, analog	1	
575 24	Screened cable, BNC/4 mm	1	
313 07	Hand-held stop watch I, mechanical	1	
300 02	Stand base, V-shaped, small	1	2
300 11	Saddle base	1	1
300 41	Stand rod, 25 cm, 12 mm Ø	1	1
300 43	Stand rod, 75 cm, 12 mm diam.	1	1
301 01	Leybold multiclamp	1	2
608 100	Stand ring with clamp, 70 mm diam.	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>		<b>1</b>
524 220	CASSY Lab 2		1
524 034	Timer box		1
524 073	Laser motion sensor S		1
300 40	Stand rod, 10 cm, 12 mm diam.		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

The change in the observed frequency for a relative motion of the transmitter and receiver with respect to the propagation medium is called the acoustic Doppler effect. If the transmitter with the frequency  $f_0$  moves at a velocity  $v$  relative to a receiver at rest, the receiver measures the frequency

$$f = \frac{f_0}{1 - \frac{v}{c}}$$

$c$ : velocity of sound


If, on the other hand, the receiver moves at a velocity  $v$  relative to a transmitter at rest, we can say

$$f = f_0 \cdot \left(1 + \frac{v}{c}\right)$$

The change in the frequency  $f - f_0$  is proportional to the frequency  $f_0$ . Investigation of the acoustic Doppler effect on ultrasonic waves thus suggests itself. In the experiment P1.7.6.1, two identical ultrasonic transducers are used as the transmitter and the receiver, and differ only in their connection. One transducer is mounted on a measuring trolley with electric drive, while the other transducer is at rest on the lab bench. The frequency of the received signal is measured using a digital counter. To determine the speed of the transducer in motion, the time  $\Delta t$  which the measuring trolley requires to traverse the measuring distance is measured using a stopclock. In the experiment P1.7.6.2, two identical ultrasonic transducers are used as transmitter and receiver. One transducer is mounted on a measuring trolley with electric drive, while the other transducer is at rest on the lab bench. The frequency of the received signal is measured using a high-resolution digital counter inside the CASSY. The Sensor-CASSY 2 in conjunction with the laser motion sensor S measures the speed of the transducer in motion.



Propagation of sound with the sound source and the observer at rest



## P1.7.7 FOURIER ANALYSIS

P1.7.7.1  
Investigating fast Fourier transforms: simulation of Fourier analysis and Fourier synthesis

P1.7.7.2  
Fourier analysis of the periodic signals of a function generator

P1.7.7.3  
Fourier analysis of an electric oscillator circuit

P1.7.7.4  
Fourier analysis of sounds

Fourier analysis of an electric oscillator circuit (P1.7.7.3)

Cat. No.	Description	P1.7.7.1	P1.7.7.2	P1.7.7.3	P1.7.7.4
524 220	CASSY Lab 2	1	1	1	1
522 621	Function generator S 12		1		
524 013	Sensor-CASSY 2		1	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair		1	4	
562 14	Coil, 500 turns			2	
578 15	Capacitor, 1 μF, STE 2/19			2	
579 10	Push button (NO), STE 2/19			1	
577 19	Resistor, 1 Ω, STE 2/19			1	
577 20	Resistor, 10 Ω, STE 2/19			1	
577 21	Resistor, 5.1 Ω, STE 2/19			1	
577 23	Resistor, 20 Ω, STE 2/19			1	
577 32	Resistor, 100 Ω, STE 2/19			1	
576 74	Plug-in board, DIN A4, STE			1	
524 059	Microphone S				1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1

Fourier analysis and synthesis of sound waves are important tools in acoustics. Thus, for example, knowing the harmonics of a sound is important for artificial generation of sounds or speech.


The experiments P1.7.7.1 and 1.7.7.2 investigate Fourier transforms of periodic signals which are either numerically simulated or generated using a function generator.

In the experiment P1.7.7.3, the frequency spectrum of coupled electric oscillator circuits is compared with the spectrum of an uncoupled oscillator circuit. The Fourier transform of the uncoupled, damped oscillation is a Lorentz curve

$$L(f) = L_0 \cdot \frac{\gamma^2}{(f - f_0)^2 + \gamma^2}$$

in which the width increases with the ohmic resistance of the oscillator circuit. The Fourier-transformed signal of the coupled oscillator circuits shows the split into two distributions lying symmetrically around the uncoupled signal, with their spacing depending on the coupling of the oscillator circuits.

The aim of the experiment P1.7.7.4 is to conduct Fourier analysis of sounds having different tone colors and pitches. As examples, the vowels of the human voice and the sounds of musical instruments are analyzed. The various vowels of a language differ mainly in the amplitudes of the harmonics. The fundamental frequency  $f_0$  depends on the pitch of the voice. This is approx. 200 Hz for high-pitched voices and approx. 80 Hz for low-pitched voices. The vocal tone color is determined by variations in the excitation of the harmonics. The audible tones of musical instruments are also determined by the excitation of harmonics.




Coupled oscillator circuit with its Fourier transform

P1.7.8

ULTRASOUND IN MEDIA

P1.7.8.1  
Optical determination of the velocity of sound in liquids

P1.7.8.2  
Laser diffraction at an ultrasonic wave in liquids (Debye-Sears-Effect)



Laser diffraction at an ultrasonic wave in liquids (Debye-Sears-Effect) (P1.7.8.2)

Cat. No.	Description	P1.7.8.1	P1.7.8.2
417 11	Ultrasound generator, 4 MHz	1	1
460 32	Optical bench with standardised profile, 1 m	1	1
460 374	Optics rider, 90/50	5	4
471 791	Diode laser, 635 nm, 1 mW	1	1
460 02	Lens in frame, f=50 mm	1	
460 25	Prism table	1	1
477 02	Plate glass cell (cuvette), 50 x 100 x 50 mm	1	1
460 380	Cantilever arm	1	1
382 35	Thermometer, -10...+50 °C/0.1 K	1	1
300 41	Stand rod, 25 cm, 12 mm Ø	1	1
301 01	Leybold multiclamp	1	1
441 531	Screen	1	1
675 3410	Water, pure, 5 l	1	1
672 1210	Glycerin, 99 %, 250 ml	1	
671 9740	Ethanol, denaturated, 250 ml	1	
673 5700	Sodium chloride, 250 g	1	

Today's acousto-optic modulators are important building parts for telecommunication and rely on the interaction of sound and light in media. Density variations created by ultrasound are used as diffraction gratings. Experiment P1.7.8.1 measures the wavelength of a standing ultrasound wave in different liquids. The local variation of density in the liquid is made visible on screen by geometrical projection.

Experiment P1.7.8.2 demonstrates the classic Debye-Sears-Effect, i.e. the diffraction of laser light by a phase grating created by ultrasound in a liquid. This is the basic principle of an acousto-optic modulator.



Debye-Sears Effect, Diffraction at an ultrasound grating (P1.7.8.2)



**P1.8.1**  
**BAROMETRIC**  
**MEASUREMENTS**

P1.8.1.1  
Definition of pressure

P1.8.1.2  
Hydrostatic pressure as a  
non-directional quantity

Definition of pressure (P1.8.1.1)

Cat. No.	Description	P1.8.1.1	P1.8.1.2
361 30	Gas syringes with holder, set of 2	1	
315 456	Slotted weight, 100 g	6	
300 02	Stand base, V-shaped, small	1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	
311 77	Steel tape measure, 2 m		1
361 57	Liquid pressure gauge with U-tube manometer	1	
361 575	Glass vessel for liquid pressure gauge		1

In a gas or liquid at rest, the same pressure applies at all points:

$$p = \frac{F}{A}$$

It is measurable as the distributed force  $F$  acting perpendicularly on an area  $A$ .

The experiment P1.8.1.1. aims to illustrate the definition of pressure as the ratio of force and area by experimental means using two gas syringes of different diameters. The pressure is the same in both gas syringes. Thus, we can say for the forces  $F_1$  and  $F_2$  acting on the gas syringes

$$\frac{F_1}{F_2} = \frac{A_1}{A_2}$$


$A_1, A_2$ : cross-section areas

The experiment P1.8.1.2 explores the hydrostatic pressure

$$p = \rho \cdot g \cdot h$$

$\rho$ : density,  $g$ : gravitational acceleration

in a water column subject to gravity. The pressure is measured as a function of the immersion depth  $h$  using a liquid pressure gauge. The displayed pressure remains constant when the gauge is turned in all directions at a constant depth. The pressure is thus a non-directional quantity.



Pressure-gauge reading as a function of the immersion depth (P1.8.1.2)

## P1.8.2

### BOUYANCY

P1.8.2.1  
Confirming Archimedes' principle

P1.8.2.2  
Measuring the buoyancy as a function of the immersion depth



Confirming Archimedes' principle (P1.8.2.1)

Cat. No.	Description	P1.8.2.1	P1.8.2.2
362 02	Archimedes' cylinder	1	1
315 011	Hydrostatic balance	1	
315 31	Set of weights, 10 mg to 200 g	1	
664 111	Beaker, DURAN, 100 ml, tall	1	
664 113	Beaker, DURAN, 250 ml, tall	1	1
672 1210	Glycerin, 99 %, 250 ml	1	1
671 9720	Ethanol, denaturated, 1 l	1	1
314 141	Precision dynamometer, 1 N		1
311 77	Steel tape measure, 2 m		1

Archimedes' principle states that the buoyancy force  $F$  acting on any immersed body corresponds to the weight  $G$  of the displaced liquid.

The experiment P1.8.2.1 verifies Archimedes' principle. In this experiment, a hollow cylinder and a solid cylinder which fits snugly inside it are suspended one beneath the other on the beam of a balance. The deflection of the balance is compensated to zero. When the solid cylinder is immersed in a liquid, the balance shows the reduction in weight due to the buoyancy of the body in the liquid. When the same liquid is filled in the hollow cylinder the deflection of the balance is once again compensated to zero, as the weight of the filled liquid compensates the buoyancy.

In the experiment P1.8.2.2, the solid cylinder is immersed in various liquids to the depth  $h$  and the weight

$$G = \rho \cdot g \cdot A \cdot h$$

$\rho$ : density,  $g$ : gravitational acceleration,  $A$ : cross-section


of the displaced liquid is measured as the buoyancy  $F$  using a precision dynamometer. The experiment confirms the relationship

$$F \sim \rho$$

As long as the immersion depth is less than the height of the cylinder, we can say:

$$F \sim h$$

At greater immersion depths the buoyancy remains constant.



Measuring the buoyancy as a function of the immersion depth (P1.8.2.2)





## P1.8.3

### VISCOSITY

#### P1.8.3.1

Assembling a falling-ball viscosimeter to determine the viscosity of viscous fluids

#### P1.8.3.2

Falling-ball viscosimeter: measuring the viscosity of sugar solutions as a function of the concentration

#### P1.8.3.3

Falling-ball viscosimeter: measuring the viscosity of Newtonian liquids as a function of the temperature

Assembling a falling-ball viscosimeter to determine the viscosity of viscous fluids (P1.8.3.1)

Cat. No.	Description	P1.8.3.1	P1.8.3.2	P1.8.3.3
379 001	Guinea-and-feather apparatus	1		
336 21	Holding magnet	1		
352 54	Steel ball, 16 mm	1		
336 25	Holding magnet adapter with a release mechanism	1		
575 471	Counter S	1		
510 48	Magnets, 35 mm Ø, pair	1		
300 01	Stand base, V-shaped, large	1		
300 41	Stand rod, 25 cm, 12 mm Ø	1		
300 44	Stand rod, 100 cm, 12 mm diam.	1		
301 01	Leybold multiclamp	1		
301 11	Clamp with jaw clamp	1		
311 77	Steel tape measure, 2 m	1		
672 1210	Glycerin, 99 %, 250 ml	6		
590 08ET2	Measuring cylinders, 100 ml, set of 2	1*		
311 54	Precision vernier callipers	1*		
OHC S-200E	Electronic balance, CS200E	1*		
665 906	Höppler falling ball viscometer		1	1
313 07	Hand-held stop watch I, mechanical		1	1
666 7681	Circulation thermostat SC 100-S5P			1
667 194	Silicone tubing, 7 mm diam., 1 m			2
675 3410	Water, pure, 5 l			2

\* additionally recommended

The falling-ball viscometer is used to determine the viscosity of liquids by measuring the falling time of a ball. The substance under investigation is filled in the vertical tube of the viscosimeter, in which a ball falls through a calibrated distance of 100 mm. The resulting falling time  $t$  is a measure of the dynamic viscosity  $\eta$  of the liquid according to the equation

$$\eta = K \cdot (\rho_1 - \rho_2) \cdot t$$

$\rho_2$ : density of the liquid under study

whereby the constant  $K$  and the ball density  $\rho_1$  may be read from the test certificate of the viscosimeter.

The object of the experiment P1.8.3.1 is to set up a falling-ball viscosimeter and to study this measuring method, using the viscosity of glycerine as an example.

The experiment P1.8.3.2 investigates the relationship between viscosity and concentration using concentrated sugar solutions at room temperature.

In the experiment P1.8.3.3, the temperature regulation chamber of the viscosimeter is connected to a circulation thermostat to measure the dependency of the viscosity of a Newtonian fluid (e. g. olive oil) on the temperature.

P1.8.3

VISCOSITY

P1.8.3.4

Hagen-Poiseuille Law



Hagen-Poiseuille Law (P1.8.3.4)

Cat. No.	Description	P1.8.3.4
665 195	Capillary tube, 300 mm x 5 mm, 0.8 mm diam.	1
665 196	Capillary tube, 300 mm x 8 mm, 2 mm diam.	1
665 205	Capillary tube, 300 mm x 8 mm, 1 mm diam.	1
667 510	Capillary tube	1
460 21	Holder for plug-in elements	1
590 02ET2	Clip plugs, small, set of 2	1
LDS 00001	Stopwatch, digital	1
382 21	Stirring thermometer, -30...+110 °C	1
311 02	Metal rule, 1 m	1
362 05	Discharge vessel	1
665 752	Measuring cylinder, 25 ml, with plastic base	1
608 160	Mohr clip, 50 mm	1
665 227	Connector, straight, 4 ... 8 mm diam.	1
667 194	Silicone tubing, 7 mm diam., 1 m	1
667 197	Silicone tubing, 4 mm diam., 1 m	1
300 01	Stand base, V-shaped, large	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
300 44	Stand rod, 100 cm, 12 mm diam.	1
301 09	Bosshead S	2
648 01	Storage tray S8-FN	1
675 3400	Water, pure, 1 l	1
460 135	Ocular with scale	1*
460 317	Optical bench, S1 profile, 0.5 m	1*
460 312	Clamp rider with clamp, 45/35	2*

\* additionally recommended

For a Newtonian liquid (i.e. laminar case) the volume flow rate  $J$  and the pressure difference  $\Delta p$  are related by the Hagen-Poiseuille law:

$$J = \frac{\pi \cdot \Delta p \cdot r^4}{8 \cdot \eta \cdot L}$$

$\eta$ : Viscosity of the liquid,

$L$ : Length of the capillary tube,

$r$ : radius of the capillary tube

In the experiment P1.8.3.4 the Hagen-Poiseuille law can be verified by using several kind of capillary tubes and various effective level differences to get variant pressure difference.



## P1.8.4 SURFACE TENSION

P1.8.4.1  
Measuring the surface tension using the „break-away“ method

P1.8.4.2  
Measuring the surface tension using the „break-away“ method - Recording and evaluating with CASSY

Measuring the surface tension using the „break-away“ method (P1.8.4.1)

Cat. No.	Description	P1.8.4.1	P1.8.4.2
367 46	Surface tension determination device	1	1
664 175	Crystallisation dish, 95 mm diam., 300 ml	1	1
314 111	Precision dynamometer , 0.1 N	1	
311 53	Vernier callipers	1	1
300 76	Laboratory stand II	1	1
300 02	Stand base, V-shaped, small	1	1
300 43	Stand rod, 75 cm, 12 mm diam.	1	
301 08	Clamp with hook	1	
671 9740	Ethanol, denaturated, 250 ml	1	1
675 3400	Water, pure, 1 l	1	1
524 060	Force sensor S, ±1N		1
524 013	Sensor-CASSY 2		1
524 220	CASSY Lab 2		1
300 42	Stand rod, 47 cm, 12 mm diam.		1
301 01	Leybold multiclamp		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

To determine the surface tension  $\sigma$  of a liquid, a metal ring is suspended horizontally from a precision dynamometer or a force sensor. The metal ring is completely immersed in the liquid, so that the entire surface is wetted. The ring is then slowly pulled out of the liquid, drawing a thin sheet of liquid behind it. The liquid sheet tears when the tensile force exceeds a limit value

$$F = \sigma \cdot 4\pi \cdot R$$

$R$ : edge radius

The experiments P1.8.4.1 and P1.8.4.2 determines the surface tension of water and ethanol. It is shown that water has a particularly high surface tension in comparison to other liquids (literature value for water:  $0.073 \text{ Nm}^{-1}$ , for ethanol:  $0.022 \text{ Nm}^{-1}$ ).

### P1.8.5

#### INTRODUCTORY EXPERIMENTS ON AERODYNAMICS

##### P1.8.5.1

Static pressure in a reduced cross-section  
- Measuring the pressure with the precision manometer

##### P1.8.5.2

Determining the volume flow with a Venturi tube - Measuring the pressure with the precision manometer

##### P1.8.5.3

Determining the wind speed with a Prandtl pressure probe - Measuring the pressure with the precision manometer

##### P1.8.5.4

Static pressure in a reduced cross-section - Measuring the pressure with a pressure sensor and CASSY

##### P1.8.5.5

Determining the volume flow with a Venturi tube - Measuring the pressure with a pressure sensor and CASSY

##### P1.8.5.6

Determining the wind speed with a Prandtl pressure probe - Measuring the pressure with a pressure sensor and CASSY



Static pressure in a reduced cross-section - Measuring the pressure with the precision manometer (P1.8.5.1)

Cat. No.	Description	P1.8.5.1-2	P1.8.5.3	P1.8.5.4-5	P1.8.5.6
373 041	Suction and pressure fan	1	1	1	1
373 091	Venturi tube with multimanoscope	1		1	
373 10	Precision manometer	1	1		
300 02	Stand base, V-shaped, small	2	1	1	
300 41	Stand rod, 25 cm, 12 mm Ø	1		1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	1		
301 01	Leybold multiclamp	1		1	
373 13	Prandtl pressure probe		1		1
524 005	Mobile-CASSY 2			1	1
524 066	Pressure sensor S, ±70 hPa			1	1

The study of aerodynamics relies on describing the flow of air through a tube using the continuity equation and the Bernoulli equation. These state that regardless of the cross-section  $A$  of the tube, the volume flow

$$\dot{V} = v \cdot A$$

$v$ : flow speed

and the total pressure

$$p_0 = p + p_s \quad \text{where} \quad p_s = \frac{\rho}{2} \cdot v^2$$

$p$ : static pressure,  $p_s$ : dynamic pressure,  $\rho$ : density of air

remain constant as long as the flow speed remains below the speed of sound.

*Note:* In the experiments P1.8.5.1 – P1.8.5.3, the precision manometer is used to measure pressures. In addition to a pressure scale, it is provided with a further scale which indicates the flow speed directly when measuring with the pressure head sensor. In the experiments P1.8.5.4 – P1.8.5.6 the pressure is measured with a pressure sensor and recorded using the universal measuring instrument Mobile-CASSY.

In order to verify these two equations, the static pressure in a Venturi tube is measured for different cross-sections in the experiments P1.8.5.1 and P1.8.5.4. The static pressure decreases in the reduced cross-section, as the flow speed increases here.

The experiments P1.8.5.2 and P1.8.5.5 uses the Venturi tube to measure the volume flow. Using the pressure difference  $\Delta p = p_2 - p_1$  between two points with known cross-sections  $A_1$  and  $A_2$ , we obtain

$$v_1 \cdot A_1 = \sqrt{\frac{2 \cdot \Delta p \cdot A_2^2}{\rho \cdot (A_2^2 - A_1^2)}}$$

The experiments P1.8.5.3 and P1.8.5.6 aims to determine flow speeds. Here, dynamic pressure (also called the "pressure head") is measured using the pressure head sensor after *Prandtl* as the difference between the total pressure and the static pressure, and this value is used to calculate the speed at a known density  $\rho$ .



## P1.8.6 MEASURING AIR RESISTANCE

- P1.8.6.1  
Measuring the air resistance as a function of the wind speed - Measuring the wind speed with the precision manometer
- P1.8.6.2  
Drag coefficient  $c_w$ : relationship between air resistance and body shape - Measuring the wind speed with the precision manometer
- P1.8.6.3  
Pressure curve on an airfoil profile - Measuring the pressure with the precision manometer
- P1.8.6.4  
Measuring the air resistance as a function of the wind speed - Measuring the pressure with a pressure sensor and CASSY
- P1.8.6.5  
Drag coefficient  $c_w$ : relationship between air resistance and body shape - Measuring the wind speed with a pressure sensor and CASSY
- P1.8.6.6  
Pressure curve on an airfoil profile - Measuring the pressure with a pressure sensor and CASSY

Measuring the air resistance as a function of the wind speed - Measuring the wind speed with the precision manometer (P1.8.6.1)

Cat. No.	Description	P1.8.6.1-2	P1.8.6.3	P1.8.6.4-5	P1.8.6.6
373 041	Suction and pressure fan	1	1	1	1
373 06	Open aerodynamics working section	1	1	1	1
373 071	Aerodynamics accessories 1	1		1	
373 075	Measurement trolley for wind tunnel	1		1	
373 14	Sector dynamometer, 0.65 N	1		1	
373 13	Prandtl pressure probe	1	1	1	1
373 10	Precision manometer	1	1		
300 02	Stand base, V-shaped, small	1	2	1	1
300 11	Saddle base	1		1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	1	1	
373 70	Air foil model		1		1
<b>524 005</b>	<b>Mobile-CASSY 2</b>			<b>1</b>	<b>1</b>
524 066	Pressure sensor S, $\pm 70$ hPa			1	1

A flow of air exercises a force  $F_w$  on a body in the flow which is parallel to the direction of the flow; this force is called the air resistance. This force depends on the flow speed  $v$ , the cross-section  $A$  of the body perpendicular to the flow direction and the shape of the body. The influence of the body shape is described using the so-called drag coefficient  $c_w$ , whereby the air resistance is determined as:

$$F_w = c_w \cdot \frac{\rho}{2} \cdot v^2 \cdot A$$

Note: In the experiments P1.8.6.1 - P1.8.6.3, the precision manometer is used to measure pressures. In addition to a pressure scale, it is provided with a further scale which indicates the flow speed directly when measuring with the pressure head sensor. In the experiments P1.8.6.4 - P1.8.6.6 the pressure is measured with a pressure sensor and recorded using the universal measuring instrument Mobile-CASSY.

The experiments P1.8.6.1 and P1.8.6.4 examines the relationship between the air resistance and the flow speed using a circular disk. The flow speed is measured using a pressure head sensor and the air resistance with a dynamometer.

The experiments P1.8.6.2 and P1.8.6.5 determines the drag coefficient  $c_w$  for various flow bodies with equal cross-sections. The flow speed is measured using a pressure head sensor and the air resistance with a dynamometer.

The aim of the experiments P1.8.6.3 and P1.8.6.6 is to measure the static pressure  $p$  at various points on the underside of an airfoil profile. The measured curve not only illustrates the air resistance, but also explains the lift acting on the airfoil.

## P1.8.7

### MEASUREMENTS IN A WIND TUNNEL

P1.8.7.1  
Recording the airfoil profile polars in a wind tunnel

P1.8.7.2  
Measuring students' own airfoils and panels in the wind tunnel

P1.8.7.3  
Verifying the Bernoulli equation -  
Measuring the pressure with the precision manometer

P1.8.7.4  
Verifying the Bernoulli equation -  
Measuring the pressure with a pressure sensor and CASSY



Recording the airfoil profile polars in a wind tunnel (P1.8.7.1)

Cat. No.	Description	P1.8.7.1-2	P1.8.7.3	P1.8.7.4
373 12	Wind tunnel	1	1	1
373 041	Suction and pressure fan	1	1	1
373 075	Measurement trolley for wind tunnel	1	1	1
373 08	Aerodynamics accessories 2	1		
373 14	Sector dynamometer, 0.65 N	1		
373 13	Prandtl pressure probe		1	1
373 10	Precision manometer		1	
524 005	Mobile-CASSY 2			1
524 066	Pressure sensor S, ±70 hPa			1

The wind tunnel provides a measuring configuration for quantitative experiments on aerodynamics that ensures an airflow which has a constant speed distribution with respect to both time and space. Among other applications, it is ideal for measurements on the physics of flight.

In the experiment P1.8.7.1, the air resistance  $f_W$  and the lift  $F_A$  of an airfoil are measured as a function of the angle of attack  $\alpha$  of the airfoil against the direction of flow. In a polar diagram,  $F_W$  is graphed as a function of  $F_A$  with the angle of attack  $\alpha$  as the parameter. From this polar diagram, we can read e. g. the optimum angle of attack.

In the experiment P1.8.7.2, the students perform comparable measurements on airfoils of their own design. The aim is to determine what form an airfoil must have to obtain the smallest possible quotient  $F_W / F_A$  at a given angle of attack  $\alpha$ .

The experiments P1.8.7.3 and P1.8.7.4 verify the Bernoulli equation. The difference between the total pressure and the static pressure is measured as a function of the cross-section, whereby the cross-section of the wind tunnel is gradually reduced by means of a built-in ramp. If we assume that the continuity equation applies, the cross-section  $A$  provides a measure of the flow speed  $v$  due to

$$v = \frac{v_0 \cdot A_0}{A}$$

$v_0$ : flow speed at cross-section  $A_0$

The experiment verifies the following relationship, which follows from the Bernoulli equation:

$$\varphi \sim \frac{1}{A^2}$$

# P2 HEAT



P2.1 THERMAL EXPANSION	65
P2.2 HEAT TRANSFER	68
P2.3 HEAT AS A FORM OF ENERGY	70
P2.4 PHASE TRANSITIONS	74
P2.5 KINETIC THEORY OF GASES	77
P2.6 THERMODYNAMIC CYCLE	81

# P2 HEAT



## P2.1 THERMAL EXPANSION

- P2.1.1 Thermal expansion of solids
- P2.1.2 Thermal expansion of liquids
- P2.1.3 Thermal anomaly of water

65  
66  
67

## P2.2 HEAT TRANSFER

- P2.2.1 Thermal conductivity
- P2.2.2 Solar collector

68  
69

## P2.3 HEAT AS A FORM OF ENERGY

- P2.3.1 Mixing temperatures
- P2.3.2 Heat capacities
- P2.3.3 Converting mechanical energy into heat
- P2.3.4 Converting electrical energy into heat

70  
71  
72  
73

## P2.4 PHASE TRANSITIONS

- P2.4.1 Latent heat and vaporization heat
- P2.4.2 Measuring vapor pressure
- P2.4.3 Critical temperature

74  
75  
76

## P2.5 KINETIC THEORY OF GASES

- P2.5.1 Brownian motion of molecules
- P2.5.2 Gas laws
- P2.5.3 Specific heat of gases
- P2.5.4 Real gases

77  
78  
79  
80

## P2.6 THERMODYNAMIC CYCLE

- P2.6.1 Hot-air engine: qualitative experiments
- P2.6.2 Hot-air engine: quantitative experiments
- P2.6.3 Heat pump

81-82  
83-84  
85





## P2.1.1

### THERMAL EXPANSION OF SOLIDS

P2.1.1.1  
Thermal expansion of solids -  
Measuring using STM equipment

P2.1.1.2  
Thermal expansion of solids -  
Measuring using the expansion  
apparatus

P2.1.1.3  
Measuring the linear expansion  
of solids as a function of temperature

P2.1.1.4  
Thermal expansion of solids -  
Recording and evaluating with CASSY

Thermal expansion of solids - Measuring using the expansion apparatus (P2.1.1.2)

Cat. No.	Description	P2.1.1.1	P2.1.1.2	P2.1.1.3	P2.1.1.4
301 21	Stand base MF	2			2
301 27	Stand rod, 50 cm, 10 mm diam.	2			2
301 26	Stand rod, 25 cm, 10 mm diam.	1			
301 25	Support block	2			
301 09	Bosshead S	2			2
666 555	Universal clamp, 0...80 mm	1			1
664 248	Erlenmeyer flask, Boro 3.3, 50 ml, narrow neck	1			1
667 2545	Rubber stopper, one 7-mm hole, 17...23 mm diam.	1			1
665 226	Connector, straight, 6 ... 8 mm diam.	1			1
667 194	Silicone tubing, 7 mm diam., 1 m	1	1	2	1
664 183	Petri dishes	1			
314 04ET5	Support clip, for plugging in, set of 5	1			
340 82	Double scale	1			
381 331	Pointer for linear expansion	1			
381 332	Aluminium tube, 44 cm x 8 mm diam.	1			1
381 333	Iron tube, 44 cm x 8 mm diam.	1			1
311 77	Steel tape measure, 2 m	1			1
303 22	Alcohol burner, metal	1			
381 341	Longitudinal expansion apparatus D		1	1	
361 152	Dial gauge with holder		1	1	
382 34	Thermometer, -10...+110 °C/0.2 K		1	1	1
303 28	Steam generator		1		
664 185	Petri dishes		1		
666 7681	Circulation thermostat SC 100-S5P			1	
675 3410	Water, pure, 5 l			2	
686 53ET5	Round tins with cap, set of 5				1
524 013	Sensor-CASSY 2				1
524 220	CASSY Lab 2				1

Cat. No.	Description	P2.1.1.1	P2.1.1.2	P2.1.1.3	P2.1.1.4
524 082	Rotary motion sensor S				1
666 615	Universal bosshead				1
302 68	Stand ring with stem, 13 cm diam.				1
666 685	Wire gauze, 160 mm x 160 mm				1
666 711	Butane gas burner				1
666 712ET3	Butane cartridge, 190 g, set of 3				1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)				1

The relationship between the length  $s$  and the temperature  $\vartheta$  of a liquid is approximately linear:

$$s = s_0 \cdot (1 + \alpha \cdot \vartheta)$$

$$s_0: \text{length at } 0^\circ\text{C}, \vartheta: \text{temperature in } ^\circ\text{C}$$

The linear expansion coefficient  $\alpha$  is determined by the material of the solid body. We can conduct measurements on this topic using e.g. thin tubes through which hot water or steam flows.

In the experiment P2.1.1.1, steam is channeled through different tube samples. The thermal expansion is measured in a simple arrangement, and the dependency on the material is demonstrated.

The experiment P2.1.1.2 measures the increase in length of various tube samples between room temperature and steam temperature using the expansion apparatus. The effective length  $s_0$  of each tube can be defined as 200, 400 or 600 mm.

In the experiment P2.1.1.3, a circulation thermostat is used to heat the water, which flows through various tube samples. The expansion apparatus measures the change in the lengths of the tubes as a function of the temperature  $\vartheta$ .

In the experiment P2.1.1.4, steam is channeled through different tube samples. The thermal expansion is measured with a rotary motion sensor.

## P2.1.2

### THERMAL EXPANSION OF LIQUIDS

#### P2.1.2.1

Determining the volumetric expansion coefficient of liquids



Determining the volumetric expansion coefficient of liquids (P2.1.2.1\_c)

Cat. No.	Description	P2.1.2.1 (G)
382 15	Dilatometer	1
524 005	Mobile-CASSY 2	1
666 767	Hotplate, 1500 W, 180 mm diam.	1
664 104	Beaker, DURAN, 400 ml, squat	1
315 05	Single-pan suspension balance 311	1
300 02	Stand base, V-shaped, small	1
300 42	Stand rod, 47 cm, 12 mm diam.	1
301 01	Leybold multiclamp	2
666 555	Universal clamp, 0...80 mm	2
671 9720	Ethanol, denaturated, 1 l	1

In general, liquids expands more than solids when heated. The relationship between the Volume  $V$  and the temperature  $\vartheta$  of a liquid is approximately linear here:

$$V = V_0 \cdot (1 + \gamma \cdot \vartheta)$$


$V_0$ : volume at 0 °C,  $\vartheta$ : temperature in °C

When determining the volumetric expansion coefficient  $\gamma$ , it must be remembered that the vessel in which the liquid is heated also expands.

In the experiment P2.1.2.1, the volumetric expansion coefficients of water and methanol are determined using a volume dilatometer made of glass. An attached riser tube with a known cross-section is used to measure the change in volume. i.e. the change in volume is determined from the rise height of the liquid.

P2.1.3  
THERMAL ANOMALY  
OF WATER

P2.1.3.1  
Investigating the density  
maximum of water



Investigating the density maximum of water (P2.1.3.1\_c)


Cat. No.	Description	P2.1.3.1 (G)
667 505	Anomaly of water apparatus	1
666 8451	Magnetic stirrer	1
664 195	Glass tanks	1
665 009	Funnel, PP, 75 mm diam.	1
307 66	Tubing (rubber)	1
300 42	Stand rod, 47 cm, 12 mm diam.	1
666 555	Universal clamp, 0...80 mm	1
301 01	Leybold multiclamp	1
300 02	Stand base, V-shaped, small	1
608 100	Stand ring with clamp, 70 mm diam.	1
524 005	Mobile-CASSY 2	1

When heated from a starting temperature of 0 °C, water demonstrates a critical anomaly: it has a negative volumetric expansion coefficient up to 4 °C, i.e. it contracts when heated. After reaching zero at 4 °C, the volumetric expansion coefficient takes on a positive value. As the density corresponds to the reciprocal of the volume of a quantity of matter, water has a density maximum at 4 °C.

The experiment P2.1.3.1 verifies the density maximum of water by measuring the expansion in a vessel with riser tube. Starting at room temperature, the complete setup is cooled in a constantly stirred water bath to about 1 °C, or alternatively allowed to gradually reach the ambient temperature after cooling in an ice chest or refrigerator. The rise height  $h$  is measured as a function of the temperature  $\vartheta$ . As the change in volume is very slight in relation to the total volume  $V_0$ , we obtain the density

$$\rho(\vartheta) = \rho(0^\circ\text{C}) \cdot \left( 1 - \frac{A}{V_0} \cdot h(\vartheta) \right)$$

A: cross-section of riser tube



Relative density of water as a function of the temperature

## P2.2.1

### THERMAL CONDUCTIVITY

#### P2.2.1.1

Determining the thermal conductivity of building materials using the single-plate method

#### P2.2.1.2

Determining the thermal conductivity of building materials using the heat flux plate principle

#### P2.2.1.3

Damping temperature fluctuations using multiple-layered walls



Determining the thermal conductivity of building materials using the single-plate method (P2.2.1.1)


Cat. No.	Description	P2.2.1.1	P2.2.1.2	P2.2.1.3
389 29	Calorimetric chamber	1	1	1
389 30	Building materials for calorimetric chamber	1	1	1
521 25	Transformer, 2...12 V, 120 W	1	1	1
524 013	Sensor-CASSY 2	1	1	1
524 220	CASSY Lab 2	1	1	1
524 0673	NiCr-Ni adapter S, type K	1	2	2
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K	2	3	3
501 451	Connecting leads, 19 A, 50 cm, black, pair	1		
501 33	Connecting lead, 32 A, 100 cm, black	4	2	2
450 64	Halogen lamp, 12 V, 50/100 W			1
450 63	Halogen bulb, 12 V/100 W, G6.35			1
300 11	Saddle base			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

In the equilibrium state, the heat flow through a plate with the cross-section area  $A$  and the thickness  $d$  depends on the temperature difference  $\vartheta_2 - \vartheta_1$  between the front and rear sides and on the thermal conductivity  $\lambda$  of the plate material:

$$\frac{\Delta Q}{\Delta t} = \lambda \cdot A \cdot \frac{\vartheta_2 - \vartheta_1}{d}$$

The object of the experiments P2.2.1.1 und P2.2.1.2 is to determine the thermal conductivity of building materials. In these experiments, sheets of building materials are placed in the heating chamber and their front surfaces are heated. The temperatures  $\vartheta_1$  and  $\vartheta_2$  are measured using measuring sensors. The heat flow is determined either from the electrical power of the hot plate or by measuring the temperature using a reference material with known thermal conductivity  $\lambda_0$  which is pressed against the sheet of the respective building material from behind.

The experiment P2.2.1.3 demonstrates the damping of temperature variations by means of two-layer walls. The temperature changes between day and night are simulated by repeatedly switching a lamp directed at the outside surface of the wall on and off. This produces a temperature "wave" which penetrates the wall; the wall in turn damps the amplitude of this wave. This experiment measures the temperatures  $\vartheta_A$  on the outer surface,  $\vartheta_Z$  between the two layers and  $\vartheta_I$  on the inside as a function of time.



Temperature variations in multi-layer walls (P2.2.1.3)



## P2.2.2 SOLAR COLLECTOR

P2.2.2.1  
Determining the efficiency of a solar collector as a function of the throughput volume of water

P2.2.2.2  
Determining the efficiency of a solar collector as a function of the heat insulation

Determining the efficiency of a solar collector as a function of the throughput volume of water (P2.2.2.1\_b)

Cat. No.	Description	P2.2.2.1-2 (b)
389 50	Solar collector	1
579 220	Water pump, STE 2/50	1
450 72	Floodlight lamp, 1000 W, with light shades	1
521 35	Variable extra-low voltage transformer S	1
524 005	Mobile-CASSY 2	1
524 0673	NiCr-Ni adapter S, type K	1
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K	2
311 77	Steel tape measure, 2 m	1
313 17	Hand-held stop-watch II, mechanical	1
300 02	Stand base, V-shaped, small	2
300 41	Stand rod, 25 cm, 12 mm Ø	1
300 42	Stand rod, 47 cm, 12 mm diam.	1
300 43	Stand rod, 75 cm, 12 mm diam.	1
301 01	Leybold multiclamp	3
666 555	Universal clamp, 0...80 mm	1
590 06	Plastic beaker	1
604 431	Silicone tubing, 5 mm diam., 1 m	1
604 432	Silicone tubing, 6 mm diam., 1 m	1
604 434	Silicone tubing, 8 mm diam., 1 m	1
665 226	Connector, straight, 6 ... 8 mm diam.	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1

A solar collector absorbs radiant energy to heat the water flowing through it. When the collector is warmer than its surroundings, it loses heat to its surroundings through radiation, convection and heat conductivity. These losses reduce the efficiency

$$\eta = \frac{\Delta Q}{\Delta E}$$

i. e. the ratio of the emitted heat quantity  $\Delta Q$  to the absorbed radiant energy  $\Delta E$ .

In the experiments P2.2.2.1 and P2.2.2.2, the heat quantity  $\Delta Q$  emitted per unit of time is determined from the increase in the temperature of the water flowing through the apparatus, and the radiant energy absorbed per unit of time is estimated on the basis of the power of the lamp and its distance from the absorber. The throughput volume of the water and the heat insulation of the solar collector are varied in the course of the experiment.

### P2.3.1

#### MIXING TEMPERATURES

##### P2.3.1.1

##### Mixing temperature of water



Mixing temperature of water (P2.3.1.1\_a)

Cat. No.	Description	P2.3.1.1 (a)
384 161	Cover for Dewar vessel	1
386 48	Dewar vessel calorimeter, 250 ml	1
382 34	Thermometer, -10...+110 °C/0.2 K	1
315 23	Single-pan suspension balance 610 Tara	1
313 07	Hand-held stop watch I, mechanical	1
666 767	Hotplate, 1500 W, 180 mm diam.	1
664 104	Beaker, DURAN, 400 ml, squat	2

When cold water with the temperature  $\vartheta_1$  is mixed with warm or hot water having the temperature  $\vartheta_2$ , an exchange of heat takes place until all the water reaches the same temperature. If no heat is lost to the surroundings, we can formulate the following for the mixing temperature:

$$\vartheta_m = \frac{m_1}{m_1 + m_2} \vartheta_1 + \frac{m_2}{m_1 + m_2} \vartheta_2$$

$m_1$ ,  $m_2$ : mass of cold and warm water respectively

Thus the mixing temperature  $\vartheta_m$  is equivalent to a weighted mean value of the two temperatures  $\vartheta_1$  and  $\vartheta_2$ .

The use of the Dewar flask in the experiment P2.3.1.1 essentially prevents the loss of heat to the surroundings. This vessel has a double wall; the intermediate space is evacuated and the interior surface is mirrored. The water is stirred thoroughly to ensure a complete exchange of heat. This experiment measures the mixing temperature  $\vartheta_m$  for different values for  $\vartheta_1$ ,  $\vartheta_2$ ,  $m_1$  and  $m_2$ .

### P2.3.2

#### HEAT CAPACITIES

##### P2.3.2.1

Determining the specific heat of solids



Determining the specific heat of solids (P2.3.2.1\_c)

Cat. No.	Description	P2.3.2.1 (c)
384 161	Cover for Dewar vessel	1
386 48	Dewar vessel calorimeter, 250 ml	1
<b>524 005</b>	<b>Mobile-CASSY 2</b>	<b>1</b>
384 34	Heating apparatus	1
384 35	Copper shot, 200 g	1
384 36	Glass shot, 100 g	1
315 76	Lead shot, 200 g	1
315 23	Single-pan suspension balance 610 Tara	1
303 28	Steam generator	1
664 104	Beaker, DURAN, 400 ml, squat	1
667 194	Silicone tubing, 7 mm diam., 1 m	1
300 02	Stand base, V-shaped, small	1
300 42	Stand rod, 47 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
666 555	Universal clamp, 0...80 mm	1
667 614	Heat protective gloves	1

When a body is heated or cooled, the absorbed heat capacity  $\Delta Q$  is proportional to the change in temperature  $\Delta\vartheta$  and to the mass  $m$  of the body:

$$\Delta Q = c \cdot m \cdot \Delta\vartheta$$

The proportionality factor  $c$ , the specific heat capacity of the body, is a quantity which depends on the respective material.

To determine the specific heat capacity in experiment P2.3.2.1, various materials in particle form are weighed, heated in steam to the temperature  $\vartheta_1$  and poured into a weighed-out quantity of water with the temperature  $\vartheta_2$ . After careful stirring, heat exchange ensures that the particles and the water have the same temperature  $\vartheta_m$ . The heat quantity released by the particles:

$$\Delta Q_1 = c_1 \cdot m_1 \cdot (\vartheta_1 - \vartheta_m)$$

$m_1$ : mass of particles

$c_1$ : specific heat capacity of particles

is equal to the quantity absorbed by the water

$$\Delta Q_2 = c_2 \cdot m_2 \cdot (\vartheta_m - \vartheta_2)$$

$m_2$ : mass of water

The specific heat capacity of water  $c_2$  is assumed as a given. The temperature  $\vartheta_1$  corresponds to the temperature of the steam. Therefore, the specific heat quantity  $c_1$  can be calculated from the measurement quantities  $\vartheta_2$ ,  $\vartheta_m$ ,  $m_1$  and  $m_2$ .

### P2.3.3

#### CONVERTING MECHANICAL ENERGY INTO HEAT

##### P2.3.3.1

Converting mechanical energy into heat energy - Recording and evaluating measured values manually

##### P2.3.3.2

Converting mechanical energy into heat energy - Recording and evaluating with CASSY



Converting mechanical energy into heat energy - Recording and evaluating measured values manually (P2.3.3.1\_a)

Cat. No.	Description	P2.3.3.1 (a)	P2.3.3.2
388 00	Equivalent of heat, basic apparatus	1	1
388 01	Water calorimeter	1	1
388 02	Copper-block calorimeter	1	1
388 03	Aluminium-block calorimeter	1	1
388 04	Aluminium-block calorimeter, large	1	1
388 05	Thermometer for calorimeters, +15...35 °C/0.2 K	1	
388 24	Weight, 5 kg	1	1
524 013	Sensor-CASSY 2		1
524 220	CASSY Lab 2		1
524 074	Timer S		1
524 0673	NiCr-Ni adapter S, type K		1
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K		1
337 46	Fork-type light barrier		1
501 16	Multi-core cable, 6-pole, 1.5 m		1
300 02	Stand base, V-shaped, small		1
301 11	Clamp with jaw clamp		1
300 40	Stand rod, 10 cm, 12 mm diam.		1
300 41	Stand rod, 25 cm, 12 mm Ø		1
301 07	Simple bench clamp		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

Energy is a fundamental quantity of physics. This is because the various forms of energy can be converted from one to another and are thus equivalent to each other, and because the total energy is conserved in the case of conversion in a closed system.

These experiments P2.3.3.1 und P2.3.3.2 show the equivalence of mechanical and heat energy. A hand crank is used to turn various calorimeter vessels on their own axes, and friction on a nylon belt causes them to become warmer. The friction force is equivalent to the weight  $G$  of a suspended weight. For  $n$  turns of the calorimeter, the mechanical work is thus

$$W_n = G \cdot n \cdot \pi \cdot d$$

$d$ : diameter of calorimeter

This results in an increase in the temperature of the calorimeter which corresponds to the specific heat capacity

$$Q_n = m \cdot c \cdot (\vartheta_n - \vartheta_0)$$

$c$ : specific heat capacity,  $m$ : mass,  
 $\vartheta_n$ : temperature after  $n$  turns

To confirm the relationship

$$Q_n = W_n$$

the two quantities are plotted together in a diagram. In the experiment P2.3.3.1, the measurement is conducted and evaluated manually point by point. The experiment P2.3.3.2 takes advantage of the computer-assisted measuring system CASSY.





**P2.3.4**  
**CONVERTING ELECTRICAL ENERGY INTO HEAT**

P2.3.4.1  
Converting electrical energy into heat energy - Measuring with a voltmeter and an ammeter

P2.3.4.2  
Converting electrical energy into heat energy - Measuring with the joule and wattmeter

P2.3.4.3  
Converting electrical energy into heat energy - Measuring with CASSY

P2.3.4.4  
Converting electrical energy into heat energy - Measuring with Mobile-CASSY

Converting electrical energy into heat energy - Measuring with the joule and wattmeter (P2.3.4.2\_a)

Cat. No.	Description	P2.3.4.1 (c)	P2.3.4.2 (a)	P2.3.4.3	P2.3.4.4
384 20	Lid for 250 ml Dewar vessel with heater	1			
386 48	Dewar vessel calorimeter, 250 ml	1			
524 005	Mobile-CASSY 2	1			1
524 0431	30-A-Box	1			
313 07	Hand-held stop watch I, mechanical	1			
664 103	Beaker, DURAN, 250 ml, squat	1			
665 755	Measuring cylinder, 250 ml, with plastic base	1			
521 35	Variable extra-low voltage transformer S	1	1	1	1
501 28	Connecting lead, 32 A, 50 cm, black	1			
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	1	1	1
388 02	Copper-block calorimeter		1	1	1
388 03	Aluminium-block calorimeter		1	1	1
388 04	Aluminium-block calorimeter, large		1	1	1
388 05	Thermometer for calorimeters, +15...35 °C/0.2 K		1		
388 06	Connecting leads, pair		1	1	1
531 831	Joule and wattmeter		1		
524 013	Sensor-CASSY 2			1	
524 220	CASSY Lab 2			1	
524 0673	NiCr-Ni adapter S, type K			1	
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K			1	
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1	

Just like mechanical energy, electrical energy can also be converted into heat. We can use e.g. a calorimeter vessel with a wire winding to which a voltage is connected to demonstrate this fact. When a current flows through the wire, Joule heat is generated and heats the calorimeter.

The supplied electrical energy

$$W(t) = U \cdot I \cdot t$$

is determined in the experiment P2.3.4.1 by measuring the voltage  $U$ , the current  $I$  and the time  $t$ , in the experiment P2.3.4.4 by measuring with the Mobile-CASSY and in the experiment P2.3.4.2 measured directly using the Joule and Wattmeter. This results in a change in the temperature of the calorimeter which corresponds to the specific heat capacity

$$Q(t) = m \cdot c \cdot (\vartheta(t) - \vartheta(0))$$

$c$ : specific heat capacity

$m$ : mass

$\vartheta(t)$ : temperature at time  $t$

To confirm the equivalence

$$Q(t) = W(t)$$

the two quantities are plotted together in a diagram.

In the experiment P2.3.4.3, the equivalence of electrical energy  $E_{el}$  and thermal energy  $E_{th}$  is established experimentally. The supplied electrical energy  $E_{el}$  is converted into heat  $E_{th}$  in the heating coil (or heating spiral). This leads to a temperature rise in the calorimeter (or water, in which the heating spiral is immersed). As the current  $I$  and the temperature  $\vartheta$  are measured simultaneously as functions of the time  $t$ , the constant voltage  $U$  being known, the two energy forms can be registered quantitatively in units of wattsecond (Ws) and Joule (J) so that their numerical equivalence can be demonstrated experimentally:  $E_{el} = E_{th}$ .

## P2.4.1

### LATENT HEAT AND VAPORIZATION HEAT

#### P2.4.1.1

Determining the specific vaporization heat of water

#### P2.4.1.2

Determining the specific latent heat of ice



Determining the specific vaporization heat of water (P2.4.1.1\_c)

Cat. No.	Description	P2.4.1.1 (c)	P2.4.1.2 (c)
386 48	Dewar vessel calorimeter, 250 ml	1	1
384 17	Water separator	1	
524 005	Mobile-CASSY 2	1	1
315 23	Single-pan suspension balance 610 Tara	1	1
303 28	Steam generator	1	
667 194	Silicone tubing, 7 mm diam., 1 m	1	
664 104	Beaker, DURAN, 400 ml, squat	1	1
300 02	Stand base, V-shaped, small	1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	
301 01	Leybold multiclamp	2	
666 555	Universal clamp, 0...80 mm	2	
303 25	Immersion heater		1
590 06	Plastic beaker		1

When a substance is heated at a constant pressure, its temperature generally increases. When that substance undergoes a phase transition, however, the temperature does not increase even when more heat is added, as the heat is required for the phase transition. As soon as the phase transition is complete, the temperature once more increases with the additional heat supplied. Thus, for example, the specific evaporation heat  $Q_V$  per unit of mass is required for evaporating water, and the specific melting heat  $Q_S$  per unit of mass is required for melting ice.

To determine the specific evaporation heat  $Q_V$  of water, pure steam is fed into the calorimeter in the experiment P2.4.1.1, in which cold water is heated to the mixing temperature  $\vartheta_m$ . The steam condenses to water and gives off heat in the process; the condensed water is cooled to the mixing temperature. The experiment measures the starting temperature  $\vartheta_2$  and the mass  $m_2$  of the cold water, the mixing temperature  $\vartheta_m$  and the total mass

$$m = m_1 + m_2$$

By comparing the amount of heat given off and absorbed, we can derive the equation

$$Q_V = \frac{m_1 \cdot c \cdot (\vartheta_m - \vartheta_1) + m_2 \cdot c \cdot (\vartheta_m - \vartheta_2)}{m_1}$$

$$\vartheta_1 \approx 100 \text{ °C}, c: \text{ specific heat capacity of water}$$

In the experiment P2.4.1.2, pure ice is filled in a calorimeter, where it cools water to the mixing temperature  $\vartheta_m$ , in order to determine the specific melting heat. The ice absorbs the melting heat and melts into water, which warms to the mixing temperature. Analogously to the experiment P2.4.1.1, we can say for the specific melting heat:

$$Q_S = \frac{m_1 \cdot c \cdot (\vartheta_m - \vartheta_1) + m_2 \cdot c \cdot (\vartheta_m - \vartheta_2)}{m_1}$$

$$\vartheta_1 = 0 \text{ °C}$$



## P2.4.2

### MEASURING VAPOR PRESSURE

#### P2.4.2.1

Recording the vapor-pressure curve of water - Pressures up to 1 bar

#### P2.4.2.2

Recording the vapor-pressure curve of water - Pressures up to 50 bar

Recording the vapor-pressure curve of water - Pressures up to 50 bar (P2.4.2.2)

Cat. No.	Description	P2.4.2.1	P2.4.2.2
664 315	Double-necked, round-bottom flask	1	
665 305	Adapter	1	
667 186	Vacuum rubber tubing, 8 mm diam.	1	
665 255	Three-way valve, T-shaped, ST nozzles	1	
378 031	Hose nozzle, DN 16 KF	1	
378 045ET2	Centring rings, DN 16 KF, set of 2	1	
378 050	Clamping ring, DN 10/16 KF	1	
378 701	Leybold high-vacuum grease	1	
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	
524 220	CASSY Lab 2	1	
524 065	Absolute pressure sensor S, 0...1500 hPa	1	
501 11	Extension cable, 15 pin	1	
688 808	Stand rod, 10 x 223 mm, with thread M6	1	
524 045	Temperature box, NiCr-Ni/NTC	1	
666 216	Temperature probe, NiCr-Ni, fast	1	
300 02	Stand base, V-shaped, small	1	
300 43	Stand rod, 75 cm, 12 mm diam.	1	
666 555	Universal clamp, 0...80 mm	1	
301 01	Leybold multiclamp	3	
302 68	Stand ring with stem, 13 cm diam.	1	
666 685	Wire gauze, 160 mm x 160 mm	1	
666 711	Butane gas burner	1	1
666 712ET3	Butane cartridge, 190 g, set of 3	1	1
667 614	Heat protective gloves	1	1
385 16	High-pressure steam boiler		1
664 109	Beaker, DURAN, 25 ml, squat		1
300 01	Stand base, V-shaped, large		1
667 6131	Safety goggles for wearing over glasses		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	

The vapour pressure  $p$  of a liquid-vapor mixture in a closed system depends on the temperature  $T$ . Above the critical temperature, the vapor pressure is undefined. The substance is gaseous and cannot be liquefied no matter how high the pressure. The increase in the vapor-pressure curve  $p(T)$  is determined by several factors, including the molar evaporation heat  $q_v$  of the substance:

$$T \cdot \frac{dp}{dT} = \frac{q_v}{v_1 - v_2} \quad (\text{Clausius-Clapeyron})$$

$T$ : absolute temperature

$v_1$ : molar volume of vapor

$v_2$ : molar volume of liquid

As we can generally ignore  $v_2$  and  $q_v$  hardly varies with  $T$ , we can derive a good approximation from the law of ideal gases:

$$\ln p = \ln p_0 - \frac{q_v}{R \cdot T}$$

In the experiment P2.4.2.1, the vapor pressure curve of water below the normal boiling point is recorded with the computer-assisted measuring system CASSY. The water is placed in a glass vessel, which was sealed beforehand while the water was boiling at standard pressure. The vapor pressure  $p$  is measured as a function of the temperature  $T$  when cooling and subsequently heating the system, respectively.

The high-pressure steam apparatus is used in the experiment P2.4.2.2 for measuring pressures of up to 50 bar. The vapor pressure can be read directly from the manometer of this device. A thermometer supplies the corresponding temperature. The measured values are recorded and evaluated manually point by point.

## P2.4.3

### CRITICAL TEMPERATURE

#### P2.4.3.1

Observing the phase transition between the liquid and the gas phase at the critical point




Observing the phase transition between the liquid and the gas phase at the critical point (P2.4.3.1\_c)

Cat. No.	Description	P2.4.3.1 (c)
371 401	Pressure chamber for demonstrating the critical temperature	1
450 60	Lamp housing with cable	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
460 20	Condenser with diaphragm holder	1
521 210	Transformer, 6/12 V	1
460 03	Lens in frame, f=100 mm	1
460 43	Small optical bench	1
300 01	Stand base, V-shaped, large	1
301 01	Leybold multiclamp	3
524 005	Mobile-CASSY 2	1
666 7681	Circulation thermostat SC 100-S5P	1
667 194	Silicone tubing, 7 mm diam., 1 m	2
675 3410	Water, pure, 5 l	2

The critical point of a real gas is defined by the critical pressure  $p_c$ , the critical density  $\rho_c$  and the critical temperature  $T_c$ . Below the critical temperature, the substance is gaseous for a sufficiently great molar volume – it is termed a vapor – and is liquid at a sufficiently small molar volume. Between these extremes, a liquid-vapor mix exists, in which the vapor component increases with the molar volume. As liquid and vapor have different densities, they are separated in a gravitational field. As the temperature rises, the density of the liquid decreases and that of the vapor increases, until finally at the critical temperature both densities have the value of the critical density. Liquid and vapor mix completely, and the phase boundary disappears. Above the critical temperature, the substance is gaseous, regardless of the molar volume.

The experiment P2.4.3.1 investigates the behavior of sulfur hexafluoride ( $\text{SF}_6$ ) close to the critical temperature. The critical temperature of this substance is  $T_c = 318.7 \text{ K}$  and the critical pressure is  $p_c = 37.6 \text{ bar}$ . The substance is enclosed in a pressure chamber designed so that hot water or steam can flow through the mantle. The dissolution of the phase boundary between liquid and gas while heating the substance, and its restoration during cooling, are observed in projection on the wall. As the system approaches the critical point, the substance scatters short-wave light particularly intensively; the entire contents of the pressure chamber appears red-brown. This critical opalescence is due to the variations in density, which increase significantly as the system approaches the critical point.

*Note:* The dissolution of the phase boundary during heating can be observed best when the pressure chamber is heated as slowly as possible using a circulation thermostat.



Contents of the pressure chamber: below, at the and above the critical temperature



## P2.5.1

### BROWNIAN MOTION OF MOLECULES


P2.5.1.1  
Brownian movement of smoke particles

Brownian movement of smoke particles (P2.5.1.1)

Cat. No.	Description	P2.5.1.1
MIK 738865	Microscope EduLED, mono 2	1
372 51	Smoke chamber	1
450 60	Lamp housing with cable	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
460 20	Condenser with diaphragm holder	1
521 210	Transformer, 6/12 V	1
300 02	Stand base, V-shaped, small	1

A particle which is suspended in a gas constantly executes a motion which changes in its speed and in all directions. *J. Perrin* first explained this molecular motion, discovered by *R. Brown*, which is caused by bombardment of the particles with the gas molecules. The smaller the particle is, the more noticeably it moves. The motion consists of a translational component and a rotation, which also constantly changes.

In the experiment P2.5.1.1, the motion of smoke particles in the air is observed using a microscope.



Schematic diagram of Brownian motion of molecules

### P2.5.2

#### GAS LAWS

##### P2.5.2.1

Pressure-dependency of the volume of a gas at a constant temperature (Boyle-Mariotte's law)

##### P2.5.2.2

Temperature-dependency of the volume of a gas at a constant pressure (Gay-Lussac's law)

##### P2.5.2.3


Temperature-dependency of the pressure of a gas at a constant volume (Amontons' law)



Temperature-dependency of the volume of a gas at a constant pressure (Gay-Lussac's law) (P2.5.2.2\_c)

Cat. No.	Description	P2.5.2.1 (a)	P2.5.2.1 (c)	P2.5.2.2-3 (c)
382 00	Gas thermometer	1	1	1
300 02	Stand base, V-shaped, small	1	1	1
300 42	Stand rod, 47 cm, 12 mm diam.	1	1	1
301 11	Clamp with jaw clamp	2	2	2
375 58	Hand vacuum pump	1	1	1
524 005	Mobile-CASSY 2		1	1
524 064	Pressure sensor S, ±2000 hPa		1	1
665 223ET10	Connector, T-shaped, 8 mm diam., 10 pieces		1	1
524 220	CASSY Lab 2		1*	1*
666 767	Hotplate, 1500 W, 180 mm diam.			1
664 103	Beaker, DURAN, 250 ml, squat			1

\* additionally recommended



Pressure-dependency of the volume of a gas at a constant temperature (P2.5.2.1)

The gas thermometer consists of a glass tube closed at the bottom end, in which a mercury stopper seals the captured air at the top. The volume of the air column is determined from its height and the cross-section of the glass tube. When the pressure at the open end is altered using a hand pump, this changes the pressure on the sealed side correspondingly. The temperature of the entire gas thermometer can be varied using a water bath.

In the experiment P2.5.2.1, the air column is maintained at a constant room temperature  $T$ . At an external pressure  $p_0$ , it has a volume of  $V_0$  bounded by the mercury stopper. The pressure  $p$  in the air column is reduced by evacuating air at the open end, and the increased volume  $V$  of the air column is determined for different pressure values  $p$ . The evaluation confirms the relationship

$$p \cdot V = p_0 \cdot V_0 \quad \text{for } T = \text{const. (Boyle-Mariotte's law)}$$

In the experiment P2.5.2.2, the gas thermometer is placed in a water bath of a specific temperature which is allowed to gradually cool. The open end is subject to the ambient air pressure, so that the pressure in the air column is constant. This experiment measures the volume  $V$  of the air column as a function of the temperature  $T$  of the water bath. The evaluation confirms the relationship

$$V \propto T \quad \text{for } p = \text{const. (Gay-Lussac's law)}$$

In the experiment P2.5.2.3, the pressure  $p$  in the air column is constantly reduced by evacuating the air at the open end so that the volume  $V$  of the air column also remains constant as the temperature drops. This experiment measures the pressure  $p$  of the air column as a function of the temperature  $T$  of the water bath. The evaluation confirms the relationship

$$p \propto T \quad \text{for } V = \text{const. (Amontons' law)}$$



### P2.5.3 SPECIFIC HEAT OF GASES

P2.5.3.1  
Determining the adiabatic exponent  $c_p/c_v$  of air after Rüchardt

P2.5.3.2  
Determining the adiabatic exponent  $c_p/c_v$  of various gases using the gas elastic resonance apparatus

Determining the adiabatic exponent  $c_p/c_v$  of air after Rüchardt (P2.5.3.1)

Cat. No.	Description	P2.5.3.1	P2.5.3.2
371 051	Oscillation tube with Mariotte's bottle	1	
313 07	Hand-held stop watch I, mechanical	1	
317 19	Aneroid barometer	1	
590 06	Plastic beaker	1	
675 3100	Vaseline, 50 g	1	
371 07	Gas elastic resonance apparatus		1
531 120	Multimeter LD analog 20		1
522 561	Function generator P		1
300 02	Stand base, V-shaped, small		1
660 980	Fine regulating valve for minican gas canisters		1
660 985	Minican pressurised gas canister, neon		1
660 999	Minican pressurised gas canister, carbon dioxide		1
665 255	Three-way valve, T-shaped, ST nozzles		1
667 194	Silicone tubing, 7 mm diam., 1 m		1
604 481	Rubber tubing, 1 m x 4 mm diam., DIN 12865		1
604 510	Tubing connector, 4...15 mm		1
500 422	Connecting lead 19 A, 50 cm, blue		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1

In the case of adiabatic changes in state, the pressure  $p$  and the volume  $V$  of a gas demonstrate the relationship

$$p \cdot V^\kappa = \text{const.}$$

whereby the adiabatic exponent is defined as

$$\kappa = \frac{c_p}{c_v}$$

i.e. the ratio of the specific heat capacities  $c_p$  and  $c_v$  of the respective gas.

The experiment P2.5.3.1 determines the adiabatic exponent of air from the oscillation period of a ball which caps and seals a gas volume in a glass tube, whereby the oscillation of the ball around the equilibrium position causes adiabatic changes in the state of the gas. In the equilibrium position, the force of gravity and the opposing force resulting from the pressure of the enclosed gas are equal. A deflection from the equilibrium position by  $\Delta x$  causes the pressure to change by

$$\Delta p = -\kappa \cdot p \cdot \frac{A \cdot \Delta x}{V}$$

$A$ : cross-section of riser tube

which returns the ball to the equilibrium position. The ball thus oscillates with the frequency

$$f_0 = \frac{1}{2\pi} \cdot \sqrt{\frac{\kappa \cdot p \cdot A^2}{m \cdot V}}$$

around its equilibrium position.

In the experiment P2.5.3.2, the adiabatic exponent is determined using the gas elastic resonance apparatus. Here, the air column is sealed by a magnetic piston which is excited to forced oscillations by means of an alternating electromagnetic field. The aim of the experiment is to find the characteristic frequency  $f_0$  of the system, i.e. the frequency at which the piston oscillates with maximum amplitude. Other gases, such as carbon dioxide and nitrogen, can alternatively be used in this experiment.

### P2.5.4

#### REAL GASES

##### P2.5.4.1

##### Joule-Thomson effect



Joule-Thomson effect (P2.5.4.1)


Cat. No.	Description	P2.5.4.1
371 56	Joule-Thomson apparatus	1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
524 045	Temperature box, NiCr-Ni/NTC	1
666 216	Temperature probe, NiCr-Ni, fast	2
524 064	Pressure sensor S, ±2000 hPa	1
667 186	Vacuum rubber tubing, 8 mm diam.	1
667 184	Tubing clamps, 10...16 mm, set of 10	1
661 017	Pressure reducing valve for carbon dioxide, helium, argon	1
661 018	Pressure reducing valve for nitrogen	1
664 569	Open jaw spanner for gas canisters	1
661 020	Trolley for two gas canisters	1
661 012	Compressed gas bottle, carbon dioxide	1
661 013	Compressed gas bottle, hydrogen	1
300 01	Stand base, V-shaped, large	1
300 42	Stand rod, 47 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
666 555	Universal clamp, 0...80 mm	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The Joule-Thomson effect occurs when a real gas or gas mixture undergoes a change in temperature by reduction of the pressure  $p$ . The strength and direction of the change in temperature is described by the Joule-Thomson coefficient  $\mu$ :

$$\mu_{JT} = \left( \frac{\partial T}{\partial p} \right)_H$$


$H$  : enthalpy = const.

In the experiment P2.5.4.1 the flowing gas is throttled through a membrane and expanded behind it. The volume occupied by the gas can increase behind the obstacle, so the average particle distance is increased. This causes a change in the temperature of the gas. The pressure and the temperature difference before and after the membrane is measured in this experiment.



Temperature change against pressure





**P2.6.1**  
**HOT-AIR ENGINE:  
QUALITATIVE EXPERIMENTS**

P2.6.1.1  
Operating a hot-air engine as a thermal engine

P2.6.1.5  
Operating the hot-air engine P as a thermal engine

Operating a hot-air engine as a thermal engine (P2.6.1.1)

Cat. No.	Description	P2.6.1.1	P2.6.1.5
388 182	Hot-air engine	1	
562 11	U-core with yoke	1	
562 121	Clamping device with spring clip	1	
562 21	Mains coil, 500 turns	1	
562 18	Coil, 50 turns, extra-low voltage	1	
501 33	Connecting lead, 32 A, 100 cm, black	2	
388 181	Immersion pump	1*	
521 231	Low-voltage power supply, 3/6/9/12 V	1*	
667 194	Silicone tubing, 7 mm diam., 1 m	2*	
604 313	Wide-mouthed can, 10 l	1*	
388 176	Hot-air engine P		1
303 22	Alcohol burner, metal		1
524 005	Mobile-CASSY 2		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1

\* additionally recommended




Diagram illustrating the principle of operation of a hot-air engine as a heat engine

The hot-air engine (invented by *R. Stirling*, 1816) is the oldest thermal engine, along with the steam engine. In greatly simplified terms, its thermodynamic cycle consists of an isothermic compression at low temperature, an isochoric application of heat, an isothermic expansion at high temperature and an isochoric emission of heat. The displacement piston and the working piston are connected to a crankshaft via tie rods, whereby the displacement piston leads the working piston by 90°. When the working piston is at top dead center (a), the displacement piston is moving downwards, displacing the air into the electrically heated zone of the cylinder. Here, the air is heated, expands and forces the working piston downward (b). The mechanical work is transferred to the flywheel. When the working piston is at bottom dead center (c), the displacement piston is moving upwards, displacing the air into the water-cooled zone of the cylinder. The air cools and is compressed by the working cylinder (d). The flywheel delivers the mechanical work required to execute this process.

The experiment P2.6.1.1 qualitatively investigates the operation of the hot-air engine as a thermal engine. Mechanical power is derived from the engine by braking at the brake hub. The voltage of the heating filament is varied in order to demonstrate the relationship between the thermal power supplied and the mechanical power removed from the system. The no-load speed of the motor for each case is used as a measure of the mechanical power produced in the system.

The experiment P2.6.1.5 qualitatively investigates the operation of the hot-air engine P as a thermal engine. The hot-air engine P is a transparent model to show the basics of energy transformation. The hot-air engine P is driven by an alcohol burner and transforms via a small generator the mechanical energy into a voltage which can be measured.

## P2.6.1


### HOT-AIR ENGINE: QUALITATIVE EXPERIMENTS

#### P2.6.1.3

Operating the hot-air engine as a heat pump and a refrigerator

#### P2.6.1.6

Operating the hot-air engine P as a heat pump and a refrigerator



Operating the hot-air engine P as a heat pump and a refrigerator (P2.6.1.6)


Cat. No.	Description	P2.6.1.3	P2.6.1.6
388 182	Hot-air engine	1	
388 19	Thermometer for hot-air engine	1	
347 35	Experiment motor	1	
347 36	Control unit for experiment motor	1	
388 181	Immersion pump	1*	
521 231	Low-voltage power supply, 3/6/9/12 V	1*	1
667 194	Silicone tubing, 7 mm diam., 1 m	2*	
604 313	Wide-mouthed can, 10 l	1*	
388 176	Hot-air engine P		1
524 005	Mobile-CASSY 2		1
524 0673	NiCr-Ni adapter S, type K		1
666 1261	Temperature probe, NiCr-Ni, fast, type K		2
300 11	Saddle base		1
301 01	Leybold multiclamp		1
300 41	Stand rod, 25 cm, 12 mm Ø		1
590 13	Stand rod with bore holes		1
340 89ET5	Coupling plugs, 4 mm, set of 5		1
501 861	Crocodile-clips, polished, set of 6		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1

\* additionally recommended

Depending on the direction of rotation of the crankshaft, the hot-air engine operates as either a heat pump or a refrigerating machine when its flywheel is externally driven. When the displacement piston is moving upwards while the working piston is at bottom dead center, it displaces the air in the top part of the cylinder. The air is then compressed by the working piston and transfers its heat to the cylinder head, i.e. the hot-air motor operates as a heat pump. When run in the opposite direction, the working piston causes the air to expand when it is in the top part of the cylinder, so that the air draws heat from the cylinder head; in this case the hot-air engine operates as a refrigerating machine.

The experiment P2.6.1.3 qualitatively investigates the operation of the hot-air engine as a heat pump and a refrigerating machine. In order to demonstrate the relationship between the externally supplied mechanical power and the heating or refrigerating power, respectively, the speed of the electric motor is varied and the change in temperature observed.

In the experiment P2.6.1.6 the operation of the hot-air engine P as a heat pump and a refrigerating machine is shown. By varying the speed of the electric motor, i.e. the mechanical power supplied to the hot-air engine P, the impact on heating or refrigerating power is observed.



## P2.6.2 HOT-AIR ENGINE: QUANTITATIVE EXPERIMENTS

P2.6.2.1  
Frictional losses in the hot-air engine (calorific determination)

P2.6.2.2  
Determining the efficiency of the hot-air engine as a heat engine

P2.6.2.3  
Determining the efficiency of the hot-air engine as a refrigerator

Frictional losses in the hot-air engine (calorific determination) (P2.6.2.1)

Cat. No.	Description	P2.6.2.1	P2.6.2.2	P2.6.2.3
388 182	Hot-air engine	1	1	1
388 221	Accessories for hot air engine for power measurement	1	1	1
347 35	Experiment motor	1		1
347 36	Control unit for experiment motor	1		1
575 471	Counter S	1	1	1
337 46	Fork-type light barrier	1	1	1
501 16	Multi-core cable, 6-pole, 1.5 m	1	1	1
313 17	Hand-held stop-watch II, mechanical	1	1	1
382 35	Thermometer, -10...+50 °C/0.1 K	1	1	1
300 02	Stand base, V-shaped, small	1	2	1
300 41	Stand rod, 25 cm, 12 mm Ø	1	1	1
590 06	Plastic beaker	1	1	1
388 181	Immersion pump	1*	1*	1*
521 231	Low-voltage power supply, 3/6/9/12 V	1*	1*	1*
667 194	Silicone tubing, 7 mm diam., 1 m	2*	2*	2*
604 313	Wide-mouthed can, 10 l	1*	1*	1*
562 11	U-core with yoke		1	
562 121	Clamping device with spring clip		1	
562 21	Mains coil, 500 turns		1	
562 18	Coil, 50 turns, extra-low voltage		1	
531 120	Multimeter LDanalog 20		1	1
531 130	Multimeter LDanalog 30		1	1
314 141	Precision dynamometer, 1 N		1	
300 42	Stand rod, 47 cm, 12 mm diam.		1	
300 51	Stand rod, right-angled		1	
301 01	Leybold multiclamp		2	
342 61	Weights, 50 g, set of 12		1	
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair		1	1

Cat. No.	Description	P2.6.2.1	P2.6.2.2	P2.6.2.3
501 33	Connecting lead, 32 A, 100 cm, black		3	3
521 35	Variable extra-low voltage transformer S			1

\* additionally recommended

When the hot-air engine is operated as a heat engine, each engine cycle withdraws the amount of heat  $Q_1$  from reservoir 1, generates the mechanical work  $W$  and transfers the difference  $Q_2 = Q_1 - W$  to reservoir 2. The hot-air engine can also be made to function as a refrigerating machine while operated in the same rotational direction by externally applying the mechanical work  $W$ . In both cases, the work  $W_F$  converted into heat in each cycle through the friction of the piston in the cylinder must be taken into consideration.

In order to determine the work of friction  $W_F$  in the experiment P2.6.2.1, the temperature increase  $\Delta T_F$  in the cooling water is measured while the hot-air engine is driven using an electric motor and the cylinder head is open.

The experiment P2.6.2.2 determines the efficiency

$$\eta = \frac{W}{W + Q_2}$$

of the hot-air engine as a heat engine. The mechanical work  $W$  exerted on the axle in each cycle can be calculated using the external torque  $N$  of a dynamometrical brake which brakes the hot-air engine to a speed  $f$ . The amount of heat  $Q_2$  given off corresponds to a temperature increase  $\Delta T$  in the cooling water.

The experiment P2.6.2.3 determines the efficiency

$$\eta = \frac{Q_2}{Q_1 - Q_2}$$

of the hot-air engine as a refrigerating machine. Here, the hot-air engine with closed cylinder head is driven using an electric motor and  $Q_1$  is determined as the electrical heating energy required to maintain the cylinder head at the ambient temperature.

## P2.6.2

### HOT-AIR ENGINE: QUANTITATIVE EXPERIMENTS

#### P2.6.2.4

pV diagram of the hot-air engine  
as a heat engine - Recording and  
evaluating with CASSY



pV diagram of the hot-air engine as a heat engine - Recording and evaluating with CASSY (P2.6.2.4)

Cat. No.	Description	P2.6.2.4
388 182	Hot-air engine	1
562 11	U-core with yoke	1
562 121	Clamping device with spring clip	1
562 21	Mains coil, 500 turns	1
562 18	Coil, 50 turns, extra-low voltage	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>
524 220	CASSY Lab 2	1
524 082	Rotary motion sensor S	1
524 064	Pressure sensor S, $\pm 2000$ hPa	1
309 48ET2	Fishing line, set of 2	1
352 08ET2	Helical spring, 25 N/m, set of 2	1
501 33	Connecting lead, 32 A, 100 cm, black	2
388 181	Immersion pump	1*
521 231	Low-voltage power supply, 3/6/9/12 V	1*
667 194	Silicone tubing, 7 mm diam., 1 m	2*
604 313	Wide-mouthed can, 10 l	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

\* additionally recommended

Thermodynamic cycles are often described as a closed curve in a  $pV$  diagram ( $p$ : pressure,  $V$ : volume). The work added to or withdrawn from the system (depending on the direction of rotation) corresponds to the area enclosed by the curve.

In the experiment P2.6.2.4, the  $pV$  diagram of the hot air engine as a heat engine is recorded using the computer-assisted measured value recording system CASSY. The pressure sensor measures the pressure  $p$  in the cylinder and a displacement sensor measures the position  $s$ , from which the volume is calculated, as a function of the time  $t$ . The measured values are displayed on the screen directly in a  $pV$  diagram. In the further evaluation, the mechanical work performed as piston friction per cycle


$$W = -\int p \cdot dV$$

and from this the mechanical power

$$P = W \cdot f$$

$f$ : no-load speed

are calculated and plotted in a graph as a function of the no-load speed.



pV diagram of the hot air engine

P2.6.3

HEAT PUMP

P2.6.3.1

Determining the efficiency of the heat pump as a function of the temperature differential

P2.6.3.2

Investigating the function of the expansion valve of the heat pump

P2.6.3.3

Analyzing the cyclical process of the heat pump with the Mollier diagram



Determining the efficiency of the heat pump as a function of the temperature differential (P2.6.3.1\_a)

Cat. No.	Description	P2.6.3.1 (a)	P2.6.3.2 (a)	P2.6.3.3 (a)
389 521	Heat pump	1	1	1
531 831	Joule and wattmeter	1		1
666 209	Digital thermometer with 4 inputs	1	1	1
666 193	Temperature probe, NiCr-Ni, 1.5 mm	2	2	3
313 12	Digital stop-watch	1	1	1
729 769	RS-232 cable, 9-pole	1*	1*	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

\* additionally recommended


The heat pump extracts heat from a reservoir with the temperature  $T_1$  through vaporization of a coolant and transfers heat to a reservoir with the temperature  $T_2$  through condensation of the coolant. In the process, compression in the compressor (a-b) greatly heats the gaseous coolant. It condenses in the liquefier (c-d) and gives up the released condensation heat  $\Delta Q_2$  to the reservoir  $T_2$ . The liquefied coolant is filtered and fed to the expansion valve (e-f) free of bubbles. This regulates the supply of coolant to the vaporizer (g-h). In the vaporizer, the coolant once again becomes a gas, withdrawing the necessary evaporation heat  $\Delta Q_1$  from the reservoir  $T_1$ .

The aim of the experiment P2.6.3.1 is to determine the efficiency

$$\varepsilon = \frac{\Delta Q_2}{\Delta W}$$

of the heat pump as a function of the temperature differential  $\Delta T = T_2 - T_1$ . The heat quantity  $\Delta Q_2$  released is determined from the heating of water reservoir  $T_2$ , while the applied electrical energy  $\Delta W$  is measured using the joule and wattmeter.

In the experiment P2.6.3.2, the temperatures  $T_f$  and  $T_h$  are recorded at the outputs of the expansion valve and the vaporizer. If the difference between these two temperatures falls below a specific limit value, the expansion valve chokes off the supply of coolant to the vaporizer. This ensures that the coolant in the vaporizer is always vaporized completely. In the experiment P2.6.3.3, a Mollier diagram, in which the pressure  $p$  is graphed as a function of the specific enthalpy  $h$  of the coolant, is used to trace the energy transformations of the heat pump. The pressures  $p_1$  and  $p_2$  in the vaporizer and liquefier, as well as the temperatures  $T_a$ ,  $T_b$ ,  $T_c$  and  $T_f$  of the coolant are used to determine the corresponding enthalpy values  $h_a$ ,  $h_b$ ,  $h_c$  and  $h_f$ . This experiment also measures the heat quantities  $\Delta Q_2$  and  $\Delta Q_1$  released and absorbed per unit of time. This in turn is used to determine the amount of coolant  $\Delta m$  circulated per unit of time.



Heat pump pT (389 521) with schematic diagram of all functional components

P2.2.2.1  
Determining the efficiency of a solar collector  
as a function of the throughput volume of water



# P3 ELECTRICITY



P3.1	ELECTROSTATICS	89
P3.2	FUNDAMENTALS OF ELECTRICITY	104
P3.3	MAGNETOSTATICS	111
P3.4	ELECTROMAGNETIC INDUCTION	115
P3.5	ELECTRICAL MACHINES	122
P3.6	DC AND AC CIRCUITS	126
P3.7	ELECTROMAGNETIC OSCILLATIONS AND WAVES	134
P3.8	FREE CHARGE CARRIERS IN A VACUUM	140
P3.9	ELECTRICAL CONDUCTION IN GASES	145

# P3 ELECTRICITY



## P3.1 ELECTROSTATICS

P3.1.1 Basic experiments on electrostatics	89-90
P3.1.2 Coulomb's law	91-93
P3.1.3 Field lines and equipotential lines	94-96
P3.1.4 Effects of force in an electric field	97-98
P3.1.5 Charge distributions on electrical conductors	99
P3.1.6 Definition of capacitance	100
P3.1.7 Plate capacitor	101-103

## P3.2 FUNDAMENTALS OF ELECTRICITY

P3.2.1 Charge transfer with drops of water	104
P3.2.2 Ohm's law	105
P3.2.3 Kirchhoff's laws	106-107
P3.2.4 Circuits with electrical measuring instruments	108
P3.2.5 Conducting electricity by means of electrolysis	109
P3.2.6 Experiments on electrochemistry	110

## P3.3 MAGNETOSTATICS

P3.3.1 Basic experiments on magnetostatics	111
P3.3.2 Magnetic dipole moment	112
P3.3.3 Effects of force in a magnetic field	113
P3.3.4 Biot-Savart's law	114

## P3.4 ELECTROMAGNETIC INDUCTION

P3.4.1 Voltage impulse	115
P3.4.2 Induction in a moving conductor loop	116
P3.4.3 Induction by means of a variable magnetic field	117
P3.4.4 Eddy currents	118
P3.4.5 Transformer	119-120
P3.4.6 Measuring the earth's magnetic field	121

## P3.5 ELECTRICAL MACHINES

P3.5.1 Basic experiments on electrical machines	122
P3.5.2 Electric generators	123
P3.5.3 Electric motors	124
P3.5.4 Three-phase machines	125

## P3.6 DC AND AC CIRCUITS

P3.6.1 Circuit with capacitor	126
P3.6.2 Circuit with coil	127
P3.6.3 Impedances	128
P3.6.4 Measuring-bridge circuits	129
P3.6.5 Measuring AC voltages and AC currents	130
P3.6.6 Electrical work and power	131-132
P3.6.7 Electromechanical devices	133

## P3.7 ELECTROMAGNETIC OSCILLATIONS AND WAVES

P3.7.1 Electromagnetic oscillator circuit	134
P3.7.2 Decimeter-range waves	135
P3.7.3 Propagation of decimeter-range waves along lines	136
P3.7.4 Microwaves	137
P3.7.5 Propagation of microwaves along lines	138
P3.7.6 Directional characteristic of dipole radiation	139

## P3.8 FREE CHARGE CARRIERS IN A VACUUM

P3.8.1 Tube diode	140
P3.8.2 Tube triode	141
P3.8.3 Maltese-cross tube	142
P3.8.4 Perrin tube	143
P3.8.5 Thomson tube	144

## P3.9 ELECTRICAL CONDUCTION IN GASES

P3.9.1 Spontaneous and non-spontaneous discharge	145
P3.9.2 Gas discharge at reduced pressure	146
P3.9.3 Cathode rays and canal rays	147



P3.1.1

BASIC EXPERIMENTS  
ON ELECTROSTATICS

P3.1.1.1

Basic electrostatics experiments  
with the field electrometer




Basic electrostatics experiments with the field electrometer (P3.1.1.1)

Cat. No.	Description	P3.1.1.1
540 10	Field electrometer	1
540 11	Electrostatics set 1	1
540 12	Electrostatics set 2	1
300 02	Stand base, V-shaped, small	1
300 43	Stand rod, 75 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
501 861	Crocodile-clips, polished, set of 6	1
501 20	Connecting lead, 32 A, 25 cm, red	1

The field electrometer is a classic apparatus for demonstrating electrical charges. Its metallized pointer, mounted on needle bearings, is conductively connected to a fixed metal support. When an electrical charge is transferred to the metal support via a pluggable metal plate or a Faraday's cup, part of the charge flows onto the pointer. The pointer is thus repelled, indicating the charge.

In the experiment P3.1.1.1, the electrical charges are generated by rubbing two materials together (more precisely, by intensive contact followed by separation), and demonstrated using the field electrometer. This experiment proves that charges can be transferred between different bodies. Additional topics include charging of an electrometer via induction, screening induction via a metal screen and discharge in ionized air.



### P3.1.1

#### BASIC EXPERIMENTS ON ELECTROSTATICS

##### P3.1.1.2


Basic electrostatics experiments  
with the electrometer amplifier



Basic electrostatics experiments with the electrometer amplifier (P3.1.1.2\_a)

Cat. No.	Description	P3.1.1.2 (a)
532 14	Electrometer amplifier	1
562 791	Plug-in power supply, 12 V AC	1
578 25	Capacitor, 1 nF, STE 2/19	1
578 10	Capacitor, 10 nF, STE 2/19	1
532 16	Connecting rod	1
531 120	Multimeter LDanalog 20	1
541 00	Friction rods, PVC and acrylic	1
541 21	Leather	1
686 63	Polyethylene friction foils, set of 10	1
546 12	Faraday's cup	1
590 011	Clamping plug	1
542 51	Induction plate, 8 cm x 4 cm	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
500 424	Connecting lead 19 A, 50 cm, black	1
666 711	Butane gas burner	1*
666 712ET3	Butane cartridge, 190 g, set of 3	1*

\* additionally recommended



Measuring charges with the electrometer amplifier

The electrometer amplifier is an impedance converter with an extremely high-ohm voltage input ( $\geq 10^{13} \Omega$ ) and a low-ohm voltage output ( $\leq 1 \Omega$ ). By means of capacitive connection of the input and using a Faraday's cup to collect charges, this device is ideal for measuring extremely small charges. Experiments on contact and friction electricity can be conducted with a high degree of reliability.

The experiment P3.1.1.2 investigates how charges can be separated through rubbing two materials together. It shows that one of the materials carries positive charges, and the other negative charges, and that the absolute values of the charges are equal. If we measure the charges of both materials at the same time, they cancel each other out. The sign of the charge of a material does not depend on the material alone, but also on the properties of the other material.

## P3.1.2

### COULOMB'S LAW

#### P3.1.2.1

Confirming Coulomb's law -  
Measuring with the torsion  
balance, Schürholz design



Confirming Coulomb's law - Measuring with the torsion balance, Schürholz design (P3.1.2.1)

Cat. No.	Description	P3.1.2.1
516 01	Torsion balance, Schürholz design	1
516 20	Electrostatics, accessories	1
516 04	Scale on stand	1
521 721	High-voltage power supply, 25 kV	1
501 051	Cable for high voltages, 1.5 m	1
590 13	Stand rod with bore holes	1
300 11	Saddle base	1
532 14	Electrometer amplifier	1
562 791	Plug-in power supply, 12 V AC	1
578 25	Capacitor, 1 nF, STE 2/19	1
578 10	Capacitor, 10 nF, STE 2/19	1
531 120	Multimeter LDanalog 20	1
546 12	Faraday's cup	1
590 011	Clamping plug	1
532 16	Connecting rod	1
471 830	He-Ne Laser, linearly polarised	1
300 02	Stand base, V-shaped, small	2
300 42	Stand rod, 47 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
313 07	Hand-held stop watch I, mechanical	1
311 02	Metal rule, 1 m	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1
500 414	Connecting lead, 19 A, 25 cm, black	1
500 424	Connecting lead 19 A, 50 cm, black	1
500 444	Connecting lead 19 A, 100 cm, black	2
501 43	Connecting lead, 32 A, 200 cm, yellow/green	1

According to Coulomb's law, the force acting between two point-shaped electrical charges  $Q_1$  and  $Q_2$  at a distance  $r$  from each other can be determined using the formula

$$F = \frac{1}{4\pi\epsilon_0} \cdot \frac{Q_1 \cdot Q_2}{r^2}$$

$$\text{where } \epsilon_0 = 8.85 \cdot 10^{-12} \frac{\text{As}}{\text{Vm}} \text{ (permittivity)}$$

The same force acts between two charged fields when the distance  $r$  between the sphere midpoints is significantly greater than the sphere diameter, so that the uniform charge distributions of the spheres is undisturbed. In other words, the spheres in this geometry may be treated as points.

In the experiment P3.1.2.1, the coulomb force between two charged spheres is measured using the torsion balance. The heart of this extremely sensitive measuring instrument is a rotating body elastically mounted between two torsion wires, to which one of the two spheres is attached. When the second sphere is brought into close proximity with the first, the force acting between the two charged spheres produces torsion of the wires; this can be indicated and measured using a light pointer. The balance must be calibrated if the force is to be measured in absolute terms. The coulomb force is measured as a function of the distance  $r$ . For this purpose, the second sphere, mounted on a stand, is brought close to the first one. Then, at a fixed distance, the charge of one sphere is reduced by half. The measurement can also be carried out using spheres with opposing charges. The charges are measured using an electrometer amplifier connected as a coulomb meter. The aim of the evaluation is to verify the proportionalities

$$F \propto \frac{1}{r^2} \text{ and } F \propto Q_1 \cdot Q_2$$

and to calculate the permittivity  $\epsilon_0$ .

### P3.1.2

#### COULOMB'S LAW

##### P3.1.2.2

Confirming Coulomb's law -  
Measuring with the force sensor



Confirming Coulomb's law - Measuring with the force sensor (P3.1.2.2\_b)

Cat. No.	Description	P3.1.2.2 (b)
314 263	Bodies for electric charge, set	1
337 00	Trolley	1
460 82	Precision metal rail, 50 cm	1
460 95ET5	Clamp riders, set of 5	1
524 005	Mobile-CASSY 2	1
524 060	Force sensor S, ±1N	1
521 721	High-voltage power supply, 25 kV	1
501 051	Cable for high voltages, 1.5 m	1
590 13	Stand rod with bore holes	1
300 11	Saddle base	1
590 02ET2	Clip plugs, small, set of 2	1
532 14	Electrometer amplifier	1
562 791	Plug-in power supply, 12 V AC	1
578 25	Capacitor, 1 nF, STE 2/19	1
578 10	Capacitor, 10 nF, STE 2/19	1
546 12	Faraday's cup	1
590 011	Clamping plug	1
532 16	Connecting rod	1
300 02	Stand base, V-shaped, small	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
301 01	Leybold multiclamp	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1
500 414	Connecting lead, 19 A, 25 cm, black	1
500 424	Connecting lead 19 A, 50 cm, black	1
500 444	Connecting lead 19 A, 100 cm, black	1
501 43	Connecting lead, 32 A, 200 cm, yellow/green	1

As an alternative to measuring with the torsion balance, the coulomb force between two spheres can also be determined using the force sensor. This device consists of two bending elements connected in parallel with four strain gauges in a bridge configuration; their electrical resistance changes when a load is applied. The change in resistance is proportional to the force acting on the instrument.

In the experiment P3.1.2.2, the force sensor is connected to a measuring instrument, which displays the measured force directly. No calibration is necessary. The coulomb force is measured as a function of the distance  $r$  between the sphere midpoints, the charge  $Q_1$  of the first sphere and the charge  $Q_2$  of the second sphere. The charges of the spheres are measured using an electrometer amplifier connected as a coulomb meter. The aim of the evaluation is to verify the proportionalities

$$F \propto \frac{1}{r^2}, F \propto Q_1 \text{ and } F \propto Q_2$$

and to calculate the permittivity  $\epsilon_0$ .



## P3.1.2

### COULOMB'S LAW

#### P3.1.2.3

Confirming Coulomb's law -  
Recording and evaluating with CASSY

Confirming Coulomb's law - Recording and evaluating with CASSY (P3.1.2.3)

Cat. No.	Description	P3.1.2.3
314 263	Bodies for electric charge, set	1
337 00	Trolley	1
460 82	Precision metal rail, 50 cm	1
460 95ET5	Clamp riders, set of 5	1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
524 060	Force sensor S, ±1N	1
524 082	Rotary motion sensor S	1
521 721	High-voltage power supply, 25 kV	1
501 051	Cable for high voltages, 1.5 m	1
590 13	Stand rod with bore holes	1
300 11	Saddle base	1
590 02ET2	Clip plugs, small, set of 2	1
532 14	Electrometer amplifier	1
562 791	Plug-in power supply, 12 V AC	1
578 25	Capacitor, 1 nF, STE 2/19	1
578 10	Capacitor, 10 nF, STE 2/19	1
531 120	Multimeter LDanalog 20	1
546 12	Faraday's cup	1
590 011	Clamping plug	1
532 16	Connecting rod	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
300 02	Stand base, V-shaped, small	1
301 01	Leybold multiclamp	1
337 04	Driving weights, set	1
301 07	Simple bench clamp	1

Cat. No.	Description	P3.1.2.3
309 48ET2	Fishing line, set of 2	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1
500 414	Connecting lead, 19 A, 25 cm, black	1
500 424	Connecting lead 19 A, 50 cm, black	1
501 43	Connecting lead, 32 A, 200 cm, yellow/green	2
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

For computer-assisted measuring of the coulomb force between two charged spheres, we can also connect the force sensor to the CASSY interface. A displacement sensor (Rotary motion sensor S) is additionally required to measure the distance between the charged spheres.

The experiment P3.1.2.3 utilizes the software CASSY Lab to record the values and evaluate them. The coulomb force is measured for different charges  $Q_1$  and  $Q_2$  as a function of the distance  $r$ . The charges of the spheres are measured using an electrometer amplifier connected as a coulomb meter. The aim of the evaluation is to verify the proportionality

$$F \propto \frac{1}{r^2}$$

and to calculate of the permittivity  $\epsilon_0$ .

### P3.1.3

#### FIELD LINES AND EQUIPOTENTIAL LINES

##### P3.1.3.1

##### Displaying lines of electric flux




Displaying lines of electric flux (P3.1.3.1)

Cat. No.	Description	P3.1.3.1
541 06	Electric field lines equipment set	1
501 051	Cable for high voltages, 1.5 m	2
521 70	High-voltage power supply, 10 kV	1
MIK 74702	BMS EcoCam 5.5 mega pixel, WiFi	1*
	additionally required: PC with Windows 7 or higher	1

\* additionally recommended

The space which surrounds an electric charge is in a state which we describe as an electric field. The electric field is also present even when it cannot be demonstrated through a force acting on a sample charge. A field is best described in terms of lines of electric flux, which follow the direction of electric field strength. The orientation of these lines of electric flux is determined by the spatial arrangement of the charges generating the field.


In the experiment P3.1.3.1, small particles in an oil-filled cuvette are used to illustrate the lines of electric flux. The particles align themselves in the electric field to form chains which run along the lines of electric flux. Four different pairs of electrodes are provided to enable electric fields with different spatial distributions to be generated; these electrode pairs are mounted beneath the cuvette, and connected to a high voltage source of up to 10 kV. The resulting patterns can be interpreted as the cross-sections of two spheres, one sphere in front of a plate, a plate capacitor and a spherical capacitor.



Equipment set E-field lines (541 06)

P3.1.3  
FIELD LINES AND  
EQUIPOTENTIAL LINES

P3.1.3.2  
Displaying the equipotential lines  
of electric fields




Displaying the equipotential lines of electric fields (P3.1.3.2)

Cat. No.	Description	P3.1.3.2
545 09	Electrolytic tank	1
501 861	Crocodile-clips, polished, set of 6	1
521 231	Low-voltage power supply, 3/6/9/12 V	1
531 120	Multimeter LDanalog 20	1
686 64ET5	Metal needles, set of 5	1
590 011	Clamping plug	1
590 13	Stand rod with bore holes	1
300 41	Stand rod, 25 cm, 12 mm Ø	1
301 01	Leybold multiclamp	1
300 11	Saddle base	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2

In a two-dimensional cross-section of an electric field, points of equal potential form a line. The direction of these isoelectric lines, just like the lines of electric flux, are determined by the spatial arrangement of the charges generating the field.

The experiment P3.1.3.2 measures the isoelectric lines for bodies with different charges. To do this, a voltage is applied to a pair of electrodes placed in an electrolytic tray filled with distilled water. An AC voltage is used to avoid potential shifts due to electrolysis at the electrodes. A voltmeter measures the potential difference between the 0 V electrode and a steel needle immersed in the water. To display the isoelectric lines, the points of equal potential difference are localized and drawn on graph paper. In this way, it is possible to observe and study two-dimensional sections through the electric field in a plate capacitor, a Faraday's cup, a dipole, an image charge and a slight curve.




Measurement example: equipotential lines around a needle tip

### P3.1.3

#### FIELD LINES AND EQUIPOTENTIAL LINES

P3.1.3.3  
Measuring the potential inside a plate capacitor

P3.1.3.4  
Measuring the potential around a charged sphere



Measuring the potential around a charged sphere (P3.1.3.4\_b)

Cat. No.	Description	P3.1.3.3 (b)	P3.1.3.4 (b)
524 080	Electric field meter S	1	1
540 540	Accessories for electric field meter S	1	1
524 005	<b>Mobile-CASSY 2</b>	1	1
311 02	Metal rule, 1 m	1	1
521 70	High-voltage power supply, 10 kV	1	1
460 317	Optical bench, S1 profile, 0.5 m	1	
460 312	Clamp rider with clamp, 45/35	2	
300 11	Saddle base	2	3
300 41	Stand rod, 25 cm, 12 mm Ø	2	
301 01	Leybold multiclamp	1	
500 600	Safety connection lead, 10 cm, yellow/green	1	1
500 621	Safety connecting lead, 50 cm, red	1	1
500 622	Safety connecting lead, 50 cm, blue	1	
500 641	Safety connecting lead, 100 cm, red	1	1
500 642	Safety connecting lead, 100 cm, blue	1	1
667 193	PVC tubing, 7 mm diam., 1 m	1	1
666 716	Valve for gas cartridge	1	1
666 715	Cartridge	1	1
543 021	Sphere on insulated stand rod		1
500 95	Safety adapter sockets, red, set of 6		1

Using a flame probe, the electric potential around a charged object can be investigated in all three dimensions and the equipotential surfaces can be determined.

In the experiment P3.1.3.3, the electric potential of a plate capacitor is investigated. The equipotential surfaces parallel to the capacitor plates are identified by measuring the electrical potential at different positions but with constant distance to the capacitor plates. In addition, the dependence of the variation of the electric potential on the distance to the capacitor plates is determined and used to calculate the electric field strength.

The aim of the experiment P3.1.3.4 is to investigate the electric potential around a charged sphere. The equipotential surfaces are concentric spherical shells around the charged sphere. They are identified by measuring the electrical potential at different positions but with constant distance to the surface of the sphere. In addition, the dependence of the variation of the electric potential on the distance to the surface of the sphere is determined and used to calculate the electric field strength.



P3.1.4  
EFFECTS OF FORCE IN  
AN ELECTRIC FIELD

P3.1.4.1  
Measuring the force of an  
electric charge in a homogeneous  
electric field



Measuring the force of an electric charge in a homogeneous electric field (P3.1.4.1)

Cat. No.	Description	P3.1.4.1
516 32	Current balance	1
314 081	Precision dynamometer, 0.01 N	1
314 263	Bodies for electric charge, set	1
541 00	Friction rods, PVC and acrylic	1
541 21	Leather	1
544 22	Parallel plate capacitor	1
300 75	Laboratory stand I	1
521 70	High-voltage power supply, 10 kV	1
501 051	Cable for high voltages, 1.5 m	2
471 830	He-Ne Laser, linearly polarised	1
441 53	Screen, translucent	1
300 01	Stand base, V-shaped, large	1
300 02	Stand base, V-shaped, small	1
300 11	Saddle base	1
300 42	Stand rod, 47 cm, 12 mm diam.	2
301 01	Leybold multiclamp	1
500 414	Connecting lead, 19 A, 25 cm, black	1

In a homogeneous electric field, the force  $F$  acting on an elongated charged body is proportional to the total charge  $Q$  and the electric field strength  $E$ . Thus, the formula

$$F = Q \cdot E$$

applies.

In the experiment P3.1.4.1, the greatest possible charge  $Q$  is transferred to an electrostatic spoon from a plastic rod. The electrostatic spoon is within the electric field of a plate capacitor and is aligned parallel to the plates. To verify the proportional relationship between the force and the field strength, the force  $F$  acting on the electrostatic spoon is measured at a known plate distance  $d$  as a function of the capacitor voltage  $U$ . The electric field  $E$  is determined using the equation

$$E = \frac{U}{d}$$

The measuring instrument in this experiment is a current balance, a differential balance with light-pointer read-out, in which the force to be measured is compensated by the spring force of a precision dynamometer.

### P3.1.4

#### EFFECTS OF FORCE IN AN ELECTRIC FIELD

##### P3.1.4.2

Kirchhoff's voltage balance:  
Measuring the force between two charged plates of a plate capacitor

##### P3.1.4.3

Measuring the force between a charged sphere and a metal plate



Measuring the force between a charged sphere and a metal plate (P3.1.4.3\_b)

Cat. No.	Description	P3.1.4.2 (b)	P3.1.4.3 (b)
516 37	Electrostatics with current balance, accessories	1	1
516 31	Vertically adjustable stand	1	1
524 005	Mobile-CASSY 2	1	1
524 060	Force sensor S, ±1N	1	1
314 265	Support for conductor loops	1	1
521 70	High-voltage power supply, 10 kV	1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	1
300 02	Stand base, V-shaped, small	1	1
301 01	Leybold multiclamp	1	1
500 410	Connecting lead, 19 A, 25 cm, yellow/green	1	
500 420	Connecting lead, 19 A, 50 cm, yellow/green	2	
500 444	Connecting lead 19 A, 100 cm, black	1	
541 00	Friction rods, PVC and acrylic		1
541 21	Leather		1
500 440	Connecting lead, 19 A, 100 cm, yellow/green		1

The force in an electric field is measured using a force sensor connected to a measuring instrument. The force sensor consists of two bending elements connected in parallel with four strain gauges in a bridge configuration; their electrical resistance changes when a load is applied. The change in resistance is proportional to the force acting on the sensor. The measuring instrument displays the measured force directly.

In the experiment P3.1.4.2 a balance is set up in order to measure the force

$$F = \frac{1}{2} \cdot \epsilon_0 \cdot \frac{U^2}{d^2} \cdot A$$

$$\text{where } \epsilon_0 = 8.85 \cdot 10^{-12} \frac{\text{As}}{\text{Vm}} \text{ (permittivity)}$$

acting between the two charged plates of a plate capacitor. At a given area  $A$ , the measurement is conducted as a function of the plate distance  $d$  and the voltage  $U$ . The aim of the evaluation is to confirm the proportionalities

$$F \propto \frac{1}{d^2} \text{ and } F \propto U^2$$

and to determine the permittivity  $\epsilon_0$ .

The experiment P3.1.4.3 consists of a practical investigation of the principle of the image charge. Here, the attractive force acting on a charged sphere in front of a metal plate is measured. This force is equivalent to the force of an equal, opposite charge at twice the distance  $2d$ . Thus, it is described by the formula

$$F = \frac{1}{4\pi\epsilon_0} \cdot \frac{Q^2}{(2d)^2}$$

First, the force for a given charge  $Q$  is measured as a function of the distance  $d$ . The measurement is then repeated with half the charge. The aim of the evaluation is to confirm the proportionalities

$$F \propto \frac{1}{d^2} \text{ and } F \propto Q^2$$

## P3.1.5

### CHARGE DISTRIBUTIONS ON ELECTRICAL CONDUCTORS

#### P3.1.5.1

Investigating the charge distribution on the surface of electrical conductors

#### P3.1.5.2

Electrostatic induction with the hemispheres after Cavendish



Electrostatic induction with the hemispheres after Cavendish (P3.1.5.2)

Cat. No.	Description	P3.1.5.1	P3.1.5.2
543 071	Conical conductor on insulating stand	1	
546 12	Faraday's cup	2	
542 52	Test plate, 4 cm x 4 cm	1	
521 70	High-voltage power supply, 10 kV	1	1
501 051	Cable for high voltages, 1.5 m	1	1
532 14	Electrometer amplifier	1	1
562 791	Plug-in power supply, 12 V AC	1	1
578 25	Capacitor, 1 nF, STE 2/19	1	1
578 10	Capacitor, 10 nF, STE 2/19	1	1
531 120	Multimeter LDanalog 20	1	1
590 011	Clamping plug	1	
532 16	Connecting rod	1	1
540 52	Demonstration insulator	1	
501 861	Crocodile-clips, polished, set of 6	1	
300 11	Saddle base	1	3
500 610	Safety connecting lead, 25 cm, yellow/green	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	1
500 424	Connecting lead 19 A, 50 cm, black	1	2
500 444	Connecting lead 19 A, 100 cm, black	1	
543 021	Sphere on insulated stand rod		1
543 05	Cavendish hemispheres, pair		1
340 89ET5	Coupling plugs, 4 mm, set of 5		1
300 41	Stand rod, 25 cm, 12 mm Ø		2
301 01	Leybold multiclamp		2
590 13	Stand rod with bore holes		1

In static equilibrium, the interior of a metal conductor or a hollow body contains neither electric fields nor free electron charges. On the outer surface of the conductor, the free charges are distributed in such a way that the electric field strength is perpendicular to the surface at all points, and all points have equal potential.

In the experiment P3.1.5.1, an electric charge is collected from a charged hollow metal sphere using a charge spoon, and measured using a coulomb meter. It becomes apparent that the charge density is greater, the smaller the bending radius of the surface is. This experiment also shows that no charge can be taken from the interior of the hollow body.

The experiment P3.1.5.2 reconstructs a historic experiment first performed by Cavendish. A metal sphere is mounted on an insulated base. Two hollow hemispheres surround the sphere completely, but without touching it. When one of the hemispheres is charged, the charge distributes itself uniformly over both hemispheres, while the inside sphere remains uncharged. If the inside sphere is charged and then surrounded by the hemispheres, the two hemispheres again show equal charges, and the inside sphere is uncharged.

### P3.1.6

#### DEFINITION OF CAPACITANCE

P3.1.6.1  
Determining the capacitance  
of a sphere in free space

P3.1.6.2  
Determining the capacitance  
of a sphere in front of a metal plate



Determining the capacitance of a sphere in free space (P3.1.6.1)

Cat. No.	Description	P3.1.6.1	P3.1.6.2
543 00	Set of 3 conducting spheres	1	1
521 70	High-voltage power supply, 10 kV	1	1
501 051	Cable for high voltages, 1.5 m	1	1
532 14	Electrometer amplifier	1	1
562 791	Plug-in power supply, 12 V AC	1	1
578 25	Capacitor, 1 nF, STE 2/19	1	1
578 10	Capacitor, 10 nF, STE 2/19	1	1
531 120	Multimeter LDanalog 20	1	1
546 12	Faraday's cup	1	1
590 011	Clamping plug	1	1
532 16	Connecting rod	1	1
590 13	Stand rod with bore holes	1	1
300 11	Saddle base	2	3
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	1
500 414	Connecting lead, 19 A, 25 cm, black	1	1
500 424	Connecting lead 19 A, 50 cm, black	1	1
500 444	Connecting lead 19 A, 100 cm, black	1	2
501 43	Connecting lead, 32 A, 200 cm, yellow/green	1	1
587 66	Reflection plate		1
501 861	Crocodile-clips, polished, set of 6		1
311 77	Steel tape measure, 2 m		1
300 42	Stand rod, 47 cm, 12 mm diam.		1

The potential difference  $U$  of a charged conductor in an insulated mounting in free space with reference to an infinitely distant reference point is proportional to the charge  $Q$  of the body. We can express this using the relationship

$$Q = C \cdot U$$

and call  $C$  the capacitance of the body. Thus, for example, the capacitance of a sphere with the radius  $r$  in a free space is

$$C = 4\pi\epsilon_0 \cdot r$$

because the potential difference of the charged sphere with respect to an infinitely distant reference point is

$$U = \frac{1}{4\pi\epsilon_0} \cdot \frac{Q}{r}$$

$$\text{where } \epsilon_0 = 8.85 \cdot 10^{-12} \frac{\text{As}}{\text{Vm}} \text{ (permittivity)}$$


The experiment P3.1.6.1 determines the capacitance of a sphere in a free space by charging the sphere with a known high voltage  $U$  and measuring its charge  $Q$  using an electrometer amplifier connected as a coulomb meter. The measurement is conducted for different sphere radii  $r$ . The aim of the evaluation is to verify the proportionalities

$$Q \propto U \text{ and } C \propto r$$

The experiment P3.1.6.2 shows that the capacitance of a body also depends on its environment, e.g. the distance to other earthed conductors. In this experiment, spheres with the radii  $r$  are arranged at a distance  $s$  from an earthed metal plate and charged using a high voltage  $U$ . The capacitance of the arrangement is now

$$C = 4\pi\epsilon_0 \cdot r \cdot \left(1 + \frac{r}{2s}\right)$$

The aim of the evaluation is to confirm the proportionality between the charge  $Q$  and the potential difference  $U$  at any given distance  $s$  between the sphere and the metal plate.



## P3.1.7

### PLATE CAPACITOR

#### P3.1.7.1

Determining the capacitance of a plate capacitor - Measuring the charge with the electrometer amplifier

#### P3.1.7.2

Parallel and series connection of capacitors - Measuring the charge with the electrometer amplifier

Determining the capacitance of a plate capacitor - Measuring the charge with the electrometer amplifier (P3.1.7.1)

Cat. No.	Description	P3.1.7.1	P3.1.7.2
544 23	Capacitor assembly kit	1	1
522 27	Power supply, 450 V	1	1
504 48	Two-way switch	1	1
531 120	Multimeter LDanalog 20	2	2
532 14	Electrometer amplifier	1	1
578 10	Capacitor, 10 nF, STE 2/19	1	1
578 31	Capacitor, 0.1 μF, STE 2/19	1	1
532 16	Connecting rod	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	4	5
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	2

A plate capacitor is the simplest form of a capacitor. Its capacitance depends on the plate area  $A$  and the plate spacing  $d$ . The capacitance increases when an insulator with the dielectric constant  $\epsilon_r$  is placed between the two plates. The total capacitance is

$$C = \epsilon_r \epsilon_0 \cdot \frac{A}{d}$$

$$\text{where } \epsilon_0 = 8.85 \cdot 10^{-12} \frac{\text{As}}{\text{Vm}} \text{ (permittivity)}$$

In the experiment P3.1.7.1, this relationship is investigated using a demountable capacitor assembly with variable geometry. Capacitor plates with the areas  $A = 40 \text{ cm}^2$  and  $A = 80 \text{ cm}^2$  can be used, as well as various plate-type dielectrics. The distance can be varied in steps of one millimeter.

The experiment P3.1.7.2 determines the total capacitance  $C$  of the demountable capacitor with the two plate pairs arranged at a fixed distance and connected first in parallel and then in series, compares these with the individual capacitances  $C_1$  and  $C_2$  of the two plate pairs. The evaluation confirms the relationship

$$C = C_1 + C_2$$

for parallel connection and

$$\frac{1}{C} = \frac{1}{C_1} + \frac{1}{C_2}$$

for serial connection.

P3.1.7

PLATE CAPACITOR

P3.1.7.3

Determining the capacitance of a plate capacitor – Measuring the charge with the I-measuring amplifier D



Determining the capacitance of a plate capacitor – Measuring the charge with the I-measuring amplifier D (P3.1.7.3)

Cat. No.	Description	P3.1.7.3
544 22	Parallel plate capacitor	1
521 65	Tube power supply, 0...500 V	1
504 48	Two-way switch	1
532 00	I-measuring amplifier D	1
531 120	Multimeter LDanalog 20	1
531 130	Multimeter LDanalog 30	1
536 221	Measuring resistor, 100 MΩ	1
500 421	Connecting lead 19 A, 50 cm, red	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	3
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1

Calculation of the capacitance of a plate capacitor using the formula

$$C = \epsilon_0 \cdot \frac{A}{d}$$

A: plate area

d: plate spacing

where  $\epsilon_0 = 8.85 \cdot 10^{-12} \frac{\text{As}}{\text{Vm}}$  (permittivity)

ignores the fact that part of the electric field of the capacitor extends beyond the edge of the plate capacitor, and that consequently a greater charge is stored for a specific potential difference between the two capacitors. For example, for a plate capacitor grounded on one side and having the area

$$A = \pi \cdot r^2$$


the capacitance is given by the formula

$$C = \epsilon_0 \left( \frac{\pi \cdot r^2}{d} + 3.7724 \cdot r + r \cdot \ln \left( \frac{\pi r}{d} \right) + \dots \right)$$

In the experiment P3.1.7.3, the capacitance C of a plate capacitor is measured as a function of the plate spacing d with the greatest possible accuracy. This experiment uses a plate capacitor with a plate radius of 13 cm and a plate spacing which can be continuously varied between 0 and 70 mm. The aim of the evaluation is to plot the measured values in the form

$$C = f \left( \frac{1}{d} \right)$$

and compare them with the values to be expected according to theory.



**P3.1.7**  
**PLATE CAPACITOR**

P3.1.7.4  
Measuring the electric field strength inside a plate capacitor

P3.1.7.5  
Measuring the electric field strength inside a plate capacitor as a function of the dielectrics

P3.1.7.6  
Measuring the electric field strength of a charged sphere in front of a conductive plate (image charge)

Measuring the electric field strength inside a plate capacitor (P3.1.7.4\_c)

Cat. No.	Description	P3.1.7.4 (c)	P3.1.7.5 (b)	P3.1.7.6 (c)
524 080	Electric field meter S	1	1	1
540 540	Accessories for electric field meter S	1	1	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>		<b>1</b>
524 220	CASSY Lab 2	1		1
521 70	High-voltage power supply, 10 kV	1		1
460 317	Optical bench, S1 profile, 0.5 m	1	1	
460 312	Clamp rider with clamp, 45/35	2	2	
500 600	Safety connection lead, 10 cm, yellow/green	1		1
500 641	Safety connecting lead, 100 cm, red	1		1
500 642	Safety connecting lead, 100 cm, blue	1		1
<b>524 005</b>	<b>Mobile-CASSY 2</b>		<b>1</b>	
531 120	Multimeter LDanalog 20		1	
522 27	Power supply, 450 V		1	
504 45	On-off switch, single pole		1	
500 421	Connecting lead 19 A, 50 cm, red		3	
500 422	Connecting lead 19 A, 50 cm, blue		1	
500 442	Connecting lead 19 A, 100 cm, blue		1	
543 021	Sphere on insulated stand rod			1
311 02	Metal rule, 1 m			1
300 11	Saddle base			2
500 95	Safety adapter sockets, red, set of 6			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1		1

Using the electric field meter S the electric field strength  $E$  in a plate capacitor can be measured. The electric field strength depends on the applied voltage  $U$  and the distance  $d$  of the capacitor plates:

$$E = \frac{U}{d}$$

Alternatively, the electrical field strength  $E$  can be calculated from the charge  $Q$  on the capacitor plates:

$$E = \frac{Q}{\epsilon_0 \cdot \epsilon_r \cdot A}$$

Here,  $E$  depends on the area of the plates  $A$  and the permittivity  $\epsilon_r$  of the material between the capacitor plates as well.

In the experiment P3.1.7.4 the dependance of the electric field strength  $E$  on the applied voltage  $U$  and the plate spacing  $d$  is determined. First, keeping the distance of the plates constant, the value of the applied voltage  $U$  is varied and the electric field strength is measured. Then, the voltage  $U$  is kept constant and the dependance of the electric field strength  $E$  on the plate spacing  $d$  is determined.

The aim of the experiment P3.1.7.5 is to investigate the influence of the permittivity  $\epsilon_r$  on the field strength  $E$ . First, keeping the applied voltage  $U$  constant a dielectric (glass, plastics) is placed between the capacitor plates and the electric field strength is measured. Second, the charged capacitor is disconnected from the power supply. Then, the dielectric is removed and the field strength measured again.


In the experiment P3.1.7.6, the electric field strength on the surface of a conductive plate with distance  $r$  to a charged sphere is measured. The field gradient in front of the plate is equivalent to the case where instead of the plate a sphere with opposite charge is situated in twice the distance to the sphere (mirror or image charge). This leads to a doubling in field strength compared to a free-standing sphere.

### P3.2.1

#### CHARGE TRANSFER WITH DROPS OF WATER

##### P3.2.1.1

Generating an electric current through the motion of charged drops of water



Generating an electric current through the motion of charged drops of water (P3.2.1.1)

Cat. No.	Description	P3.2.1.1
665 843	Burette, clear glass, 10 ml	1
522 27	Power supply, 450 V	1
532 14	Electrometer amplifier	1
532 16	Connecting rod	1
546 12	Faraday's cup	1
578 25	Capacitor, 1 nF, STE 2/19	1
578 26	Capacitor, 2.2 nF, STE 2/19	1
578 10	Capacitor, 10 nF, STE 2/19	1
578 22	Capacitor, 100 pF, STE 2/19	1
531 120	Multimeter LD analog 20	1
501 641	Two-way adapters, red, set of 6	1
550 41	Constantan resistance wire, 0.25 mm diameter, 100 m	1
501 861	Crocodile-clips, polished, set of 6	1
664 120	Beaker, PP, 50 ml, squat	1
301 21	Stand base MF	2
301 27	Stand rod, 50 cm, 10 mm diam.	1
301 26	Stand rod, 25 cm, 10 mm diam.	1
301 01	Leybold multiclamp	1
666 555	Universal clamp, 0...80 mm	1
500 412	Connecting lead, 19 A, 25 cm, blue	1
500 424	Connecting lead 19 A, 50 cm, black	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2
500 444	Connecting lead 19 A, 100 cm, black	2
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1*</b>
524 220	CASSY Lab 2	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

\* additionally recommended

Each charge transport is an electric current. The electrical current strength (or more simply the "current")

$$I = \frac{\Delta Q}{\Delta t}$$

is the charge  $\Delta Q$  transported per unit of time  $\Delta t$ . For example, in a metal conductor,  $\Delta Q$  is given by the number  $\Delta N$  of free electrons which flow through a specific conductor cross-section per unit of time  $\Delta t$ . We can illustrate this relationship using charged water droplets.

In the experiment P3.2.1.1, charged water drops drip out of a burette at a constant rate

$$\dot{N} = \frac{\Delta N}{\Delta t}$$


$N$ : number of water drops

into a Faraday's cup, and gradually charge the latter. Each individual drop of water transports approximately the same charge  $q$ . The total charge  $Q$  in the Faraday's cup is measured using an electrometer amplifier connected as a coulomb meter. This charge shows a step-like curve as a function of the time  $t$ , as can be recorded using CASSY. At a high drip rate  $N$ , a very good approximation is

$$Q = \dot{N} \cdot q \cdot t$$


The current is then

$$I = \dot{N} \cdot q$$



Charge transfer by water drops





### P3.2.2

#### OHM'S LAW

##### P3.2.2.1

Verifying Ohm's law and measuring specific resistances

##### P3.2.2.2

Verifying Ohm's law - Measuring with CASSY

Verifying Ohm's law and measuring specific resistances (P3.2.2.1)

Cat. No.	Description	P3.2.2.1	P3.2.2.2
550 57	Apparatus for resistance measurements	1	
521 49	AC/DC power supply, 0...12 V	1	
531 120	Multimeter LDanalog 20	2	
501 23	Connecting lead, 32 A, 25 cm, black	1	
501 33	Connecting lead, 32 A, 100 cm, black	3	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	
<b>524 006</b>	<b>Pocket-CASSY</b>		<b>1</b>
524 220	CASSY Lab 2		1
524 0621	UIP sensor S		1
576 74	Plug-in board, DIN A4, STE		1
501 48	Bridging plugs, STE 2/19, set of 10		1
567 18	Wrapping plate for wires		1
550 46	Chrome-nickel resistance wire, 0.25 mm diameter, 100 m		1
579 331	Plug-in holder, STE		2
579 13	Toggle switch, STE 2/19		1
577 32	Resistor, 100 Ω, STE 2/19		1
521 231	Low-voltage power supply, 3/6/9/12 V		1
501 44	Connecting leads, 19 A, 25 cm, red/blue, pair		3
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

In circuits consisting of metal conductors, Ohm's law

$$U = R \cdot I$$

represents a very close approximation of the actual circumstances. In other words, the voltage drop  $U$  in a conductor is proportional to the current  $I$  through the conductor. The proportionality constant  $R$  is called the resistance of the conductor. For the resistance, we can say

$$R = \rho \cdot \frac{s}{A}$$

$\rho$ : resistivity of the conductor material

$s$ : length of wire

$A$ : cross-section of wire

The experiment P3.2.2.1 verifies the proportionality between the current and voltage for metal wires of different materials, thicknesses and lengths, and calculates the resistivity of each material.

The experiment P3.2.2.2 verifies the proportionality between the current and voltage for a metal wire using the CASSY system.

### P3.2.3

#### KIRCHHOFF'S LAWS

##### P3.2.3.1


Measuring current and voltage at resistors connected in parallel and in series

##### P3.2.3.2

Voltage division with a potentiometer

##### P3.2.3.3

Principle of a Wheatstone bridge



Measuring current and voltage at resistors connected in parallel and in series (P3.2.3.1)

Cat. No.	Description	P3.2.3.1	P3.2.3.2	P3.2.3.3
576 74	Plug-in board, DIN A4, STE	1	1	1
577 36	Resistor, 220 $\Omega$ , STE 2/19	1		1
577 38	Resistor, 330 $\Omega$ , STE 2/19	1		2
577 40	Resistor, 470 $\Omega$ , STE 2/19	1	1	1
577 44	Resistor, 1 k $\Omega$ , STE 2/19	1		1
577 53	Resistor, 5.6 k $\Omega$ , STE 2/19	1		
577 56	Resistor, 10 k $\Omega$ , STE 2/19	1		
577 68	Resistor, 100 k $\Omega$ , STE 2/19	1		
501 48	Bridging plugs, STE 2/19, set of 10	1	1	1
521 45	DC power supply 0... $\pm$ 15 V	1	1	1
531 120	Multimeter LDanalog 20	2	2	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	3	3	2
577 28	Resistor, 47 $\Omega$ , STE 2/19		1	
577 32	Resistor, 100 $\Omega$ , STE 2/19		2	
577 34	Resistor, 150 $\Omega$ , STE 2/19		1	
577 90	Potentiometer, 220 $\Omega$ , STE 4/50		1	
577 92	Potentiometer, 1 k $\Omega$ , STE 4/50			1

Kirchhoff's laws are of fundamental importance in calculating the component currents and voltages in branching circuits. The so-called "node rule" states that the sum of all currents flowing into a particular junction point in a circuit is equal to the sum of all currents flowing away from this junction point. The "mesh rule" states that in a closed path the sum of all voltages through the loop in any arbitrary direction of flow is zero. Kirchhoff's laws are used to derive a system of linear equations which can be solved for the unknown current and voltage components.

The experiment P3.2.3.1 examines the validity of Kirchhoff's laws in circuits with resistors connected in parallel and in series. The result demonstrates that two resistors connected in series have a total resistance  $R$

$$R = R_1 + R_2$$

while for parallel connection of resistors, the total resistance  $R$  is

$$\frac{1}{R} = \frac{1}{R_1} + \frac{1}{R_2}$$

In the experiment P3.2.3.2, a potentiometer is used as a voltage divider in order to tap a lower voltage component  $U_1$  from a voltage  $U$ .  $U$  is present at the total resistance of the potentiometer. In a no-load, zero-current state, the voltage component

$$U_1 = \frac{R_1}{R} \cdot U$$

can be tapped at the variable component resistor  $R_1$ . The relationship between  $U_1$  and  $R_1$  at the potentiometer under load is no longer linear.


The experiment P3.2.3.3 examines the principle of a Wheatstone bridge, in which "unknown" resistances can be measured through comparison with "known" resistances.

### P3.2.3

#### KIRCHHOFF'S LAWS

##### P3.2.3.4

Determining resistances using a Wheatstone bridge



Determining resistances using a Wheatstone bridge (P3.2.3.4)


Cat. No.	Description	P3.2.3.4
536 02	Demonstration bridge	1
536 121	Measuring resistor, 10 Ω	1
536 131	Measuring resistor, 100 Ω	1
536 141	Measuring resistor, 1 kΩ	1
536 776	Decade resistor, 0...1 kΩ	1
536 777	Decade resistor, 0...100 Ω	1
536 778	Decade resistor, 0...10 Ω	1
536 779	Decade resistor, 0...1 Ω	1
521 45	DC power supply 0...±15 V	1
531 13	Galvanometer C.A 403	1
501 28	Connecting lead, 32 A, 50 cm, black	3
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1

In modern measuring practice, the bridge configuration developed in 1843 by *Ch. Wheatstone* is used almost exclusively.

In the experiment P3.2.3.4, a voltage  $U$  is applied to a 1 m long measuring wire with a constant cross-section. The ends of the wire are connected to an unknown resistor  $R_x$  and a variable resistor  $R$  arranged behind it, whose value is known precisely. A sliding contact divides the measuring wire into two parts with the lengths  $s_1$  and  $s_2$ . The slide contact is connected to the node between  $R_x$  and  $R$  via an ammeter which is used as a zero indicator. Once the current has been regulated to zero, the relationship

$$R_x = \frac{s_1}{s_2} \cdot R$$

applies. Maximum accuracy is achieved by using a symmetrical experiment setup, i. e. when the slide contact over the measuring wire is set in the middle position so that the two sections  $s_1$  and  $s_2$  are the same length.



Circuit diagram of Wheatstone bridge

### P3.2.4


#### CIRCUITS WITH ELECTRICAL MEASURING INSTRUMENTS

##### P3.2.4.1

The ammeter as an ohmic resistor in a circuit

##### P3.2.4.2

The voltmeter as an ohmic resistor in a circuit



*The ammeter as an ohmic resistor in a circuit (P3.2.4.1)*

Cat. No.	Description	P3.2.4.1	P3.2.4.2
521 45	DC power supply 0...±15 V	1	1
576 74	Plug-in board, DIN A4, STE	1	1
577 33	Resistor, 82 Ω, STE 2/19	3	
577 52	Resistor, 4.7 kΩ, STE 2/19	1	1
531 110	Multimeter LDanalog 10	2	2
501 48	Bridging plugs, STE 2/19, set of 10	1	
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	3	3
577 75	Resistor, 680 kΩ, STE 2/19		1
577 71	Resistor, 220 kΩ, STE 2/19		1

One important consequence of Kirchhoff's laws is that the internal resistance of an electrical measuring instrument affects the respective current or voltage measurement. Thus, an ammeter increases the overall resistance of a circuit by the amount of its own internal resistance and thus measures a current value which is too low whenever the internal resistance is above a negligible level. A voltmeter measures a voltage value which is too low when its internal resistance is not great enough with respect to the resistance at which the voltage drop is to be measured.

In the experiment P3.2.4.1, the internal resistance of an ammeter is determined by measuring the voltage which drops at the ammeter during current measurement. It is subsequently shown that the deflection of the ammeter pointer is reduced by half, or that the current measuring range is correspondingly doubled, by connecting a second resistor equal to the internal resistance in parallel to the ammeter.

The experiment P3.2.4.2 determines the internal resistance of a voltmeter by measuring the current flowing through it. In this experiment, the measuring range is extended by connecting a second resistor with a value equal to the internal resistance to the voltmeter in series.

P3.2.5

CONDUCTING ELECTRICITY  
BY MEANS OF ELECTROLYSIS

P3.2.5.1

Determining the Faraday constant



Determining the Faraday constant (P3.2.5.1)

Cat. No.	Description	P3.2.5.1
664 350	Water electrolysis apparatus	1
382 35	Thermometer, -10...+50 °C/0.1 K	1
531 832	Digital multimeter P	1
521 45	DC power supply 0...±15 V	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
649 45	Tray, 552 mm x 459 mm x 48 mm	1
674 7920	Sulfuric acid, diluted, approx. 2 N, 500 ml	1

In electrolysis, the processes of electrical conduction entails liberation of material. The quantity of liberated material is proportional to the transported charge  $Q$  flowing through the electrolyte. This charge can be calculated using the Faraday constant  $F$ , a universal constant which is related to the unit charge  $e$  by means of Avogadro's number  $N_A$ .

$$F = N_A \cdot e$$

When we insert the molar mass  $n$  for the material quantity and take the valence  $z$  of the separated ions into consideration, we obtain the relationship

$$Q = n \cdot F \cdot z$$

In the experiment P3.2.5.1, a specific quantity of hydrogen is produced in an electrolysis apparatus after Hofmann to determine the Faraday constant. The valence of the hydrogen ions is  $z = 1$ . The molar mass  $n$  of the liberated hydrogen atoms is calculated using the laws of ideal gas on the basis of the volume  $V$  of the hydrogen collected at an external pressure  $p$  and room temperature  $T$ :

$$n = 2 \cdot \frac{pV}{RT}$$

$$\text{where } R = 8.314 \frac{\text{J}}{\text{mol} \cdot \text{K}} \text{ (universal gas constant)}$$

At the same time, the electric work  $W$  is measured which is expended for electrolysis at a constant voltage  $U_0$ . The transported charge quantity is then

$$Q = \frac{W}{U_0}$$

### P3.2.6

#### EXPERIMENTS ON ELECTROCHEMISTRY

##### P3.2.6.1

Generating electric current with a Daniell cell

##### P3.2.6.2

Measuring the voltage at simple galvanic elements

##### P3.2.6.3

Determining the standard potentials of corresponding redox pairs



Measuring the voltage at simple galvanic elements (P3.2.6.2)

Cat. No.	Description	P3.2.6.1-3
667 4041	Measuring unit S	1
664 3951	Electrochemistry workstation	1
661 125	Electrochemistry chemicals set	1

In galvanic cells, electrical energy is generated using an electrochemical process. The electrochemistry workplace enables you to investigate the physical principles which underlie such processes.

In the experiment P3.2.6.1, a total of four Daniell cells are assembled. These consist of one half-cell containing a zinc electrode in a  $ZnSO_4$  solution and one half-cell containing a copper electrode in a  $CuSO_4$  solution. The voltage produced by multiple cells connected in series is measured and compared with the voltage from a single cell. The current of a single cell is used to drive an electric motor.

The experiment P3.2.6.2 combines half-cells of corresponding redox pairs of the type metal/metal cation to create simple galvanic cells. For each pair, the object is to determine which metal represents the positive and which one the negative pole, and to measure the voltage between the half-cells. From this, a voltage series for the corresponding redox pairs can be developed.

The experiment P3.2.6.3 uses a platinum electrode in 1-mol hydrochloric acid as a simple standard hydrogen electrode in order to permit direct measurement of the standard potentials of corresponding redox pairs of the type metal/metal cation and nonmetallic anion/non-metallic substance directly.

## P3.3.1

### BASIC EXPERIMENTS ON MAGNETOSTATICS

#### P3.3.1.1

Displaying lines of magnetic flux

#### P3.3.1.2

Basics of electromagnetism



Basics of electromagnetism (P3.3.1.2)

Cat. No.	Description	P3.3.1.1	P3.3.1.2
560 701	Magnetic field demonstration set	1	
521 55	High current power supply	1	1
501 26	Connecting lead, 32 A, 50 cm, blue	1	1
501 30	Connecting lead, 32 A, 100 cm, red	1	
501 31	Connecting lead, 32 A, 100 cm, blue	1	
MIK 74702	BMS EcoCam 5.5 mega pixel, WiFi	1*	
560 15	Equipment for electromagnetism		1
513 511	Magnetic needle on base, needle bearing		1
510 21	Horseshoe magnet with yoke		1
510 12	Cylindrical bar magnets, pair		1
514 72ET5	Shaker for iron filings, set of 5		1
514 73	Iron powder		1
314 111	Precision dynamometer , 0.1 N		1
300 02	Stand base, V-shaped, small		1
300 43	Stand rod, 75 cm, 12 mm diam.		1
301 01	Leybold multiclamp		3
666 555	Universal clamp, 0...80 mm		1
540 52	Demonstration insulator		2
300 11	Saddle base		2
501 35	Connecting lead, 32 A, 200 cm, red		1
501 36	Connecting lead, 32 A, 200 cm, blue		1
	additionally required: PC with Windows 7 or higher	1	

\* additionally recommended

Magnetostatics studies the spatial distribution of magnetic fields in the vicinity of permanent magnets and stationary currents as well as the force exerted by a magnetic field on magnets and currents. Basic experiments on this topic can be carried out without complex experiment setups.

In the experiment P3.3.1.1, magnetic fields are observed by spreading iron filings over a smooth surface so that they align themselves with the lines of magnetic flux. By this means it becomes possible to display the magnetic field of a straight conductor, the magnetic field of a conductor loop and the magnetic field of a coil.

The experiment P3.3.1.2 combines a number of fundamental experiments on electromagnetic phenomena. First, the magnetic field surrounding a current-carrying conductor is illustrated. Then the force exerted by two current-carrying coils on each other and the deflection of a current-carrying coil in the magnetic field of a second coil are demonstrated.



Displaying lines of magnetic flux

## P3.3.2

### MAGNETIC DIPOLE MOMENT

#### P3.3.2.1

Measuring the magnetic dipole moments of long magnetic needles



Measuring the magnetic dipole moments of long magnetic needles (P3.3.2.1)

Cat. No.	Description	P3.3.2.1
516 01	Torsion balance, Schürholz design	1
516 21	Magnetostatics, accessories	1
516 04	Scale on stand	1
510 50ET2	Bar magnets, 60 x 13 x 5 mm, set of 2	1
450 60	Lamp housing with cable	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
460 20	Condenser with diaphragm holder	1
521 210	Transformer, 6/12 V	1
300 02	Stand base, V-shaped, small	1
300 42	Stand rod, 47 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1

Although only magnetic dipoles occur in nature, it is useful in some cases to work with the concept of highly localized "magnetic charges". Thus, we can assign pole strengths or "magnetic charges"  $q_m$  to the pole ends of elongated magnetic needles on the basis of their length  $d$  and their magnetic moment  $m$ :

$$q_m = \frac{m}{d}$$

The pole strength is proportional to the magnetic flux  $\Phi$ :

$$\Phi = \mu_0 \cdot q_m$$

$$\text{where } \mu_0 = 4\pi \cdot 10^{-7} \frac{\text{Vs}}{\text{Am}} \text{ (permeability)}$$

Thus, for the spherical surface with a small radius  $r$  around the pole (assumed as a point source), the magnetic field is

$$B = \frac{1}{4\pi\mu_0} \cdot \frac{q_m}{r^2}$$

At the end of a second magnetic needle with the pole strength  $q'_m$ , the magnetic field exerts a force

$$F = q'_m \cdot B$$

and consequently

$$F = \frac{1}{4\pi\mu_0} \cdot \frac{q_m \cdot q'_m}{r^2}$$

In formal terms, this relationship is equivalent to Coulomb's law governing the force between two electrical charges.

The experiment P3.3.2.1 measures the force  $F$  between the pole ends of two magnetized steel needles using the torsion balance. The experiment setup is similar to the one used to verify Coulomb's law. The measurement is initially carried out as a function of the distance  $r$  of the pole ends. To vary the pole strength  $q_m$ , the pole ends are exchanged, and multiple steel needles are mounted next to each other in the holder.



## P3.3.3 EFFECTS OF FORCE IN A MAGNETIC FIELD


P3.3.3.1  
Measuring the force acting on current-carrying conductors in the field of a horseshoe magnet

P3.3.3.2  
Measuring the force acting on current-carrying conductors in a homogeneous magnetic field – Recording with CASSY

P3.3.3.3  
Measuring the force acting on current-carrying conductors in the magnetic field of an air coil – Recording with CASSY

P3.3.3.4  
Basic measurements for the electrodynamic definition of the ampere

P3.3.3.5  
Demonstration of the force between a current-carrying conductor and a permanent magnet



Measuring the force acting on current-carrying conductors in the field of a horseshoe magnet (P3.3.3.1\_b)

Cat. No.	Description	P3.3.3.1 (b)	P3.3.3.2	P3.3.3.3	P3.3.3.4 (b)	P3.3.3.5
510 22	Large horseshoe magnet with yoke	1				
314 265	Support for conductor loops	1	1	1	1	
516 34	Conductor loops for force measurement	1	1	1		
521 55	High current power supply	1	1	1	1	
524 005	Mobile-CASSY 2	1			1	
524 060	Force sensor S, ±1N	1	1	1	1	
300 02	Stand base, V-shaped, small	1	1	1	1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	1	1	1	
301 01	Leybold multiclamp	1	1	1	1	
501 30	Connecting lead, 32 A, 100 cm, red	1	2	2	1	
501 31	Connecting lead, 32 A, 100 cm, blue	1	2	2	1	
562 11	U-core with yoke		1			
562 14	Coil, 500 turns		2			
562 25	Pole-shoe yoke		1			
524 013	Sensor-CASSY 2		1	1		
524 220	CASSY Lab 2		1	1		
524 0431	30-A-Box		1	1		
521 501	AC/DC power supply, 0...15 V/0...5 A		1	1		
501 26	Connecting lead, 32 A, 50 cm, blue		2	1	1	
516 244	Field coil, d = 120 mm			1		
516 249	Stand for tubes and coils			1		
516 33	Conductors for Ampere definition				1	
516 31	Vertically adjustable stand				1	
510 20	Horseshoe magnet, small					1
315 233	Electronic balance MAULTec S					1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A					1
500 444	Connecting lead 19 A, 100 cm, black					1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1	1		

To measure the force acting on a current-carrying conductor in a magnetic field, conductor loops are attached to a force sensor. The force sensor contains two bending elements arranged in parallel with four strain gauges connected in a bridge configuration; their resistance changes in proportion to the force when a strain is applied. The force sensor is connected to a measuring instrument, or alternatively to the CASSY computer interface device. When using CASSY a 30 ampere box is recommended for current measurement.

In the experiment P3.3.3.1, the conductor loops are placed in the magnetic field of a horseshoe magnet. This experiment measures the force  $F$  as a function of the current  $I$ , the conductor length  $s$  and the angle  $\alpha$  between the magnetic field and the conductor, and reveals the relationship

$$F = I \cdot s \cdot B \cdot \sin \alpha$$

In the experiment P3.3.3.2, a homogeneous magnetic field is generated using an electromagnet with U-core and pole-piece attachment. This experiment measures the force  $F$  as a function of the current  $I$ . The measurement results for various conductor lengths  $s$  are compiled and evaluated in a graph.

The experiment P3.3.3.3 uses an air coil to generate the magnetic field. The magnetic field is calculated from the coil parameters and compared with the values obtained from the force measurement.

The object of the experiment P3.3.3.4 is the electrodynamic definition of the ampere. Here, the current is defined on the basis of the force exerted between two parallel conductors of infinite length which carry an identical current. When  $r$  represents the distance between the two conductors, the force per unit of length of the conductor is:

$$\frac{F}{s} = \mu_0 \cdot \frac{I_1 I_2}{2\pi \cdot r}$$

This experiment uses two conductors approx. 30 cm long, placed just a few millimeters apart. The forces  $F$  are measured as a function of the different current levels  $I$  and distances  $r$ .

In the experiment P3.3.3.5 the interaction forces between a current-carrying conductor and a permanent magnet are recorded using a balance. Depending on the current direction, the „weight“ of the magnet is increased or decreased by the current related forces.

## P3.3.4

### BIOT-SAVART'S LAW

#### P3.3.4.1

Measuring the magnetic field for a straight conductor and on circular conductor loops

#### P3.3.4.2

Measuring the magnetic field of an air coil

#### P3.3.4.3

Measuring the magnetic field of a pair of coils in the Helmholtz configuration



Measuring the magnetic field for a straight conductor and on circular conductor loops (P3.3.4.1\_b)

Cat. No.	Description	P3.3.4.1 (b)	P3.3.4.1 (c)	P3.3.4.2 (b)	P3.3.4.2 (c)	P3.3.4.3 (b)	P3.3.4.3 (c)
516 235	Current conductors, set of 4	1	1				
524 005	Mobile-CASSY 2	1	1	1	1	1	1
524 0381	Combi B sensor S	1					
501 11	Extension cable, 15 pin	1	1	1	1	1	1
521 55	High current power supply	1		1		1	
460 21	Holder for plug-in elements	1	1				
460 43	Small optical bench	1	1			1	1
301 01	Leybold multiclamp	2	3			3	3
300 01	Stand base, V-shaped, large	1	1			1	1
501 644	Two-way adapters, black, set of 6	1	1				
501 30	Connecting lead, 32 A, 100 cm, red	1	1	1	1	1	1
501 31	Connecting lead, 32 A, 100 cm, blue	1	1	1	1	1	2
524 0383	Axial B sensor S, ±0.3 mT			1		1	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A			1		1	1
300 43	Stand rod, 75 cm, 12 mm diam.			1			
516 242	Coil with variable number of turns per unit length			1	1		
516 249	Stand for tubes and coils			1	1		
524 0382	Axial B sensor S, ±1000 mT			1		1	
300 11	Saddle base			1	1		
555 604	Pair of Helmholtz coils					1	1
501 26	Connecting lead, 32 A, 50 cm, blue					1	1

In principle, it is possible to calculate the magnetic field of any current-carrying conductor using Biot and Savart's law. However, analytical solutions can only be derived for conductors with certain symmetries, e.g. for an infinitely long straight wire, a circular conductor loop and a cylindrical coil. Biot and Savart's law can be verified easily using these types of conductors.

In the experiment P3.3.4.1, the magnetic field of a long, straight conductor is measured for various currents  $I$  as a function of the distance  $r$  from the conductor. The result is a quantitative confirmation of the relationship

$$B = \frac{\mu_0}{2\pi} \cdot \frac{I}{r}$$

In addition, the magnetic fields of circular coils with different radii  $R$  are measured as a function of the distance  $x$  from the axis through the center of the coil. The measured values are compared with the values which are calculated using the equation

$$B = \frac{\mu_0}{2} \cdot \frac{I \cdot R^2}{(R^2 + x^2)^{3/2}}$$

The measurements can be carried out using the combi B sensor. This device contains two Hall sensors which one is extremely sensitive to fields parallel to the probe axis and the second one is sensitive perpendicular to the probe axis. Using a sensitive magnetic field sensor, the experiment is performed at low current levels.

The experiment P3.3.4.2 investigates the magnetic field of an air coil in which the length  $L$  can be varied for a constant number of turns  $N$ . For the magnetic field the relationship


$$B = \mu_0 \cdot I \cdot \frac{N}{L}$$

applies.

The experiment P3.3.4.3 examines the homogeneity of the magnetic field in a pair of Helmholtz coils. The magnetic field along the axis through the coil centers is recorded in several measurement series; the spacing  $a$  between the coils is varied from measurement series to measurement series. When  $a$  is equal to the coil radius, the magnetic field is essentially independent of the location  $x$  on the coil axis.

P3.4.1  
VOLTAGE IMPULSE

P3.4.1.1  
Generating a voltage surge in a conductor loop with a moving permanent magnet



Generating a voltage surge in a conductor loop with a moving permanent magnet (P3.4.1.1)

Cat. No.	Description	P3.4.1.1
510 11	Cylindrical bar magnet	2
562 13	Coil, 250 turns	1
562 14	Coil, 500 turns	1
562 15	Coil, 1 000 turns	1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

Each change in the magnetic flux  $\Phi$  through a conductor loop induces a voltage  $U$ , which has a level proportional to the change in the flux. Such a change in the flux is caused e. g. when a permanent magnet is moved inside a fixed conductor loop. In this case, it is common to consider not only the time-dependent voltage


$$U = -\frac{d\Phi}{dt}$$

but also the voltage surge

$$\int_{t_1}^{t_2} U(t) dt = \Phi(t_1) - \Phi(t_2)$$

This corresponds to the difference in the magnetic flux densities before and after the change.

In the experiment P3.4.1.1, the voltage surge is generated by manually inserting a bar magnet into an air coil, or pulling it out of a coil. The curve of the voltage  $U$  over time is measured and the area inside the curve is evaluated. This is always equal to the flux  $\Phi$  of the permanent magnet inside the air coil independent of the speed at which the magnet is moved, i. e. proportional to the number of turns of the coil for equal coil areas.



Induced voltages of a moving magnet

### P3.4.2

#### INDUCTION IN A MOVING CONDUCTOR LOOP

##### P3.4.2.1

Measuring the induction voltage in a conductor loop moved through a magnetic field



Measuring the induction voltage in a conductor loop moved through a magnetic field (P3.4.2.1\_b)

Cat. No.	Description	P3.4.2.1 (b)
516 40	Induction apparatus with wire loops	1
510 48	Magnets, 35 mm Ø, pair	6
347 35	Experiment motor	1
347 36	Control unit for experiment motor	1
<b>524 005</b>	<b>Mobile-CASSY 2</b>	<b>1</b>
524 040	µV box	1

When a conductor loop with the constant width  $b$  is withdrawn from a homogeneous magnetic field  $B$  with the speed

$$v = \frac{dx}{dt}$$

the magnetic flux changes over the time  $dt$  by the value

$$d\Phi = -B \cdot b \cdot dx$$


This change in flux induces the voltage

$$U = B \cdot b \cdot v$$


in the conductor loop.

In the experiment P3.4.2.1, a slide on which induction loops of various widths are mounted is moved between the two pole pieces of a magnet. The object is to measure the induction voltage  $U$  as a function of the magnetic flux density  $B$ , the width  $b$  and the speed  $v$  of the induction loops. The aim of the evaluation is to verify the proportionalities

$$U \propto B, U \propto b, U \propto v$$



Induction voltage in a moved conductor loop



### P3.4.3

#### INDUCTION BY MEANS OF A VARIABLE MAGNETIC FIELD

P3.4.3.1  
Measuring the induction voltage in a conductor loop for a variable magnetic field - with triangular wave-form power supply

P3.4.3.2  
Measuring the induction voltage in a conductor loop for a variable magnetic field - with Power-CASSY as variable source of current

Measuring the induction voltage in a conductor loop for a variable magnetic field - with triangular wave-form power supply (P3.4.3.1\_a)

Cat. No.	Description	P3.4.3.1 (a)	P3.4.3.1 (b)	P3.4.3.2
516 249	Stand for tubes and coils	1	1	1
516 244	Field coil, d = 120 mm	1	1	1
516 241	Induction coils, set of 3	1	1	1
521 56	Delta current power supply	1	1	
524 013	Sensor-CASSY 2	1		1
524 220	CASSY Lab 2	1		1
524 040	µV box	1	1	1
524 0431	30-A-Box	1	1	
500 422	Connecting lead 19 A, 50 cm, blue	1	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	2	2
524 005	Mobile-CASSY 2		1	
524 011USB	Power-CASSY USB			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1		1

A change in the homogeneous magnetic field  $B$  inside a coil with  $N_1$  windings and the area  $A_1$  over time induces the voltage

$$U = N_1 \cdot A_1 \cdot \frac{dB}{dt}$$


in the coil.

In the experiments P3.4.3.1 and P3.4.3.2, induction coils with different areas and numbers of turns are arranged in a cylindrical field coil through which alternating currents of various frequencies, amplitudes and signal forms flow. In the field coil, the currents generate the magnetic field

$$B = \mu_0 \cdot \frac{N_2}{L_2} \cdot I$$

$$\text{where } \mu_0 = 4\pi \cdot 10^{-7} \frac{\text{Vs}}{\text{Am}} \text{ (permeability)}$$

and  $I(t)$  is the time-dependent current level,  $N_2$  the number of turns and  $L_2$  the overall length of the coil. The curve over time  $U(t)$  of the voltages induced in the induction coils is recorded using the computer-based CASSY measuring system. This experiment explores how the voltage is dependent on the area and the number of turns of the induction coils, as well as on the frequency, amplitude and signal form of the exciter current.



Induction in a conduction loop for a variable magnetic field

### P3.4.4

#### EDDY CURRENTS

##### P3.4.4.1

Waltenhofen's pendulum:  
demonstration of an  
eddy-current brake

##### P3.4.4.2

Demonstrating the operating  
principle of an AC power meter



Waltenhofen's pendulum: demonstration of an eddy-current brake (P3.4.4.1)

Cat. No.	Description	P3.4.4.1	P3.4.4.2
560 34	Waltenhofen's pendulum	1	
342 07	Clamp with knife-edge bearings	1	
562 11	U-core with yoke	1	1
562 13	Coil, 250 turns	2	1
560 31	Bored pole pieces, pair	1	
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1	
300 02	Stand base, V-shaped, small	1	
301 01	Leybold multiclamp	1	2
300 51	Stand rod, right-angled	1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	
501 28	Connecting lead, 32 A, 50 cm, black	3	
560 32	Rotatable aluminium disc		1
562 10	Yoke		1
562 15	Coil, 1 000 turns		1
562 18	Coil, 50 turns, extra-low voltage		2
562 34	Coil holder		1
510 22	Large horseshoe magnet with yoke		1
521 39	Variable extra-low voltage transformer		1
537 32	Rheostat, 10 ohms		1
531 120	Multimeter LDanalog 20		2
313 07	Hand-held stop watch I, mechanical		1
300 01	Stand base, V-shaped, large		1
300 41	Stand rod, 25 cm, 12 mm Ø		1
501 33	Connecting lead, 32 A, 100 cm, black		9

When a metal disk is moved into a magnetic field, eddy currents are produced in the disk. The eddy currents generate a magnetic field which interacts with the inducing field to resist the motion of the disk. The energy of the eddy currents, which is liberated by the Joule effect, results from the mechanical work which must be performed to overcome the magnetic force.

In the experiment P3.4.4.1, the occurrence and suppression of eddy currents is demonstrated using Waltenhofen's pendulum. The aluminum plate swings between the pole pieces of a strong electromagnet. As soon as the magnetic field is switched on, the pendulum is arrested when it enters the field. The pendulum oscillations of a slitted plate, on the other hand, are only slightly attenuated, as only weak eddy currents can form.

The experiment P3.4.4.2 examines the workings of an alternating current meter. In principle, the AC meter functions much like an asynchronous motor with squirrel-cage rotor. A rotating aluminium disk is mounted in the air gap between the poles of two magnet systems. The current to be measured flows through the bottom magnet system, and the voltage to be measured is applied to the top magnet system. A moving magnetic field is formed which generates eddy currents in the aluminum disk. The moving magnetic field and the eddy currents produce an asynchronous angular momentum

$$N_1 \propto P$$

proportional to the electrical power  $P$  to be measured. The angular momentum accelerates the aluminum disk until it attains equilibrium with its counter-torque

$$N_2 \propto \omega$$

$\omega$ : angular velocity of disk

generated by an additional permanent magnet embedded in the turning disk. Consequently, at equilibrium

$$N_1 = N_2$$

the angular velocity of the disk is proportional to the electrical power  $P$ .



**P3.4.5**  
**TRANSFORMER**

- P3.4.5.1  
Voltage and current transformation with a transformer
- P3.4.5.2  
Voltage transformation with a transformer under load
- P3.4.5.3  
Recording the voltage and current of a transformer under load as a function of time

Voltage and current transformation with a transformer (P3.4.5.1)

Cat. No.	Description	P3.4.5.1	P3.4.5.2	P3.4.5.3 (b)
562 801	Transformer for students' experiments	1	1	1
531 120	Multimeter LDanalog 20	2	2	
521 35	Variable extra-low voltage transformer S	1	1	
500 444	Connecting lead 19 A, 100 cm, black	6	7	6
537 34	Rheostat, 100 ohms	1	1	
459 23	Acrylic glass screen on rod	1		
514 72ET5	Shaker for iron filings, set of 5	1		
514 73	Iron filings	1		
524 013	Sensor-CASSY 2			1
524 011USB	Power-CASSY USB			1
524 220	CASSY Lab 2			1
500 414	Connecting lead, 19 A, 25 cm, black			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1

Regardless of the physical design of the transformer, the voltage transformation of a transformer without load is determined by the ratio of the respective number of turns

$$\frac{U_2}{U_1} = \frac{N_2}{N_1} \text{ (when } I_2 = 0 \text{)}$$

The current transformation in short-circuit operation is inversely proportional to the ratio of the number of turns

$$\frac{I_2}{I_1} = \frac{N_1}{N_2} \text{ (when } U_2 = 0 \text{)}$$

The behavior of the transformer under load, on the other hand, depends on its particular physical design. This fact can be demonstrated using the transformer for students' experiments.

The aim of the experiment P3.4.5.1 is to measure the voltage transformation of a transformer without load and the current transformation of a transformer in short-circuit mode. At the same time, the difference between an isolating transformer and an autotransformer is demonstrated.

The experiment P3.4.5.2 examines the ratio between primary and secondary voltage in a "hard" and a "soft" transformer under load. In both cases, the lines of magnetic flux of the transformer are revealed using iron filings on a glass plate placed on top of the transformer.

In the experiment P3.4.5.3, the primary and secondary voltages and the primary and secondary currents of a transformer under load are recorded as time-dependent quantities using the computer-based CASSY measuring system. The CASSY software determines the phase relationships between the four quantities directly and additionally calculates the time-dependent power values of the primary and secondary circuits.


### P3.4.5

#### TRANSFORMER

P3.4.5.4  
Power transmission of a transformer

P3.4.5.5  
Experiments with high currents

P3.4.5.6  
High-voltage experiments with  
a two-pronged lightning rod



Power transmission of a transformer (P3.4.5.4\_b)

Cat. No.	Description	P3.4.5.4 (a)	P3.4.5.4 (b)	P3.4.5.5	P3.4.5.6
562 11	U-core with yoke	1	1	1	1
562 121	Clamping device with spring clip	1	1	1	1
562 13	Coil, 250 turns	2	2		
524 013	Sensor-CASSY 2	2	1		
524 220	CASSY Lab 2	1	1		
521 35	Variable extra-low voltage transformer S	1			
537 34	Rheostat, 100 ohms	1	1		
500 414	Connecting lead, 19 A, 25 cm, black	2	1		2
500 444	Connecting lead 19 A, 100 cm, black	8	6		
524 011USB	Power-CASSY USB		1		
562 21	Mains coil, 500 turns			1	1
562 20	Ring-shaped melting ladle			1	
562 32	Melting ring			1	
562 19	Coil, 5 turns			1	
562 31	Sheet-metal strips, set of 5			1	
562 17	Coil, 23 000 turns				1
540 52	Demonstration insulator				2
300 11	Saddle base				2
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1		

As an alternative to the transformer for students' experiments, the demountable transformer with a full range of coils is available which simply slide over the arms of the U-core, making them easily interchangeable. The experiments described for the transformer for students' experiments (P3.4.5.1-3) can of course be performed just as effectively using the demountable transformer, as well as a number of additional experiments.

The experiment P3.4.5.4 examines the power transmission of a transformer. Here, the RMS values of the primary and secondary voltage and the primary and secondary current are measured on a variable load resistor  $R = 0 - 100 \Omega$  using the computer-based CASSY measuring system. The phase shift between the voltage and current on the primary and secondary sides is determined at the same time. In the evaluation, the primary power  $P_1$ , the secondary power  $P_2$  and the efficiency

$$\eta = \frac{P_2}{P_1}$$

are calculated and displayed in a graph as a function of the load resistance  $R$ .

In the experiment P3.4.5.5, a transformer is assembled in which the primary side with 500 turns is connected directly to the mains voltage. In a melting ring with one turn or a welding coil with five turns on the secondary side, extremely high currents of up to 100 A can flow, sufficient to melt metals or spot-weld wires.

In the experiment P3.4.5.6, a transformer is assembled in which the primary side with 500 turns is connected directly to the mains voltage. Using a secondary coil with 23,000 turns, high voltages of up to 10 kV are generated, which can be used to produce electric arcs in horn-shaped spark electrodes.





### P3.4.6

#### MEASURING THE EARTH'S MAGNETIC FIELD

##### P3.4.6.1

Measuring the earth's magnetic field with a rotating induction coil (earth inductor)

##### P3.4.6.2

Measuring the earth's magnetic field with a rotating induction coil and wireless CASSY

Measuring the earth's magnetic field with a rotating induction coil and wireless CASSY (P3.4.6.2)

Cat. No.	Description	P3.4.6.1	P3.4.6.2
555 604	Pair of Helmholtz coils	1	1
524 013	Sensor-CASSY 2	1	
524 220	CASSY Lab 2	1	1
524 040	$\mu$ V box	1	
501 35	Connecting lead, 32 A, 200 cm, red	1	
501 36	Connecting lead, 32 A, 200 cm, blue	1	
347 35	Experiment motor	1*	
347 36	Control unit for experiment motor	1*	
348 22	Gyroscope axle		1
524 018	Pocket-CASSY 2 Bluetooth		1
524 019	Rechargeable battery for Pocket-CASSY 2 Bluetooth		1
524 0031	Bluetooth dongle		1
524 0621	UIP sensor S		1
666 615	Universal bosshead		1
501 44	Connecting leads, 19 A, 25 cm, red/blue, pair		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

\* additionally recommended

When a circular induction loop with  $N$  turns and a radius  $R$  rotates in a homogeneous magnetic field  $B$  around its diameter as its axis, it is permeated by a magnetic flux of

$$\Phi(t) = N \cdot \pi \cdot R^2 \cdot n(t) \cdot B$$

$n(t)$ : normal vector of a rotating loop

If the angular velocity  $\omega$  is constant, we can say that

$$\Phi(t) = N \cdot \pi \cdot R^2 \cdot B_{\perp} \cdot \cos \omega t$$


Where  $B_{\perp}$  is the effective component of the magnetic field perpendicular to the axis of rotation. We can determine the magnetic field from the amplitude of the induced voltage

$$U_0 = N \cdot \pi \cdot R^2 \cdot B_{\perp} \cdot \omega$$

To achieve the maximum measuring accuracy, we need to use the largest possible coil.

In the experiment P3.4.6.1 the voltage  $U(t)$  induced in the earth's magnetic field for various axes of rotation is measured using the computer-based CASSY measuring system. The amplitude and frequency of the recorded signals and the respective active component  $B_{\perp}$  are used to calculate the earth's magnetic field. The aim of the evaluation is to determine the total value, the horizontal component and the angle of inclination of the earth's magnetic field.

In the experiment P3.4.6.2 the voltage  $U(t)$  induced in the earth's magnetic field is measured in the rotating system and transmitted wirelessly using the computer-based Pocket-CASSY measuring system. The amplitude and frequency of the recorded signals are used to calculate the horizontal component of the earth's magnetic field.



Induced voltage in the rotating coil

### P3.5.1

#### BASIC EXPERIMENTS ON ELECTRICAL MACHINES

P3.5.1.1  
Investigating the interactions of forces of rotors and stators

P3.5.1.2  
Simple induction experiments with electromagnetic rotors and stators



Simple induction experiments with electromagnetic rotors and stators (P3.5.1.2)

Cat. No.	Description	P3.5.1.1	P3.5.1.2
563 480	ELM basic set	1	1
727 81	Basic machine unit	1	1
560 61	Magnet model, cubical	1	
521 485	AC/DC power supply, 0...12 V/3 A	1	
500 422	Connecting lead 19 A, 50 cm, blue	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1
531 120	Multimeter LDanalog 20		1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair		1

The term "electrical machines" is used to refer to both motors and generators. Both devices consist of a stationary stator and a rotating armature or rotor. The function of the motors is due to the interaction of the forces arising through the presence of a current-carrying conductor in a magnetic field, and that of the generators is based on induction in a conductor loop moving within a magnetic field.

The action of forces between the magnetic field and the conductor is demonstrated in the experiment P3.5.1.1 using permanent and electromagnetic rotors and stators. A magnet model is used to represent the magnetic fields.

The object of the experiment P3.5.1.2 is to carry out qualitative measurements on electromagnetic induction in electromagnetic rotors and stators.



### P3.5.2

#### ELECTRIC GENERATORS

- P3.5.2.1  
Generating AC voltage using a revolving-field generator and a stationary-field generator
- P3.5.2.2  
Generating DC voltage using a stationary-field generator
- P3.5.2.3  
Generating AC voltage using a generator with electromagnetic rotating pole (power-plant generator)
- P3.5.2.4  
Generating voltage with an AC-DC generator (generator with electromagnetic stationary pole)
- P3.5.2.5  
Generating voltage using self-exciting generators

Generating AC voltage using a revolving-field generator and a stationary-field generator (P3.5.2.1\_b)

Cat. No.	Description	P3.5.2.1 (b)	P3.5.2.2 (b)	P3.5.2.3-4 (b)	P3.5.2.5 (b)
563 480	ELM basic set	1	1	1	1
727 81	Basic machine unit	1	1	1	1
563 303	ELM hand-cranked gear	1	1	1	1
301 300	Demonstration panel frame	1	1	1	1
531 120	Multimeter LDanalog 20	1	1	2	1
531 282	Multimeter Metrahit Pro	1			
537 36	Rheostat, 1000 ohms	1			
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1		1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	1	2	2
563 23	ELM three-pole rotor		1*		1
575 214	Oscilloscope 30 MHz, two-channel, analogous		1*		
575 24	Screened cable, BNC/4 mm		1*		
521 485	AC/DC power supply, 0...12 V/3 A			1	1
500 422	Connecting lead 19 A, 50 cm, blue				1

\* additionally recommended

Electric generators exploit the principle of electromagnetic induction discovered by Faraday to convert mechanical into electrical energy. We distinguish between revolving-armature generators (excitation of the magnetic field in the stator, induction in the rotor) and revolving-field generators (excitation of the magnetic field in the rotor, induction in the stator).

Both types of generators are assembled in the experiment P3.5.2.1 using permanent magnets. The induced AC voltage  $U$  is measured as a function of the speed  $f$  of the rotor. Also, the electrical power  $P$  produced at a fixed speed is determined as a function of the load resistance  $R$ .

The experiment P3.5.2.2 demonstrates the use of a commutator to rectify the AC voltage generated in the rotor of a rotating-armature generator. The number of rectified half-waves per rotor revolution increases when the two-pole rotor is replaced with a three-pole rotor.

The experiments P3.5.2.3 and P3.5.2.4 investigate generators which use electromagnets instead of permanent magnets. Here, the induced voltage depends on the excitation current of the magnetic field. The excitation current can be used to vary the generated power without changing the speed of the rotor or the frequency of the AC voltage. This principle is used in power-plant generators. In the AC/DC generator, the voltage can also be tapped via the commutator in rectified form.

The experiment P3.5.2.5 examines generators in which the magnetic field of the stator is amplified by the generator current by means of self-excitation. The stator and rotor windings are conductively connected with each other. We distinguish between serieswound generators, in which the rotor, stator and load are all connected in series, and shuntwound generators, in which the stator and the load are connected in parallel to the rotor.

### P3.5.3

#### ELECTRIC MOTORS

##### P3.5.3.1

Experiments on DC motor with two-pole rotor

##### P3.5.3.2

Experiments on DC motor with three-pole rotor

##### P3.5.3.3

Experiments with a universal motor in series and shunt connection

##### P3.5.3.4

Assembling an AC synchronous motor



Experiments on DC motor with two-pole rotor (P3.5.3.1\_b)

Cat. No.	Description	P3.5.3.1 (b)	P3.5.3.2 (b)	P3.5.3.3 (b)	P3.5.3.4 (b)
563 480	ELM basic set	1	1	1	1
727 81	Basic machine unit	1	1	1	1
301 300	Demonstration panel frame	1	1	1	1
531 120	Multimeter LDanalog 20	2	2	2	
521 35	Variable extra-low voltage transformer S	1	1	1	1
451 281	Stroboscope	1	1	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	1	2	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	2	2	2
563 23	ELM three-pole rotor		1	1*	
314 151	Precision dynamometer, 2 N		1	1	
314 161	Precision dynamometer, 5 N		1	1	
309 50	Demonstration cord		1	1	
666 470	Holder with clamp, height-adjustable, CPS		1	1	
300 41	Stand rod, 25 cm, 12 mm Ø		1	1	
563 303	ELM hand-cranked gear				1
726 50	Plug-in board, 297 mm x 300 mm, STE				1
579 13	Toggle switch, STE 2/19				1
579 06	Lamp holder, E10, top, STE 2/19				1
505 181	Bulbs, 24 V/3 W, E10, set of 5				1

\* additionally recommended

Electric motors exploit the force acting on current-carrying conductors in magnetic fields to convert electrical energy into mechanical energy. We distinguish between asynchronous motors, in which the rotor is supplied with AC or DC voltage via a commutator, and synchronous motors, which have no commutator, and whose frequencies are synchronized with the frequency of the applied voltage.

The experiment P3.5.3.1 investigates the basic function of an electric motor with commutator. The motor is assembled using a permanent magnet as stator and a two-pole rotor. The polarity of the rotor current determines the direction in which the rotor turns. This experiment measures the relationship between the applied voltage  $U$  and the no-load speed  $f_0$  as well as, at a fixed speed, the current  $I$  consumed as a function of the load-dependent speed  $f$ .

The use of the three-pole rotor is the object of the experiment P3.5.3.2. The rotor starts turning automatically, as an angular momentum (torque) acts on the rotor for any position in the magnetic field. To record the torque curve  $M(f)$ , the speed  $f$  of the rotor is recorded as a function of a counter-torque  $M$ . In addition, the mechanical power produced is compared with the electrical power consumed.

The experiment P3.5.3.3 takes a look at the so-called universal motor, in which the stator and rotor fields are electrically excited. The stator and rotor coils are connected in series ("serieswound") or in parallel ("shunt-wound") to a common voltage source. This motor can be driven both with DC and AC voltage, as the torque acting on the rotor remains unchanged when the polarity is reversed. The torque curve  $M(f)$  is recorded for both circuits. The experiment shows that the speed of the shuntwound motor is less dependent on the load than that of the series-wound motor.

In the experiment P3.5.3.4, the rotor coil of the AC synchronous motor is synchronized with the frequency of the applied voltage using a hand crank, so that the rotor subsequently continues running by itself.



### P3.5.4

#### THREE-PHASE MACHINES

##### P3.5.4.1

Experiments with a three-phase revolving-armature generator

##### P3.5.4.2

Experiments with a three-phase revolving-field generator

##### P3.5.4.3

Comparing star and delta connections on a three-phase generator

##### P3.5.4.4

Assembling synchronous and asynchronous three-phase motors

Experiments with a three-phase revolving-armature generator (P3.5.4.1\_b)

Cat. No.	Description	P3.5.4.1 (b)	P3.5.4.2 (b)	P3.5.4.3 (b)	P3.5.4.4 (b)
563 480	ELM basic set	1	1	1	1
563 481	ELM supplementary set	1	1	1	1
727 81	Basic machine unit	1	1	1	1
563 303	ELM hand-cranked gear	1	1	1	
301 300	Demonstration panel frame	1	1	1	1
531 120	Multimeter LDanalog 20	3	3	2	1
501 451	Connecting leads, 19 A, 50 cm, black, pair	3	4	6	2
575 214	Oscilloscope 30 MHz, two-channel, analogous	1*	1*		
575 24	Screened cable, BNC/4 mm	2*	2*		
313 07	Hand-held stop watch I, mechanical	1*	1*		
521 485	AC/DC power supply, 0...12 V/3 A		1	1	
726 50	Plug-in board, 297 mm x 300 mm, STE			1	
579 06	Lamp holder, E10, top, STE 2/19			3	
505 14	Bulbs, 6 V/3 W, E10, set of 10			1	
501 48	Bridging plugs, STE 2/19, set of 10			1	
500 414	Connecting lead, 19 A, 25 cm, black			3	3
563 12	ELM squirrel-cage rotor				1
521 291	Three-phase extra-low voltage transformer				1

\* additionally recommended

In the real world, power is supplied mainly through the generation of three-phase AC, usually referred to simply as "three-phase current". Consequently, three-phase generators and motors are extremely significant in actual practice. In principle, their function is analogous to that of AC machines. As with AC machines, we differentiate between revolving-armature and revolving-field generators, and between asynchronous and synchronous motors.

The simplest configuration for generating three-phase current, a revolving-armature generator which rotates in a permanent magnetic field, is assembled in the experiment P3.5.4.1 using a threepole rotor.

The experiment P3.5.4.2 examines the more common revolving-field generator, in which the magnetic field of the rotor in the stator coils is induced by phase-shifted AC voltages. In both cases, instruments for measuring current and voltage, and for observing the phase shift for a slowly turning rotor, are connected between two taps. For faster rotor speeds, the phase shift is measured using an oscilloscope.

In the experiment P3.5.4.3, loads are connected to the three-phase generator in star and delta configuration. In the star configuration, the relationship

$$\frac{U_{aa}}{U_{a0}} = \sqrt{3}$$

is verified for the voltages  $U_{aa}$  between any two outer conductors as well as  $U_{a0}$  between the outer and neutral conductors. For the currents  $I_1$  flowing to the loads and the currents  $I_2$  flowing through the generator coils in delta configuration, the result is

$$\frac{I_1}{I_2} = \sqrt{3}$$

The experiment P3.5.4.4 examines the behavior of asynchronous and synchronous machines when the direction of rotation is reversed.

### P3.6.1

#### CIRCUIT WITH CAPACITOR

##### P3.6.1.1


Charging and discharging a capacitor when switching DC on and off

##### P3.6.1.2

Determining the capacitive reactance of a capacitor in an AC circuit


##### P3.6.1.3

Charging and discharging a capacitor when switching DC on and off - Measuring with a mutlimeter



Charging and discharging a capacitor when switching DC on and off (P3.6.1.1\_a)

Cat. No.	Description	P3.6.1.1 (a)	P3.6.1.1 (b)	P3.6.1.2 (a)	P3.6.1.2 (b)	P3.6.1.3
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1
578 15	Capacitor, 1 $\mu$ F, STE 2/19	3	3	3	3	
577 40	Resistor, 470 $\Omega$ , STE 2/19	1	1			
577 44	Resistor, 1 k $\Omega$ , STE 2/19	1	1			
577 48	Resistor, 2.2 k $\Omega$ , STE 2/19	1	1			
522 621	Function generator S 12	1	1	1	1	
575 214	Oscilloscope 30 MHz, two-channel, analogous	1		1		
575 24	Screened cable, BNC/4 mm	2		2		
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	3	1	3	
524 005	Mobile-CASSY 2		1		1	
577 19	Resistor, 1 $\Omega$ , STE 2/19			1		
577 20	Resistor, 10 $\Omega$ , STE 2/19			1	1	
578 16	Capacitor, 4.7 $\mu$ F, STE 2/19					1
578 12	Capacitor, 10 $\mu$ F, STE 2/50					1
577 76	Resistor, 1 M $\Omega$ , STE 2/19					1
582 81	Change-over switch, STE 4/50					1
501 48	Bridging plugs, STE 2/19, set of 10					1
521 485	AC/DC power supply, 0...12 V/3 A					1
531 120	Multimeter LDanalog 20					1
313 07	Hand-held stop watch I, mechanical					1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair					2



Schematic circuit diagram

To investigate the behavior of capacitors in DC and AC circuits, the voltage  $U_C$  at a capacitor is measured using a two-channel oscilloscope, and the current  $I_C$  through the capacitor is additionally calculated from the voltage drop across a resistor  $R$  connected in series. Alternatively a Mobile-CASSY 2 is used. The circuits for conducting these measurements are assembled on a plug-in board using the STE plug-in system. A function generator is used as a voltage source with variable amplitude and variable frequency.

In the experiment P3.6.1.1, the function generator generates periodic square-wave signals which simulate switching a DC voltage on and off. The square-wave signals are displayed on channel I of the oscilloscope, and the capacitor voltage or capacitor current is displayed on oscilloscope channel II. Alternatively current and voltage are recorded by Mobile-CASSY. The aim of the experiment is to determine the time constant

$$\tau = R \cdot C$$

for various capacitances  $C$  from the exponential curve of the respective charging or discharge current  $I_C$ .

In the experiment P3.6.1.2, an AC voltage with the amplitude  $U_0$  and the frequency  $f$  is applied to a capacitor. The voltage  $U_C(t)$  and the current  $I_C(t)$  are displayed simultaneously on the oscilloscope. Alternatively current and voltage are recorded by Mobile-CASSY. The experiment shows that in this circuit the current leads the voltage by  $90^\circ$ . In addition, the proportionality between the voltage amplitude  $U_0$  and the current amplitude  $I_0$  is confirmed, and for the proportionality constant


$$Z_C = \frac{U_0}{I_0}$$

the relationship

$$Z_C = -\frac{1}{2\pi f \cdot C}$$

is revealed.

In experiment P3.6.1.3 a large capacitor is charged and discharged with a long time constant and the process can be recorded with a multimeter and a stopwatch.



### P3.6.2

#### CIRCUIT WITH COIL

P3.6.2.1  
Measuring the current in a coil when switching DC on and off

P3.6.2.2  
Determining the inductive reactance of a coil in an AC circuit

Measuring the current in a coil when switching DC on and off (P3.6.2.1\_a)

Cat. No.	Description	P3.6.2.1 (a)	P3.6.2.1 (b)	P3.6.2.2 (a)	P3.6.2.2 (b)
576 74	Plug-in board, DIN A4, STE	1	1	1	1
590 84	Coil, 1000 turns, STE 2/50	2	2	2	2
577 19	Resistor, 1 Ω, STE 2/19	1		1	
577 20	Resistor, 10 Ω, STE 2/19	1	1	1	1
577 24	Resistor, 22 Ω, STE 2/19	1	1		
577 28	Resistor, 47 Ω, STE 2/19	1	1		
501 48	Bridging plugs, STE 2/19, set of 10	1	1	1	1
522 621	Function generator S 12	1	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1		1	
575 24	Screened cable, BNC/4 mm	2		2	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	3	1	3
524 005	Mobile-CASSY 2		1		1

To investigate the behavior of coils in DC and AC circuits, the voltage  $U_L$  at a coil is measured using a two-channel oscilloscope, and the current  $I_L$  through the coil is additionally calculated from the voltage drop across a resistor  $R$  connected in series. Alternatively a Mobile-CASSY 2 can be used. The circuits for conducting these measurements are assembled on a plug-in board using the STE plug-in system for electricity/electronics. A function generator is used as a voltage source with variable amplitude and variable frequency.

In the experiment P3.6.2.1, the function generator generates periodic square-wave signals which simulate switching a DC voltage on and off. The square-wave signals are displayed on channel I of the oscilloscope, and the coil voltage or coil current is displayed on oscilloscope channel II. Alternatively current and voltage are recorded by Mobile-CASSY. The aim of the experiment is to determine the time constant

$$\tau = \frac{L}{R}$$

for different inductances  $L$  from the exponential curve of the coil voltage  $U_L$ .


In the experiment P3.6.2.2, an AC voltage with the amplitude  $U_0$  and the frequency  $f$  is applied to a coil. The voltage  $U_L(t)$  and the current  $I_L(t)$  are displayed simultaneously on the oscilloscope. Alternatively current and voltage are recorded by Mobile-CASSY. The experiment shows that in this circuit the current lags behind the voltage by  $90^\circ$ . In addition, the proportionality between the voltage amplitude  $U_0$  and the current amplitude  $I_0$  is confirmed, and, for the proportionality constant

$$Z_L = \frac{U_0}{I_0}$$

the relationship

$$Z_L = 2\pi f \cdot L$$

is revealed.



Schematic circuit diagram


### P3.6.3

#### IMPEDANCES

**P3.6.3.1**  
Determining the impedance in circuits with capacitors and ohmic resistors

**P3.6.3.2**  
Determining the impedance in circuits with coils and ohmic resistors

**P3.6.3.3**  
Determining the impedance in circuits with capacitors and coils



Determining the impedance in circuits with capacitors and coils (P3.6.3.3\_a)

Cat. No.	Description	P3.6.3.1 (a)	P3.6.3.1 (b)	P3.6.3.2 (a)	P3.6.3.2 (b)	P3.6.3.3 (a)	P3.6.3.3 (b)
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1	1
577 19	Resistor, 1 Ω, STE 2/19	1		1			
577 32	Resistor, 100 Ω, STE 2/19	1	1	1	1	1	1
578 12	Capacitor, 10 μF, STE 2/50	1	1				
578 15	Capacitor, 1 μF, STE 2/19	1	1			1	1
578 31	Capacitor, 0.1 μF, STE 2/19	1	1				
522 621	Function generator S 12	1	1	1	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1		1		1	
575 24	Screened cable, BNC/4 mm	2		2		2	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	3	1	3	1	3
<b>524 005</b>	<b>Mobile-CASSY 2</b>		<b>1</b>		<b>1</b>		<b>1</b>
590 83	Coil, 500 turns, STE 2/50			1	1	1	1
590 84	Coil, 1000 turns, STE 2/50			1	1	1	1
577 20	Resistor, 10 Ω, STE 2/19					1	1
578 16	Capacitor, 4.7 μF, STE 2/19					1	1

The current  $I(t)$  and the voltage  $U(t)$  in an AC circuit are measured as time-dependent quantities using a dual-channel oscilloscope or a Mobile-CASSY 2. A function generator is used as a voltage source with variable amplitude  $U_0$  and variable frequency  $f$ . The measured quantities are then used to determine the absolute value of the total impedance

$$Z = \frac{U_0}{I_0}$$

and the phase shift  $\varphi$  between the current and the voltage.

A resistor  $R$  is combined with a capacitor  $C$  in the experiment P3.6.3.1, and an inductor  $L$  in the experiment P3.6.3.2. These experiments confirm the relationship

$$Z_s = \sqrt{R^2 + Z_l^2} \text{ and } \tan \varphi_s = \frac{Z_l}{R}$$

$$\text{with } Z_l = -\frac{1}{2\pi f \cdot C} \text{ resp. } Z_l = 2\pi f \cdot L$$

for series connection and

$$\frac{1}{Z_p} = \sqrt{\frac{1}{R^2} + \frac{1}{Z_l^2}} \text{ and } \tan \varphi_p = \frac{R}{Z_l}$$

for parallel connection.

The experiment P3.6.3.3 examines the oscillator circuit as the series and parallel connection of capacitance and inductance. The total impedance of the series circuit

$$Z_s = 2\pi f \cdot L - \frac{1}{2\pi f \cdot C}$$

disappears at the resonance frequency


$$f_r = \frac{1}{2\pi \cdot \sqrt{LC}}$$

i.e. at a given current  $I$  the total voltage  $U$  at the capacitor and the coil is zero, because the individual voltages  $U_C$  and  $U_L$  are equal and opposite. For parallel connection, we can say

$$\frac{1}{Z_p} = \frac{1}{2\pi f \cdot L} - 2\pi f \cdot C$$

At the resonance frequency, the impedance of this circuit is infinitely great; in other words, at a given voltage  $U$  the total current  $I$  in the incumingine is zero, as the two individual currents  $I_C$  and  $I_L$  are equal and opposed.





### P3.6.4 MEASURING-BRIDGE CIRCUITS

P3.6.4.1  
Determining capacitive reactance  
with a Wien measuring bridge

P3.6.4.2  
Determining inductive reactance  
with a Maxwell measuring bridge

Determining capacitive reactance with a Wien measuring bridge (P3.6.4.1\_b)

Cat. No.	Description	P3.6.4.1 (b)	P3.6.4.1 (c)	P3.6.4.2 (b)	P3.6.4.2 (c)
576 74	Plug-in board, DIN A4, STE	1	1	1	1
577 32	Resistor, 100 Ω, STE 2/19	1	1	1	1
577 93	Potentiometer, 1 kΩ, 10-turn, STE 4/50	1	1	2	2
578 15	Capacitor, 1 μF, STE 2/19	1	1		
578 16	Capacitor, 4.7 μF, STE 2/19	1	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analog	1		1	
575 24	Screened cable, BNC/4 mm	1		1	
522 621	Function generator S 12	1	1	1	1
501 48	Bridging plugs, STE 2/19, set of 10	1	1	2	2
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	2	1	2
<b>524 005</b>	<b>Mobile-CASSY 2</b>		<b>1</b>		<b>1</b>
590 83	Coil, 500 turns, STE 2/50			1	1
590 84	Coil, 1000 turns, STE 2/50			1	1

The Wheatstone measuring bridge is one of the most effective means of measuring ohmic resistance in DC and AC circuits. Capacitive and inductive reactance can also be determined by means of analogous circuits. These measuring bridges consist of four passive bridge arms which are connected to form a rectangle, an indicator arm with a null indicator and a supply arm with the voltage source. Inserting variable elements in the bridge arm compensates the current in the indicator arm to zero. Then, for the component resistance values, the fundamental compensation condition

$$Z_1 = Z_2 \cdot \frac{Z_3}{Z_4}$$

applies, from which the measurement quantity  $Z_1$  is calculated.

The experiment P3.6.4.1 investigates the principle of a Wien measuring bridge for measuring a capacitive reactance  $Z_1$ . In this configuration,  $Z_2$  is a fixed capacitive reactance,  $Z_3$  is a fixed ohmic resistance and  $Z_4$  is a variable ohmic resistance. For zero compensation, the following applies regardless of the frequency of the AC voltage:


$$\frac{1}{C_1} = \frac{1}{C_2} \cdot \frac{R_3}{R_4}$$

An oscilloscope or Mobile-CASSY can be used as a zero indicator.

In the experiment P3.6.4.2, a Maxwell measuring bridge is assembled to determine the inductive reactance  $Z_1$ . As the resistive component of  $Z_1$  is also to be compensated, this circuit is somewhat more complicated. Here,  $Z_2$  is a variable ohmic resistance,  $Z_3$  is a fixed ohmic resistance and  $Z_4$  is a parallel connection consisting of a capacitive reactance and a variable ohmic resistor. For the purely inductive component, the following applies with respect to zero compensation:

$$2\pi f \cdot L_1 = R_2 \cdot R_3 \cdot 2\pi f \cdot C_4$$

$f$ : AC voltage frequency



### P3.6.5

#### MEASURING AC VOLTAGES AND AC CURRENTS

##### P3.6.5.1

Frequency response and form factor  
of a multimeter



Frequency response and form factor of a multimeter (P3.6.5.1)

Cat. No.	Description	P3.6.5.1
531 120	Multimeter LDanalog 20	2
536 131	Measuring resistor, 100 $\Omega$	1
536 211	Measuring resistor, 10 M $\Omega$	1
522 621	Function generator S 12	1
500 424	Connecting lead 19 A, 50 cm, black	5
575 214	Oscilloscope 30 MHz, two-channel, analogous	1*
575 24	Screened cable, BNC/4 mm	1*

\* additionally recommended

When measuring voltages and currents in AC circuits at higher frequencies, the indicator of the meter no longer responds in proportion to the voltage or current amplitude. The ratio of the reading value to the true value as a function of frequency is referred to as the "frequency response". When measuring AC voltages or currents in which the shape of the signal deviates from the sinusoidal oscillation, a further problem occurs. Depending on the signal form, the meter will display different current and voltage values at the same frequency and amplitude. This phenomenon is described by the wave form factor.

The experiment P3.6.5.1 determines the frequency response and wave form factor of a multimeter. Signals of a fixed amplitude and varying frequencies are generated using a function generator and measured using the multimeter.



### P3.6.6

#### ELECTRICAL WORK AND POWER

##### P3.6.6.1

Determining the heating power of an ohmic load in an AC circuit as a function of the applied voltage

##### P3.6.6.2

Determining the electric work of an immersion heater using an AC power meter

Determining the electric work of an immersion heater using an AC power meter (P3.6.6.2)

Cat. No.	Description	P3.6.6.1	P3.6.6.2
590 50	Lid with heater	1	
384 52	Aluminium calorimeter	1	
313 07	Hand-held stop watch I, mechanical	1	1
382 34	Thermometer, -10...+110 °C/0.2 K	1	1
531 120	Multimeter LDanalog 20	1	1
531 130	Multimeter LDanalog 30	1	1
521 35	Variable extra-low voltage transformer S	1	
590 06	Plastic beaker	1	1
501 23	Connecting lead, 32 A, 25 cm, black	4	
501 28	Connecting lead, 32 A, 50 cm, black	2	
560 331	Alternating current meter		1
301 339	Stand bases, pair		1
303 25	Immersion heater		1
500 624	Safety connecting lead, 50 cm, black		4

The relationship between the power  $P$  at an ohmic resistance  $R$  and the applied voltage  $U$  can be expressed with the relationship

$$P = \frac{U^2}{R}$$

The same applies for AC voltage when  $P$  is the power averaged over time and  $U$  is replaced by the RMS value

$$U_{\text{rms}} = \frac{U_0}{\sqrt{2}}$$

$U_0$ : amplitude of AC voltage

The relationship

$$P = U \cdot I$$

can also be applied to ohmic resistors in AC circuits when the direct current  $I$  is replaced by the RMS value of the AC

$$I_{\text{rms}} = \frac{I_0}{\sqrt{2}}$$

$I_0$ : amplitude of AC

In the experiment P3.6.6.1, the electrical power of an immersion heater for extra-low voltage is determined from the Joule heat emitted per unit of time and compared with the applied voltage  $U_{\text{rms}}$ . This experiment confirms the relationship

$$P \propto U_{\text{rms}}^2$$

In the experiment P3.6.6.2, an AC power meter is used to determine the electrical work  $W$  which must be performed to produce one liter of hot water using an immersion heater. For comparison purposes, the voltage  $U_{\text{rms}}$ , the current  $I_{\text{rms}}$  and the heating time  $t$  are measured and the relationship

$$W = U_{\text{rms}} \cdot I_{\text{rms}} \cdot t$$

is verified.

### P3.6.6

#### ELECTRICAL WORK AND POWER

##### P3.6.6.3


Quantitative comparison of DC power and AC power in an incandescent lamp

##### P3.6.6.4

Determining the crest factors of various AC signal forms

##### P3.6.6.5

Determining the active and reactive power in AC circuits



Determining the active and reactive power in AC circuits (P3.6.6.5)

Cat. No.	Description	P3.6.6.3	P3.6.6.4 (a)	P3.6.6.5
531 831	Joule and wattmeter	1	1	1
505 14	Bulbs, 6 V/3 W, E10, set of 10	1		
579 06	Lamp holder, E10, top, STE 2/19	2		
576 71	Plug-in board section, STE	2		
521 485	AC/DC power supply, 0...12 V/3 A	1		
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	2	2
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2		2
522 621	Function generator S 12		1	
536 131	Measuring resistor, 100 Ω		1	
575 214	Oscilloscope 30 MHz, two-channel, analogous		1	1
575 24	Screened cable, BNC/4 mm		1	1
521 35	Variable extra-low voltage transformer S			1
537 35	Rheostat, 330 ohms			1
517 021	Capacitor, 40 μF			1
562 11	U-core with yoke			1
562 121	Clamping device with spring clip			1
562 15	Coil, 1 000 turns			1
575 35	Adapter, BNC/4 mm, 2-pole			1
504 45	On-off switch, single pole			1
500 421	Connecting lead 19 A, 50 cm, red			1

The electrical power of a time-dependent voltage  $U(t)$  at any load resistance is also a function of time:

$$P(t) = U(t) \cdot I(t)$$

$I(t)$ : time-dependent current through the load resistor

Thus, for periodic currents and voltages, we generally consider the power averaged over one period  $T$ . This quantity is often referred to as the active power  $P_W$ . It can be measured electronically for any DC or AC voltages using the joule and wattmeter.

In the experiment P3.6.6.3, two identical incandescent light bulbs are operated with the same electrical power. One bulb is operated with DC voltage, the other with AC voltage. The equality of the power values is determined directly using the joule and wattmeter, and additionally by comparing the lamp brightness levels. This equality is reached when the DC voltage equals the RMS value of the AC voltage.

The object of the experiment P3.6.6.4 is to determine the crest factors, i. e. the quotients of the amplitude  $U_0$  and the RMS value  $U_{rms}$  for different AC voltage signal forms generated using a function generator by experimental means. The amplitude is measured using the joule and wattmeter connected to a PC. The RMS value is calculated from the power  $P$  measured at an ohmic resistor  $R$  using the joule and wattmeter according to the formula

$$U_{eff} = \sqrt{P \cdot R}$$


The experiment P3.6.6.5 measures the current  $I_{rms}$  through a given load and the active power  $P_W$  for a fixed AC voltage  $U_{rms}$ . To verify the relationship

$$P_W = U_{rms} \cdot I_{rms} \cdot \cos \varphi$$

the phase shift  $\varphi$  between the voltage and the current is additionally determined using an oscilloscope. This experiment also shows that the active power for a purely inductive or capacitive load is zero, because the phase shift is  $\varphi = 90^\circ$ . The apparent power

$$P_s = U_{rms} \cdot I_{rms}$$

is also referred to as reactive power in this case.



### P3.6.7

#### ELECTROMECHANICAL DEVICES

P3.6.7.1  
Demonstrating the function of a bell


P3.6.7.2  
Demonstrating the function of a relay

Demonstrating the function of a relay (P3.6.7.2\_b)

Cat. No.	Description	P3.6.7.1	P3.6.7.2 (b)
561 071	Bell/relay set	1	1
301 339	Stand bases, pair	1	1
521 210	Transformer, 6/12 V	1	1
579 10	Push button (NO), STE 2/19	1	
500 444	Connecting lead 19 A, 100 cm, black	2	7
579 30	Adjustable contact, STE 2/19		1
579 13	Toggle switch, STE 2/19		1
576 71	Plug-in board section, STE		2
579 06	Lamp holder, E10, top, STE 2/19		2
505 131	Bulbs, 6 V/5 W, E10, set of 10		1

In the experiment P3.6.7.1, an electric bell is assembled using a hammer interrupter (Wagner interrupter). The hammer interrupter consists of an electromagnet and an oscillating armature. In the resting state, the oscillating armature touches a contact, thus switching the electromagnet on. The electromagnet attracts the oscillating armature, which strikes a bell. At the same time, this action interrupts the circuit, and the oscillating armature returns to the resting position.

The experiment P3.6.7.2 demonstrates how a relay functions. A control circuit operates an electromagnet which attracts the armature of the relay. When the electromagnet is switched off, the armature returns to the resting position. When the armature touches a contact, a second circuit is closed, which e.g. supplies power to a lamp. When the contact is configured so that the armature touches it in the resting state, we call this a break contact; the opposite case is termed a make contact.



Function principle of a relay (P3.6.7.2)

### P3.7.1

#### ELECTROMAGNETIC OSCILLATOR CIRCUIT

##### P3.7.1.1

Free electromagnetic oscillations

##### P3.7.1.2

De-damping of electromagnetic oscillations through inductive three-point coupling after Hartley



Free electromagnetic oscillations (P3.7.1.1\_a)

Cat. No.	Description	P3.7.1.1 (a)	P3.7.1.2 (a)
517 011	Coil with high inductance	1	1
517 021	Capacitor, 40 $\mu\text{F}$	1	1
301 339	Stand bases, pair	2	2
501 48	Bridging plugs, STE 2/19, set of 10	1	1
521 45	DC power supply 0... $\pm$ 15 V	1	
531 120	Multimeter LDanalog 20	1	1
313 07	Hand-held stop watch I, mechanical	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	1
576 74	Plug-in board, DIN A4, STE		1
578 76	Transistor, BC 140, NPN, emitter bottom, STE 4/50		1
577 68	Resistor, 100 k $\Omega$ , STE 2/19		1
576 86	Monocell holder STE 2/50		1
685 48ET5	Batteries 1.5 V (D, mono), set of 5		1
579 13	Toggle switch, STE 2/19		1
500 424	Connecting lead 19 A, 50 cm, black		3

Electromagnetic oscillation usually occurs in a frequency range in which the individual oscillations cannot be seen by the naked eye. However, this is not the case in an oscillator circuit consisting of a high-capacity capacitor ( $C = 40 \mu\text{F}$ ) and a high-inductance coil ( $L = 500 \text{ H}$ ). Here, the oscillation period is about 1 s, so that the voltage and current oscillations can be observed directly on a pointer instrument or CASSY.

The experiment P3.7.1.1 investigates the phenomenon of free electromagnetic oscillations. The damping is so low that multiple oscillation periods can be observed and their duration measured with e. g. a stop-clock. In the process, the deviations between the observed oscillation periods and those calculated using Thomson's equation

$$T = 2\pi \cdot \sqrt{L \cdot C}$$

are observed. These deviations can be explained by the current-dependency of the inductance, as the permeability of the iron core of the coil depends on the magnetic field strength.

In the experiment P3.7.1.2, an oscillator circuit after Hartley is used to "de-damp" the electromagnetic oscillations in the circuit, or in other words to compensate the ohmic energy losses in a feedback loop by supplying energy externally. Oscillator circuits of this type are essential components in transmitter and receiver circuits used in radio and television technology. A coil with center tap is used, in which the connection points are connected with the emitter, base and collector of a transistor via AC. The base current controls the collector current synchronously with the oscillation to compensate for energy losses.



## P3.7.2 DECIMETER-RANGE WAVES

P3.7.2.1  
Radiation characteristic and polarization of decimeter waves

P3.7.2.2  
Amplitude modulation of decimeter waves

P3.7.2.4  
Estimating the dielectric constant of water in the decimeter-wave range

Estimating the dielectric constant of water in the decimeter-wave range (P3.7.2.4)

Cat. No.	Description	P3.7.2.1	P3.7.2.2	P3.7.2.4
587 551	UHF wave generator	1	1	1
531 110	Multimeter LDanalog 10	1		
300 11	Saddle base	2	3	1
501 38	Connecting lead, 32 A, 200 cm, black	2		
522 621	Function generator S 12		1	
532 20	AC/DC amplifier 30 W		1	
587 08	Broad-band speaker		1	
575 24	Screened cable, BNC/4 mm		1	
501 33	Connecting lead, 32 A, 100 cm, black		4	
587 54	Dipoles in water tank			1

It is possible to excite electromagnetic oscillations in a straight conductor in a manner analogous to an oscillator circuit. An oscillator of this type emits electromagnetic waves, and their radiated intensity is greatest when the conductor length is equivalent to exactly one half the wavelength (we call this a  $\lambda/2$  dipole). Experiments on this topic are particularly successful with wavelengths in the decimeter range. We can best demonstrate the existence of such decimeter waves using a second dipole which also has the length  $\lambda/2$ , and from which the voltage is applied to an incandescent lamp or (via a high-frequency rectifier) to a measuring instrument.

The experiment P3.7.2.1 investigates the radiation characteristic of a  $\lambda/2$  dipole for decimeter waves. Here, the receiver is aligned parallel to the transmitter and moved around the transmitter. In a second step, the receiver is rotated with respect to the transmitter in order to demonstrate the polarization of the emitted decimeter waves.

The experiment P3.7.2.2 deal with the transmission of audio-frequency signals using amplitude-modulated decimeter waves. In amplitude modulation a decimeter-wave signal

$$E(t) = E_0 \cdot \cos(2\pi \cdot f \cdot t)$$

is modulated through superposing of an audio-frequency signal  $u(t)$  in the form

$$E_{AM}(t) = E_0 \cdot (1 + k_{AM} \cdot u(t)) \cdot \cos(2\pi \cdot f \cdot t)$$

$k_{AM}$ : coupling coefficient

The experiment P3.7.2.4 demonstrates the dielectric nature of water. In water, decimeter waves of the same frequency propagate with a shorter wavelength than in air. Therefore, a receiver dipole tuned for reception of the wavelength in air is no longer adequately tuned when placed in water.

### P3.7.3

#### PROPAGATION OF DECIMETER-RANGE WAVES ALONG LINES

P3.7.3.1  
Determining the current and voltage  
maxima on a Lecher line

P3.7.3.2  
Investigating the current and voltage  
on a Lecher line with a loop dipole



Determining the current and voltage maxima on a Lecher line (P3.7.3.1)


Cat. No.	Description	P3.7.3.1-2
587 551	UHF wave generator	1
587 56	Lecher system with accessories	1
311 77	Steel tape measure, 2 m	1
300 11	Saddle base	3

*E. Lecher* (1890) was the first to suggest using two parallel wires for directional transmission of electromagnetic waves. Using such Lecher lines, as they are known today, electromagnetic waves can be transmitted to any point in space. They are measured along the line as a voltage  $U(x,t)$  propagating as a wave, or as a current  $I(x,t)$ .

In the experiment P3.7.3.1, a Lecher line open at the wire ends and a shorted Lecher line are investigated. The waves are reflected at the ends of the wires, so that standing waves are formed. The current is zero at the open end, while the voltage is zero at the shorted end. The current and voltage are shifted by  $\lambda/4$  with respect to each other, i. e. the wave antinodes of the voltage coincide with the wave nodes of the current. The voltage maxima are located using a probe with an attached incandescent lamp. An induction loop with connected incandescent lamp is used to detect the current maxima. The wavelength  $\lambda$  is determined from the intervals  $d$  between the current maxima or voltage maxima. We can say

$$d = \frac{\lambda}{2}$$

In the experiment P3.7.3.2, a transmitting dipole ( $\lambda/2$  folded dipole) is attached to the end of the Lecher line. Subsequently, it is no longer possible to detect any voltage or current maxima on the Lecher line itself. A current maximum is detectable in the middle of the dipole, and voltage maxima at the dipole ends.



Current and voltage maxima on a Lecher line



### P3.7.4 MICROWAVES

P3.7.4.1  
Directional characteristic and polarization of microwaves in front of a horn antenna

P3.7.4.2  
Absorption of microwaves

P3.7.4.3  
Interference of microwaves

P3.7.4.4  
Diffraction of microwaves

P3.7.4.5  
Refraction of microwaves

P3.7.4.6  
Total reflection of microwaves



Diffraction of microwaves (P3.7.4.4)

Cat. No.	Description	P3.7.4.1-2	P3.7.4.3	P3.7.4.4	P3.7.4.5	P3.7.4.6
737 01	Gunn oscillator	1	1	1	1	1
737 020	Gunn power supply with amplifier	1	1	1	1	1
737 21	Large horn antenna	1	1	1	1	1
737 35	Electric field probe	1	1	1	1	1
688 809	Stand rod 245 mm long, with thread M6	1	1	1	1	1
737 27	Physics microwave accessories I	1	1	1		
531 120	Multimeter LDanalog 20	1	1	1	1	1
300 11	Saddle base	2	2	4	2	1
501 022	BNC cable, 2 m	2	2	2	2	2
501 461	Connecting leads, 19 A, 100 cm, black, pair	1	1	1	1	1
737 390	Set of Microwave Absorbers	1*	1*	1*	1*	1*
737 275	Physics microwave accessories II		1	1	1	1
311 77	Steel tape measure, 2 m		1			
300 02	Stand base, V-shaped, small					1

\* additionally recommended

Microwaves are electromagnetic waves in the wavelength range between 0.1 mm and 100 mm. They are generated e.g. in a cavity resonator, where the frequency is determined by the volume of the cavity resonator. An E-field probe is used to detect the microwaves; this device measures the parallel component of the electric field. The output signal of the probe is proportional to the square of the field strength, and thus to the intensity.

The experiment P3.7.4.1 investigates the orientation and polarization of the microwave field in front of a radiating horn antenna. Here, the field in front of the horn antenna is measured point by point in both the longitudinal and transverse directions using the E-field probe. To determine the polarization, a rotating polarization grating made of thin metal strips is used; in this apparatus, the electric field can only form perpendicular to the metal strips. The polarization grating is set up between the horn antenna and the E-field probe. This experiment shows that the electric field vector of the radiated microwaves is perpendicular to the long side of the horn radiator.

The experiment P3.7.4.2 deals with the absorption of microwaves. Working on the assumption that reflections may be ignored, the absorption in different materials is calculated using both the incident and the transmitted intensity. This experiment reveals a fact which has had a profound impact on modern cooking: microwaves are absorbed particularly intensively by water.

In the experiment P3.7.4.3, standing microwaves are generated by reflection at a metal plate. The intensity, measured at a fixed point between the horn antenna and the metal plate, changes when the metal plate is shifted longitudinally. The distance between two intensity maxima corresponds to one half the wavelength. Inserting a dielectric in the beam path shortens the wavelength.

The experiments P3.7.4.4 and P3.7.4.5 show that many of the properties of microwaves are comparable to those of visible light. The diffraction of microwaves at an edge, a single slit, a double slit and an obstacle are investigated. Additionally, the refraction of microwaves is demonstrated and the validity of Snell's law of refraction is confirmed.

The experiment P3.7.4.6 investigates total reflection of microwaves at media with lower refractive indices. We know from wave mechanics that the reflected wave penetrates about three to four wavelengths deep into the medium with the lower refractive index, before traveling along the boundary surface in the form of surface waves. This is verified in an experiment by placing an absorber (e.g. a hand) on the side of the medium with the lower refractive index close to the boundary surface and observing the decrease in the reflected intensity.

### P3.7.5

#### PROPAGATION OF MICROWAVES ALONG LINES

##### P3.7.5.1

Guiding of microwaves along a Lecher line

##### P3.7.5.2

Qualitative demonstration of guiding of microwaves along a metal waveguide

##### P3.7.5.3

Determining the standing-wave ratio of a rectangular wave-guide for a variable reflection factor



Guiding of microwaves along a Lecher line (P3.7.5.1)

Cat. No.	Description	P3.7.5.1	P3.7.5.2	P3.7.5.3 (a)
737 01	Gunn oscillator	1	1	1
737 020	Gunn power supply with amplifier	1	1	
737 21	Large horn antenna	1	1	
737 35	Electric field probe	1	1	
688 809	Stand rod 245 mm long, with thread M6	1	1	
737 275	Physics microwave accessories II	1		
531 120	Multimeter LDanalog 20	1	1	
300 11	Saddle base	2	1	
501 022	BNC cable, 2 m	2	2	
501 461	Connecting leads, 19 A, 100 cm, black, pair	1	1	
737 390	Set of Microwave Absorbers	1*		
737 27	Physics microwave accessories I		1	
737 021	Gunn Power Supply with SWR Meter			1
737 095	Fixed Attenuator			1
737 111	Slotted Measuring Line			1
737 03	Coax Detector			1
737 09	Variable Attenuator			1
737 14	Waveguide Termination			1
737 10	Moveable Short			1
737 399	Thumb screw M4, set of 10			1
737 15	Support for waveguide components			1
301 21	Stand base MF			2
501 01	BNC cable, 0.25 m			1
501 02	BNC cable, 1 m			2


\* additionally recommended

To minimize transmission losses over long distances, microwaves can also be transmitted along lines. For this application, metal waveguides are most commonly used; Lecher lines, consisting of two parallel wires, are less common.

Despite this, the experiment P3.7.5.1 investigates the guiding of microwaves along a Lecher line. The voltage along the line is measured using the E-field probe. The wavelengths are determined from the spacing of the maxima.

The experiment P3.7.5.2 demonstrates the guiding of microwaves along a hollow metal waveguide. First, the E-field probe is used to verify that the radiated intensity at a position beside the horn antenna is very low. Next, a flexible metal waveguide is set up and bent so that the microwaves are guided to the E-field probe, where they are measured at a greater intensity.

Quantitative investigations on guiding microwaves in a rectangular waveguide are conducted in the experiment P3.7.5.3. Here, standing microwaves are generated by reflection at a shorting plate in a waveguide, and the intensity of these standing waves is measured as a function of the location in a measuring line with movable measuring probe. The wavelength in the waveguide is calculated from the distance between two intensity maxima or minima. A variable attenuator is set up between the measuring line and the short which can be used to attenuate the intensity of the returning wave by a specific factor, and thus vary the standing-wave ratio.



**P3.7.6**  
**DIRECTIONAL CHARACTERISTIC OF DIPOLE RADIATION**

P3.7.6.1  
Directional characteristic of a helix antenna - Recording measured values manually

P3.7.6.2  
Directional characteristic of a Yagi antenna - Recording measured values manually

P3.7.6.3  
Directional characteristic of a helix antenna - Recording measured values with computer

P3.7.6.4  
Directional characteristic of a Yagi antenna - Recording measured values with computer

P3.7.6.5  
Directional characteristic of a Patch antenna - Recording measured values with computer

Directional characteristic of a helix antenna - Recording measured values manually (P3.7.6.1)

Cat. No.	Description	P3.7.6.1	P3.7.6.2	P3.7.6.3	P3.7.6.4	P3.7.6.5
737 440	Helical Antenna Kit	1		1		
737 03	Coax Detector	1		1		1
737 407	Antenna Stand with Amplifier	1	1			
737 020	Gunn power supply with amplifier	1	1			
737 01	Gunn oscillator	1	1	1	1	1
737 21	Large horn antenna	1	1	1	1	1
688 809	Stand rod 245 mm long, with thread M6	2	2			
737 390	Set of Microwave Absorbers	1	1	1	1	1
531 120	Multimeter LDanalog 20	1	1			
300 11	Saddle base	1	1			
501 022	BNC cable, 2 m	1	1			
575 24	Screened cable, BNC/4 mm	1	1			
501 461	Connecting leads, 19 A, 100 cm, black, pair	2	2			
737 415	Set of wire antennas		1		1	
737 405	Rotating Antenna Platform			1	1	1
737 05	PIN Modulator			1	1	1
737 06	Unidirectional Line			1	1	1
737 15	Support for waveguide components			1		1
301 21	Stand base MF			2	2	2
301 26	Stand rod, 25 cm, 10 mm diam.			2	2	2
501 02	BNC cable, 1 m			1		1
737 428	Set Mikrostrip Antennas					1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1	1	1

Directional antennas radiate the greater part of their electromagnetic energy in a particular direction and/or are most sensitive to reception from this direction. All directional antennas require dimensions which are equivalent to multiple wavelengths. In the microwave range, this requirement can be fulfilled with an extremely modest amount of cost and effort. Thus, microwaves are particularly suitable for experiments on the directional characteristics of antennas.

In the experiment P3.7.6.1, the directional characteristic of a helical antenna is recorded. As the microwave signal is excited with a linearly polarizing horn antenna, the rotational orientation of the helical antenna (clockwise or counterclockwise) is irrelevant. The measurement results are represented in the form of a polar diagram, from which the unmistakable directional characteristic of the helical antenna can be clearly seen.

In the experiment P3.7.6.2, a dipole antenna is expanded using parasitic elements to create a Yagi antenna, to improve the directional properties of the dipole arrangement. Here, a total of four shorter elements are placed in front of the dipole as directors, and a slightly longer element placed behind the dipole serves as a reflector. The directional factor of this arrangement is determined from the polar diagram.

In the experiments P3.7.6.3 and P3.7.6.4, the antennas are placed on a turntable which is driven by an electric motor; the angular turntable position is transmitted to a computer. The antennas receive the amplitude-modulated microwave signals, and frequency-selective and phase-selective detection are applied to suppress noise. The received signals are preamplified in the turntable. After filtering and amplification, they are passed on to the computer. For each measurement, the included software displays the receiving power logarithmically in a polar diagram.

In the experiment P3.7.6.5 the directional characteristic of patch antennas is investigated. Depending on the number of patches, the directional characteristics change. The measurement is done on a computer controlled turntable.

# ELECTRICITY

## FREE CHARGE CARRIERS IN A VACUUM

### P3.8.1

#### TUBE DIODE

P3.8.1.1  
Recording the characteristic  
of a tube diode

P3.8.1.2  
Half-wave rectification  
using a tube diode



Recording the characteristic of a tube diode (P3.8.1.1)

Cat. No.	Description	P3.8.1.1	P3.8.1.2
555 610	Demonstration diode	1	1
555 600	Tube stand	1	1
521 65	Tube power supply, 0...500 V	1	1
531 120	Multimeter LDanalog 20	2	
531 130	Multimeter LDanalog 30	1	
500 641	Safety connecting lead, 100 cm, red	5	2
500 642	Safety connecting lead, 100 cm, blue	4	3
536 191	Measuring resistor, 10 kΩ		1
521 40	Variable transformer, 0...250 V		1
575 214	Oscilloscope 30 MHz, two-channel, analogous		1
575 231	Probe 100 MHz, 1:1 / 10:1		1
575 24	Screened cable, BNC/4 mm		1


A tube diode contains two electrodes: a heated cathode, which emits electrons due to thermionic emission, and an anode. A positive potential between the anode and the cathode generates an emission current to the anode, carried by the free electrons. If this potential is too low, the emission current is prevented by the space charge of the emitted electrons, which screen out the electrical field in front of the cathode. When the potential between the anode and the cathode is increased, the isoelectric lines penetrate deeper into the space in front of the cathode, and the emission current increases. This increase of the current with the potential is described by the Schottky-Langmuir law:

$$I \propto U^{\frac{3}{2}}$$

This current increases until the space charge in front of the cathode has been overcome and the saturation value of the emission current has been reached. On the other hand, if the negative potential applied to the anode is sufficient, the electrons cannot flow to the anode and the emission current is zero.

In the experiment P3.8.1.1, the characteristic of a tube diode is recorded, i.e. the emission current is measured as the function of the anode potential. By varying the heating voltage, it can be demonstrated that the saturation current depends on the temperature of the cathode.

The experiment P3.8.1.2 demonstrates half-wave rectification of the AC voltage signal using a tube diode. For this experiment, an AC voltage is applied between the cathode and the anode via an isolating transformer, and the voltage drop is measured at a resistor connected in series. This experiment reveals that the diode blocks when the voltage is reversed.



Anode current  $I_A$  as a function of the anode voltage  $U_A$



### P3.8.2 TUBE TRIODE

P3.8.2.1  
Recording the characteristic field of a tube triode

P3.8.2.2  
Amplifying voltages with a tube triode

Recording the characteristic field of a tube triode (P3.8.2.1)

Cat. No.	Description	P3.8.2.1	P3.8.2.2
555 612	Demonstration triode	1	1
555 600	Tube stand	1	1
521 65	Tube power supply, 0...500 V	1	1
531 120	Multimeter LDanalog 20	2	
531 130	Multimeter LDanalog 30	1	
500 622	Safety connecting lead, 50 cm, blue	1	2
500 641	Safety connecting lead, 100 cm, red	6	3
500 642	Safety connecting lead, 100 cm, blue	4	3
536 251	Measuring resistor, 100 kΩ		1
522 621	Function generator S 12		1
575 214	Oscilloscope 30 MHz, two-channel, analogous		1
575 231	Probe 100 MHz, 1:1 / 10:1		1
575 24	Screened cable, BNC/4 mm		1


In a tube triode, the electrons pass through the mesh of a grid on their way from the cathode to the anode. When a negative voltage  $U_G$  is applied to the grid, the emission current  $I_A$  to the anode is reduced; a positive grid voltage increases the anode current. In other words, the anode current can be controlled by the grid voltage.

The experiment P3.8.2.1 records the family of characteristics of the triode, i.e. the anode current  $I_A$  as a function of the grid voltage  $U_G$  and the anode voltage  $U_A$

The experiment P3.8.2.2 demonstrates how a tube triode can be used as an amplifier. A suitable negative voltage  $U_G$  is used to set the working point of the triode on the characteristic curve  $I_A(U_A)$  so that the characteristic is as linear as possible in the vicinity of the working point. Once this has been set, small changes in the grid voltage  $\delta U_G$  cause a change in the anode voltage  $\delta U_A$  by means of a proportional change in the anode current  $\delta I_A$ . The ratio:

$$V = \left| \frac{\delta U_A}{\delta U_G} \right|$$

is known as the gain.



Characteristic field of a tube triode

# ELECTRICITY

## FREE CHARGE CARRIERS IN A VACUUM

### P3.8.3

#### MALTESE-CROSS TUBE

P3.8.3.1  
Demonstrating the linear propagation of electrons in a field-free space

P3.8.3.2  
Deflection of electrons in an axial magnetic field



*Deflection of electrons in an axial magnetic field (P3.8.3.2)*

Cat. No.	Description	P3.8.3.1	P3.8.3.2
555 620	Maltese cross tube	1	1
555 600	Tube stand	1	1
521 70	High-voltage power supply, 10 kV	1	1
510 48	Magnets, 35 mm Ø, pair	1	
500 611	Safety connecting lead, 25 cm, red	1	1
500 621	Safety connecting lead, 50 cm, red	1	2
500 641	Safety connecting lead, 100 cm, red	1	2
500 642	Safety connecting lead, 100 cm, blue	1	2
500 644	Safety connecting lead, 100 cm, black	2	2
555 604	Pair of Helmholtz coils		1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A		1
500 622	Safety connecting lead, 50 cm, blue		1


In the Maltese cross tube, the electrons are accelerated by the anode to a fluorescent screen, where they can be observed as luminescent phenomena. A Maltese cross is arranged between the anode and the fluorescent screen, and its shadow can be seen on the screen. The Maltese cross has its own separate lead, so that it can be connected to any desired potential.

The experiment P3.8.3.1 confirms the linear propagation of electrons in a field-free space. In this experiment, the Maltese cross is connected to the anode potential and the shadow of the Maltese cross in the electron beam is compared with the light shadow. We can conclude from the observed coincidence of the shadows that electrons propagate in a straight line. The Maltese cross is then disconnected from any potential. The resulting space charges around the Maltese cross give rise to a repulsive potential, so that the image on the fluorescent screen becomes larger.

In the experiment P3.8.3.2 an axial magnetic field is applied using an electromagnet. The shadow cross turns and shrinks as a function of the coil current. When a suitable relationship between the high voltage and the coil current is set, the cross is focused almost to a point, and becomes larger again when the current is increased further. The explanation for this magnetic focusing may be found in the helical path of the electrons in the magnetic field.



*Shadow of the maltese cross on the fluorescent screen*



### P3.8.4

#### PERRIN TUBE

##### P3.8.4.1

Hot-cathode emission in a vacuum: determining the polarity and estimating the specific charge of the emitted charge carriers

##### P3.8.4.2

Generating Lissajou figures through electron deflection in crossed alternating magnetic fields

##### P3.8.4.3

Generating Lissajou figures through electron deflection in parallel alternating electrical and magnetic field

*Hot-cathode emission in a vacuum: determining the polarity and estimating the specific charge of the emitted charge carriers (P3.8.4.1)*

Cat. No.	Description	P3.8.4.1	P3.8.4.2	P3.8.4.3
555 622	Perrin tube	1	1	1
555 600	Tube stand	1	1	1
555 604	Pair of Helmholtz coils	1	1	1
521 70	High-voltage power supply, 10 kV	1	1	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1		
540 091	Electroscope	1		
300 11	Saddle base	1		
501 051	Cable for high voltages, 1.5 m	1		
500 611	Safety connecting lead, 25 cm, red	1	1	1
500 621	Safety connecting lead, 50 cm, red	2	2	2
500 622	Safety connecting lead, 50 cm, blue	1	1	1
500 641	Safety connecting lead, 100 cm, red	4	3	3
500 642	Safety connecting lead, 100 cm, blue	2	3	3
500 644	Safety connecting lead, 100 cm, black	2	2	2
562 14	Coil, 500 turns	1		
521 35	Variable extra-low voltage transformer S		1	
522 621	Function generator S 12		1	1
300 761	Support blocks, set of 6		1	
521 40	Variable transformer, 0...250 V			1

In the Perrin tube, the electrons are accelerated through an anode with pin-hole diaphragm onto a fluorescent screen. Deflection plates are mounted at the opening of the pin-hole diaphragm for horizontal electrostatic deflection of the electron beam. A Faraday's cup, which is set up at an angle of 45° to the electron beam, can be charged by the electrons deflected vertically upward. The charge current can be measured using a separate connection.

In the experiment P3.8.4.1, the current through a pair of Helmholtz coils is set so that the electron beam is incident on the Faraday's cup of the Perrin tube. The Faraday's cup is connected to an electroscopes which has been pre-charged with a known polarity. The polarity of the electron charge can be recognized by the direction of electroscopes deflection when the Faraday's cup is struck by the electron beam. At the same time, the specific electron charge can be estimated. The following relationship applies:

$$\frac{e}{m} = \frac{2U_A}{(B \cdot r)^2} \quad U_A: \text{anode voltage}$$

The bending radius  $r$  of the orbit is predetermined by the geometry of the tube. The magnetic field  $B$  is calculated from the current  $I$  through the Helmholtz coils.

In the experiment P3.8.4.2, the deflection of electrons in crossed alternating magnetic fields is used to produce Lissajou figures on the fluorescent screen. This experiment demonstrates that the electrons respond to a change in the electromagnetic fields with virtually no lag.

In the experiment P3.8.4.3, the deflection of electrons in parallel electric and magnetic alternating fields is used to produce Lissajou figures on the fluorescent screen.

# ELECTRICITY

## FREE CHARGE CARRIERS IN A VACUUM

### P3.8.5


#### THOMSON TUBE

##### P3.8.5.1

Investigating the deflection of electrons in electrical and magnetic fields

##### P3.8.5.2

Assembling a velocity filter (Wien filter) to determine the specific electron charge



Investigating the deflection of electrons in electrical fields (P3.8.5.1)

Cat. No.	Description	P3.8.5.1-2
555 624	Electron deflection tube	1
555 600	Tube stand	1
555 604	Pair of Helmholtz coils	1
521 70	High-voltage power supply, 10 kV	2
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1
500 611	Safety connecting lead, 25 cm, red	2
500 621	Safety connecting lead, 50 cm, red	1
500 622	Safety connecting lead, 50 cm, blue	1
500 641	Safety connecting lead, 100 cm, red	2
500 642	Safety connecting lead, 100 cm, blue	4
500 644	Safety connecting lead, 100 cm, black	2

In the Thomson tube, the electrons pass through a slit behind the anode and fall glancingly on a fluorescent screen placed in the beam path at an angle. A plate capacitor is mounted at the opening of the slit diaphragm which can electrostatically deflect the electron beam vertically. In addition, Helmholtz coils can be used to generate an external magnetic field which can also deflect the electron beam.

The experiment P3.8.5.1 investigates the deflection of electrons in electric and magnetic fields. For different anode voltages  $U_A$ , the beam path of the electrons is observed when the deflection voltage  $U_p$  at the plate capacitor is varied. Additionally, the electrons are deflected in the magnetic field of the Helmholtz coils by varying the coil current  $I$ . The point at which the electron beam emerges from the fluorescent screen gives us the radius  $R$  of the orbit. When we insert the anode voltage in the following equation, we can obtain an experimental value for the specific electron charge

$$\frac{e}{m} = \frac{2U_A}{(B \cdot r)^2}$$

whereby the magnetic field  $B$  is calculated from the current  $I$ .

In the experiment P3.8.5.2, a velocity filter (Wien filter) is constructed using crossed electrical and magnetic fields. Among other things, this configuration permits a more precise determination of the specific electron charge. At a fixed anode voltage  $U_A$ , the current  $I$  of the Helmholtz coils and the deflection voltage  $U_p$  are set so that the effects of the electric field and the magnetic field just compensate each other. The path of the beam is then virtually linear, and we can say:

$$\frac{e}{m} = \frac{1}{2U_A} \cdot \left( \frac{U_p}{B \cdot d} \right)^2$$

$d$ : plate spacing of the plate capacitor



Investigating the deflection of electrons in magnetic fields





### P3.9.1

#### SPONTANEOUS AND NON-SPONTANEOUS DISCHARGE

##### P3.9.1.1

Non-spontaneous gas discharge: comparison between the charge transport in a gas triode and a high-vacuum triode

##### P3.9.1.2

Ignition and extinction of spontaneous gas discharge

*Non-spontaneous gas discharge: comparison between the charge transport in a gas triode and a high-vacuum triode (P3.9.1.1)*

Cat. No.	Description	P3.9.1.1	P3.9.1.2
555 614	Gas triode	1	1
555 612	Demonstration triode	1	
555 600	Tube stand	1	1
521 65	Tube power supply, 0...500 V	1	1
531 130	Multimeter LDanalog 30	1	1
531 120	Multimeter LDanalog 20	2	1
500 641	Safety connecting lead, 100 cm, red	6	4
500 642	Safety connecting lead, 100 cm, blue	4	2

A gas becomes electrically conductive, i. e. gas discharge occurs, when a sufficient number of ions or free electrons as charge carriers are present in the gas. As the charge carriers recombine with each other, new ones must be produced constantly. We speak of self-maintained gas discharge when the existing charge carriers produce a sufficient number of new charge carriers through the process of collision ionization. In non-self-maintained gas discharge, free charge carriers are produced by external effects, e. g. by the emission of electrons from a hot cathode.

The experiment P3.9.1.1 looks at non-self-maintained gas discharge. The comparison of the current-voltage characteristics of a high-vacuum triode and a He gas triode shows that additional charge carriers are created in a gas triode. Some of the charge carriers travel to the grid of the gas triode, where they are measured using a sensitive ammeter to determine their polarity.

The experiment P3.9.1.2 investigates self-maintained discharge in a He gas triode. Without cathode heating, gas discharge occurs at an ignition voltage  $U_2$ . This gas discharge also maintains itself at lower voltages, and only goes out when the voltage falls below the extinction voltage  $U_1$ . Below the ignition voltage  $U_2$ , non-self-maintained discharge can be triggered, e. g. by switching on the cathode heating.

# ELECTRICITY

## ELECTRICAL CONDUCTION IN GASES

### P3.9.2

#### GAS DISCHARGE AT REDUCED PRESSURE

##### P3.9.2.1

Investigating spontaneous gas discharge in air as a function of pressure



Investigating spontaneous gas discharge in air as a function of pressure (P3.9.2.1)

Cat. No.	Description	P3.9.2.1
554 161	Discharge tube, canal rays	1
378 752	Rotary-vane vacuum pump D 2.5 E	1
378 023	Male ground joint, ST 19/26, DN 16 KF	1
378 015	Cross piece, DN 16 KF	1
378 050	Clamping ring, DN 10/16 KF	5
378 045ET2	Centring rings, DN 16 KF, set of 2	3
378 777	Fine vacuum ball valve, DN 16 KF	1
378 776	Variable leak valve, DN 16 KF	1
378 5131	Pirani vacuum gauge with display	1
378 701	Leybold high-vacuum grease	1
521 70	High-voltage power supply, 10 kV	1
501 051	Cable for high voltages, 1.5 m	2
378 764	Exhaust filter AF 8	1*

\* additionally recommended

Glow discharge is a special form of gas discharge. It maintains itself at low pressures with a relatively low current density, and is connected with spectacular luminous phenomena. Research into these phenomena provided fundamental insights into the structure of the atom.

In the experiment P3.9.2.1, a cylindrical glass tube is connected to a vacuum pump and slowly evacuated. A high voltage is applied to the electrodes at the end of the glass tube. No discharge occurs at standard pressure. However, when the pressure is reduced to a certain level, current flows, and a luminosity is visible. When the gas pressure is further reduced, multiple phases can be observed: First, a luminous "thread" joins the anode and the cathode. Then, a column of light extends from the anode until it occupies almost the entire space. A glowing layer forms on the cathode. The column gradually becomes shorter and breaks down into multiple layers, while the glowing layer becomes larger. The layering of the luminous zone occurs because after collision excitation, the exciting electrons must traverse an acceleration distance in order to acquire enough energy to re-excite the atoms. The spacing of the layers thus illustrates the free path length.



### P3.9.3

#### CATHODE RAYS AND CANAL RAYS

##### P3.9.3.1

Magnetic deflection of cathode and canal rays

Magnetic deflection of cathode and canal rays (P3.9.3.1)

Cat. No.	Description	P3.9.3.1
554 161	Discharge tube, canal rays	1
378 752	Rotary-vane vacuum pump D 2.5 E	1
378 023	Male ground joint, ST 19/26, DN 16 KF	1
378 015	Cross piece, DN 16 KF	1
378 050	Clamping ring, DN 10/16 KF	5
378 045ET2	Centring rings, DN 16 KF, set of 2	3
378 777	Fine vacuum ball valve, DN 16 KF	1
378 776	Variable leak valve, DN 16 KF	1
378 5131	Pirani vacuum gauge with display	1
378 701	Leybold high-vacuum grease	1
521 70	High-voltage power supply, 10 kV	1
501 051	Cable for high voltages, 1.5 m	2
510 48	Magnets, 35 mm Ø, pair	1
378 764	Exhaust filter AF 8	1*


\* additionally recommended

Cathode and canal rays can be observed in a gas discharge tube which contains only a residual pressure of less than 0.1 mbar. When a high voltage is applied, more and more electrons are liberated from the residual gas on collision with the cathode. The electrons travel to the anode virtually unhindered, and some of them manage to pass through a hole to the glass wall behind it. Here they are observed as fluorescence phenomena. The luminosity also appears behind the cathode, which is also provided with a hole. A tightly restricted canal ray consisting of positive ions passes straight through the hole until it hits the glass wall. In the experiment P3.9.3.1, the cathode rays, i. e. the electrons, and the canal rays are deflected using a magnet. From the observation that the deflection of the canal rays is significantly less, we can conclude that the ions have a lower specific charge

P3.9.3.1 Magnetic deflection of cathode and canal rays (close up)






# P4 ELECTRONICS



P4.1	COMPONENTS AND BASIC CIRCUITS	151
P4.2	OPERATIONAL AMPLIFIER	159
P4.3	OPEN- AND CLOSED-LOOP CONTROL	161
P4.5	DIGITAL ELECTRONICS	162

# P4 ELECTRONICS



## P4.1 COMPONENTS AND BASIC CIRCUITS

P4.1.1 Current and voltage sources	151-152
P4.1.2 Special resistors	153
P4.1.3 Diodes	154
P4.1.4 Diode circuits	155
P4.1.5 Transistors	156
P4.1.6 Transistor circuits	157
P4.1.7 Optoelectronics	158

## P4.2 OPERATIONAL AMPLIFIER

P4.2.1 Internal design of an operational amplifier	159
P4.2.2 Operational amplifier circuits	160

## P4.3 OPEN- AND CLOSED-LOOP CONTROL

P4.3.2 Closed-loop control	161
----------------------------	-----

## P4.5 DIGITAL ELECTRONICS

P4.5.1 Simple combinations	162
P4.5.2 Logic circuits	163
P4.5.3 Analog inputs and outputs	164



## P4.1.1

### CURRENT AND VOLTAGE SOURCES

#### P4.1.1.1

Determining the internal resistance of a battery

#### P4.1.1.2

Operating a DC power supply as constant-current and constant-voltage source

Determining the internal resistance of a battery (P4.1.1.1)

Cat. No.	Description	P4.1.1.1	P4.1.1.2
576 86	Monocell holder STE 2/50	1	
576 71	Plug-in board section, STE	1	
685 48ET5	Batteries 1.5 V (D, mono), set of 5	1	
531 120	Multimeter LDanalog 20	2	
537 32	Rheostat, 10 ohms	1	1
501 23	Connecting lead, 32 A, 25 cm, black	5	
521 501	AC/DC power supply, 0...15 V/0...5 A		1
501 30	Connecting lead, 32 A, 100 cm, red		1
501 31	Connecting lead, 32 A, 100 cm, blue		1
531 130	Multimeter LDanalog 30		1*
501 25	Connecting lead, 32 A, 50 cm, red		1*
501 26	Connecting lead, 32 A, 50 cm, blue		1*

\* additionally recommended

The voltage  $U_0$  generated in a voltage source generally differs from the terminal voltage  $U$  measured at the connections as soon as a current  $I$  is drawn from the voltage source. A resistance  $R_i$  must therefore exist within the voltage source, across which a part of the generated voltage drops. This resistance is called the internal resistance of the voltage source.

In the experiment P4.1.1.1, a rheostat as an ohmic load is connected to a battery to determine the internal resistance. The terminal voltage  $U$  of the battery is measured for different loads, and the voltage values are plotted over the current  $I$  through the rheostat. The internal resistance  $R_i$  is determined using the formula

$$U = U_0 - R_i \cdot I$$

by drawing a best-fit straight line through the measured values. A second diagram illustrates the power

$$P = U \cdot I$$

as a function of the load resistance. The power is greatest when the load resistance has the value of the internal resistance  $R_i$ .


The experiment P4.1.1.2 demonstrates the difference between a constant-voltage source and a constant-current source using a DC power supply in which both modes are implemented. The voltage and current of the power supply are limited to the respective values  $U_0$  and  $I_0$ . The terminal voltage  $U$  and the current  $I$  consumed are measured for various load resistances  $R$ . When the load resistance  $R$  is reduced, the terminal voltage retains a constant value  $U_0$  as long as the current  $I$  remains below the set limit value  $I_0$ . The DC power supply operates as a constant-voltage source with an internal resistance of zero. When the load resistance  $R$  is increased, the current consumed remains constant at  $I_0$  as long as the terminal voltage does not exceed the limit value  $U_0$ . The DC power supply operates as a constant-current source with infinite internal resistance.

### P4.1.1

#### CURRENT AND VOLTAGE SOURCES

##### P4.1.1.3

Recording the current-voltage characteristics of a solar battery as a function of the irradiance




Recording the current-voltage characteristics of a solar battery as a function of the irradiance (P4.1.1.3)

Cat. No.	Description	P4.1.1.3
578 63	Solar module, 2 V/0.3 A, STE 4/100	1
576 74	Plug-in board, DIN A4, STE	1
576 77	Board holders, STE, pair	1
577 90	Potentiometer, 220 $\Omega$ , STE 4/50	1
501 48	Bridging plugs, STE 2/19, set of 10	1
531 120	Multimeter LDanalog 20	2
450 64	Halogen lamp, 12 V, 50/100 W	1
450 63	Halogen bulb, 12 V/100 W, G6.35	1
521 25	Transformer, 2...12 V, 120 W	1
300 11	Saddle base	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2
501 461	Connecting leads, 19 A, 100 cm, black, pair	1

The solar cell is a semiconductor photoelement in which irradiance is converted directly to electrical energy at the p-n junction. Often, multiple solar cells are combined to create a solar battery.

In the experiment P4.1.1.3 the current-voltage characteristics of a solar battery are recorded for different irradiance levels. The irradiance is varied by changing the distance of the light source. The characteristic curves reveal the characteristic behavior. At a low load resistance, the solar battery supplies an approximately constant current. When it exceeds a critical voltage (which depends on the irradiance), the solar battery functions increasingly as a constant-voltage source.



Current-voltage characteristics for different illuminance levels





## P4.1.2 SPECIAL RESISTORS

P4.1.2.1  
Recording the current-voltage characteristic of an incandescent lamp

P4.1.2.2  
Recording the current-voltage characteristic of a varistor

P4.1.2.3  
Measuring the temperature-dependency of PTC and NTC resistors

P4.1.2.4  
Measuring the light-dependency of photoresistors

Recording the current-voltage characteristic of an incandescent lamp (P4.1.2.1)

Cat. No.	Description	P4.1.2.1	P4.1.2.2	P4.1.2.3	P4.1.2.4
505 08	Bulbs, 12 V/3 W, E10, set of 10	1			
579 06	Lamp holder, E10, top, STE 2/19	1			
524 011USB	Power-CASSY USB	1			
524 220	CASSY Lab 2	1			
578 00	VDR resistor, STE 2/19		1		
576 71	Plug-in board section, STE	1	1	2	
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1	1	1	
531 120	Multimeter LDanalog 20	2	2	2	
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	2	2	
500 441	Connecting lead 19 A, 100 cm, red	1	1	1	
578 06	PTC probe, 30 $\Omega$ , STE 2/19			1	
578 04	NTC probe, 4.7 k $\Omega$ , STE 2/19			1	
666 767	Hotplate, 1500 W, 180 mm diam.			1	
382 34	Thermometer, -10...+110 $^{\circ}$ C/0.2 K			1	
664 104	Beaker, DURAN, 400 ml, squat			1	
578 02	Photoresistor LDR 05, STE 2/19				1
579 05	Lamp holder, E10, lateral, STE 2/19				1
505 131	Bulbs, 6 V/5 W, E10, set of 10				1
521 210	Transformer, 6/12 V				1
311 77	Steel tape measure, 2 m				1
501 461	Connecting leads, 19 A, 100 cm, black, pair				1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1			


Many materials do not conduct voltage and current in proportion to one another. Their resistance depends on the current level. In technical applications, elements in which the resistance depends significantly on the temperature, the luminous intensity or another physical quantity are increasingly important.

In the experiment P4.1.2.1, the computer-assisted measured-value recording system CASSY is used to record the current-voltage characteristic of an incandescent lamp. As the incandescent filament heats up when current is applied, and its resistance depends on the temperature, different characteristic curves are generated when the current is switched on and off. The characteristic also depends on the rate of increase  $dU/dt$  of the voltage.

The experiment P4.1.2.2 records the current-voltage characteristic of a varistor (VDR voltage dependent resistor). Its characteristic is non-linear in its operating range. At higher currents, it enters the so-called "rise range", in which the ohmic component of the total resistance increases.

The aim of the experiment P4.1.2.3 is to measure the temperature characteristics of an NTC thermistor resistor and a PTC thermistor resistor. The respective measured values can be described using empirical equations in which only the rated value  $R_0$ , the reference temperature  $T_0$  and a material constant appear as parameters.

The subject of the experiment P4.1.2.4 is the characteristic of a CdS light-dependent resistor (LDR). Its resistance varies from approx. 100  $\Omega$  to approx. 10 M $\Omega$ , depending on the brightness. The resistance is measured as a function of the distance from an incandescent lamp which illuminates the light-dependent resistor.



Current-voltage characteristic of the light bulb

### P4.1.3

#### DIODES

##### P4.1.3.1

Recording the current-voltage characteristics of diodes

##### P4.1.3.2

Recording the current-voltage characteristics of Zener diodes (Z-diodes)

##### P4.1.3.3

Recording the current-voltage characteristics of light-emitting diodes (LED)

##### P4.1.3.4


Recording the current-voltage characteristics of diodes with CASSY

##### P4.1.3.5

Recording the current-voltage characteristics of Zener diodes (Z-diodes) with CASSY

##### P4.1.3.6

Recording the current-voltage characteristics of light-emitting diodes (LED) with CASSY



Recording the current-voltage characteristics of diodes (P4.1.3.1\_a)


Cat. No.	Description	P4.1.3.1 (a)	P4.1.3.2 (a)	P4.1.3.3 (a)	P4.1.3.4	P4.1.3.5	P4.1.3.6
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1	1
578 51	Diode, 1N 4007, STE 2/19	1			1		
578 50	Diode, AA 118, germanium, STE 2/19	1			1		
577 32	Resistor, 100 Ω, STE 2/19	1	1	1	1	1	1
521 485	AC/DC power supply, 0...12 V/3 A	1	1	1			
531 120	Multimeter LDanalog 20	2	2	2			
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	2	2	2	2	2
500 441	Connecting lead 19 A, 100 cm, red	1	1	1	1	1	1
578 55	Zener diode, 6.2 V, STE 2/19		1			1	
578 54	Zener diode, 9.1 V, STE 2/19		1			1	
578 57	Light emitting diode, green, STE 2/19			1			1
578 47	Light emitting diode, yellow, STE 2/19			1			1
578 48	Light emitting diode, red, STE 2/19			1			1
578 49	Infrared diode, lateral, STE 2/19			1			1
524 013	Sensor-CASSY 2				1	1	1
524 220	CASSY Lab 2				1	1	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)				1	1	1

Virtually all aspects of electronic circuit technology rely on semiconductor components. The semiconductor diodes are among the simplest of these. They consist of a semiconductor crystal in which an n-conducting zone is adjacent to a p-conducting zone. Capture of the charge carriers, i.e. the electrons in the n-conducting and the "holes" in the p-conducting zones, forms a low-conductivity zone at the junction called the depletion layer. The size of this zone is increased when electrons or holes are removed from the depletion layer by an external electric field with a certain orientation. The direction of this electric field is called the reverse direction. Reversing the electric field drives the respective charge carriers into the depletion layer, allowing current to flow more easily through the diode.


In the experiments P4.1.3.1 and P4.1.3.4, the current-voltage characteristics of an Si-diode (silicon diode) and a Ge-diode (germanium diode) are measured and graphed manually point by point or with CASSY. The aim is to compare the current in the reverse direction and the threshold voltage as the most important specifications of the two diodes.

The objective of the experiments P4.1.3.2 and P4.1.3.5 is to measure the current-voltage characteristic of a zener or Z-diode. Here, special attention is paid to the breakdown voltage in the reverse direction, as when this voltage level is reached the current rises abruptly. The current is due to charge carriers in the depletion layer, which, when accelerated by the applied voltage, ionize additional atoms of the semiconductor through collision.

The experiments P4.1.3.3 and P4.1.3.6 compares the characteristics of infrared, red, yellow and green light-emitting diodes.



Recording the current-voltage characteristics of light-emitting diodes (LED)



**P4.1.4**  
**DIODE CIRCUITS**

P4.1.4.1  
Rectifying AC voltage using diodes

P4.1.4.2  
Voltage-limiting with a Z-diode

P4.1.4.3  
Testing polarity with light-emitting diodes

P4.1.4.4  
Rectifying AC voltage using diodes with CASSY

Rectifying AC voltage using diodes (P4.1.4.1\_a)

Cat. No.	Description	P4.1.4.1 (a)	P4.1.4.1 (b)	P4.1.4.2	P4.1.4.3	P4.1.4.4
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1
578 51	Diode, 1N 4007, STE 2/19	4	4			4
579 06	Lamp holder, E10, top, STE 2/19	1	1	1		1
505 08	Bulbs, 12 V/3 W, E10, set of 10	1	1	1		1
501 48	Bridging plugs, STE 2/19, set of 10	1	1			1
521 485	AC/DC power supply, 0...12 V/3 A	1	1	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1				
575 24	Screened cable, BNC/4 mm	1				
531 120	Multimeter LDanalog 20	1		2		
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	3	3	1	3
<b>524 005</b>	<b>Mobile-CASSY 2</b>		<b>1</b>			
578 55	Zener diode, 6.2 V, STE 2/19			1		
577 42	Resistor, 680 Ω, STE 2/19			1	1	
578 57	Light emitting diode, green, STE 2/19				1	
578 48	Light emitting diode, red, STE 2/19				1	
<b>524 013</b>	<b>Sensor-CASSY 2</b>					<b>1</b>
524 220	CASSY Lab 2					1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)					1



Diodes, zener diodes (or Z-diodes) and light-emitting diodes are used today in virtually every electronic circuit.

The experiment P4.1.4.1 explores the function of half-wave and full-wave rectifiers in the rectification of AC voltages. The half-wave rectifier assembled using a single diode blocks the first half-wave of every AC cycle and conducts only the second half-wave (assuming the diode is connected with the corresponding polarity). The full-wave rectifier, assembled using four diodes in a bridge configuration, uses both half-waves of the AC voltage.

The experiment P4.1.4.2 demonstrates how a Z-diode can be used to protect against voltage surges. As long as the applied voltage is below the breakdown voltage  $U_Z$  of the Z-diode, the Z-diode acts as an insulator and the voltage  $U$  is unchanged. At voltages above  $U_Z$ , the current flowing through the Z-diode is so high that  $U$  is limited to  $U_Z$ .

The aim of the experiment P4.1.4.3 is to assemble a circuit for testing the polarity of a voltage using a green and a red light emitting diode (LED). The circuit is tested with both DC and AC voltage.

The experiment P4.1.4.4 explores the function of half-wave and full-wave rectifiers in the rectification of AC voltages. The half-wave rectifier assembled using a single diode blocks the first half-wave of every AC cycle and conducts only the second half-wave (assuming the diode is connected with the corresponding polarity). The full-wave rectifier, assembled using four diodes in a bridge configuration, uses both half-waves of the AC voltage. The computer-based system CASSY is used to record the voltages.

### P4.1.5

#### TRANSISTORS

P4.1.5.1  
Investigating the diode properties of transistor junctions

P4.1.5.2  
Recording the characteristics of a transistor

P4.1.5.3  
Recording the characteristics of a field-effect transistor

P4.1.5.4  
Recording the characteristics of a transistor with CASSY

P4.1.5.5  
Recording the characteristics of a field-effect transistor with CASSY



Recording the characteristics of a transistor (P4.1.5.2\_a)


Cat. No.	Description	P4.1.5.1	P4.1.5.2 (a)	P4.1.5.3 (a)	P4.1.5.4	P4.1.5.5
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1
578 67	Transistor, BD 137, NPN, emitter bottom, STE 4/50	1	1		1	
578 68	Transistor, BD 138, PNP, e.b., STE 4/50	1				
577 32	Resistor, 100 $\Omega$ , STE 2/19	1		1		
521 485	AC/DC power supply, 0...12 V/3 A	1	1			
531 120	Multimeter LDanalog 20	2	3	2		
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	3	4	3	2	2
577 44	Resistor, 1 k $\Omega$ , STE 2/19		1	1		
577 64	Resistor, 47 k $\Omega$ , STE 2/19		1	1		
577 90	Potentiometer, 220 $\Omega$ , STE 4/50		1	1		
577 92	Potentiometer, 1 k $\Omega$ , STE 4/50		1	1		
501 48	Bridging plugs, STE 2/19, set of 10		1	1		
578 77	Transistor (field effect), BF244, STE 4/50			1		1
578 51	Diode, 1N 4007, STE 2/19			1		
521 45	DC power supply 0... $\pm$ 15 V			1		
521 210	Transformer, 6/12 V			1		
575 214	Oscilloscope 30 MHz, two-channel, analogous			1		
575 24	Screened cable, BNC/4 mm			2		
500 422	Connecting lead 19 A, 50 cm, blue			1		
577 56	Resistor, 10 k $\Omega$ , STE 2/19				1	1
578 31	Capacitor, 0.1 $\mu$ F, STE 2/19				1	1
524 013	Sensor-CASSY 2				1	1
524 011USB	Power-CASSY USB				1	1
524 220	CASSY Lab 2				1	1
501 44	Connecting leads, 19 A, 25 cm, red/blue, pair				1	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)				1	1

Transistors are among the most important semiconductor components in electronic circuit technology. We distinguish between bipolar transistors, in which the electrons and holes are both involved in conducting current, and field-effect transistors, in which the current is carried solely by electrons. The electrodes of a bipolar transistor are called the emitter, the base and the collector. The transistor consists of a total of three n-conducting and p-conducting layers, in the order npn or pnp. The base layer, located in the middle, is so thin that charge carriers originating at one junction can cross to the other junction. In field-effect transistors, the conductivity of the current-carrying channel is changed using an electrical field, without applying power. The element which generates this field is called the gate. The input electrode of a field-effect transistor is known as the source, and the output electrode is called the drain.

The experiment P4.1.5.1 examines the principle of the bipolar transistor and compares it with a diode. Here, the difference between an npn and a pnp transistor is explicitly investigated.

The experiments P4.1.5.2 and P4.1.5.4 examine the properties of an npn transistor on the basis of its characteristics. This experiment measures the input characteristic, i.e. the base current  $I_B$  as a function of the base-emitter voltage  $U_{BE}$ , the output characteristic, i.e. the collector current  $I_C$  as a function of the collector-emitter voltage  $U_{CE}$  at a constant base current  $I_B$  and the collector current  $I_C$  as a function of the base current  $I_B$  at a constant collector-emitter voltage  $U_{CE}$ .

In the experiments P4.1.5.3 and P4.1.5.5, the characteristic of a field-effect transistor, i.e. the drain current  $I_D$ , is recorded and diagrammed as a function of the voltage  $U_{DS}$  between the drain and source at a constant gate voltage  $U_G$ .



## P4.1.6 TRANSISTOR CIRCUITS

- P4.1.6.1  
The transistor as an amplifier
- P4.1.6.2  
The transistor as a switch
- P4.1.6.3  
The transistor as a sine-wave generator (oscillator)
- P4.1.6.4  
The transistor as a function generator
- P4.1.6.5  
The field-effect transistor as an amplifier
- P4.1.6.6  
The field-effect transistor as a switch

The transistor as an amplifier (P4.1.6.1\_a)

Cat. No.	Description	P4.1.6.1 (a)	P4.1.6.2	P4.1.6.3	P4.1.6.4	P4.1.6.5	P4.1.6.6
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1	1
578 67	Transistor, BD 137, NPN, emitter bottom, STE 4/50	1	1				
577 44	Resistor, 1 kΩ, STE 2/19	1	1	2			
577 56	Resistor, 10 kΩ, STE 2/19	1		3		1	1
577 64	Resistor, 47 kΩ, STE 2/19	1		2			1
577 80	Variable resistor, 10 kΩ, STE 2/19	1	1				
577 82	Variable resistor, 47 kΩ, STE 2/19	1					
578 38	Capacitor (electrolytic), 47 μF, STE 2/19	1				1	
578 39	Capacitor (electrolytic), 100 μF, STE 2/19	1					
578 40	Capacitor (electrolytic), 470 μF, STE 2/19	1			1		
501 48	Bridging plugs, STE 2/19, set of 10	1	1	1	1	1	1
522 621	Function generator S 12	1				1	1
521 485	AC/DC power supply, 0...12 V/3 A	1					
575 214	Oscilloscope 30 MHz, two-channel, analogous	1		1	1	1	1
575 24	Screened cable, BNC/4 mm	2		2	2	2	2
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	4	2		2	2
501 451	Connecting leads, 19 A, 50 cm, black, pair	1					1
578 02	Photoresistor LDR 05, STE 2/19		1				
578 06	PTC probe, 30 Ω, STE 2/19		1				
579 06	Lamp holder, E10, top, STE 2/19		1		2		
505 08	Bulbs, 12 V/3 W, E10, set of 10		1				
579 13	Toggle switch, STE 2/19		1				
579 38	Heating element 100 Ohm, 2 W, STE 2/50		1				
521 45	DC power supply 0...±15 V	1	1	1	1	1	1
531 120	Multimeter LDanalog 20	2	1			1	1
578 76	Transistor, BC 140, NPN, emitter bottom, STE 4/50			2	2		

Cat. No.	Description	P4.1.6.1 (a)	P4.1.6.2	P4.1.6.3	P4.1.6.4	P4.1.6.5	P4.1.6.6
577 58	Resistor, 15 kΩ, STE 2/19			2	2	1	
577 68	Resistor, 100 kΩ, STE 2/19			2		1	
577 81	Variable resistor, 4.7 kΩ, STE 2/19			2			
578 22	Capacitor, 100 pF, STE 2/19			2			
578 23	Capacitor, 220 pF, STE 2/19			2			
578 35	Capacitor, 1 μF, STE 2/19			2		2	
578 16	Capacitor, 4.7 μF, STE 2/19			2			
501 28	Connecting lead, 32 A, 50 cm, black			1	3	1	
577 46	Resistor, 1.5 kΩ, STE 2/19				2		
578 41	Capacitor (electrolytic), 220 μF, bipolar, STE 2/19				1		
578 13	Capacitor, 0.22 μF, STE 2/19				1		
578 33	Capacitor, 0.47 μF, STE 2/19				1		
578 51	Diode, 1N 4007, STE 2/19				2		
505 191	Bulbs, 15 V/2 W, E10, set of 5				1		
578 77	Transistor (field effect), BF244, STE 4/50					1	1
577 61	Resistor, 33 kΩ, STE 2/19					1	
577 657	Resistor, 68 kΩ, STE 2/19					1	
577 76	Resistor, 1 MΩ, STE 2/19					1	
578 36	Capacitor, 2.2 μF, STE 2/19					1	
577 92	Potentiometer, 1 kΩ, STE 4/50						1

Transistor circuits are investigated on the basis of a number of examples. These include the basic connections of a transistor as an amplifier, the transistor as a light-dependent or temperature-dependent electronic switch, the Wien bridge oscillator as an example of a sine-wave generator, the astable multivibrator, basic circuits with field-effect transistors as amplifiers as well as the field-effect transistor as a low-frequency switch.

### P4.1.7


#### OPTOELECTRONICS

##### P4.1.7.1

Recording the characteristics of a phototransistor connected as a photodiode

##### P4.1.7.2

Assembling a purely optical transmission line



Recording the characteristics of a phototransistor connected as a photodiode (P4.1.7.1\_a)


Cat. No.	Description	P4.1.7.1 (a)	P4.1.7.2
576 74	Plug-in board, DIN A4, STE	1	1
578 61	Phototransistor, lateral, STE 2/19	1	1
577 32	Resistor, 100 $\Omega$ , STE 2/19	1	
577 56	Resistor, 10 k $\Omega$ , STE 2/19	1	3
579 05	Lamp holder, E10, lateral, STE 2/19	1	
505 08	Bulbs, 12 V/3 W, E10, set of 10	1	
501 48	Bridging plugs, STE 2/19, set of 10	1	1
521 485	AC/DC power supply, 0...12 V/3 A	1	
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	
575 24	Screened cable, BNC/4 mm	2	
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	2
578 57	Light emitting diode, green, STE 2/19		1
578 58	Light emitting diode, red, lateral, STE 2/19		1
578 68	Transistor, BD 138, PNP, e.b., STE 4/50		1
578 85	Operational amplifier, LM 741, STE 4/50		1
577 28	Resistor, 47 $\Omega$ , STE 2/19		1
577 40	Resistor, 470 $\Omega$ , STE 2/19		1
577 44	Resistor, 1 k $\Omega$ , STE 2/19		1
577 48	Resistor, 2.2 k $\Omega$ , STE 2/19		1
577 64	Resistor, 47 k $\Omega$ , STE 2/19		1
578 16	Capacitor, 4.7 $\mu$ F, STE 2/19		2
578 39	Capacitor (electrolytic), 100 $\mu$ F, STE 2/19		1
578 40	Capacitor (electrolytic), 470 $\mu$ F, STE 2/19		1
521 45	DC power supply 0... $\pm$ 15 V		1
522 621	Function generator S 12		1
579 29	Earphone		1

Cat. No.	Description	P4.1.7.1 (a)	P4.1.7.2
500 414	Connecting lead, 19 A, 25 cm, black		3
500 424	Connecting lead 19 A, 50 cm, black		1

Optoelectronics deals with the application of the interactions between light and electrical charge carriers in optical and electronic devices. Optoelectronic arrangements consist of a light-emitting, a light-transmitting and a light-sensitive element. The light beam is controlled electrically.

The subject of the experiment P4.1.7.1 is a phototransistor without base terminal connection used as a photodiode. The current-voltage characteristics are displayed on an oscilloscope for the unilluminated, weakly illuminated and fully illuminated states. It is revealed that the characteristic of the fully illuminated photodiode is comparable with that of a Z-diode, while no conducting-state behavior can be observed in the unilluminated state.

The experiment P4.1.7.2 demonstrates optical transmission of the electrical signals of a function generator to a loudspeaker. The signals modulate the light intensity of an LED by varying the on-state current; the light is transmitted to the base of a phototransistor via a flexible light waveguide. The phototransistor is connected in series to the speaker, so that the signals are transmitted to the loudspeaker.



## P4.2.1

### INTERNAL DESIGN OF AN OPERATIONAL AMPLIFIER

#### P4.2.1.1

Discrete assembly of an operational amplifier as a transistor circuit

Discrete assembly of an operational amplifier as a transistor circuit (P4.2.1.1)


Cat. No.	Description	P4.2.1.1
576 75	Plug-in board, DIN A3, STE	2
577 20	Resistor, 10 $\Omega$ , STE 2/19	2
577 36	Resistor, 220 $\Omega$ , STE 2/19	1
577 38	Resistor, 330 $\Omega$ , STE 2/19	1
577 40	Resistor, 470 $\Omega$ , STE 2/19	1
577 44	Resistor, 1 k $\Omega$ , STE 2/19	8
577 52	Resistor, 4.7 k $\Omega$ , STE 2/19	2
577 56	Resistor, 10 k $\Omega$ , STE 2/19	4
577 68	Resistor, 100 k $\Omega$ , STE 2/19	1
577 93	Potentiometer, 1 k $\Omega$ , 10-turn, STE 4/50	1
578 31	Capacitor, 0.1 $\mu$ F, STE 2/19	2
578 39	Capacitor (electrolytic), 100 $\mu$ F, STE 2/19	1
578 51	Diode, 1N 4007, STE 2/19	4
578 55	Zener diode, 6.2 V, STE 2/19	1
578 69	Transistor, BC 550, NPN, emitter bottom, STE 4/50	3
578 71	Transistor, BC 550, NPN, emitter top, STE 4/50	1
578 72	Transistor, BC 560, PNP, emitter top, STE 4/50	1
501 48	Bridging plugs, STE 2/19, set of 10	5
522 621	Function generator S 12	1
521 45	DC power supply 0... $\pm$ 15 V	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1
575 24	Screened cable, BNC/4 mm	2
500 414	Connecting lead, 19 A, 25 cm, black	5
500 424	Connecting lead 19 A, 50 cm, black	2
500 444	Connecting lead 19 A, 100 cm, black	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1*

Cat. No.	Description	P4.2.1.1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
531 183	Digital multimeter 3340	1*

\* additionally recommended

Many electronics applications place great demands on the amplifier. The ideal characteristics include an infinite input resistance, an infinitely high voltage gain and an output voltage which is independent of load and temperature. These requirements can be satisfactorily met using an operational amplifier.

In the experiment P4.2.1.1, an operational amplifier is assembled from discrete elements as a transistor circuit. The key components of the circuit are a difference amplifier on the input side and an emitter-follower stage on the output side. The gain and the phase relation of the output signals are determined with respect to the input signals in inverting and non-inverting operation. This experiment additionally investigates the frequency characteristic of the circuit.




Setup of the operational amplifier for the non-inverting case

### P4.2.2

#### OPERATIONAL AMPLIFIER CIRCUITS

- P4.2.2.1  
Unconnected operational amplifier (comparator)
- P4.2.2.2  
Inverting operational amplifier
- P4.2.2.3  
Non-inverting operational amplifier
- P4.2.2.4  
Adder and subtractor
- P4.2.2.5  
Differentiator and integrator



Unconnected operational amplifier (comparator) (P4.2.2.1)

Cat. No.	Description	P4.2.2.1	P4.2.2.2	P4.2.2.3	P4.2.2.4	P4.2.2.5
576 74	Plug-in board, DIN A4, STE	1	1	1	1	1
578 85	Operational amplifier, LM 741, STE 4/50	1	1	1	1	1
577 56	Resistor, 10 kΩ, STE 2/19	1	2	2	2	1
577 61	Resistor, 33 kΩ, STE 2/19	2	1		1	
577 62	Resistor, 39 kΩ, STE 2/19	1				
577 68	Resistor, 100 kΩ, STE 2/19	1	1		4	1
577 74	Resistor, 470 kΩ, STE 2/19	1				
577 96	Potentiometer, 100 kΩ, STE 4/50	2	1			1
578 26	Capacitor, 2.2 nF, STE 2/19	2				1
578 28	Capacitor, 10 nF, STE 2/19	1				1
578 51	Diode, 1N 4007, STE 2/19	1				
501 48	Bridging plugs, STE 2/19, set of 10	1	1	1	1	1
522 621	Function generator S 12	1	1	1		1
521 45	DC power supply 0...±15 V	1	1	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	1	1		1
575 24	Screened cable, BNC/4 mm	2	2	2		2
500 424	Connecting lead 19 A, 50 cm, black	4	4	5	3	3
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	2	2	2	2
577 44	Resistor, 1 kΩ, STE 2/19		1		1	1
577 50	Resistor, 3.3 kΩ, STE 2/19		1			
577 52	Resistor, 4.7 kΩ, STE 2/19		1	1		1
577 64	Resistor, 47 kΩ, STE 2/19		2			
577 80	Variable resistor, 10 kΩ, STE 2/19		1	1		
531 120	Multimeter LAnalog 20		1	1	1	
577 32	Resistor, 100 Ω, STE 2/19			1		
577 40	Resistor, 470 Ω, STE 2/19			1	1	
577 46	Resistor, 1.5 kΩ, STE 2/19			1	1	
577 48	Resistor, 2.2 kΩ, STE 2/19			1		
577 58	Resistor, 15 kΩ, STE 2/19			1		

Cat. No.	Description	P4.2.2.1	P4.2.2.2	P4.2.2.3	P4.2.2.4	P4.2.2.5
577 38	Resistor, 330 Ω, STE 2/19				1	
577 60	Resistor, 22 kΩ, STE 2/19				1	
500 421	Connecting lead 19 A, 50 cm, red				1	
577 76	Resistor, 1 MΩ, STE 2/19					1
578 15	Capacitor, 1 μF, STE 2/19					1
578 16	Capacitor, 4.7 μF, STE 2/19					1
578 76	Transistor, BC 140, NPN, emitter bottom, STE 4/50					1

The operational amplifier is an important analogue component in modern electronics. Originally designed as a calculating component for analogue computers, it has been introduced into an extremely wide range of applications as an amplifier. The experiment P4.2.2.1 shows that the unconnected operational amplifier overdrives for even the slightest voltage differential at the inputs. It generates a maximum output signal with a sign corresponding to that of the input-voltage differential.

In the experiments P4.2.2.2 and 4.2.2.3, the output of the operational amplifier is fed back to the inverting and non-inverting inputs via resistor  $R_2$ . The initial input signal applied via resistor  $R_1$  is amplified in the inverting operational amplifier by the factor

$$V = -\frac{R_2}{R_1}$$

and in the non-inverting module by the factor

$$V = \frac{R_2}{R_1} + 1$$

The experiment P4.2.2.4 demonstrates the addition of multiple input signals and the subtraction of input signals.


The aim of the experiment P4.2.2.5 is to use the operational amplifier as a differentiator and an integrator. For this purpose, a capacitor is connected to the input resp. the feedback loop of the operational amplifier. The output signals of the differentiator are proportional to the change in the input signals, and those of the integrator are proportional to the integral of the input signals.



P4.3.2  
CLOSED-LOOP CONTROL

P4.3.2.2  
Brightness control with CASSY

P4.3.2.3  
Voltage control with CASSY




Voltage control with CASSY (P4.3.2.3)

Cat. No.	Description	P4.3.2.2	P4.3.2.3
576 74	Plug-in board, DIN A4, STE	1	1
579 05	Lamp holder, E10, lateral, STE 2/19	1	
505 10	Bulbs, 3.8 V/0.27 W, E10, set of 10	1	1
579 13	Toggle switch, STE 2/19	1	1
578 02	Photoresistor LDR 05, STE 2/19	1	
577 20	Resistor, 10 $\Omega$ , STE 2/19	1	
577 23	Resistor, 20 $\Omega$ , STE 2/19	1	
577 28	Resistor, 47 $\Omega$ , STE 2/19	1	
577 32	Resistor, 100 $\Omega$ , STE 2/19	1	
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
524 031	Current source box	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2	2
579 43	Motor and tachogenerator, STE 4/19/50		2
307 641ET5	PVC tubing, 6 mm diam., 5 m		1
579 06	Lamp holder, E10, top, STE 2/19		3
501 48	Bridging plugs, STE 2/19, set of 10		1
524 011USB	Power-CASSY USB		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

Modern technology without control engineering cannot be imagined. Practical examples such as a heating control or voltage control are familiar to everybody. In the following experiments, various controls from the two-point regulator to the PID controller are presented and investigated.

The aim of the experiments P4.3.2.2 and P4.3.2.3 is the computer-aided realization of closed control loops. In the one case, a PID controller is assembled and used to control an incandescent lamp whose brightness is measured using a photoresistor. The other configuration controls a generator which supplies a constant voltage independently of the load.



Voltage regulation of the motor-generator setup


## P4.5.1

### SIMPLE COMBINATIONS

P4.5.1.1  
AND, OR, XOR, and NAND  
operations with two variables

P4.5.1.2  
De Morgan's laws

P4.5.1.3  
Operations with three variables



De Morgan's laws (P4.5.1.2)

Cat. No.	Description	P4.5.1.1-3
571 401	Logic board 1	1

Digital devices are built on the simple concept of repeated application of just a few basic circuits. Operations using these circuits are governed by the rules of Boolean algebra, sometimes also called "logic algebra" when applied to digital circuit technology.

The experiment P4.5.1.1 introduces all operations with one or two variables used in digital technology. The aim is to verify the laws which apply in Boolean algebra, i.e. those describing commutation, idempotents, absorption and negation.

The experiment P4.5.1.2 demonstrates de Morgan's laws in practical application.  $\text{Not}(\text{And}(A,B)) = \text{Or}(\text{Not}(A);\text{Not}(B))$

The object of the experiment P4.5.1.3 is to verify the associative and distributive laws through experiment when operating three variables.

## P4.5.2

### LOGIC CIRCUITS

P4.5.2.1  
AND, NAND, OR and XOR  
operations with four variables


P4.5.2.3  
Multiplexers and demultiplexers

P4.5.2.4  
Adders

P4.5.2.5  
Flipflops

P4.5.2.6  
Counters

P4.5.2.7  
Shift registers



Flipflops (P4.5.2.5)

Cat. No.	Description	P4.5.2.1	P4.5.2.3-7
571 401	Logic board 1	1	
571 402	Logic board 2		1

A combinatorial circuit performs operations such that the output variables are only determined by the input variables, but not by past states. A sequential circuit is additionally able to store the states of individual variables. The output variables also depend on the result of preceding events, which is represented by the switching state of flipflops.

As an approach to the structure of complex combinatorial circuits, the experiment P4.5.2.1 applies the understanding of basic operations previously learned to the logical operation of four inputs.

The experiment P4.5.2.3 demonstrates how a multiplexer is used to switch multiple inputs onto a single output and a demultiplexer distributes the signals of a single input line to multiple output lines.

The experiment P4.5.2.4 investigates discrete and complex adders as key components of a computer.

The aim of the experiment P4.5.2.5 is to study the function of flipflops. It deals with the various demands on the behavior of these fundamental components of sequential circuits, which are required for assembling RS, D and JK flipflops.


The experiment P4.5.2.6 is showing the properties of counters in different configurations.

The experiment P4.5.2.7 investigates the shift register as a working block of a deserializer.

### P4.5.3

#### ANALOG INPUTS AND OUTPUTS

##### P4.5.3.1 DA and AD Converter



DA and AD Converter (P4.5.3.1)

Cat. No.	Description	P4.5.3.1
571 402	Logic board 2	1
531 120	Multimeter LDanalog 20	1

Interfacing digital circuitry to the often analog world requires the use of A/D and D/A converters. Transfer characteristics are investigated. The experiment P4.5.3.1 shows the function of analog to digital conversion and digital to analog conversion.

# P5 OPTICS



P5.1 GEOMETRICAL OPTICS	167
P5.2 DISPERSION AND CHROMATICS	171
P5.3 WAVE OPTICS	176
P5.4 POLARIZATION	189
P5.5 LIGHT INTENSITY	195
P5.6 VELOCITY OF LIGHT	198
P5.7 SPECTROMETER	202
P5.8 PHOTONICS	206

# P5 OPTICS



## P5.1 GEOMETRICAL OPTICS

P5.1.1 Reflection, refraction	167
P5.1.2 Laws of imaging	168
P5.1.3 Image distortion	169
P5.1.4 Optical instruments	170

## P5.2 DISPERSION AND CHROMATICS

P5.2.1 Refractive index and dispersion	171
P5.2.3 Color mixing	172
P5.2.4 Absorption spectra	173-174
P5.2.5 Reflection spectra	175

## P5.3 WAVE OPTICS

P5.3.1 Diffraction	176-179
P5.3.2 Two-beam interference	180
P5.3.3 Newton's Rings	181-182
P5.3.4 Michelson interferometer	183-184
P5.3.5 Other types of interferometers	185-186
P5.3.6 White-light Reflection Holography	187
P5.3.7 Transmission Holography	188

## P5.4 POLARIZATION

P5.4.1 Basic experiments	189
P5.4.2 Birefringence	190
P5.4.3 Optical activity, polarimetry	191
P5.4.4 Kerr effect	192
P5.4.5 Pockels effect	193
P5.4.6 Faraday effect	194

## P5.5 LIGHT INTENSITY

P5.5.1 Quantities and measuring methods of lighting engineering	195
P5.5.2 Laws of radiation	196-197

## P5.6 VELOCITY OF LIGHT


P5.6.1 Measurement according to Foucault/Michelson	198
P5.6.2 Measuring with short light pulses	199
P5.6.3 Measuring with an electronically modulated signal	200-201

## P5.7 SPECTROMETER

P5.7.1 Prism spectrometer	202
P5.7.2 Grating spectrometer	203-205

## P5.8 PHOTONICS

P5.8.2 Basic Optics	206-208
P5.8.3 Optical applications	209-211
P5.8.4 Optical imaging and colour	212
P5.8.5 Laser basics	213-216
P5.8.6 Solid state laser	217-219
P5.8.7 Optical fibres	220-222
P5.8.8 Technical applications	223-227



**P5.1.1**  
**REFLECTION, REFRACTION**

P5.1.1.1  
Reflection of light at straight and curved mirrors

P5.1.1.2  
Refraction of light at straight surfaces and investigation of ray paths in prisms and lenses

Reflection of light at straight and curved mirrors (P5.1.1.1)

Cat. No.	Description	P5.1.1.1-2
463 52	Optical disc	1
450 60	Lamp housing with cable	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
521 210	Transformer, 6/12 V	1
460 43	Small optical bench	1
463 51	Diaphragm with 5 slits	1
460 08	Lens in frame, f=150 mm	1
300 01	Stand base, V-shaped, large	1
301 01	Leybold multiclamp	4
300 41	Stand rod, 25 cm, 12 mm Ø	1

Frequently, the propagation of light can be adequately described simply by defining the ray path. Examples of this are the ray paths of light in mirrors, in lenses and in prisms using sectional models.

The experiment P5.1.1.1 examines how a mirror image is formed by reflection at a plane mirror and demonstrates the reversibility of the ray path. The law of reflection is experimentally validated:

$$\alpha = \beta$$

$\alpha$ : angle of incidence,  $\beta$ : angle of reflection

Further experiment objectives deal with the reflection of a parallel light beam in the focal point of a concave mirror, the existence of a virtual focal point for reflection in a convex mirror, the relationship between focal length and bending radius of the curved mirror and the creation of real and virtual images for reflection at a curved mirror

The experiment P5.1.1.2 deals with the change of direction when light passes from one medium into another. The law of refraction discovered by W. Snell is quantitatively verified:

$$\frac{\sin \alpha}{\sin \beta} = \frac{n_2}{n_1}$$

$\alpha$ : angle of incidence,  $\beta$ : angle of refraction,

$n_1$ : refractive index of medium 1 (here air),

$n_2$ : refractive index of medium 2 (here glass)

This experiment topic also studies total reflection at the transition from a medium with a greater refractive index to one with a lesser refractive index, the concentration of a parallel light beam at the focal point of a collecting lens, the existence of a virtual focal point when a parallel light beam passes through a dispersing lens, the creation of real and virtual images when imaging with lenses and the ray path through a prism.



P5.1.2

LAWS OF IMAGING

P5.1.2.1  
Determining the focal lengths at collecting and dispersing lenses using collimated light

P5.1.2.2  
Determining the focal lengths at collecting lenses through autocollimation

P5.1.2.3  
Determining the focal lengths at collecting lenses using Bessel's method

P5.1.2.4  
Verifying the imaging laws with a collecting lens



Determining the focal lengths at collecting lenses using Bessel's method (P5.1.2.3\_b)

Cat. No.	Description	P5.1.2.1 (b)	P5.1.2.2 (b)	P5.1.2.3-4 (b)
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	1	1
450 60	Lamp housing with cable	1	1	1
460 20	Condenser with diaphragm holder	1	1	1
521 210	Transformer, 6/12 V	1	1	1
460 02	Lens in frame, f=50 mm	1		1
460 03	Lens in frame, f=100 mm	1		1
460 04	Lens in frame, f=200 mm	1		
460 06	Lens in frame, f=-100 mm	1		
441 53	Screen, translucent	1		1
460 310	Optical bench, S1 profile, 1 m	1	1	1
460 311	Clamp rider with clamp	2	2	2
460 312	Clamp rider with clamp, 45/35	1	1	1
311 77	Steel tape measure, 2 m	1	1	1
460 08	Lens in frame, f=150 mm		1	
460 09	Lens in frame, f=300 mm		1	
461 66	Objects for investigating images, pair		1	1
460 28	Plane mirror, 14 cm x 9 cm, with ball joint		1	

The focal lengths of lenses are determined by a variety of means. The basis for these are the laws of imaging.

In the experiment P5.1.2.1, an observation screen is set up parallel to the optical axis so that the path of a parallel light beam can be observed on the screen after passing through a collecting or dispersing lens. The focal length is determined directly as the distance between the lens and the focal point

In autocollimation, experiment P5.1.2.2 a parallel light beam is reflected by a mirror behind a lens so that the image of an object is viewed right next to that object. The distance  $d$  between the object and the lens is varied until the object and its image are exactly the same size. At this point, the focal length is

$$f = d$$

In the Bessel method, experiment P5.1.2.3 the object and the observation screen are set up at a fixed overall distance  $s$  apart. Between these points there are two lens positions  $x_1$  and  $x_2$  at which a sharply focused image of the object is produced on the observation screen. From the lens laws, we can derive the following relationship for the focal length

$$f = \frac{1}{4} \cdot \left( s - \frac{(x_1 - x_2)^2}{s} \right)$$

In the experiment P5.1.2.4, the object height  $G$ , the object width  $g$ , the image height  $B$  and the image width  $b$  are measured directly for a collecting lens in order to confirm the lens laws. The focal length can be calculated using the formula:

$$f = \frac{g \cdot b}{g + b}$$






## P5.1.3 IMAGE DISTORTION

- P5.1.3.1  
Spherical aberration in lens imaging
- P5.1.3.2  
Astigmatism and curvature of image field in lens imaging
- P5.1.3.3  
Lens imaging distortions (barrel and cushion) and coma
- P5.1.3.4  
Chromatic aberration in lens imaging

Spherical aberration in lens imaging (P5.1.3.1\_b)

Cat. No.	Description	P5.1.3.1 (b)	P5.1.3.2 (b)	P5.1.3.3 (b)	P5.1.3.4 (b)
450 60	Lamp housing with cable	1	1	1	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	1	1	1
460 20	Condenser with diaphragm holder	1	1	1	1
521 210	Transformer, 6/12 V	1	1	1	1
461 61	Diaphragms for spherical aberration, pair	1			
461 66	Objects for investigating images, pair	1	1	1	
460 08	Lens in frame, $f=150$ mm	1	1	1	1
460 26	Iris diaphragm	1	1		1
441 53	Screen, translucent	1	1	1	1
460 310	Optical bench, S1 profile, 1 m	1	1	1	1
460 311	Clamp rider with clamp	2	2	2	2
460 312	Clamp rider with clamp, 45/35	2	2	2	2
460 02	Lens in frame, $f=50$ mm			1	
467 95	Colour filter set, primary				1



Intersections of paraxial and abaxial rays

A spherical lens only images a point in an ideal point when the imaging ray traces intersect the optical axis at small angles, and the angle of incidence and angle of refraction are also small when the ray passes through the lens. As this condition is only fulfilled to a limited extent in practice, aberrations (image defects) are unavoidable.

The experiments P5.1.3.1 and P5.1.3.2 deal with aberrations of image sharpness. In a ray path parallel to the optical axis, paraxial rays are united at a different distance from abaxial rays. This effect, known as "spherical aberration", is particularly apparent in lenses with sharp curvatures. Astigmatism and curvature of field may be observed when imaging long objects with narrow light beams. The focal plane is in reality a curved surface, so that the image on the observation screen becomes increasingly fuzzy toward the edges when the middle is sharply focused. Astigmatism is the phenomenon whereby a tightly restricted light beam does not produce a point-type image, but rather two lines which are perpendicular to each other with a finite spacing with respect to the axis.

The experiment P5.1.3.3 explores aberrations of scale. Blocking light rays in front of the lens causes a barrel-shaped distortion, i. e. a reduction in the imaging scale with increasing object size. Screening behind the lens results in cushion-type aberrations. "Coma" is the term for one-sided, plume-like or blob-like distortion of the image when imaged by a beam of light passing through the lens at an oblique angle. The experiment P5.1.3.4 examines chromatic aberrations. These are caused by a change in the refractive index with the wavelength, and are thus unavoidable when not working with non-monochromatic light.

P5.1.4

OPTICAL INSTRUMENTS

P5.1.4.1  
Magnifier and microscope

P5.1.4.2  
Kepler's telescope and  
Galileo's telescope



Magnifier and microscope (P5.1.4.1\_b)

Cat. No.	Description	P5.1.4.1 (b)	P5.1.4.2 (b)
450 60	Lamp housing with cable	1	
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	
460 20	Condenser with diaphragm holder	1	
521 210	Transformer, 6/12 V	1	
460 22	Holder with spring clips	1	
311 09	Glass scale, 5 cm	1	
460 02	Lens in frame, f=50 mm	1	1
460 03	Lens in frame, f=100 mm	1	1
460 08	Lens in frame, f=150 mm	1	
460 04	Lens in frame, f=200 mm	1	1
460 310	Optical bench, S1 profile, 1 m	1	1
460 311	Clamp rider with clamp	2	
460 312	Clamp rider with clamp, 45/35	4	2
441 53	Screen, translucent	1	
311 77	Steel tape measure, 2 m	1	
460 05	Lens in frame, f=500 mm		1
460 06	Lens in frame, f=-100 mm		1
311 22	Vertical rule		1
300 11	Saddle base		1



Ray path through the Kepler telescope

The magnifier, the microscope and the telescope are introduced as optical instruments which primarily increase the angle of vision. The design principle of each of these instruments is reproduced on the optical bench. For quantitative conclusions, the common definition of magnification is used:

$$V = \frac{\tan \psi}{\tan \varphi}$$

$\psi$ : angle of vision with instrument

$\varphi$ : angle of vision without instrument

In the experiment P5.1.4.1, small objects are observed from a short distance. First, a collecting lens is used as a magnifier. Then, a microscope in its simplest form is assembled using two collecting lenses. The first lens, the objective, produces a real, magnified and inverted intermediate image. The second lens, the ocular (or eyepiece) is used as a magnifier to view this intermediate image. For the total magnification of the microscope, the following applies:

$$V_M = V_{ob} \cdot V_{oc}$$

$V_{ob}$ : imaging scale of objective

$V_{oc}$ : imaging scale of ocular

Here,  $V_{oc}$  corresponds to the magnification of the magnifier.

$$V_{oc} = \frac{s_0}{f_{oc}}$$

$s_0$ : clear field of vision

$f_{oc}$ : focal length of ocular

The aim of the experiment P5.1.4.2 is to observe distant objects using a telescope. The objective and the ocular of a telescope are arranged so that the back focal point of the objective coincides with the front focal point of the ocular. A distinction is made between the Galilean telescope, which uses a dispersing lens as an ocular and produces an erect image, and the Kepler telescope, which produces an inverted image because its ocular is a collecting lens. In both cases, the total magnification can be determined as:

$$V_T = \frac{f_{ob}}{|f_{oc}|}$$

$f_{ob}$ : focal length of objective

$f_{oc}$ : focal length of ocular



## P5.2.1 REFRACTIVE INDEX AND DISPERSION

P5.2.1.1  
Determining the refractive index and dispersion of flint glass and crown glass

P5.2.1.2  
Determining the refractive index and dispersion of liquids

Determining the refractive index and dispersion of liquids (P5.2.1.2\_b)

Cat. No.	Description	P5.2.1.1 (b)	P5.2.1.2 (b)
465 22	Prism, crown glass	1	
465 32	Prism, flint glass	1	
460 25	Prism table	1	1
460 22	Holder with spring clips	1	1
450 60	Lamp housing with cable	1	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	1
460 20	Condenser with diaphragm holder	1	1
521 210	Transformer, 6/12 V	1	1
468 03	Light filter, red	1	1
468 07	Light filter, yellow-green	1	1
468 11	Light filter, blue with violet	1	1
460 08	Lens in frame, $f=150$ mm	1	1
460 310	Optical bench, S1 profile, 1 m	1	1
460 311	Clamp rider with clamp	1	1
460 312	Clamp rider with clamp, 45/35	3	3
311 77	Steel tape measure, 2 m	1	1
465 51	Hollow prism		1
665 002	Funnel		1
675 2100	Toluene, 250 ml		1
675 0410	Turpentine oil, rectified, 250 ml		1
675 4760	Cinnamic ethylester, 100 ml		1


Dispersion is the term for the fact that the refractive index  $n$  is different for different-colored light. Often, dispersion also refers to the quantity  $dn/d\lambda$ , i.e. the quotient of the change in the refractive index  $dn$  and the change in the wavelength  $d\lambda$ .

In the experiment P5.2.1.1, the angle of minimum deviation  $\varphi$  is determined for a flint glass and a crown glass prism at the same refracting angle  $\varepsilon$ . This enables determination of the refractive index of the respective prism material according to the formula

$$n = \frac{\sin \frac{1}{2}(\varepsilon + \varphi)}{\sin \frac{1}{2}\varepsilon}$$

The measurement is conducted for several different wavelengths, so that the dispersion can also be quantitatively measured.

In the experiment P5.2.1.2, an analogous setup is used to investigate dispersion in liquids. Toluol, turpentine oil, cinnamic ether, alcohol and water are each filled into a hollow prism in turn, and the differences in the refractive index and dispersion are observed.



Ray path through a prism

### P5.2.3


#### COLOR MIXING

##### P5.2.3.2

Demonstration of additive color mixing

##### P5.2.3.3

Demonstration of subtractive color mixing



*Demonstration of additive color mixing (P5.2.3.2)*

Cat. No.	Description	P5.2.3.2	P5.2.3.3
459 046	Triple lamp, 12 V/3 x 6 W	1	
562 791	Plug-in power supply, 12 V AC	1	1
460 03	Lens in frame, f=100 mm	1	1
460 22	Holder with spring clips	1	1
688 045	Sliding diaphragms, set of 6	1	
467 95	Colour filter set, primary	1	1
467 96	Colour filter set, secondary	1	1
467 97	Triple colour filter	1	
441 53	Screen, translucent	1	1
460 310	Optical bench, S1 profile, 1 m	1	1
460 311	Clamp rider with clamp	1	1
460 312	Clamp rider with clamp, 45/35	2	3
460 313	Clamp rider with fixing column	1	1
459 032	Halogen lamp, 12 V/20 W		1
460 02	Lens in frame, f=50 mm		1
311 77	Steel tape measure, 2 m		1

The colour recognition of the human eye is determined by three types of light receptor cones in the retina. Comparison of the different colours (wavelength ranges) of the visible spectrum with the sensitivity of the different types of cone reveals division into the primary colours: red, green and blue. Combinations of two primary colours result in the secondary colours: cyan, magenta and yellow. Therefore secondary colour filters absorb the third primary colour. A combination of all three primary colours results in white.

In experiment P5.2.3.2 several colour filters (red, green, blue) are placed next to each other in front of a lamp with 3 light sources and in the image on a screen the overlap of these primary colours show additive colour mixing.

In experiment P5.2.3.3 subtractive colour mixing is shown by placing colour filters (yellow, magenta, cyan) partially overlapping in the light beam of a lamp.



## P5.2.4 ABSORPTION SPECTRA

P5.2.4.1  
Absorption spectra of tinted glass samples

P5.2.4.2  
Absorption spectra of colored liquids


Absorption spectra of tinted glass samples (P5.2.4.1\_b)

Cat. No.	Description	P5.2.4.1 (b)	P5.2.4.2 (b)
466 05	Direct vision prism	1	1
467 96	Colour filter set, secondary	1	
468 01	Light filter, dark red	1	
468 09	Light filter, blue-green	1	
468 11	Light filter, blue with violet	1	
460 22	Holder with spring clips	1	
460 25	Prism table	1	2
450 60	Lamp housing with cable	1	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	1
460 20	Condenser with diaphragm holder	1	1
521 210	Transformer, 6/12 V	1	1
441 53	Screen, translucent	1	1
460 03	Lens in frame, f=100 mm	1	1
460 310	Optical bench, S1 profile, 1 m	1	1
460 311	Clamp rider with clamp	2	2
460 312	Clamp rider with clamp, 45/35	3	3
477 14	Plate glass cell (cuvette), 50 x 50 x 20 mm		1
672 0810	Fuchsine, 25 g		1

The colors we perceive when looking through colored glass or liquids are created by the transmitted component of the spectral colors.

In the experiment P5.2.4.1, the light passing through colored pieces of glass from an incandescent light bulb is viewed through a direct-vision prism and compared with the continuous spectrum of the lamp light. The original, continuous spectrum with the continuum of spectral colors disappears. All that remains is a band with the color components of the filter.

In the experiment P5.2.4.2, the light passing through colored liquids from an incandescent light bulb is viewed through a direct-vision prism and compared with the continuous spectrum of the lamp light. The original, continuous spectrum with the continuum of spectral colors disappears. All that remains is a band with the color components of the liquid.



Absorption spectra of tinted glass samples (without filter set, magenta, yellow, cyan)

### P5.2.4

#### ABSORPTION SPECTRA

##### P5.2.4.3


Absorption spectra of tinted glass samples - Recording and evaluating with a spectrophotometer

##### P5.2.4.4

Absorption and fluorescence spectra of coloured liquids - Recording and evaluating with a spectrophotometer

##### P5.2.4.5

Absorption spectra of PMMA optical waveguide - Recording and evaluating with a spectrophotometer



Absorption spectra of tinted glass samples - Recording and evaluating with a spectrophotometer (P5.2.4.3)

Cat. No.	Description	P5.2.4.3	P5.2.4.4	P5.2.4.5
467 96	Colour filter set, secondary	1		
468 01	Light filter, dark red	1		
468 09	Light filter, blue-green	1		
468 11	Light filter, blue with violet	1	1	
460 22	Holder with spring clips	1	1	
450 60	Lamp housing with cable	1	1	
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	1	
460 20	Condenser with diaphragm holder	1	1	
521 210	Transformer, 6/12 V	1	1	
467 251	Compact spectrometer, physics (spectral photometer)	1	1	1
460 251	Fibre holder	1	1	1
460 310	Optical bench, S1 profile, 1 m	1	1	
460 311	Clamp rider with clamp	3	4	
477 14	Plate glass cell (cuvette), 50 x 50 x 20 mm		1	
460 25	Prism table		1	
300 11	Saddle base		1	2
300 40	Stand rod, 10 cm, 12 mm diam.		1	
301 01	Leybold multiclamp		1	
604 5672	Double microspatula, steel, 150 mm		1	
672 0110	Fluoresceine, 25 g		1	
451 17	Lamp socket, E27, Euro plug			1
505 302	Halogen Bulb 230 V/ 46 W, E27			1
579 44	Light waveguide, set of 2			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

In the experiment P5.2.4.3, the light from an incandescent light bulb passing through coloured pieces of glass is recorded with a spectrometer and compared with the continuous spectrum of the lamp light. The original, continuous spectrum with the continuum of spectral colors disappears. All that remains is a band with the colour components of the filter. The transmission coefficient and the optical density of the coloured pieces of glass are calculated.

In the experiment P5.2.4.4, the light from an incandescent light bulb passing through a coloured liquid is recorded using a spectrometer. The fluorescence of the coloured liquid is recorded under a right angle. A blue filter is used to clearly separate fluorescence and light scattering. Both, absorption and fluorescence spectra are compared with the continuous spectrum of the lamp light.

In the experiment P5.2.4.5, light passing through an optical fibre is recorded by a compact spectrometer. The higher order overtones of molecular oscillations create spectral ranges of high absorption, leaving ranges of high transmission in between, the so called „optical windows“.



Absorption and fluorescence spectra of coloured liquids (P5.2.4.4)



**P5.2.5**  
**REFLECTION SPECTRA**


P5.2.5.1  
Reflection spectra of different materials – Recording and evaluating with a spectrophotometer

Reflection spectra of different materials – Recording and evaluating with a spectrophotometer (P5.2.5.1)

Cat. No.	Description	P5.2.5.1
567 06	Conductors/insulators, set of 6	1
460 22	Holder with spring clips	1
450 60	Lamp housing with cable	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
460 20	Condenser with diaphragm holder	1
521 210	Transformer, 6/12 V	1
467 251	Compact spectrometer, physics (spectral photometer)	1
460 251	Fibre holder	1
460 310	Optical bench, S1 profile, 1 m	1
460 311	Clamp rider with clamp	3
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The colors we perceive of opaque objects are induced by the reflected component of the spectral colors.

In the experiment P5.2.5.1, the light from an incandescent light bulb reflected by different materials is recorded using a spectrometer. The reflection coefficients are calculated and compared.



Reflectivity of different materials against the wavelength

P5.3.1  
DIFFRACTION

P5.3.1.1  
Diffraction at a slit, at a post  
and at a circular iris diaphragm

P5.3.1.2  
Diffraction at a double slit and  
multiple slits

P5.3.1.3  
Diffraction at one- and  
two-dimensional gratings



Diffraction at a double slit and multiple slits (P5.3.1.2)

Cat. No.	Description	P5.3.1.1	P5.3.1.2	P5.3.1.3
469 91	Diaphragm with 3 single slits	1		
469 96	Diaphragm with 3 diffraction holes	1		
469 97	Diaphragm with 3 diffraction objects	1		
460 22	Holder with spring clips	1	1	1
471 830	He-Ne Laser, linearly polarised	1	1	1
460 01	Lens in frame, f=5 mm	1	1	1
460 02	Lens in frame, f=50 mm	1	1	1
460 32	Optical bench with standardised profile, 1 m	1	1	1
460 370	Optics rider, 60/34	4	4	4
441 53	Screen, translucent	1	1	1
300 11	Saddle base	1	1	1
469 84	Diaphragm with 3 double slits of different slit widths		1	
469 85	Diaphragm with 4 double slits of different slit spacing		1	
469 86	Diaphragm with 5 multiple slits with different slit numbers		1	
469 87	Diaphragm with 3 gratings			1
469 88	Diaphragm with 2 wire-mesh gratings			1



Interference pattern of experiment P5.3.1.2 on the screen

The experiment P5.3.1.1 looks at the intensity minima for diffraction at a slit. Their angles  $\varphi_k$  with respect to the optical axis for a slit of the width  $b$  is given by the relationship

$$\sin \varphi_k = k \cdot \frac{\lambda}{b} \quad (k = 1; 2; 3; \dots)$$

$\lambda$ : wavelength of the light

In accordance with Babinet's theorem, diffraction at a post produces similar results. In the case of diffraction at a circular iris diaphragm with the radius  $r$ , concentric diffraction rings may be observed; their intensity minima can be found at the angles  $\varphi_k$  using the relationship

$$\sin \varphi_k = k \cdot \frac{\lambda}{r} \quad (k = 0.610; 1.116; 1.619; \dots)$$

The experiment P5.3.1.2 explores diffraction at a double slit. The constructive interference of secondary waves from the first slit with secondary waves from the second slit produces intensity maxima; at a given distance  $d$  between slit midpoints, the angles  $\varphi_n$  of these maxima are specified by

$$\sin \varphi_n = n \cdot \frac{\lambda}{d} \quad (n = 0; 1; 2; \dots)$$

The intensities of the various maxima are not constant, as the effect of diffraction at a single slit is superimposed on the diffraction at a double slit. In the case of diffraction at more than two slits with equal spacings  $d$ , the positions of the interference maxima remain the same. Between any two maxima, we can also detect  $N-2$  secondary maxima; their intensities decrease for a fixed slit width  $b$  and increasing number of slits  $N$ .

The experiment P5.3.1.3 investigates diffraction at a line grating and a crossed grating. We can consider the crossed grating as consisting of two line gratings arranged at right angles to each other. The diffraction maxima are points at the "nodes" of a straight, square matrix pattern.





P5.3.1  
DIFFRACTION

P5.3.1.4  
Diffraction at a single slit -  
Recording and evaluating  
with CASSY

P5.3.1.5  
Diffraction at a double slit and  
multiple slits - Recording and  
evaluating with CASSY

Diffraction at a single slit - Recording and evaluating with CASSY (P5.3.1.4)

Cat. No.	Description	P5.3.1.4	P5.3.1.5
460 14	Adjustable slit	1	
471 830	He-Ne Laser, linearly polarised	1	1
578 62	Solar cell, STE 2/19	1	1
460 21	Holder for plug-in elements	1	1
460 01	Lens in frame, f=5 mm	1	1
460 02	Lens in frame, f=50 mm	1	1
460 33	Optical bench with standardised profile, 2 m	1	1
460 374	Optics rider, 90/50	4	4
460 383	Sliding rider, 90/50	1	1
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
524 040	µV box	1	1
524 082	Rotary motion sensor S	1	1
301 07	Simple bench clamp	1	1
309 48ET2	Fishing line, set of 2	1	1
342 61	Weights, 50 g, set of 12	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1
469 84	Diaphragm with 3 double slits of different slit widths		1
469 85	Diaphragm with 4 double slits of different slit spacing		1
469 86	Diaphragm with 5 multiple slits with different slit numbers		1
460 22	Holder with spring clips		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

A photoelement with a narrow light opening is used to measure the diffraction intensities; this sensor can be moved perpendicularly to the optical axis on the optical bench, and its lateral position can be measured using a displacement transducer. The measured values are recorded and evaluated using the software CASSY Lab.

The experiment P5.3.1.4 investigates diffraction at slit of variable width. The recorded measured values for the intensity  $I$  are compared with the results of a model calculation for small diffraction angles  $\varphi$  which uses the slit width  $b$  as a parameter:

$$I \propto \left( \frac{\sin\left(\frac{\pi b}{\lambda} \varphi\right)}{\frac{\pi b}{\lambda} \varphi} \right)^2 \text{ where } \varphi = \frac{s}{L}$$

$\lambda$ : wavelength of the light

$s$ : lateral shift of photoelement

$L$ : distance between object and photoelement

The experiment P5.3.1.5 explores diffraction at multiple slits. In the model calculation performed for comparison purposes, the slit width  $b$  and the slit spacing  $d$  are both used as parameters.

$$I \propto \left( \frac{\sin\left(\frac{\pi b}{\lambda} \varphi\right)}{\frac{\pi b}{\lambda} \varphi} \right)^2 \cdot \left( \frac{\sin\left(\frac{N\pi d}{\lambda} \varphi\right)}{\sin\left(\frac{\pi d}{\lambda} \varphi\right)} \right)^2$$

$N$ : number of illuminated slits

P5.3.1

DIFFRACTION

P5.3.1.6

Diffraction at a single slit -  
Recording and evaluating  
with VideoCom

P5.3.1.7

Diffraction at a double slit and  
multiple slits - Recording and  
evaluating with VideoCom

P5.3.1.8

Diffraction at a half-plane -  
Recording and evaluating  
with VideoCom



Diffraction at a half-plane - Recording and evaluating with VideoCom (P5.3.1.8)

Cat. No.	Description	P5.3.1.6	P5.3.1.7	P5.3.1.8
460 14	Adjustable slit	1		
471 830	He-Ne Laser, linearly polarised	1	1	1
472 401	Polarisation filter	1	1	1
337 47USB	VideoCom	1	1	1
460 01	Lens in frame, f=5 mm	1	1	1
460 02	Lens in frame, f=50 mm	1	1	
460 11	Lens in frame, f=500 mm	1	1	1
460 32	Optical bench with standardised profile, 1 m	1	1	1
460 373	Optics rider, 60/50	7	7	6
469 84	Diaphragm with 3 double slits of different slit widths		1	
469 85	Diaphragm with 4 double slits of different slit spacing		1	
469 86	Diaphragm with 5 multiple slits with different slit numbers		1	
460 22	Holder with spring clips		1	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

Diffraction at a single slit P5.3.1.6 or double slit and multiple slits P5.3.1.7 can also be measured as a one-dimensional spatial intensity distribution using the single-line CCD camera VideoCom (here used without the camera lens). The VideoCom software enables fast, direct comparison of the measured intensity distributions with model calculations in which the wavelength  $\lambda$ , the focal length  $f$  of the imaging lens, the slit width  $b$  and the slit spacing  $d$  are all used as parameters. These parameters agree closely with the values arrived at through experiment.

It is also possible to investigate diffraction at a half-plane P5.3.1.8. Thanks to the high-resolution CCD camera, it becomes easy to follow the intensity distribution over more than 20 maxima and minima and compare it with the result of a model calculation. The model calculation is based on Kirchhoff's formulation of Huygens' principle. The intensity  $I$  at point  $x$  in the plane of observation is calculated from the amplitude of the electric field strength  $E$  at this point using the formula

$$I(x) = |E(x)|^2$$


The field strength is obtained through the phase-correct addition of all secondary waves originating from various points  $x'$  in the diffraction plane, from the half-plane boundary  $x' = 0$  to  $x' = \infty$ :

$$E(x) \sim \int_0^{\infty} \exp(i \cdot \varphi(x, x')) \cdot dx'$$

Here,

$$\varphi(x, x') = \frac{2\pi}{\lambda} \cdot \frac{(x - x')^2}{2L}$$

In the phase shift of the secondary wave which travels from point  $x'$  in the diffraction plane to point  $x$  in the observation plane as a function of the direct wave. The parameters in the model calculation are the wavelength  $\lambda$  and the distance  $L$  between the diffraction plane and the observation plane. Here too, the agreement with the values obtained in the experiment is close.



Measured (black) and calculated (red) intensity distributions (P5.3.1.6, P5.3.1.8)



## P5.3.1

### DIFFRACTION

#### P5.3.1.9

Investigation of the spatial coherence of an extended light source

Investigation of the spatial coherence of an extended light source (P5.3.1.9)

Cat. No.	Description	P5.3.1.9
451 062	Spectral lamp, Hg 100	1
451 16	Housing for spectral lamps	1
451 30	Universal choke, 230 V, 50 Hz	1
460 32	Optical bench with standardised profile, 1 m	1
460 370	Optics rider, 60/34	2
460 373	Optics rider, 60/50	1
460 374	Optics rider, 90/50	3
468 07	Light filter, yellow-green	1
460 22	Holder with spring clips	2
688 045	Sliding diaphragms, set of 6	1
460 14	Adjustable slit	1
469 85	Diaphragm with 4 double slits of different slit spacing	1
460 02	Lens in frame, f=50 mm	1
460 135	Ocular with scale	1

Coherence is the property of waves that enables them to exhibit stationary interference patterns. The spatial coherence of a light source can be examined in a Young's double-slit interferometer. A light source illuminates a double slit with slit width  $b$  and distance  $g$ . If the partial beams emitted by the light source are coherent at the position of the two slits an interference pattern can be observed after the double slit. The condition for coherent illumination of the two slits is

$$\Delta s = a \cdot \sin \alpha = \frac{1}{2} \cdot \frac{a}{L} \cdot (g + b) < \frac{\lambda}{2}$$

The experiment P5.3.1.9 explores the condition for spatial coherence. The light source is a single slit of variable width illuminated by a Hg spectral lamp. Combined with a filter this results in a monochromatic light source with variable width  $a$ . At a distance  $L$  double slits of different distances of the slits  $g$  (and fixed slit width  $b$ ) are illuminated. For each distance  $g$  the width  $a$  of the adjustable single slit is determined where the interference pattern after the double slit vanishes. Then, the coherence condition is no longer fulfilled.

## P5.3.2

### TWO-BEAM INTERFERENCE

#### P5.3.2.1

Interference at a Fresnel's mirror with an He-Ne laser

#### P5.3.2.2

Lloyd's mirror experiment with an He-Ne laser

#### P5.3.2.3

Interference at Fresnel's biprism with an He-Ne laser



Interference at a Fresnel's mirror with an He-Ne laser (P5.3.2.1)

Cat. No.	Description	P5.3.2.1-2	P5.3.2.3
471 830	He-Ne Laser, linearly polarised	1	1
471 05	Fresnel's mirror, adjustable	1	
460 01	Lens in frame, f=5 mm	1	1
460 04	Lens in frame, f=200 mm	1	1
460 32	Optical bench with standardised profile, 1 m	1	1
460 370	Optics rider, 60/34	3	3
460 373	Optics rider, 60/50	1	1
441 53	Screen, translucent	1	1
300 11	Saddle base	1	1
311 53	Vernier callipers	1	1
311 77	Steel tape measure, 2 m	1	1
471 09	Fresnel biprism		1
460 25	Prism table		1

In these experiments, two coherent light sources are generated by recreating three experiments of great historical significance. In each of these experiments, the respective wavelength  $\lambda$  of the light used is determined by the distance  $d$  between two interference lines and the distance  $a$  of the (virtual) light sources. At a sufficiently great distance  $L$  between the (virtual) light sources and the projection screen, the relationship


$$\lambda = a \cdot \frac{d}{L}$$

obtains. The determination of the quantity  $a$  depends on the respective experiment setup.

In 1821, A. Fresnel used two mirrors inclined with respect to one another to create two virtual light sources positioned close together, which, being coherent, interfered with each other - P5.3.2.1.

In 1839, H. Lloyd demonstrated that a second, virtual light source coherent with the first can be created by reflection in a mirror. He observed interference phenomena between direct and reflected light - P5.3.2.2.

Coherent light sources can also be produced using a Fresnel biprism, first demonstrated in 1826 (P5.3.2.3). Refraction in both halves of the prism results in two virtual images, which are closer together the smaller the prism angle is.



P5.3.2.1 P5.3.2.2 P5.3.2.3

### P5.3.3

#### NEWTON'S RINGS

##### P5.3.3.1

Newton's Rings in transmitted monochromatic light

##### P5.3.3.2

Newton's rings in transmitted and reflected white light



Newton's rings in transmitted and reflected white light (P5.3.3.2)

Cat. No.	Description	P5.3.3.1	P5.3.3.2
471 111	Glass plates for Newton's rings	1	1
460 03	Lens in frame, f=100 mm	2	
460 26	Iris diaphragm	1	
460 22	Holder with spring clips	1	
460 32	Optical bench with standardised profile, 1 m	1	1
460 370	Optics rider, 60/34	6	5
451 111	Spectral lamp, Na	1	
451 062	Spectral lamp, Hg 100	1	
451 16	Housing for spectral lamps	1	
451 30	Universal choke, 230 V, 50 Hz	1	
468 30	Light filter, yellow	1	
468 31	Light filter, green	1	
468 32	Light filter, blue	1	
441 53	Screen, translucent	1	
300 11	Saddle base	1	
460 04	Lens in frame, f=200 mm		2
460 373	Optics rider, 60/50		1
460 380	Cantilever arm		1
471 88	Beam splitter		2
450 64	Halogen lamp, 12 V, 50/100 W		1
450 63	Halogen bulb, 12 V/100 W, G6.35		1
521 25	Transformer, 2...12 V, 120 W		1
501 33	Connecting lead, 32 A, 100 cm, black		2

Newton's rings are produced using an arrangement in which a convex lens with an extremely slight curvature is touching a glass plate, so that an air wedge with a spherically curved boundary surface is formed. When this configuration is illuminated with a vertically incident, parallel light beam, concentric interference rings (the Newton's rings) are formed around the point of contact between the two glass surfaces both in reflection and in transmitted light. For the path difference of the interfering partial beams, the thickness  $d$  of the air wedge is the defining factor; this distance is not in a linear relation to the distance  $r$  from the point of contact:

$$d = \frac{r^2}{2R}$$

$R$ : bending radius of convex lens

In the experiment P5.3.3.1, the Newton's rings are investigated with monochromatic, transmitted light. At a known wavelength  $\lambda$ , the bending radius  $R$  is determined from the radii  $r_n$  of the interference rings. Here, the relationship for constructive interference is:

$$d = n \cdot \frac{\lambda}{2} \quad \text{where } n = 0, 1, 2, \dots$$

Thus, for the radii of the bright interference rings, we can say:

$$r_n^2 = n \cdot R \cdot \lambda \quad \text{where } n = 0, 1, 2, \dots$$

In the experiment P5.3.3.2, the Newton's rings are studied both in reflection and in transmitted light. As the partial beams in the air wedge are shifted in phase by  $\lambda/2$  for each reflection at the glass surfaces, the interference conditions for reflection and transmitted light are complementary. The radii  $r_n$  of the bright interference lines calculated for transmitted light using the equations above correspond precisely to the radii of the dark rings in reflection. In particular, the center of the Newton's rings is bright in transmitted light and dark in reflection. As white light is used, the interference rings are bordered by colored fringes.

P5.3.3

NEWTON'S RINGS

P5.3.3.3  
Newton's rings in reflected  
monochromatic light -  
Recording and evaluating  
with VideoCom



Newton's rings in reflected monochromatic light - Recording and evaluating with VideoCom (P5.3.3.3)

Cat. No.	Description	P5.3.3.3
471 111	Glass plates for Newton's rings	1
460 32	Optical bench with standardised profile, 1 m	1
460 370	Optics rider, 60/34	4
460 373	Optics rider, 60/50	1
460 380	Cantilever arm	1
450 60	Lamp housing with cable	1
450 521	Bulbs, 12 V/30 W, E14, set of 2	1
460 20	Condenser with diaphragm holder	1
521 485	AC/DC power supply, 0...12 V/3 A	1
471 88	Beam splitter	1
460 04	Lens in frame, f=200 mm	1
468 41	Holder for interference filters	1
468 401	Interference filter, 578 nm	1
337 47USB	VideoCom	1
300 02	Stand base, V-shaped, small	1
300 40	Stand rod, 10 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
468 403	Interference filter, 436 nm	1*
468 405	Interference filter, 633 nm	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

\* additionally recommended

In the experiment P5.3.3.3 Newton's rings in reflection are investigated. The interference rings are displayed directly on the photo-transistor line of the VideoCom. This is used for identifying the diameters or radii  $r_n$  for the different wavelengths  $\lambda$  of the interference filters. The condition for destructive interference is in this case

$$r_n^2 = n \cdot \lambda$$

with  $n = 1, 2, 3, \dots$ . Therefore the radii of the dark interference rings are given by

$$d = (n-1) \cdot \frac{\lambda}{2}$$

with  $n = 1, 2, 3, \dots$



## P5.3.4 MICHELSON INTERFEROMETER


P5.3.4.1  
Setting up a Michelson interferometer on the laser optics base plate

P5.3.4.2  
Determining the wavelength of the light of an He-Ne laser using a Michelson interferometer

P5.3.4.7  
Magnetostriction using a Michelson interferometer

Setting up a Michelson interferometer on the laser optics base plate (P5.3.4.1)

Cat. No.	Description	P5.3.4.1	P5.3.4.2	P5.3.4.7
473 40	Laser optics base plate	1	1	1
471 830	He-Ne Laser, linearly polarised	1	1	1
473 411	Laser mount	1	1	1
473 421	Optics base	4	5	3
473 432	Beam divider, 50 %	1	1	1
473 431	Holder for beam divider	1	1	1
473 461	Planar mirror with fine adjustment	2	2	1
473 471	Spherical lens, $f = 2.7 \text{ mm}$	1	1	1
441 53	Screen, translucent	1	1	1
300 11	Saddle base	1	1	1
311 02	Metal rule, 1 m	1	1	1
473 48	Fine adjustment mechanism		1	
471 931	Set magnetostriction			1
562 15	Coil, 1 000 turns			1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A			1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair			1



Setup of the Michelson interferometer

In a Michelson interferometer, an optical element divides a coherent light beam into two parts. The component beams travel different paths, are reflected into each other and finally recombined. As the two component beams have a fixed phase relationship with respect to each other, interference patterns can occur when they are superposed on each other. A change in the optical path length of one component beam alters the phase relation, and thus the interference pattern as well. Thus, given a constant refractive index, a change in the interference pattern can be used to determine a change in the geometric path, e.g. changes in length due to heat expansion or the effects of electric or magnetic fields. When the geometric path is unchanged, then this configuration can be used to investigate changes in the refractive index due to variations e.g. in pressure, temperature and density

In the experiment P5.3.4.1, the Michelson interferometer is assembled on the vibration-proof laser optics base plate. This setup is ideal for demonstrating the effects of mechanical shocks and air streaking.

In the experiment P5.3.4.2, the wavelength of an He-Ne laser is determined from the change in the interference pattern when moving an interferometer mirror using the shifting distance  $\Delta s$  of the mirror. During this shift, the interference lines on the observation screen move. In evaluation, either the interference maxima or interference minima passing a fixed point on the screen while the plane mirror is shifted are counted. For the wavelength  $\lambda$ , the following equation applies:

$$\lambda = 2 \cdot \frac{\Delta s}{Z}$$

Z: number of intensity maxima or minima counted

In the experiment P5.3.4.7 a metal sample attached to one of the mirrors of a Michelson interferometer is placed in a magnetic field. Varying the current through an external, thermally isolated coil will change the length of the sample. Counting fringes, this behaviour can be quantified for different metals.

## P5.3.4

### MICHELSON INTERFEROMETER

#### P5.3.4.3

Determining the wavelength of the light of an He-Ne laser using a Michelson interferometer - Set-up on the optical bench

#### P5.3.4.4

Determination of the coherence time and the line width of spectral lines with the Michelson interferometer

#### P5.3.4.5

Investigation of the pressure induced line broadening using a Michelson interferometer

#### P5.3.4.6

Determination of the line splitting of two spectral lines using a Michelson interferometer



Determination of the coherence time and the line width of spectral lines with the Michelson interferometer (P5.3.4.4\_a)

Cat. No.	Description	P5.3.4.3	P5.3.4.4 (a)	P5.3.4.5 (a)	P5.3.4.6 (a)
471 830	He-Ne Laser, linearly polarised	1			
460 32	Optical bench with standardised profile, 1 m	1	1	1	1
460 373	Optics rider, 60/50	1	1	1	1
460 374	Optics rider, 90/50	5	7	7	7
471 88	Beam splitter	1	1	1	1
473 461	Planar mirror with fine adjustment	2	2	2	2
460 380	Cantilever arm	1	1	1	1
460 01	Lens in frame, f=5 mm	1			
473 48	Fine adjustment mechanism	1	1	1	1
441 53	Screen, translucent	1	1	1	1
300 11	Saddle base	1	1	1	1
311 02	Metal rule, 1 m	1			
451 062	Spectral lamp, Hg 100		1	1	1
451 16	Housing for spectral lamps		1	1	1
451 30	Universal choke, 230 V, 50 Hz		1	1	1
460 26	Iris diaphragm		2	2	2
468 07	Light filter, yellow-green		1	1	
460 22	Holder with spring clips		1	1	1
451 15	High pressure mercury lamp			1	
451 19	Lamp socket, E 27, multi-pin connector			1	
468 30	Light filter, yellow				1

Temporal coherence can be investigated by means of a Michelson interferometer. The maximum time difference  $\Delta t$  during which interference can be observed is called the coherence time. The coherence length is defined as the distance  $\Delta s_C$  the light travels in the coherence time. Typical coherence lengths are a few microns in incandescent lamps, some millimeters in spectral lamps and many meters in lasers. In addition, the coherence time  $\Delta t_C$  is connected to the spectral width  $\Delta \nu$  or  $\Delta \lambda$  of the light source:

$$\Delta \nu = \frac{1}{\Delta t_C} \quad \text{or} \quad \Delta \lambda = \frac{1}{c} \cdot \frac{\lambda_0^2}{\Delta t_C}$$

In the experiment P5.3.4.3, the Michelson interferometer is assembled on the optical bench. The wavelength of an He-Ne laser is determined from the change in the interference pattern when moving an interferometer mirror using the shifting distance  $\Delta s$  of the mirror.

In the experiment P5.3.4.4 the wavelength  $\lambda$  of the green spectral line of a Hg spectral lamp is determined. To measure the coherence length the positions of the movable plane mirror are measured where interference can barely be seen. From the difference in path length the coherence length  $\Delta s_C$ , the coherence time  $\Delta t_C$  and the line width  $\Delta \nu$  of the spectral line are determined.

In experiment P5.3.4.5 the coherence lengths and spectral widths of the green spectral line of a Hg spectral lamp and a high pressure mercury lamp are determined and the results are compared. The higher pressure in the high pressure mercury lamp leads to a significant broadening of the spectral line causing a shorter coherence length.

In the experiment P5.3.4.6 the mean wavelength  $\lambda$  and the line splitting  $\Delta \lambda$  of the yellow line doublet is determined. For two different proximate wavelengths  $\lambda_1$  and  $\lambda_2$  the coherent superposition of two beams leads to a beating: At distinct path length differences the contrast between bright and dark rings of the interference pattern is big while for other path length differences the contrast vanishes completely.





## P5.3.5

### OTHER TYPES OF INTERFEROMETERS

P5.3.5.1  
Setting up a Mach-Zehnder interferometer on the laser optics base plate

P5.3.5.2  
Measuring the refractive index of air with a Mach-Zehnder interferometer

Measuring the refractive index of air with a Mach-Zehnder interferometer (P5.3.5.2)


Cat. No.	Description	P5.3.5.1	P5.3.5.2
473 40	Laser optics base plate	1	1
471 830	He-Ne Laser, linearly polarised	1	1
473 411	Laser mount	1	1
473 421	Optics base	5	6
473 431	Holder for beam divider	2	2
473 432	Beam divider, 50 %	2	2
473 461	Planar mirror with fine adjustment	2	2
473 471	Spherical lens, $f = 2.7 \text{ mm}$	1	1
441 53	Screen, translucent	1	1
300 11	Saddle base	1	1
311 02	Metal rule, 1 m	1	1
473 485	Evacuatable chamber		1
375 58	Hand vacuum pump		1
667 186	Vacuum rubber tubing, 8 mm diam.		1
604 520	Connector with nipple		1
300 02	Stand base, V-shaped, small		1
666 555	Universal clamp, 0...80 mm		1

In a Mach-Zehnder interferometer, an optical element divides a coherent light beam into two parts. The component beams are deflected by mirrors and finally recombined. As the two partial beams have a fixed phase relationship with respect to each other, interference patterns can occur when they are superposed on each other. A change in the optical path length of one component beam alters the phase relation, and consequently the interference pattern as well. As the component beams are not reflected into each other, but rather travel separate paths, these experiments are easier to comprehend and didactically more effective than experiments with the Michelson interferometer. However, the Mach-Zehnder interferometer is more difficult to adjust.

In the experiment P5.3.5.1, the Mach-Zehnder interferometer is assembled on the vibration-proof laser optics base plate.

In the experiment P5.3.5.2, the refractive index of air is determined. To achieve this, an evacuatable chamber is placed in the path of one component beam of the Mach-Zehnder interferometer. Slowly evacuating the chamber alters the optical path length of the respective component beam.

*Note:* Setting up a Michelson interferometer is recommended before using a Mach-Zehnder interferometer for the first time.



Setting up a Mach-Zehnder interferometer on the laser optics base plate (P5.3.5.1)

### P5.3.5

#### OTHER TYPES OF INTERFEROMETERS

##### P5.3.5.3

Determining the wavelength of the light of an He-Ne laser using a Fabry-Perot interferometer



*Determining the wavelength of the light of an He-Ne laser using a Fabry-Perot interferometer (P5.3.5.3)*

Cat. No.	Description	P5.3.5.3
473 40	Laser optics base plate	1
471 830	He-Ne Laser, linearly polarised	1
473 411	Laser mount	1
473 421	Optics base	5
473 432	Beam divider, 50 %	2
473 431	Holder for beam divider	2
473 461	Planar mirror with fine adjustment	1
473 48	Fine adjustment mechanism	1
473 471	Spherical lens, $f = 2.7$ mm	1
441 53	Screen, translucent	1
300 11	Saddle base	1
311 02	Metal rule, 1 m	1

Besides the Michelson and Mach-Zender interferometers, more geometries can be used to split a light beam in two parts and let them overlap afterwards. One with particular practical use is the Fabry-Perot interferometer.

In the experiment P5.3.5.3 the mirrors on the laser base plate are arranged to set up a Fabry-Perot interferometer with two parallel mirrors in line. The first semi-transparent mirror splits the light beam, while the second fully reflective mirror reflects the light back to overlap with the first reflected beam. Moving the second mirror will change the length of the resonator, thus changing the interference pattern and enables the measurement of the laser wavelength.



## P5.3.6

### WHITE-LIGHT REFLECTION HOLOGRAPHY

#### P5.3.6.1

Creating white-light reflection holograms on the laser optics base plate


Creating white-light reflection holograms on the laser optics base plate (P5.3.6.1)

Cat. No.	Description	P5.3.6.1
473 40	Laser optics base plate	1
471 830	He-Ne Laser, linearly polarised	1
473 411	Laser mount	1
473 421	Optics base	3
473 441	Film holder	1
473 451	Object holder	1
473 471	Spherical lens, $f = 2.7 \text{ mm}$	1
311 02	Metal rule, 1 m	1
663 615	Socket strip for earthed plugs, 5-way	1
313 17	Hand-held stop-watch II, mechanical	1
649 11	Storage trays, 86 x 86 x 26, set of 6	1
661 234	Polyethylene bottle	3
667 016	Scissors, 200 mm, pointed	1
473 448	Holography film, 3000 lines/mm	1
473 446	Darkroom accessories	1
473 444	Photographic chemicals	1
671 8910	Iron(III) nitrate-9-hydrate, 250 g	1
672 4910	Potassium bromide, 100 g	1

In creating white-light reflection holograms, a broadened laser beam passes through a film and illuminates an object placed behind the film. Light is reflected from the surface of the object back onto the film, where it is superposed with the light waves of the original laser beam. The film consists of a light-sensitive emulsion of sufficient thickness. Interference creates standing waves within the film, i.e. a series of numerous nodes and antinodes at a distance of  $\lambda/4$  apart. The film is exposed in the planes of the anti-nodes but not in the nodes. Semi-transparent layers of metallic silver are formed at the exposed areas. To reconstruct the image, the finished hologram is illuminated with white light – the laser is not required. The light waves reflected by the semitransparent layers are superposed on each other in such a way that they have the same properties as the waves originally reflected by the object. The observer sees a three-dimensional image of the object. Light beams originating at different layers only reinforce each other when they are in phase. The in-phase condition is only fulfilled for a certain wavelength, which allows the image to be reconstructed using white light.

The object of the experiment P5.3.6.1 is to create white-light reflection holograms. This process uses a protection class 2 laser, so as to minimize the risk of eye damage for the experimenter. Both amplitude and phase holograms can be created simply by varying the photochemical processing of the exposed film.

*Recommendation:* The Michelson interferometer on the laser optics base plate is ideal for demonstrating the effects of disturbances due to mechanical shocks or air streaking in unsuitable rooms, which can prevent creation of satisfactory holograms



Experiment setup for creating white-light reflection holograms

P5.3.7

TRANSMISSION  
HOLOGRAPHY

P5.3.7.1  
Creating transmission holograms  
on the laser optics base plate




Creating transmission holograms on the laser optics base plate (P5.3.7.1)

Cat. No.	Description	P5.3.7.1
473 40	Laser optics base plate	1
471 830	He-Ne Laser, linearly polarised	1
473 411	Laser mount	1
473 421	Optics base	5
473 435	Variable beam divider	1
473 431	Holder for beam divider	1
473 441	Film holder	1
473 451	Object holder	1
473 471	Spherical lens, $f = 2.7 \text{ mm}$	2
311 02	Metal rule, 1 m	1
663 615	Socket strip for earthed plugs, 5-way	1
313 17	Hand-held stop-watch II, mechanical	1
649 11	Storage trays, 86 x 86 x 26, set of 6	1
661 234	Polyethylene bottle	3
667 016	Scissors, 200 mm, pointed	1
473 448	Holography film, 3000 lines/mm	1
473 446	Darkroom accessories	1
473 444	Photographic chemicals	1
671 8910	Iron(III) nitrate-9-hydrate, 250 g	1
672 4910	Potassium bromide, 100 g	1

In creating transmission holograms, a laser beam is split into an object beam and a reference beam, and then broadened. The object beam illuminates an object and is reflected. The reflected light is focused onto a film together with the reference beam, which is coherent with the object beam. The film records an irregular interference pattern which shows no apparent similarity with the object in question. To reconstruct the hologram, a light beam which corresponds to the reference beam is diffracted at the amplitude hologram in such a way that the diffracted waves are practically identical to the object waves. In reconstructing a phase hologram the phase shift of the reference waves is exploited. In both cases, the observer sees a three-dimensional image of the object.

The object of the experiment P5.3.7.1 is to create transmission holograms and subsequently reconstruct them. This process uses a protection class 2 laser, so as to minimize the risk of eye damage for the experimenter. Both amplitude and phase holograms can be created simply by varying the photochemical processing of the exposed film.

*Recommendation:* The Michelson interferometer on the laser optics base plate is ideal for demonstrating the effects of disturbances due to mechanical shocks or air streaking in unsuitable rooms, which can prevent creation of satisfactory holograms.



Experiment setup for creating transmission holograms



P5.4.1  
BASIC EXPERIMENTS

- P5.4.1.1  
Polarization of light through reflection at a glass plate
- P5.4.1.2  
Fresnel's laws of reflection
- P5.4.1.3  
Polarization of light through scattering in an emulsion
- P5.4.1.4  
Malus' law

Fresnel's laws of reflection (P5.4.1.2)

Cat. No.	Description	P5.4.1.1	P5.4.1.2	P5.4.1.3 (a)	P5.4.1.4 (a)
477 20	Plate glass cell (cuvette), 100 x 100 x 10 mm	1	1	1	
460 25	Prism table	1	1	1	
450 64	Halogen lamp, 12 V, 50/100 W	1	1		1
450 63	Halogen bulb, 12 V/100 W, G6.35	1	1		1
450 66	Picture slider	1	1		1
521 25	Transformer, 2...12 V, 120 W	1	1		1
460 26	Iris diaphragm	1	1	1	1
472 401	Polarisation filter	2	2	2	2
460 03	Lens in frame, f=100 mm	1	1		1
441 53	Screen, translucent	1			
460 43	Small optical bench	2	2	1	1
460 40	Swivel joint with protractor scale	1	1		
301 01	Leybold multiclamp	6	7	6	6
300 01	Stand base, V-shaped, large	2	2	1	1
501 33	Connecting lead, 32 A, 100 cm, black	2	2		2
460 08	Lens in frame, f=150 mm		1		
578 62	Solar cell, STE 2/19		1		1
460 21	Holder for plug-in elements		1		1
531 183	Digital multimeter 3340		1		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1		1
450 60	Lamp housing with cable			1	
450 511	Bulbs, 6 V/30 W, E14, set of 2			1	
460 20	Condenser with diaphragm holder			1	
521 210	Transformer, 6/12 V			1	
460 04	Lens in frame, f=200 mm			1	

The fact that light can be polarized is important evidence of the transversal nature of light waves. Natural light is unpolarized. It consists of mutually independent, unordered waves, each of which has a specific polarization state. Polarization of light is the selection of waves having a specific polarization state.

In the experiment P5.4.1.1, unpolarized light is reflected at a glass surface. When we view this through an analyzer, we see that the reflected light is at least partially polarized. The greatest polarization is observed when reflection occurs at the polarizing angle (Brewster angle)  $\alpha_p$ . The relationship

$$\tan \alpha_p = n$$

gives us the refractive index  $n$  of the glass.

Closer observation leads to Fresnel's laws of reflection, which describe the ratio of reflected to incident amplitude for different directions of polarization. These laws are quantitatively verified in the experiment P5.4.1.2.


The experiment P5.4.1.3 demonstrates that unpolarized light can also be polarized through scattering in an emulsion, e. g. diluted milk, and that polarized light is not scattered uniformly in all directions.

The aim of the experiment P5.4.1.4 is to derive Malus's law: when linearly polarized light falls on an analyzer, the intensity of the transmitted light is

$$I = I_0 \cdot \cos^2 \varphi$$

$I_0$ : intensity of incident light

$\varphi$ : angle between direction of polarization and analyzer



P5.4.2

BIREFRINGENCE

P5.4.2.1

Birefringence and polarization with calcareous spar

P5.4.2.2

Quarter-wavelength and half-wavelength plate

P5.4.2.3

Photoelasticity: Investigating the distribution of strains in mechanically stressed bodies



Quarter-wavelength and half-wavelength plate (P5.4.2.2\_b)


Cat. No.	Description	P5.4.2.1 (b)	P5.4.2.2 (b)	P5.4.2.3 (b)
472 02	Calcite crystal, birefringent	1		
460 25	Prism table	1	1	
460 26	Iris diaphragm	1	1	
472 401	Polarisation filter	1	2	2
460 02	Lens in frame, f=50 mm	1		
460 06	Lens in frame, f=-100 mm	1		
441 53	Screen, translucent	1	1	1
460 310	Optical bench, S1 profile, 1 m	1	1	1
460 311	Clamp rider with clamp	2	2	2
460 312	Clamp rider with clamp, 45/35	5	5	7
450 64	Halogen lamp, 12 V, 50/100 W	1	1	1
450 63	Halogen bulb, 12 V/100 W, G6.35	1	1	1
450 66	Picture slider	1	1	1
521 25	Transformer, 2...12 V, 120 W	1	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	2	1
472 601	Quarter-wavelength plate, 140 nm		2	2
472 59	Half-wavelength plate		1	
468 30	Light filter, yellow		1	
578 62	Solar cell, STE 2/19		1	
460 21	Holder for plug-in elements		1	
531 183	Digital multimeter 3340		1	
471 95	Photoelastic models, set of 4			1
460 08	Lens in frame, f=150 mm		2	
300 11	Saddle base			1

The validity of Snell's law of refraction is based on the premise that light propagates in the refracting medium at the same velocity in all directions. In birefringent media, this condition is only fulfilled for the ordinary component of the light beam (the ordinary ray); the law of refraction does not apply for the extraordinary ray.

The experiment P5.4.2.1 looks at birefringence of calcite (Iceland spar). We can observe that the two component rays formed in the crystal are linearly polarized, and that the directions of polarization are perpendicular to each other.

The experiment P5.4.2.2 investigates the properties of  $\lambda/4$  and  $\lambda/2$  plates and explains these in terms of their birefringence; it further demonstrates that the names for these plates refer to the path difference between the ordinary and the extraordinary rays through the plates.

In the experiment P5.4.2.3, the magnitude and direction of mechanical stresses in transparent plastic models are determined. The plastic models become optically birefringent when subjected to mechanical stress. Thus, the stresses in the models can be revealed using polarization-optical methods. For example, the plastic models are illuminated in a setup consisting of a polarizer and analyzer arranged at right angles. The stressed points in the plastic models polarize the light elliptically. Thus, the stressed points appear as bright spots in the field of view. In another configuration, the plastic models are illuminated with circularly polarized light and observed using a quarter-wavelength plate and an analyzer. Here too, the stressed points appear as bright spots in the field of view.



Photoelasticity: Investigating the distribution of strains in mechanically stressed bodies (P5.4.2.3)



Rotation of the plane of polarization with sugar solutions (P5.4.3.2\_b)

### P5.4.3

#### OPTICAL ACTIVITY, POLARIMETRY

##### P5.4.3.1

Rotation of the plane of polarization with quartz

##### P5.4.3.2

Rotation of the plane of polarization with sugar solutions

##### P5.4.3.3

Building a half-shadow polarimeter with discrete elements

##### P5.4.3.4

Determining the concentration of sugar solutions with a standard commercial polarimeter

Cat. No.	Description	P5.4.3.1 (b)	P5.4.3.2 (b)	P5.4.3.3 (b)	P5.4.3.4
472 621	Quartz, parallel	1			
472 641	Quartz, right-handed	1			
472 651	Quartz, left-handed	1			
460 22	Holder with spring clips	1	1		
450 64	Halogen lamp, 12 V, 50/100 W	1	1	1	
450 63	Halogen bulb, 12 V/100 W, G6.35	1	1	1	
450 66	Picture slider	1	1	1	
521 25	Transformer, 2...12 V, 120 W	1	1	1	
468 30	Light filter, yellow	1		1	
472 401	Polarisation filter	2	2	2	
460 03	Lens in frame, f=100 mm	1	1	1	
441 53	Screen, translucent	1	1	1	
460 310	Optical bench, S1 profile, 1 m	1	1	1	
460 311	Clamp rider with clamp	2	2	2	
460 312	Clamp rider with clamp, 45/35	4	4	5	
501 33	Connecting lead, 32 A, 100 cm, black	2	2	2	
477 20	Plate glass cell (cuvette), 100 x 100 x 10 mm		1		
460 25	Prism table		1	1	
468 03	Light filter, red		1		
468 07	Light filter, yellow-green		1		
468 11	Light filter, blue with violet		1		
666 963	Spoon-ended spatula, stainless steel, 120 mm		1	1	1
674 6050	D(+)-Saccharose, 100 g		1	1	1
688 107	Polarizing foils 38 mm Ø, set of 2			1	
688 109	Set of 100 slides cover slip 5 x 5 cm			1	
477 25	Plate glass cell (cuvette), 100 x 80 x 25 mm			1	
657 591	Polarimeter				1
664 111	Beaker, DURAN, 100 ml, tall				1
OHCS-200E	Electronic balance, CS200E				1

Optical activity is the property of some substances of rotating the plane of linearly polarized light as it passes through the material. The angle of optical rotation is measured using a device called a polarimeter.

The experiment P5.4.3.1 studies the optical activity of crystals, in this case a quartz crystal. Depending on the direction of intersection with respect to the optical axis, the quartz rotates the light clockwise ("right-handed") or counterclockwise ("left-handed"). The angle of optical rotation is closely dependent on the wavelength of the light; therefore a yellow filter is used.

The experiment P5.4.3.2 investigates the optical activity of a sugar solution. For a given cuvette length  $d$ , the angles of optical rotation  $\alpha$  of optically active solutions are proportional to the concentration  $c$  of the solution.

$$\alpha = [\alpha] \cdot c \cdot d$$

$[\alpha]$ : rotational effect of the optically active solution

The object of the experiment P5.4.3.3 is to assemble a half-shadow polarimeter from discrete components. The two main elements are a polarizer and an analyzer, between which the optically active substance is placed. Half the field of view is covered by an additional, polarizing foil, of which the direction of polarization is rotated slightly with respect to the first. This facilitates measuring the angle of optical rotation.

In the experiment P5.4.3.4, the concentrations of sugar solutions are measured using a standard commercial polarimeter and compared with the values determined by weighing.



Determining the concentration of sugar solutions with a standard commercial polarimeter (P5.4.3.4)

P5.4.4

KERR EFFECT

P5.4.4.1

Investigating the Kerr effect  
in nitrobenzene



Investigating the Kerr effect in nitrobenzene (P5.4.4.1)

Cat. No.	Description	P5.4.4.1
473 31	Kerr cell	1
450 64	Halogen lamp, 12 V, 50/100 W	1
450 63	Halogen bulb, 12 V/100 W, G6.35	1
450 66	Picture slider	1
468 03	Light filter, red	1
468 05	Light filter, yellow	1
468 07	Light filter, yellow-green	1
468 11	Light filter, blue with violet	1
472 401	Polarisation filter	2
460 03	Lens in frame, f=100 mm	1
460 25	Prism table	1
441 53	Screen, translucent	1
460 32	Optical bench with standardised profile, 1 m	1
460 373	Optics rider, 60/50	6
521 25	Transformer, 2...12 V, 120 W	1
521 70	High-voltage power supply, 10 kV	1
501 051	Cable for high voltages, 1.5 m	2
501 33	Connecting lead, 32 A, 100 cm, black	2
673 9410	Nitrobenzene, 250 ml	1

In 1875, *J. Kerr* discovered that electrical fields cause birefringence in isotropic substances. The birefringence increases quadratically with the electric field strength. For reasons of symmetry, the optical axis of birefringence lies in the direction of the electric field. The normal refractive index of the substance is changed to  $n_e$  for the direction of oscillation parallel to the applied field, and to  $n_o$  for the direction of oscillation perpendicular to it. The experiment results in the relationship

$$n_e - n_o = K \cdot \lambda \cdot E^2$$

$K$ : Kerr constant

$\lambda$ : wavelength of light used

$E$ : electric field strength

The experiment P5.4.4.1 demonstrates the Kerr effect for nitrobenzol, as the Kerr constant is particularly great for this material. The liquid is filled into a small glass vessel in which a suitable plate capacitor is mounted. The arrangement is placed between two polarization filters arranged at right angles, and illuminated with a linearly polarized light beam. The field of view is dark when no electric field is applied. When an electric field is applied, the field of view brightens, as the light beam is elliptically polarized when passing through the birefringent liquid.





## P5.4.5 POCKELS EFFECT

P5.4.5.1  
Demonstrating the Pockels effect  
in a conoscopic beam path

P5.4.5.2  
Pockels effect: transmitting  
information using modulated light

*Demonstrating the Pockels effect in a conoscopic beam path (P5.4.5.1)*

Cat. No.	Description	P5.4.5.1	P5.4.5.2
472 90	Pockels cell	1	1
521 70	High-voltage power supply, 10 kV	1	1
471 830	He-Ne Laser, linearly polarised	1	1
460 01	Lens in frame, f=5 mm	1	
460 02	Lens in frame, f=50 mm	1	
472 401	Polarisation filter	1	1
460 32	Optical bench with standardised profile, 1 m	1	1
460 370	Optics rider, 60/34	5	4
441 53	Screen, translucent	1	
300 11	Saddle base	1	1
500 604	Safety connecting lead, 10 cm, black	1	
500 641	Safety connecting lead, 100 cm, red	1	1
500 642	Safety connecting lead, 100 cm, blue	1	1
522 621	Function generator S 12		1
500 98	Safety adapter sockets, black, set of 6		1
578 62	Solar cell, STE 2/19		1
460 21	Holder for plug-in elements		1
532 20	AC/DC amplifier 30 W		1
587 08	Broad-band speaker		1
500 621	Safety connecting lead, 50 cm, red		2
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		2

The occurrence of birefringence and the alteration of existing birefringence in an electrical field as a linear function of the electric field strength is known as the Pockels effect. In terms of the visible phenomena, it is related to the Kerr effect. However, due to its linear dependency on the electric field strength, the Pockels effect can only occur in crystals without an inversion center, for reasons of symmetry.

The experiment P5.4.5.1 demonstrates the Pockels effect in a lithium niobate crystal placed in a conoscopic beam path. The crystal is illuminated with a divergent, linearly polarized light beam, and the transmitted light is viewed behind a perpendicular analyzer. The optical axis of the crystal, which is birefringent even when no electric field is applied, is parallel to the incident and exit surfaces; as a result, the interference pattern consists of two sets of hyperbolas which are rotated 90° with respect to each other. The bright lines of the interference pattern are due to light rays for which the difference  $\Delta$  between the optical paths of the extraordinary and ordinary rays is an integral multiple of the wavelength  $\lambda$ . The Pockels effect alters the difference of the main refractive indices,  $n_o - n_e$ , and consequently the position of the interference lines. When the so-called half-wave voltage  $U_{\lambda}$  is applied,  $\Delta$  changes by one half wavelength. The dark interference lines move to the position of the bright lines, and vice versa. The process is repeated each time the voltage is increased by  $U_{\lambda}$ .

The experiment P5.4.5.2 shows how the Pockels cell can be used to transmit audio-frequency signals. The output signal of a function generator with an amplitude of several volts is superposed on a DC voltage which is applied to the crystal of the Pockels cell. The intensity of the light transmitted by the Pockels cell is modulated by the superposed frequency. The received signal is output to a speaker via an amplifier and thus made audible.




*Interference pattern in the conoscopic beam path*

P5.4.6

FARADAY EFFECT

P5.4.6.1

Faraday effect: determining Verdet's constant for flint glass as a function of the wavelength



Faraday effect: determining Verdet's constant for flint glass as a function of the wavelength (P5.4.6.1\_a)

Cat. No.	Description	P5.4.6.1 (a)	P5.4.6.1 (b)
560 482	Flint glass block with holder	1	1
460 381	Rider base with thread	1	1
562 11	U-core with yoke	1	1
560 31	Bored pole pieces, pair	1	1
562 13	Coil, 250 turns	2	2
450 63	Halogen bulb, 12 V/100 W, G6.35	1	1
450 64	Halogen lamp, 12 V, 50/100 W	1	1
450 66	Picture slider	1	1
468 05	Light filter, yellow	1	1
468 09	Light filter, blue-green	1	1
468 11	Light filter, blue with violet	1	1
468 13	Light filter, violet	1	1
460 02	Lens in frame, f=50 mm	1	1
472 401	Polarisation filter	2	2
441 53	Screen, translucent	1	1
460 32	Optical bench with standardised profile, 1 m	1	1
460 373	Optics rider, 60/50	5	5
521 25	Transformer, 2...12 V, 120 W	1	1
521 39	Variable extra-low voltage transformer	1	1
531 282	Multimeter Metrahit Pro	1	
531 835	Universal measuring instrument, Physics	1	
524 0381	Combi B sensor S	1	1
501 11	Extension cable, 15 pin	1	1
300 02	Stand base, V-shaped, small	1	1
300 41	Stand rod, 25 cm, 12 mm Ø	1	1
301 01	Leybold multiclamp	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1

Cat. No.	Description	P5.4.6.1 (a)	P5.4.6.1 (b)
501 461	Connecting leads, 19 A, 100 cm, black, pair	1	1
524 005	Mobile-CASSY 2		1
524 0431	30-A-Box		1

Transparent isotropic materials become optically active in a magnetic field; in other words, the plane of polarization of linearly polarized light rotates when passing through the material. *M. Faraday* discovered this effect in 1845 while seeking a relationship between magnetic and optical phenomena. The angle of optical rotation of the plane of polarization is proportional to the illuminated length  $s$  and the magnetic field  $B$ .

$$\Delta\varphi = V \cdot B \cdot s$$

The proportionality constant  $V$  is known as Verdet's constant, and depends on the wavelength  $\lambda$  of the light and the dispersion.

$$V = \frac{e}{2mc^2} \cdot \lambda \cdot \frac{dn}{d\lambda}$$

For flint glass, the following equation approximately obtains:

$$\frac{dn}{d\lambda} = \frac{1,8 \cdot 10^{-14} \text{m}^2}{\lambda^3}$$

In the experiment P5.4.6.1, the magnetic field is initially calibrated with reference to the current through the electromagnets using a magnetic field probe, and then the Faraday effect in a flint glass square is investigated. To improve measuring accuracy, the magnetic field is reversed each time and twice the angle of optical rotation is measured. The proportionality between the angle of optical rotation and the magnetic field and the decrease of Verdet's constant with the wavelength  $\lambda$  are verified.



## P5.5.1 QUANTITIES AND MEASURING METHODS OF LIGHTING ENGINEERING

P5.5.1.1  
Determining the radiant flux density and the luminous intensity of a halogen lamp

P5.5.1.2  
Determining the luminous intensity as a function of the distance from the light source - Recording and evaluating with CASSY

P5.5.1.3  
Verifying Lambert's law of radiation

Determining the radiant flux density and the luminous intensity of a halogen lamp (P5.5.1.1\_a)

Cat. No.	Description	P5.5.1.1 (a)	P5.5.1.2 (c)	P5.5.1.3 (a)
450 64	Halogen lamp, 12 V, 50/100 W	1		1
450 63	Halogen bulb, 12 V/100 W, G6.35	1		1
450 66	Picture slider	1		
468 03	Light filter, red	1		
521 25	Transformer, 2...12 V, 120 W	1		1
557 36	Moll's thermopile	1		1
524 040	µV box	1		1
666 243	Lux sensor	1	1	
524 0511	Lux adapter S	1	1	
<b>524 005</b>	<b>Mobile-CASSY 2</b>	<b>1</b>		<b>1</b>
460 03	Lens in frame, f=100 mm	1		1
460 21	Holder for plug-in elements	1	1	
590 02ET2	Clip plugs, small, set of 2	1	1	
460 43	Small optical bench	1		2
301 01	Leybold multiclamp	3		4
300 01	Stand base, V-shaped, large	1		2
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1		1
501 33	Connecting lead, 32 A, 100 cm, black	2		2
450 60	Lamp housing with cable		1	
450 511	Bulbs, 6 V/30 W, E14, set of 2		1	
521 210	Transformer, 6/12 V		1	
<b>524 013</b>	<b>Sensor-CASSY 2</b>		<b>1</b>	
524 220	CASSY Lab 2		1	
460 310	Optical bench, S1 profile, 1 m		1	
460 311	Clamp rider with clamp		1	
460 312	Clamp rider with clamp, 45/35		1	
450 68	Halogen bulb, 12 V/50 W, G6.35			1

Cat. No.	Description	P5.5.1.1 (a)	P5.5.1.2 (c)	P5.5.1.3 (a)
460 26	Iris diaphragm			1
460 22	Holder with spring clips			1
460 40	Swivel joint with protractor scale			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1	

There are two types of physical quantities used to characterize the brightness of light sources: quantities which refer to the physics of radiation, which describe the energy radiation in terms of measurements, and quantities related to lighting engineering, which describe the subjectively perceived brightness under consideration of the spectral sensitivity of the human eye. The first group includes the irradiance  $E_e$ , which is the radiated power per unit of area  $\Phi_e$ . The corresponding unit of measure is watts per square meter. The comparable quantity in lighting engineering is illuminance  $E$ , i. e. the emitted luminous flux per unit of area  $\Phi$ , and it is measured in lumens per square meter, or lux for short.

In the experiment P5.5.1.1, the irradiance is measured using the Moll's thermopile, and the luminous flux is measured using a luxmeter. The luxmeter is matched to the spectral sensitivity of the human eye  $V(\lambda)$  by means of a filter placed in front of the photoelement. A halogen lamp serves as the light source. From its spectrum, most of the visible light is screened out using a color filter; subsequently, a heat filter is used to absorb the infrared component of the radiation.

The experiment P5.5.1.2 demonstrates that the luminous intensity is proportional to the square of the distance between a point-type light source and the illuminated surface.

The aim of the experiment P5.5.1.3 is to investigate the angular distribution of the reflected radiation from a diffusely reflecting surface, e.g. matte white paper. To the observer, the surface appears uniformly bright; however, the apparent surface area varies with the cos of the viewing angle. The dependency of the luminous intensity is described by Lambert's law of radiation:

$$E_e(\phi) = E_e(0) \cdot \cos \phi$$

P5.5.2

LAWS OF RADIATION

P5.5.2.1

Stefan-Boltzmann law: measuring the radiant intensity of a „black body“ as a function of temperature

P5.5.2.2

Stefan-Boltzmann law: measuring the radiant intensity of a „black body“ as a function of temperature - Recording and evaluating with CASSY

P5.5.2.3

Confirming the laws of radiation with Leslie's cube



Stefan-Boltzmann law: measuring the radiant intensity of a „black body“ as a function of temperature (P5.5.2.1\_b)

Cat. No.	Description	P5.5.2.1 (b)	P5.5.2.2 (b)	P5.5.2.3 (b)
555 81	Electric oven for tubes, 230 V	1	1	
389 43	Black body accessory	1	1	
502 061	Safety connecting box, with earth	1	1	
555 84	Support for electric oven	1	1	1
524 005	Mobile-CASSY 2	1		1
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K	1	1	1
524 040	µV box	1	1	1
557 36	Moll's thermopile	1	1	1
460 32	Optical bench with standardised profile, 1 m	1	1	1
460 373	Optics rider, 60/50	2	2	
460 374	Optics rider, 90/50	2	2	2
460 380	Cantilever arm	1	1	
666 555	Universal clamp, 0...80 mm	1	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1	1	1
388 181	Immersion pump	1*	1*	
521 231	Low-voltage power supply, 3/6/9/12 V	1*	1*	
667 194	Silicone tubing, 7 mm diam., 1 m	1*	1*	
604 313	Wide-mouthed can, 10 l	1*	1*	
524 013	Sensor-CASSY 2		1	
524 220	CASSY Lab 2		1	
524 0673	NiCr-Ni adapter S, type K		1	
389 261	Leslie's cube with stirrer			1
301 01	Leybold multiclamp			1
303 25	Immersion heater			1
590 06	Plastic beaker			1
665 009	Funnel, PP, 75 mm diam.			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1	

\* additionally recommended

The total radiated power  $M_B$  of a black body increases in proportion to the fourth power of its absolute temperature  $T$  (Stefan-Boltzmann's law).

$$M_B = \sigma \cdot T^4$$

$$\sigma = 5.67 \cdot 10^{-8} \text{ W m}^{-2} \text{ K}^{-4}: \text{ (Stefan-Boltzmann's constant)}$$


For all other bodies, the radiated power  $M$  is less than that of the black body, and depends on the properties of the surface of the body. The emittance of the body is described by the relationship

$$\varepsilon = \frac{M}{M_B}$$

$M$ : radiated power of body

In the two experiments P5.5.2.1 and P5.5.2.2, a cylindrical electric oven with a burnished brass cylinder is used as a "black body". The brass cylinder is heated in the oven to the desired temperature between 300 and 750 K. A thermocouple is used to measure the temperature. A water-coolable screen is positioned in front of the oven to ensure that the setup essentially measures only the temperature of the burnished brass cylinder. The measurement is conducted using a Moll's thermopile; its output voltage provides a relative measure of the radiated power  $M$ . The thermopile can be connected either to an amplifier or, via the µV box, to the CASSY computer interface device. In the former case, the measurement must be carried out manually, point by point; the latter configuration enables computer-assisted measuring and evaluation. The aim of the evaluation is to confirm Stefan-Boltzmann's law.

The experiment P5.5.2.3 uses a radiation cube after Leslie ("Leslie's cube"). This cube has four different face surfaces (metallic matte, metallic shiny, black finish and white finish), which can be heated from the inside to almost 100 °C by filling the cube with boiling water. The heat radiated by each of the surfaces is measured as a function of the falling temperature. The aim of the evaluation is to compare the emittances of the cube faces.





**P5.5.2**  
**LAWS OF RADIATION**

P5.5.2.4  
The Wien's displacement law - spectral recording of the black body radiation


The Wien's displacement law - spectral recording of the black body radiation (P5.5.2.4)

Cat. No.	Description	P5.5.2.4
467 251	Compact spectrometer, physics (spectral photometer)	1
460 251	Fibre holder	1
459 032	Halogen lamp, 12 V/20 W	1
521 485	AC/DC power supply, 0...12 V/3 A	1
460 310	Optical bench, S1 profile, 1 m	1
460 311	Clamp rider with clamp	1
460 313	Clamp rider with fixing column	1
501 451	Connecting leads, 19 A, 50 cm, black, pair	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

In the experiment P5.5.2.4 an incandescent lamp, an electric conductor is heated by an electric current causing it to glow. The emitted spectrum is continuous and can be described with Planck's radiation formula. The radiation maximum of the radiation shifts with increasing temperature  $T$  according to Wien's displacement law

$$\lambda_{\text{max}} = 2.9 \cdot 10^{-3} \text{ m}\cdot\text{K} / T$$

toward smaller wavelengths; at the same time the maximum value of the radiation increases. At the temperatures attained by normal incandescent lamps, about 2300 to 2900 K, the radiation maximum lies in the infrared spectral range. Halogen lamps reach a somewhat higher operating temperature of about 3000 K. This experiment will record the spectra of a halogen lamp at different power levels. The specification for color temperature at 12 V nominal voltage allows the wavelength dependency of spectrometer sensitivity to be determined and thus permits a corrected display for the course of spectra intensity. Subsequently voltage will be applied to a halogen lamp then reduced in steps. The lamp's color temperature at each voltage step will be determined by adapting a suitable fit function.



Spectral emission of a light bulb at different temperatures

### P5.6.1

#### MEASUREMENT ACCORDING TO FOUCAULT/MICHELSON

##### P5.6.1.1

Determining the velocity of light by means of the rotating-mirror method according to Foucault and Michelson - Measuring the image shift as a function of the rotational speed of the mirror

##### P5.6.1.2

Determining the velocity of light by means of the rotating-mirror method according to Foucault and Michelson - Measuring the image shift for the maximum rotational speed of the mirror



Determining the velocity of light by means of the rotating-mirror method according to Foucault and Michelson - Measuring the image shift as a function of the rotational speed of the mirror (P5.6.1.1)

Cat. No.	Description	P5.6.1.1	P5.6.1.2
476 40	Rotary mirror	1	1
471 830	He-Ne Laser, linearly polarised	1	1
463 20	Front-silvered mirror	1	1
460 12	Lens in frame, f=5 m	1	1
471 88	Beam splitter	1	1
460 22	Holder with spring clips	1	1
311 09	Glass scale, 5 cm	1	1
521 40	Variable transformer, 0...250 V	1	
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	
559 921	Semiconductor detector	1	
501 02	BNC cable, 1 m	1	
501 10	BNC adapter, straight	1	
300 41	Stand rod, 25 cm, 12 mm Ø	1	
300 42	Stand rod, 47 cm, 12 mm diam.	1	1
300 44	Stand rod, 100 cm, 12 mm diam.	1	1
300 01	Stand base, V-shaped, large	1	1
300 02	Stand base, V-shaped, small	4	4
300 11	Saddle base	1	
301 01	Leybold multiclamp	2	2
301 09	Bosshead S	1	
311 02	Metal rule, 1 m	1	1
537 35	Rheostat, 330 ohms		1
537 36	Rheostat, 1000 ohms		1
502 05	Measuring junction box		1
504 48	Two-way switch		1
500 644	Safety connecting lead, 100 cm, black		5

Measurement of the velocity of light by means of the rotary mirror method utilizes a concept first proposed by *L. Foucault* in 1850 and perfected by *A. A. Michelson* in 1878. In the variation utilized here, a laser beam is deviated into a fixed end mirror located next to the light source via a rotating mirror set up at a distance of  $a = 12.1$  m. The end mirror reflects the light so that it returns along the same path when the rotary mirror is at rest. Part of the returning light is imaged on a scale using a beam divider. A lens with  $f = 5$  m images the light source on the end mirror and focuses the image of the light source from the mirror on the scale. The main beam between the lens and the end mirror is parallel to the axis of the lens, as the rotary mirror is set up in the focal point of the lens.

Once the rotary mirror is turning at a high frequency  $\nu$ , the shift  $\Delta x$  of the image on the scale is observed. In the period

$$\Delta t = \frac{2a}{c}$$

which the light requires to travel to the rotary mirror and back to the end mirror, the rotary mirror turns by the angle

$$\Delta\alpha = 2\pi\nu \cdot \Delta t$$

Thus, the image shift is

$$\Delta x = 2\Delta\alpha \cdot a$$

The velocity of light can then be calculated as

$$c = 8\pi \cdot a^2 \cdot \frac{\nu}{\Delta x}$$

To determine the velocity of light, it is sufficient to measure the shift in the image at the maximum speed of the mirror, which is known (P5.6.1.2). Measuring the image shift as a function of the speed supplies more precise results (P5.6.1.1).



**P5.6.2**  
**MEASURING WITH SHORT LIGHT PULSES**

P5.6.2.1  
Determining the velocity of light in air from the path and transit time of a short light pulse

P5.6.2.2  
Determining the propagation velocity of voltage pulses in coaxial cables

Determining the velocity of light in air from the path and transit time of a short light pulse (P5.6.2.1)

Cat. No.	Description	P5.6.2.1	P5.6.2.2
476 50	Light velocity measuring instrument	1	1
460 10	Lens in frame, f=200 mm	1	
460 335	Optical bench with standardised profile, 0.5 m	1	
460 374	Optics rider, 90/50	2	
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	1
501 02	BNC cable, 1 m	3	2
311 02	Metal rule, 1 m	1	
300 01	Stand base, V-shaped, large	1	
300 44	Stand rod, 100 cm, 12 mm diam.	1	
301 01	Leybold multiclamp	1	
501 024	BNC Cable, l = 10 m		1
501 091	BNC T adapter		1
501 10	BNC adapter, straight		1
575 35	Adapter, BNC/4 mm, 2-pole		1
577 79	Variable resistor, 1 kΩ, STE 2/19		1
577 28	Resistor, 47 Ω, STE 2/19		1
300 11	Saddle base		1


The light velocity measuring instrument emits pulses of light with a pulse width of about 20 ns. After traversing a known measuring distance in both directions, the light pulses are converted into voltage pulses for observation on the oscilloscope.

In the experiment P5.6.2.1, the path of the light pulses is varied once, and the change in the transit time is measured with the oscilloscope. The velocity of light can then be calculated as quotient of the change in the transit distance and the change in the transit time. Alternatively, the total transit time of the light pulses can be measured in absolute terms using a reference pulse. In this case, the velocity of light can be calculated as quotient of the transit distance and the transit time. A quartz-controlled oscilloscope signal can be displayed on the instrument simultaneously with the measuring pulse in order to calibrate timing. Time measurement is then independent of the time base of the oscilloscope.

In the experiment P5.6.2.2, the propagation velocity of voltage pulses in coaxial cables is determined. In this configuration, the reference pulses of the light velocity measuring instrument are output to an oscilloscope and additionally fed into a 10 m long coaxial cable via a T-connector. After reflection at the cable end, the pulses return to the oscilloscope, delayed by the transit time. The propagation velocity  $v$  can be calculated from the double cable length and the time difference between the direct and reflected voltage pulses. By inserting these values in the equation

$$v = \frac{c}{\sqrt{\epsilon_r}} \quad c: \text{velocity of light in a vacuum}$$

we obtain the relative dielectricity  $\epsilon_r$  of the insulator between the inner and outer conductors of the coaxial cable. By using a variable terminating resistor  $R$  at the cable end, it becomes possible to additionally measure the reflection behaviour of voltage pulses. In particular, the special cases „open cable end“ (no phase shift at reflection), „shorted cable end“ (phase shift due to reflection) and „termination of cable end with the 50 Ω characteristic wave impedance“ (no reflection) are of special interest here.



Schematic diagram of light velocity measurement with short light pulses (P5.6.2.1)

P5.6.3

MEASURING WITH AN ELECTRONICALLY MODULATED SIGNAL

P5.6.3.1  
Determining the velocity of light using a periodical light signal at a short measuring distance

P5.6.3.2  
Determining the velocity of light in various materials in various materials



Determining the velocity of light in various materials (P5.6.3.2)

Cat. No.	Description	P5.6.3.1	P5.6.3.2
476 301	Light transmitter and receiver	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	1
460 08	Lens in frame, f=150 mm	1	1
300 11	Saddle base	2	4
311 02	Metal rule, 1 m	1	1
476 35	Tube with 2 end windows		1
476 34	Transparent plastic block		1*
477 03	Plate glass cell (cuvette), 50 x 50 x 50 mm		1*
460 25	Prism table		1*
671 9720	Ethanol, denaturated, 1 l		1*
672 1210	Glycerin, 99 %, 250 ml		1*

\* additionally recommended

In determining the velocity of light with an electronically modulated signal, a light emitting diode which pulses at a frequency of 60 MHz is used as the light transmitter. The receiver is a photodiode which converts the light signal into a 60 MHz AC voltage. A connecting lead transmits a reference signal to the receiver which is synchronized with the transmitted signal and in phase with it at the start of the measurement. The receiver is then moved by the measuring distance  $\Delta s$ , so that the received signal is phase-shifted by the additional transit time  $\Delta t$  of the light signal.

$$\Delta\varphi = 2\pi \cdot f_1 \cdot \Delta t \quad \text{where } f_1 = 60 \text{ MHz}$$


Alternatively, a medium with a greater index of refraction can be placed in the beam path. The apparent transit time to be measured is increased by means of an electronic "trick". The received signal and the reference signal are each mixed (multiplied) with a 59.9 MHz signal before being fed through a frequency filter which only passes the low frequency components with the differential frequency  $f_1 - f_2 = 0.1 \text{ MHz}$ . This mixing has no effect on the phase shift; however, this phase shift is now for a transit time  $\Delta t'$  increased by a factor of

$$\frac{f_1}{f_1 - f_2} = 600$$

In the experiment P5.6.3.1, the apparent transit time  $\Delta t'$  is measured as a function of the measuring distance  $\Delta s$ , and the velocity of light in the air is calculated according to the formula

$$c = \frac{\Delta s}{\Delta t'} \cdot \frac{f_1}{f_1 - f_2}$$

The experiment P5.6.3.2 determines the velocity of light in various propagation media. In the way of accessories, this experiment requires a tube 1 m long with two end windows, suitable for filling with water, a glass cell 5 cm wide for other liquids and an acrylic glass body 5 cm wide.



Block circuit diagram





**P5.6.3**  
**MEASURING WITH AN ELECTRONICALLY MODULATED SIGNAL**

**P5.6.3.3**  
Determining the velocity of light using a periodical light signal at a short measuring distance - measuring with the laser motion sensor S and CASSY

**P5.6.3.4**  
Determining the velocity of light for different propagation media - measuring with the laser motion sensor S and CASSY

Determining the velocity of light using a periodical light signal at a short measuring distance - measuring with the laser motion sensor S and CASSY (P5.6.3.3)

Cat. No.	Description	P5.6.3.3	P5.6.3.4
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
524 073	Laser motion sensor S	1	1
337 116	End buffers, pair	1	1
311 02	Metal rule, 1 m	1	
477 03	Plate glass cell (cuvette), 50 x 50 x 50 mm		1
476 34	Transparent plastic block		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

Modern distance meters use a periodically modulated laser beam for the measurement. They determine the phase shift between the emitted and the reflected modulated laser beam and, with the modulation frequency being known, obtain the time-of-flight  $t$  of the light on its path to and back from the reflector. Only afterwards do the distance meters calculate the distance with the aid of the known velocity of light.


In the experiment P5.6.3.3, the laser motion sensor S is used as a time-of-flight meter because it is also capable of outputting the time-of-flight  $t$  directly. The proportionality between the distance and the time-of-flight of light is confirmed, and the velocity of light is calculated.

In the experiment P5.6.3.4 water and acrylic glass of thickness  $d$  are held into the path of the beam, and then the resulting increase of the time-of-flight  $\Delta t$  is measured. With the velocity of light  $c$  in air measured in the experiment P5.6.3.3, the velocity of light  $c_M$  in matter can now be determined:

$$c_M = 2d \left( \frac{2d}{c} + \Delta t \right) = \frac{1}{\frac{1}{c} + \frac{\Delta t}{2d}}$$

Finally, the refractive index  $n$  is determined according to

$$n = \frac{c}{c_M} = c \cdot \left( \frac{1}{c} + \frac{\Delta t}{2d} \right) = 1 + \frac{c}{2d \cdot \Delta t}$$



Transit times of light at different distances

## P5.7.1

### PRISM SPECTROMETER

#### P5.7.1.1


Measuring the line spectra of inert gases and metal vapors using a prism spectrometer



Measuring the line spectra of inert gases and metal vapors using a prism spectrometer (P5.7.1.1)

Cat. No.	Description	P5.7.1.1
467 23	Spectrometer and goniometer	1
451 031	Spectral lamp, He	1
451 041	Spectral lamp, Cd	1
451 16	Housing for spectral lamps	1
451 30	Universal choke, 230 V, 50 Hz	1
521 210	Transformer, 6/12 V	1
300 02	Stand base, V-shaped, small	1
451 011	Spectral lamp, Ne	1*
451 071	Spectral lamp, Hg-Cd	1*
451 081	Spectral lamp, Tl	1*
451 111	Spectral lamp, Na	1*

\* additionally recommended



Ray path in a grating prism spectrometer (P5.7.1.1)

To assemble the prism spectrometer, a flint glass prism is placed on the prism table of a goniometer. The light of the light source to be studied passes divergently through a collimator and is incident on the prism as a parallel light beam. The arrangement exploits the wavelength-dependency of the refractive index of the prism glass: the light is refracted and each wavelength is deviated by a different angle. The deviated beams are observed using a telescope focused on infinity which is mounted on a slewable arm; this allows the position of the telescope to be determined to within a minute of arc. The refractive index is not linearly dependent on the wavelength; thus, the spectrometer must be calibrated. This is done using e.g. an He spectral lamp, as its spectral lines are known and distributed over the entire visible range.

In the experiment P5.7.1.1, the spectrometer is used to observe the spectral lines of inert gases and metal vapors which have been excited to luminance. To identify the initially "unknown" spectral lines, the angles of deviation are measured and then converted to the corresponding wavelength using the calibration curve.

Note: as an alternative to the prism spectrometer, the goniometer can also be used to set up a grating spectrometer (see P5.7.2.1)



P5.7.2  
GRATING SPECTROMETER

P5.7.2.1  
Measuring the line spectra of inert gases and metal vapors using a grating spectrometer

Measuring the line spectra of inert gases and metal vapors using a grating spectrometer (P5.7.2.1)

Cat. No.	Description	P5.7.2.1
467 23	Spectrometer and goniometer	1
471 23	Ruled grating, 6000/cm (Rowland)	1
451 031	Spectral lamp, He	1
451 111	Spectral lamp, Na	1
451 16	Housing for spectral lamps	1
451 30	Universal choke, 230 V, 50 Hz	1
521 210	Transformer, 6/12 V	1
300 02	Stand base, V-shaped, small	1
451 011	Spectral lamp, Ne	1*
451 041	Spectral lamp, Cd	1*
451 071	Spectral lamp, Hg-Cd	1*
451 081	Spectral lamp, Tl	1*

\* additionally recommended

To create a grating spectrometer, a copy of a Rowland grating is mounted on the prism table of the goniometer in place of the prism. The ray path in the grating spectrometer is essentially analogous to that of the prism spectrometer (see P 5.7.1.1). However, in this configuration the deviation of the rays by the grating is proportional to the wavelength:

$$\sin \Delta\alpha = n \cdot g \cdot \lambda$$

$n$ : diffraction order


$g$ : grating constant

$\lambda$ : wavelength

$\Delta\alpha$ : angle of deviation of  $n$ th-order spectral line

Consequently, the wavelengths of the observed spectral lines can be calculated directly from the measured angles of deviation.

In the experiment P5.7.2.1, the grating spectrometer is used to observe the spectral lines of inert gases and metal vapors which have been excited to luminance. To identify the initially "unknown" spectral lines, the angles of deviation are measured and then converted to the corresponding wavelength. The resolution of the grating spectrometer is sufficient to determine the distance between the two yellow sodium D-lines  $\lambda(D_1) - \lambda(D_2) = 0,60 \text{ nm}$  with an accuracy of 0.10 nm. However, this high resolution is achieved at the cost of a loss of intensity, as a significant part of the radiation is lost in the undiffracted zero order and the rest is distributed over multiple diffraction orders on both sides of the zero order.



Ray path in a grating spectrometer

## P5.7.2

### GRATING SPECTROMETER

#### P5.7.2.2

Assembling a grating spectrometer for measuring transmission curves

#### P5.7.2.3

Assembling a grating spectrometer for measuring spectral lines



Assembling a grating spectrometer for measuring spectral lines (P5.7.2.3)

Cat. No.	Description	P5.7.2 (b)	
		P5.7.2.2	P5.7.2.3
337 47USB	VideoCom	1	1
460 32	Optical bench with standardised profile, 1 m	1	1
460 335	Optical bench with standardised profile, 0,5 m	1	1
460 341	Swivel joint with circular scale	1	1
471 23	Ruled grating, 6000/cm (Rowland)	1	1
460 14	Adjustable slit	1	1
460 08	Lens in frame, f=150 mm	2	1
460 22	Holder with spring clips	1	1
460 373	Optics rider, 60/50	5	5
450 60	Lamp housing with cable	1	
450 511	Bulbs, 6 V/30 W, E14, set of 2	1	
460 20	Condenser with diaphragm holder	1	
521 210	Transformer, 6/12 V	1	
467 95	Colour filter set, primary	1	
467 96	Colour filter set, secondary	1	
468 03	Light filter, red	1*	
468 05	Light filter, yellow	1*	
468 07	Light filter, yellow-green	1*	
468 09	Light filter, blue-green	1*	
460 02	Lens in frame, f=50 mm		1
451 031	Spectral lamp, He		1
451 111	Spectral lamp, Na		1
451 16	Housing for spectral lamps		1
451 30	Universal choke, 230 V, 50 Hz		1
451 011	Spectral lamp, Ne		1*
451 041	Spectral lamp, Cd		1*
451 071	Spectral lamp, Hg-Cd		1*
451 081	Spectral lamp, Tl		1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1


When used in conjunction with a grating spectrometer, the single-line CCD camera VideoCom is ideal for relative measurements of spectral intensity distributions. In such measurements, each pixel of the CCD camera is assigned a wavelength

$$\lambda = d \cdot \sin \alpha$$


in the first diffraction order of the grating. The spectrometer is assembled on the optical bench using individual components. The grating in this experiment is a copy of a Rowland grating with approx. 6000 lines/cm. The diffraction pattern behind the grating is observed with VideoCom. The VideoCom software makes possible comparison of two intensity distributions, and thus recording of transmission curves of color filters or other light-permeable bodies. The spectral intensity distribution of a light source is measured both with and without filter, and the ratio of the two measurements is graphed as a function of the wavelength.

The experiment P5.7.2.2 records the transmission curves of color filters. It is revealed that simple filters are permeable for a very broad wavelength range within the visible spectrum of light, while so-called line filters have a very narrow permeability range.

In the experiment P5.7.2.3, a grating spectrometer is assembled to observe the spectral lines of inert gases and metal vapors which have been excited to luminance. The wavelength and intensity of the spectral lines are measured and compared with literature.



Transmissions curves of various color filters (P5.7.2.2)



## P5.7.2

### GRATING SPECTROMETER

#### P5.7.2.4

Determining the grating constant of a holographic grating with an He-Ne-Laser

#### P5.7.2.5

Investigating the spectrum of a xenon lamp with a holographic grating

Investigating the spectrum of a xenon lamp with a holographic grating (P5.7.2.5\_b)

Cat. No.	Description	P5.7.2.4	P5.7.2.5 (a)	P5.7.2.5 (b)
471 830	He-Ne Laser, linearly polarised	1		
460 01	Lens in frame, f=5 mm	1		
460 08	Lens in frame, f=150 mm	1	1	1
460 09	Lens in frame, f=300 mm	1	1	1
471 27	Holographic grating in frame	1	1	1
441 531	Screen	1	1	1
460 335	Optical bench with standardised profile, 0.5 m	1	1	1
460 32	Optical bench with standardised profile, 1 m	1	1	1
460 341	Swivel joint with circular scale	1	1	1
460 374	Optics rider, 90/50	5	5	6
450 80	Xenon lamp		1	1
450 83	Power supply unit for xenon lamp		1	1
460 02	Lens in frame, f=50 mm		1	1
460 14	Adjustable slit		1	1
460 382	Tilting rider, 90/50		1	1
501 25	Connecting lead, 32 A, 50 cm, red		1	1
501 26	Connecting lead, 32 A, 50 cm, blue		1	1
460 21	Holder for plug-in elements			1
460 22	Holder with spring clips			1
461 62	Set of 2 slit diaphragms			1
578 62	Solar cell, STE 2/19			1
524 013	Sensor-CASSY 2			1
524 220	CASSY Lab 2			1
524 082	Rotary motion sensor S			1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1

To assemble a grating spectrometer with very high resolution and high efficiency a holographic reflection grating with 24000 lines/cm is used. The loss of intensity is small compared to a transmission grating.

In the experiment P5.7.2.4 the grating constant of the holographic reflection grating is determined for different values of the angle of incidence. The light source used is a He-Ne-Laser with the wavelength  $\lambda = 632.8$  nm. The best value is achieved for the special case where angle of incidence and angle of diffraction are the same, the so called Littrow condition.

In the experiment P5.7.2.5 the spectrum of a xenon lamp is investigated. The diffraction pattern behind the holographic grating is recorded by varying the position of a screen or a photocell. The corresponding diffraction angle is read of the circular scale of the rail connector or measured by a rotary motion sensor. It is revealed that the spectrum of the lamp which appears white to the eye consists of a variety of different spectral lines.

## P5.8.2

### BASIC OPTICS

#### P5.8.2.1

Absorption and emission

#### P5.8.2.2

Refraction of light



Refraction of light (P5.8.2.2)

Cat. No.	Description	P5.8.2.1	P5.8.2.2
474 5216	Plano-Convex lens $f = 40$ mm, C25 mount	1	
474 5217	Plano-Convex Lens $f = 60$ mm, C25 Mount	2	
474 161	Absorption Unit	1	
474 5224	Phosphorescing Disc, C25 mount	1	
474 5261	Fluorescent Filters, Set of 3	1	
474 5302	Transmission Grating, 600 lines/mm	1	
474 107	Filter Plate Holder	2	
474 5457	Screen with rider	1	
474 306	Photodetector signal conditioning box	1	
474 108	SiPIN photodetector	1	
531 173	Digital multimeter DMM 121	1	
575 24	Screened cable, BNC/4 mm	1	
474 301	Adaptive Power Supply	1	
474 5411	LED Lamp, White	1	
474 5412	LED Lamp, Red	1	
474 5415	LED Lamp, Blue	1	
474 5442	Profile rail, 500 mm	1	1
474 209	Mounting Plate C25 with Carrier 20 mm	3	
474 2112	Adjustment holder, 4 axes, with stop ang	1	1
474 251	Transport and Storage Box #01	1	1
474 7201	Manual Emission and Absorption	1	
474 133	Optical Fibre Model		1
474 204	Collection of Mounted Models		1
474 5453	Crossed Hair Target in C25 mount		1
474 5418	Diode Laser Module, 532 nm		1
474 121	Swivel Unit with Carrier		1

Cat. No.	Description	P5.8.2.1	P5.8.2.2
474 6411	Mounting plate 40, C25		1
474 7202	Manual Brechung des Lichtes		1

The laws which are related to absorption and emission are discussed and investigated in the experiment P5.8.2.1. With an absorbing sample the Lambert-Beer law is verified by measuring the transmitted light with a photodiode. Light sources of different colours are used to excite fluorescent samples. The absorption and emission of light is visualized spectrally using an optical grating.

Within the frame of the experiment P5.8.2.2 the Snellius Law is verified quantitatively. Deflection, offset and guidance of light travelling in and through transparent materials are demonstrated. The propagation of light through an optical fibre is simulated by a wave guide model. A model for diffraction of light rays on a rain drop is also given to understand the origin of rainbows.



## P5.8.2

### BASIC OPTICS

P5.8.2.6  
Polarisation of light

P5.8.2.7  
Reflection and transmission

Reflection and transmission (P5.8.2.7)

Cat. No.	Description	P5.8.2.6	P5.8.2.7
474 5216	Plano-Convex lens $f = 40$ mm, C25 mount	1	1
474 5260	Optical Quartz Plate in C25 Mount	1	
474 5320	Quarter wave plate, C25	1	
474 5275	Half Wave Plate, C25 Mount	1	
474 1124	Polariser / Analyser with Rotator	2	2
474 306	Photodetector signal conditioning box	1	
474 321	Si PIN Photodetector	1	
531 173	Digital multimeter DMM 121	1	1
575 24	Screened cable, BNC/4 mm	1	1
474 301	Adaptive Power Supply	1	1
474 5411	LED Lamp, White	1	1
474 5412	LED Lamp, Red	1	
474 5418	Diode Laser Module, 532 nm	1	1
474 5442	Profile rail, 500 mm	1	1
474 209	Mounting Plate C25 with Carrier 20 mm	2	1
474 2112	Adjustment holder, 4 axes, with stop ang	1	1
474 251	Transport and Storage Box #01	1	1
474 7206	Manual Polarisation of light	1	
474 6413	Collimation optics in mounting plate 40		1
474 6431	Polarisation analyser 40 mm, VIS		1
474 5453	Crossed Hair Target in C25 mount		1
474 5270	Glass Plate on Rotary Disc		1
474 5271	Dichroic Mirror on Rotary Disc		1
474 5272	Front Face Mirror on Rotary Disc		1
474 5302	Transmission Grating, 600 lines/mm		1
474 6414	Photodetector for Pivot Arm		1

Cat. No.	Description	P5.8.2.6	P5.8.2.7
501 10	BNC adapter, straight		1
474 121	Swivel Unit with Carrier		1
474 238	Carrier for rotatable insert		1
474 7207	Manual Reflection and Transmission		1

Experiment P5.8.2.6 deals with the observation of polarisation of light. The law of Malus is verified and the optical activity of optical crystals are demonstrated. As application of double refractive optics the quarter and half wave plates are subject of measurements and demonstration. As light sources a green laser and a LED are used.

Although the Fresnel Laws are still valid, nowadays mirrors can be produced which seem to bypass these laws. By dielectric coating such mirrors are made either to optimize or to suppress reflection. The functionality of those coatings is based on interference which allows for instance nearly 100 % reflectivity of mirrors (e.g. for laser cavities) or anti-reflecting windows. First the experiment P5.8.2.7 demonstrates the reflection law which is verified on a metal coated mirror. The second part deals with the quantitative verification of the Fresnel Laws on a specially shaped plate and with polarized light. At last the spectral performance of dielectrically coated mirrors is investigated by means of a white light source.

P5.8.2

BASIC OPTICS

P5.8.2.8  
Diffraction of light

P5.8.2.9  
Interference of light



Diffraction of light (P5.8.2.8)


Cat. No.	Description	P5.8.2.8	P5.8.2.9
474 5216	Plano-Convex lens $f = 40$ mm, C25 mount	1	1
474 5263	Beam expander 6x	1	1
474 5298	Circular apertures	1	
474 5299	Gauze 300 mesh	1	
474 5300	Single Slit 0.06 mm	1	
474 5301	Double slit	1	
474 6417	Optical Screen with XY Scale	1	
474 6414	Photodetector for Pivot Arm	1	
531 183	Digital multimeter 3340	1	
575 24	Screened cable, BNC/4 mm	1	
501 10	BNC adapter, straight	1	
474 301	Adaptive Power Supply	1	1
474 5411	LED Lamp, White	1	1
474 5418	Diode Laser Module, 532 nm	1	1
474 5442	Profile rail, 500 mm	1	1
474 121	Swivel Unit with Carrier	1	
474 209	Mounting Plate C25 with Carrier 20 mm	1	4
474 2112	Adjustment holder, 4 axes, with stop ang	1	1
474 251	Transport and Storage Box #01	1	1
474 7208	Manual Diffraction of light	1	
474 5221	Biconcave lens $f = -20$ mm, C25 mount		1
474 5251	Fresnel Zone Plate, C25 mount		1
474 5252	Fabry Perot Insert, C25 mount		1
474 5277	Newton's Rings Optics		1

Cat. No.	Description	P5.8.2.8	P5.8.2.9
474 176	Fresnel Mirror Assembly		1
474 5457	Screen with rider		1
474 213	Adjustment Holder 1 inch, left		1
474 7209	Manual Interference of light		1

In the experiment P5.8.2.8 Fresnel and Fraunhofer types of diffraction are discussed. Investigations are performed using monochromatic laser light which will be diffracted on slits and holes of various widths and gratings. Thin wires show impressively the Babinet theorem stating that complementary masks result in the same diffraction pattern. The obtained diffraction patterns will be observed on a target screen and can be measured by a photodiode quantitatively.

Different examples of interference phenomenon are discussed and demonstrated in the experiment P5.8.2.9. Fresnel mirror, wedges and half-lenses are tools which "divide" one light source into two and superimpose their coherent portions. On a set-up proposed by Newton, interference caused by thin layers can be determined quantitatively. Since diffraction usually generates interference patterns, a Fresnel plate is used for illustrating this effect. Finally, a model of a Fabry Perot resonator demonstrates the working principle of wavelength selection in a cavity.





## P5.8.3

### OPTICAL APPLICATIONS

P5.8.3.1  
Optical interferometer

P5.8.3.2  
Refractometer

Optical interferometer (P5.8.3.1)

Cat. No.	Description	P5.8.3.1	P5.8.3.2
474 5220	Biconcave Lens $f = -10$ mm, C25 mount	1	
474 5264	Beam expander 2.7x	1	
474 169	Gas Cuvette Assembly	1	
474 171	Mach Zehnder Beam Combining Assembly	1	
474 174	Mach-Zehnder Beam Splitting Assembly	1	
474 5457	Screen with rider	1	1
474 5418	Diode Laser Module, 532 nm	1	
474 5441	Profile Rail, 300 mm	2	
474 5442	Profile rail, 500 mm	1	1
474 5449	Angle Joint, Cross Piece	1	
474 209	Mounting Plate C25 with Carrier 20 mm	3	1
474 2112	Adjustment holder, 4 axes, with stop ang	1	1
474 251	Transport and Storage Box #01	2	1
474 7210	Manual Optical Interferometer	1	
474 5225	Polariser, C25 mount		1
474 404	Prism assembly		1
474 405	Beam Bending Assembly 22.5 °deg;- 45°, left		1
474 406	Beam Bending Assembly 22.5 °deg;- 45°, right		1
474 5307	Set of Test Liquids		1
474 301	Adaptive Power Supply		1
474 5413	LED amber in C25 housing		1
474 7211	Manual Refractometer		1

While the Michelson interferometer is mainly used to determine the movement of a reflecting object in a nm scale, the Mach-Zehnder interferometer investigates transparent objects and is particularly useful for studying liquid or gas dynamics. Since the Mach-Zehnder is a unidirectional interferometer it is especially useful for measurements where the samples have to be traversed only once or in one direction. Within the frame of the experiment P5.8.3.1 both, a Michelson- and a Mach-Zehnder interferometer will be realized. The former demonstrates the principle generation of interference patterns and the use of fringe counting in metrology. The latter uses the interference pattern to visualize changes of the index of refraction as a result of changes in physical properties of gasses, like pressure or composition.

The refractometer is an essential instrument for empirical identification of pure substances, purity measurements, and quantitative analysis of solutions. It is used in the food and beverage industry as well as in the chemical and pharmaceutical industry for quality control and check-ups on compound identities. The experiments in P5.8.3.2 are performed on a model of an Abbé refractometer to show the implementation on modern, computer driven standard refractometers. This refractometer measures refraction of liquids, solutions and solid materials by observing the shadow line of the light beam on a screen.

P5.8.3

OPTICAL APPLICATIONS

P5.8.3.3  
Holography



Holography (P5.8.3.3)

Cat. No.	Description	P5.8.3.3
474 5219	Biconcave Lens $f = -5$ mm, C25 mount	2
474 5220	Biconcave Lens $f = -10$ mm, C25 mount	2
474 6418	Model for Holography	1
474 6419	Set of Tools for Holography	1
474 5303	Photographic plate, 532 nm, 63 x 63 mm	30
474 164	Set of Development Equipment	1
474 167	Holography Combining Assembly	1
474 168	Holography Beam Splitting Assembly	1
474 3024	Digital Laser Controller and Timer	1
474 5419	Diode Laser Module, 532 nm with Window	1
474 5442	Profile rail, 500 mm	2
474 209	Mounting Plate C25 with Carrier 20 mm	1
474 211	Adjustment holder, 4 axes, carrier 20 mm	1
474 165	Set of Holography Chemicals	1
675 3400	Water, pure, 1 l	1
671 0010	Pyrocatechol, 100 g	1
675 3270	Vitamin C, 50 g	1
673 7800	Sodium sulfite, 100 g	1
672 1700	Urea, 100 g	1
673 5600	Natrium carbonate, anhydrous, 100 g	1
474 251	Transport and Storage Box #01	1
474 7212	Manual Holography	1
474 5253	Photodetector, mini Si PIN, connection lead	1*
531 183	Digital multimeter 3340	1*
575 24	Screened cable, BNC/4 mm	1*
501 10	BNC adapter, straight	1*

\* additionally recommended

In the experiment P5.8.3.3 both types, transmission and reflection holograms can be recorded. In a first set-up the stability of the working place and environment is tested by an interferometer. It is extremely important for making holograms that the set-up has to be stable within the range of the optical wavelengths. The recording techniques of transmission and reflection holography will be discussed and experimentally investigated. Finally the development of the holographic plates will be performed.



Spectral analysis (P5.8.3.5)

## P5.8.3

### OPTICAL APPLICATIONS

P5.8.3.4  
Diffraction gratings

P5.8.3.5  
Spectral analysis

P5.8.3.7  
LED and laser diode

Cat. No.	Description	P5.8.3.4	P5.8.3.5	P5.8.3.7
474 5216	Plano-Convex lens $f = 40$ mm, C25 mount	1	1	1
474 5256	Biconvex lens $f = 60$ mm, C25 Mount	1		
474 5263	Beam expander 6x	1		
474 5264	Beam expander 2.7x	1		
474 5268	Transmission gratings, Set of 5	1		
474 6417	Optical Screen with XY Scale	1		1
474 306	Photodetector signal conditioning box	1	1	
474 321	Si PIN Photodetector	1		
531 173	Digital multimeter DMM 121	1	1	1
575 24	Screened cable, BNC/4 mm	1	1	1
474 5417	Spectral Lamp with Slit and Power Supply	1	1	
474 5418	Diode Laser Module, 532 nm	1		
474 5442	Profile rail, 500 mm	1	2	1
474 121	Swivel Unit with Carrier	1		1
474 6411	Mounting plate 40, C25	2		
474 209	Mounting Plate C25 with Carrier 20 mm	1	2	2
474 2112	Adjustment holder, 4 axes, with stop ang	1	1	
474 251	Transport and Storage Box #01	1	2	1
474 7213	Manual Diffraction gratings	1		
474 5218	Biconvex Lens $f = 20$ mm, C25-T Mount		1	
474 5211	Acrylic Absorption Filter		1	
474 177	Spectrometer Mirror Assembly		1	
474 178	Spectrometer Grating Assembly		1	
474 107	Filter Plate Holder		1	
474 108	SiPIN photodetector		1	
474 301	Adaptive Power Supply		1	1
474 5411	LED Lamp, White		1	1
474 7214	Manual Spectral analysis		1	
474 5222	Cylindrical Lens $f = 25$ mm, C25 Mount			1

Cat. No.	Description	P5.8.3.4	P5.8.3.5	P5.8.3.7
474 5223	Cylindrical Lens $f=80$ mm, C25 Mount			1
474 6431	Polarisation analyser 40 mm, VIS			1
474 5302	Transmission Grating, 600 lines/mm			1
474 6414	Photodetector for Pivot Arm			1
501 10	BNC adapter, straight			1
474 5412	LED Lamp, Red			1
474 5415	LED Lamp, Blue			1
474 5420	Dimo diode laser module, 630 nm (red)			1
474 2114	Adjustment holder, 4 axes, rotary insert			1
474 213	Adjustment Holder 1 inch, left			1
474 7220	Manual LED and Laser Diode			1

In the experiment P5.8.3.4, transmission gratings of different grating constants are investigated. A light source with known wavelength is used to characterize a specific grating. Finally, a two-dimensional grating is not only used to produce impressive patterns of light spots, but gives an idea about the principles of x-ray diffraction on crystal lattices or atomic layers.

The experiment P5.8.3.5 builds a model for a standard grating monochromator as well as a spectrograph. With a white light lamp a source for a continuous spectrum is provided. Here the spectrometer is used as a spectrograph and grating diffraction of first and higher orders as well as absorption spectroscopy can be demonstrated. Using the set-up in the monochromator mode in combination with a spectral lamp, line spectra are detected and features like spectral resolution or line profile are measured.

For a comprehensive study of the properties of LED and laser diodes the setup of the experiment P5.8.3.7 comes with four different light sources. One diode laser having a wavelength of 630 nm and three high brightness LED emitting white, red and blue radiation. The respective light source is plugged into the four axes adjustment holder and connected to the adaptive power supply. One lens is used for the beam collimation and cylindrical lenses for transforming the elliptical beam of the laser diode into an almost circular one.

## P5.8.4

### OPTICAL IMAGING AND COLOUR

P5.8.4.1  
Optical filter

P5.8.4.3  
Camera and imaging




Optical filter (P5.8.4.1)

Cat. No.	Description	P5.8.4.1	P5.8.4.3
474 5216	Plano-Convex lens f = 40 mm, C25 mount	1	
474 5217	Plano-Convex Lens f = 60 mm, C25 Mount	1	
474 5289	Interference Filter 550 nm, in C25 Mount	1	
474 5262	Optical filters, Set of 8	1	
474 5302	Transmission Grating, 600 lines/mm	1	
474 107	Filter Plate Holder	1	
474 6417	Optical Screen with XY Scale	1	
474 306	Photodetector signal conditioning box	1	
474 321	Si PIN Photodetector	1	
531 173	Digital multimeter DMM 121	1	
575 24	Screened cable, BNC/4 mm	1	
474 301	Adaptive Power Supply	1	1
474 5411	LED Lamp, White	1	
474 5416	LED Lamp NIR in C25 housing	1	1
474 5442	Profile rail, 500 mm	1	1
474 121	Swivel Unit with Carrier	1	
474 6411	Mounting plate 40, C25	1	
474 209	Mounting Plate C25 with Carrier 20 mm	2	1
474 2112	Adjustment holder, 4 axes, with stop ang	1	
474 213	Adjustment Holder 1 inch, left	1	
474 251	Transport and Storage Box #01	1	1
474 7216	Manual Optische Filter	1	
474 104	Focussing Optics, f = 60 mm		1
468 75	Filter, infrared barrier		1
474 281	CCD day and night camera module		1
474 9112	CCD Camera Control Software		1

Cat. No.	Description	P5.8.4.1	P5.8.4.3
524 004	Adapter, USB port/serial port		1
474 5467	Flat panel TV 19 inch		1
474 7218	Manual Camera		1
	additionally required: PC with Windows XP or higher		1

In the experiment P5.8.4.1 different kinds of filters are presented. With a set of colour filters transmission and absorption of certain spectral ranges are demonstrated. Neutral density filters are provided to dim the whole spectral range. An infrared LED and an IR filter demonstrate light filtering in the IR range.

In the experiment P5.8.4.3 a high performance industrial CCD zoom camera with computer interface is subject of a variety of investigations. The rapid development in the area of CCD sensors created a great variety of new possibilities. Most of them are introduced here and one get experienced in the manifold of parameters which needs to be set according to the requirements. The CCD camera used can be operated as day as well as night vision camera. In the latter case an IR blocking filter is switched out of the way between the objective and the CCD chip enabling the sensitivity in the near infra red region (NIR). The camera is fully controlled by a PC, the video output is connected to a TFT monitor.



## P5.8.5

### LASER BASICS

#### P5.8.5.1

Laser Doppler Anemometry  
with CASSY

Laser Doppler Anemometry with CASSY (P5.8.5.1)

Cat. No.	Description	P5.8.5.1
471 821	Head unit for 5 mW He-Ne laser	1
471 825	Power supply for 5 mW He-Ne-laser	1
470 010	Laser holder for 5 mW He-Ne laser	1
473 431	Holder for beam divider	1
473 432	Beam divider, 50 %	1
473 461	Planar mirror with fine adjustment	1
460 02	Lens in frame, f=50 mm	1
460 03	Lens in frame, f=100 mm	1
460 21	Holder for plug-in elements	1
460 22	Holder with spring clips	2
460 26	Iris diaphragm	1
461 63	Set of 4 different diaphragms	1
469 96	Diaphragm with 3 diffraction holes	1
441 53	Screen, translucent	1
460 335	Optical bench with standardised profile, 0.5 m	1
460 32	Optical bench with standardised profile, 1 m	1
460 374	Optics rider, 90/50	10
460 380	Cantilever arm	1
460 385	Extension rod	1
311 77	Steel tape measure, 2 m	1
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>
524 220	CASSY Lab 2	1
558 835	Silicon photodetector	1
532 20	AC/DC amplifier 30 W	1
577 68	Resistor, 100 kΩ, STE 2/19	1
575 24	Screened cable, BNC/4 mm	1
501 641	Two-way adapters, red, set of 6	1

Cat. No.	Description	P5.8.5.1
590 02ET2	Clip plugs, small, set of 2	1
683 70	Reflecting particles, 10 g	1
664 146	Reaction tube	1
602 404	Separation Funnel, 500 ml	1
604 433	Silicone tubing, 7 mm diam., 1 m	2
667 175	Hofmann tubing clamp, 20 mm	1
604 5672	Double microspatula, steel, 150 mm	1
602 010	Beaker, Boro 3.3, 150 ml, tall	1
604 215	Measuring beaker, clear SAN 500 ml	1
300 01	Stand base, V-shaped, large	1
300 44	Stand rod, 100 cm, 12 mm diam.	1
666 546	Stand ring with clamp, 100 mm diam.	1
500 401	Connecting lead, 19 A, 10 cm, red	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1
471 828	Adjustment goggles for He-Ne-laser	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

\* additionally recommended


Laser Doppler anemometry is a non-contact optical measurement method to obtain the velocity of a flow (fluid, gas). In the experiment P5.8.5.1 a laser Doppler anemometer is assembled. Measurements of the flow velocity of a fluid in a tube are conducted by measuring the velocity of small particles carried along in the flow. Moving through the measuring volume the particles scatter light of a laser. The scattered light is frequency shifted due to the Doppler effect. The frequency shift is determined and converted into the particle velocity, *i.e.* the flow velocity.

P5.8.5

LASER BASICS

P5.8.5.2  
Laser safety

P5.8.5.3  
Emission & absorption /  
Optical pumping



Emission & absorption / Optical pumping (P5.8.5.3)


Cat. No.	Description	P5.8.5.2	P5.8.5.3
474 5216	Plano-Convex lens f = 40 mm, C25 mount	1	
474 5220	Biconcave Lens f = -10 mm, C25 mount	1	
474 217	Scatter Probe with Holder	1	
468 77	Light filter, green	1	
474 107	Filter Plate Holder	1	1
474 306	Photodetector signal conditioning box	1	1
474 321	Si PIN Photodetector	1	
531 173	Digital multimeter DMM 121	1	
575 24	Screened cable, BNC/4 mm	1	1
474 5464	Oscilloscope, Dual Channel, Digital	1	1*
501 06	HF-Cable, BNC-BNC, 1.5 m	1	1*
474 5460	Laser Power Meter	1	
474 5462	Laser Power Sensor, 1 nW ... 50 mW	1	
474 5463	Laser Energy Sensor, 300 nJ ... 600 µJ	1	
474 309	Controller for Pulsed Laser Diode	1	
474 5428	Pulsed Diode Laser Module 908 nm	1	
474 5418	Diode Laser Module, 532 nm	1	
474 5442	Profile rail, 500 mm	1	1
474 121	Swivel Unit with Carrier	1	
474 6411	Mounting plate 40, C25	2	
474 209	Mounting Plate C25 with Carrier 20 mm	1	
474 211	Adjustment holder, 4 axes, carrier 20 mm	1	
474 2112	Adjustment holder, 4 axes, with stop ang	1	
474 122	Optics cleaning set	1	1
671 9700	Ethanol, absolute, 250 ml	1	1
474 251	Transport and Storage Box #01	2	1
474 7101	LIT: Laser Safety	1	
474 1031	Module B - Collimating optics on carrier MG-65 with Melles Griot		1
474 104	Focussing Optics, f = 60 mm		1

Cat. No.	Description	P5.8.5.2	P5.8.5.3
474 5310	Crystal in holder Nd:YAG 1064 nm		1
474 113	Laser Mirror Adjustment Holder, left		1
474 137	Spatial filter with adjustable iris		1
474 5453	Crossed Hair Target in C25 mount		1
468 74	Filter, infrared		1
474 4025	IR converter screen 0.8 - 1.6 µm		1
474 108	SiPIN photodetector		1
531 183	Digital multimeter 3340		1
474 302	Controller for Diode Laser		1
474 1021	Single Mode Diode Laser Head with Adjust		1
474 7102	LIT: Emission & Absorption/Opt. Pumping		1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m		1*

\* additionally recommended

Laser can emit dangerous radiation. To protect against any injury international safety rules like IEC 60825 or ANSI Z136 were defined. The lasers are classified into different classes with the individual maximum permissible exposure limit (MPE) which is defined as intensity, power per square centimetre (W/cm<sup>2</sup>). To classify a laser one needs to know its intensity in order to compare it with the MPE values. For pulsed laser its energy is used instead of the power of continuous wave (cw) laser. Within the experiment P5.8.5.2 one pulsed and one cw laser is classified. For this purpose the intensity needs to be determined. To calculate this value the power and the beam diameter are measured. For the pulsed laser its repetition rate and the emitted energy is measured.

The experiment P5.8.5.3 introduces to optical pumping as well as emission and absorption. Due to the pumping process spontaneous and stimulated emission is generated. The emission is temporarily as well spectroscopical measured and analyzed. The tuning of the emission wavelength of the pump diode laser due to the temperature allows the recording of the absorption spectrum. From the timely decay of the fluorescence light the lifetime of the excited state is measured and the Einstein coefficient for stimulated emission calculated.



## P5.8.5

### LASER BASICS

#### P5.8.5.4

Fabry Perot spectrum analyser

Fabry Perot spectrum analyser, basic (P5.8.5.4\_b)

Cat. No.	Description	P5.8.5.4 (b)
474 104	Focussing Optics, $f = 60$ mm	1
474 5236	Laser Mirror, $R = 100$ mm, M16 Mount	1
474 5239	Laser Mirror, $R = 100$ mm, M12 Mount	1
474 113	Laser Mirror Adjustment Holder, left	1
474 317	Piezo Actuator Controller	1
474 139	Piezo Element with Adjustment Holder	1
522 561	Function generator P	1
501 02	BNC cable, 1 m	2
501 091	BNC T adapter	1
575 24	Screened cable, BNC/4 mm	1
474 306	Photodetector signal conditioning box	1
474 108	SiPIN photodetector	1
474 5464	Oscilloscope, Dual Channel, Digital	1
474 303	HeNe Laser High Voltage supply, adjustable	1
474 5421	HeNe Pilot Laser $\varnothing 30$ mm	1
474 5445	Profile Rail, 1000 m	1
474 210	Mounting plate $\varnothing 30$ mm, carrier 20 mm	2
474 122	Optics cleaning set	1
671 9700	Ethanol, absolute, 250 ml	1
474 251	Transport and Storage Box #01	2
474 7103	Manual Fabry Perot Resonator	1
474 1404	Lens $f = 150$ mm, C25 Mount and Plate	1*
474 140	Beam expander lens, adjustment holder	1*
474 5213	Achromat $f = 20$ mm, C25 Mount	1*
474 5234	Laser Mirror, flat, M16 mount	1*

Cat. No.	Description	P5.8.5.4 (b)
474 5235	Laser Mirror, $R = 75$ mm, M16 Mount	1*
474 5237	Laser Mirror, flat, M12 Mount	1*
474 5238	Laser Mirror, $R = 75$ mm, M12 Mount	1*

\* additionally recommended

A Fabry Perot resonator is formed by two mirrors aligned parallel to each other. The resulting cavity changes transparency when the distance of the mirrors is changed by a multiple of half the wavelength. Scanning the length of the cavity creates a high resolution spectrometer. Experiment P5.8.5.4 shows an open frame scanning Fabry Perot. As scanner a Piezo element and as probe a two mode Helium Neon laser is used. The mode spectrum of the Helium Neon laser is displayed on an oscilloscope and the characteristic parameter like finesse, free spectral range, resolution and contrast are measured and discussed. Additional components for beam expansion are used to investigate the effect of technical Finesse. Additional mirrors are used to show the difference of a plane and confocal Fabry Perot arrangement.

P5.8.5


LASER BASICS

P5.8.5.5

Helium Neon laser

P5.8.5.6

Laser frequency stabilisation



Helium Neon laser, basics (P5.8.5.5\_b)

Cat. No.	Description	P5.8.5.5 (b)	P5.8.5.6
474 5243	Laser Mirror VIS 700, M16 Mount	1	
474 5246	Laser Mirror, flat, M16 Mount	1	
474 113	Laser Mirror Adjustment Holder, left	1	
474 114	Laser Mirror Adjustment Holder, right	1	
474 137	Spatial filter with adjustable iris	1	
474 1082	Module G (SiPIN) photodetector on swivel arm	1	
531 183	Digital multimeter 3340	1	
575 24	Screened cable, BNC/4 mm	1	
501 10	BNC adapter, straight	1	
474 303	HeNe Laser High Voltage supply, adjustable	1	1
474 127	Main Laser Tube with XY-Adjustment	1	
474 5422	Pilot laser 532 nm (green)	1	
474 5445	Profile Rail, 1000 m	1	
474 122	Optics cleaning set	1	
671 9700	Ethanol, absolute, 250 ml	1	
474 251	Transport and Storage Box #01	2	1
474 7104	Manual HeNe Laser	1	
474 5242	Laser Mirror OC 632, flat, M16 Mount	1*	
474 5244	Laser Mirror VIS 1000, M16 Mount	1*	
474 5245	Laser Mirror IR 713, M16 Mount	2*	
474 4025	IR converter screen 0.8 - 1.6 $\mu$ m	1*	
474 126	Littrow prism with adjustment holder	1*	
474 141	Single Mode Etalon, Adjustmend Holder	1*	
474 142	Birefringent tuner	1*	
474 409	Mode separator and photodiode		1
474 112	Polarisation Analyzer		1

Cat. No.	Description	P5.8.5.5 (b)	P5.8.5.6
474 306	Photodetector signal conditioning box		1
474 108	SiPIN photodetector		1
474 5464	Oscilloscope, Dual Channel, Digital		1
501 06	HF-Cable, BNC-BNC, 1.5 m		1
474 312	Laser frequency stabilizer LSF-01		1
474 410	Two mode HeNe-laser		1
474 5442	Profile rail, 500 mm		1
474 7111	Manual Laser frequency stabilisation		1
	additionally required: PC with Windows XP or higher		1

\* additionally recommended

The experiment P5.8.5.5 realises a Helium Neon laser from basic parts. The open frame cavity allows the variation of parameters to measure the beam profile for different cavity mirror configurations and distances. The laser tube is equipped with a Brewster's window on both sides allowing the study of polarisation and losses. Optionally, a Littrow prism selects other wavelengths than the main laser line at 632 nm, especially the orange line at 611 nm. Further line tuning is demonstrated by means of the optional birefringent tuner. 4 different lines can be demonstrated by this element. An etalon is used to obtain the single mode operation of the Helium Neon laser.

Stabilizing a gas laser emission against thermal change of the resonator is demonstrated in Experiment P5.8.5.6. Two longitudinal modes of a short Helium Neon laser are analysed and centered to the laser transition gain profile, thereby stabilizing the wavelength of emission. A PI controller with variable coefficients is used for the stabilisation.



P5.8.6  
SOLID STATE LASER

P5.8.6.1  
Diode laser



Diode laser (P5.8.6.1)

Cat. No.	Description	P5.8.6.1
474 5266	Collimating cylindrical lens $f = 20$ mm	1
474 5267	Collimating cylindrical lens $f = 80$ mm	1
474 1031	Module B - Collimating optics on carrier MG-65 with Melles Griot	1
474 5310	Crystal in holder Nd:YAG 1064 nm	1
474 113	Laser Mirror Adjustment Holder, left	1
474 112	Polarisation Analyzer	1
474 5453	Crossed Hair Target in C25 mount	1
474 4025	IR converter screen $0.8 - 1.6 \mu\text{m}$	1
474 306	Photodetector signal conditioning box	1
474 108	SiPIN photodetector	1
531 183	Digital multimeter 3340	1
575 24	Screened cable, BNC/4 mm	1
474 302	Controller for Diode Laser	1
474 1012	Diode Laser Head, Dual Axes Rotary Mount	1
474 5442	Profile rail, 500 mm	1
474 209	Mounting Plate C25 with Carrier 20 mm	2
474 122	Optics cleaning set	1
671 9700	Ethanol, absolute, 250 ml	1
474 251	Transport and Storage Box #01	1
474 7105	Manual Diode Laser	1
474 5464	Oscilloscope, Dual Channel, Digital	1*
501 06	HF-Cable, BNC-BNC, 1.5 m	1*
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	1*

\* additionally recommended

The goal of the experiment P5.8.6.1 is the study of the properties of a laser diode, i.e. the characteristic parameters like the output power and wavelength as function of the temperature. In a next step the spatial intensity distribution is measured. The more or less elliptical beam is formed by means of two cylindrical lenses into an almost circular beam.

## P5.8.6

### SOLID STATE LASER

#### P5.8.6.2

Diode laser pumped Nd:YAG laser

#### P5.8.6.3

Frequency doubling,  
1064 nm → 532 nm

#### P5.8.6.4

Frequency doubling,  
1320 nm → 660 nm

#### P5.8.6.5

Q-switch operation



Frequency doubling, 1064 nm → 532 nm (P5.8.6.3)

Cat. No.	Description	P5.8.6.2	P5.8.6.3	P5.8.6.4	P5.8.6.5
474 1031	Module B - Collimating optics on carrier MG-65 with Melles Griot	1	1	1	1
474 104	Focussing Optics, f = 60 mm	1	1	1	1
474 5310	Crystal in holder Nd:YAG 1064 nm	1	1		1
474 5311	Laser mirror in holder SHG 100	1	1		1
474 113	Laser Mirror Adjustment Holder, left	1	1	1	1
474 114	Laser Mirror Adjustment Holder, right	1	1	1	1
474 5453	Crossed Hair Target in C25 mount	1	1	1	1
474 107	Filter Plate Holder	1	1	1	1
468 74	Filter, infrared	1	1	1	1
474 4025	IR converter screen 0.8 - 1.6 μm	1	1	1	1
474 306	Photodetector signal conditioning box	1	1	1	1
474 108	SiPIN photodetector	1	1		1
474 5464	Oscilloscope, Dual Channel, Digital	1	1	1	1
501 06	HF-Cable, BNC-BNC, 1.5 m	1	1	1	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	1	1	1	1
474 302	Controller for Diode Laser	1	1	1	1
474 102	Diode Laser Head with Adjustment Holder	1	1	1	1
474 5442	Profile rail, 500 mm	1	1	1	1
474 122	Optics cleaning set	1	1	1	1
671 9700	Ethanol, absolute, 250 ml	1	1	1	1
474 251	Transport and Storage Box #01	1	1	1	1
474 7106	Manual DPSSL	1	1	1	1
474 109	KTP Crystal with Adjustment Holder		1		
468 77	Light filter, green		1		
474 137	Spatial filter with adjustable iris		1*	1*	
474 1094	Red SHG Crystal in holder			1	
474 5240	Nd:YAG rod 1.3 μm in mirror holder			1	
474 5241	Mirror SHG 1.3 μm			1	
474 5290	Coloured glass filter KG5, 50 x 50 x 3 mm			1	

Cat. No.	Description	P5.8.6.2	P5.8.6.3	P5.8.6.4	P5.8.6.5
474 1081	Modul G InGAAS photodetector			1	
474 110	Module P - Crystal for passive q-switch operation				1
474 1804	Light Chopper on Carrier with Controller				1*
474 264	Upgrade Kit for Active q-Switch				1*

\* additionally recommended

Experiment P5.8.6.2 builds a diode pumped Nd:YAG laser ground up. First the pump diode laser is characterized. Then the process of optical pumping and the emitted spontaneous fluorescence are analysed spectrally and temporally by modulation and changing the wavelength of the pump laser leading to the Einstein coefficients. In a third step the laser operation is initiated by adding the second cavity mirror. The laser threshold and efficiency are determined and by modulating the pump laser diode the so called spiking effect demonstrated. By changing the length of the laser cavity the stability criterion is verified.

Experiment P5.8.6.3 enhances the basic set-up of the diode pumped Nd:YAG laser (P5.8.6.2) by a KTP crystal module placed into the laser cavity leading to a frequency doubled green (532 nm) visible output. By adding the optional adjustable iris the transverse mode structure can be controlled for various TEM structures down to TEM00.

Experiment P5.8.6.4 is a variation of P5.8.6.3, but using differently coated mirrors and differently cut KTP, the Nd:YAG laser operates at 1320 nm and the frequency doubled visible light is red at a wavelength of 660 nm.

In the experiment P5.8.6.5 the basic set-up of the diode pumped Nd:YAG laser (P5.8.6.2) is enhanced by a passive saturable absorber module placed into the laser cavity. The initial absorption of the Cr:YAG crystal prevents continuous laser oscillation. Build-up of inversion will saturate the absorber and repeatedly create a giant and short laser puls. Optionally an active Q-switch can be used to trigger such a pulse externally.



## P5.8.6 SOLID STATE LASER

P5.8.6.6  
Pulsed diode laser

P5.8.6.7  
Diode pumped  
Nd:YVO4 micro laser

Pulsed diode laser (P5.8.6.6)

Cat. No.	Description	P5.8.6.6	P5.8.6.7
474 1031	Module B - Collimating optics on carrier MG-65 with Melles Griot	1	1
474 6412	Polarisation analyzer 40 mm	1	
474 4025	IR converter screen 0.8 - 1.6 $\mu\text{m}$	1	1
474 6414	Photodetector for Pivot Arm	1	
531 183	Digital multimeter 3340	1	
474 5464	Oscilloscope, Dual Channel, Digital	1	1
575 24	Screened cable, BNC/4 mm	1	
501 091	BNC T adapter	1	
474 341	BNC load resistor 50 Ohm	1	
474 309	Controller for Pulsed Laser Diode	1	
474 5428	Pulsed Diode Laser Module 908 nm	1	
474 5442	Profile rail, 500 mm	1	1
474 121	Swivel Unit with Carrier	1	
474 2114	Adjustment holder, 4 axes, rotary insert	1	
474 251	Transport and Storage Box #01	1	1
474 7109	LIT: Pulsed Diode Laser	1	
474 5460	Laser Power Meter	1*	
474 5463	Laser Energy Sensor, 300 nJ ... 600 $\mu\text{J}$	1*	
474 104	Focussing Optics, f = 60 mm		1
474 241	Monolithic Nd:YVO4 + KTP Core Laser		1
474 5453	Crossed Hair Target in C25 mount		1
474 107	Filter Plate Holder		1
468 74	Filter, infrared		1
474 306	Photodetector signal conditioning box		1
474 108	SiPIN photodetector		1
501 06	HF-Cable, BNC-BNC, 1.5 m		1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m		1

Cat. No.	Description	P5.8.6.6	P5.8.6.7
474 302	Controller for Diode Laser		1
474 102	Diode Laser Head with Adjustment Holder		1
474 122	Optics cleaning set		1
671 9700	Ethanol, absolute, 250 ml		1
474 7127	LIT: Micro Laser		1

\* additionally recommended

Pulsed diode lasers emit short pulses with a pulse width of 10 ... 100 nanoseconds. Similar to a flash lamp the laser can emit a very high peak power for a short time. Experiment P5.8.6.6 analyses the temporal and spatial properties of a diode laser emitting a peak power of 70 W within a pulse width of 100 ns. The electrical as well as optical pulse is monitored on an digital oscilloscope.

Green laser light is still widely produced by optical pumping and frequency doubling. Within this experiment P5.8.6.7 a so called GCL (green core laser), consisting of a Neodymium Yttrium Vanadate (Nd:YVO4) crystal which is cemented to a KTP crystal is used. Its small size of only 1.3 mm x 1.3 mm x 3 mm justifies the term „Micro Laser“. The GCL is pumped by a diode laser as known from the optical pumping of the Nd:YAG Laser (P5.8.6.2).

P5.8.7

OPTICAL FIBRES

P5.8.7.1  
Fibre laser

P5.8.7.4  
Erbium doped fibre amplifier



Fibre laser (P5.8.7.1)

Cat. No.	Description	P5.8.7.1	P5.8.7.4
474 104	Focussing Optics, $f = 60$ mm	1	
474 5308	Bandpass filter $1.5 \mu\text{m}$ in C25	1	
474 190	Fused WDM Coupler 980/1550 nm	1	
474 191	Fibre collimator with ST connector, left	1	
474 192	Fibre collimator with ST connector, right	1	
474 194	ST coupler in C25 mounted	1	
474 1898	Erbium Doped Fibre Module, 8 m	1	
474 246	Output coupling module	1	
474 5293	SM Fibre 100 m on drum, ST connector	1	
474 5296	Fibre Patch Cable ST/ST, Length 0.25 m	3	
474 5297	Fibre Patch Cable ST/ST, Length 1 m	1	
474 4025	IR converter screen $0.8 - 1.6 \mu\text{m}$	1	1
474 306	Photodetector signal conditioning box	1	1
474 108	SiPIN photodetector	1	
474 1084	InGAAS Photodetector	1	1
474 5464	Oscilloscope, Dual Channel, Digital	1	1
501 06	HF-Cable, BNC-BNC, 1.5 m	1	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	1	2
474 302	Controller for Diode Laser	1	2
474 5426	Diode laser module, ST fibre connector	1	
474 5442	Profile rail, 500 mm	1	2
474 5444	Profile Rail, 1000 mm	1	
474 209	Mounting Plate C25 with Carrier 20 mm	1	2
474 251	Transport and Storage Box #01	2	1
474 7110	LIT: Fibre Laser	1	
474 189	Erbium doped fibre 2 m module, ST connect	1*	
474 1894	Erbium Doped Fibre Module, 4 m	1*	

Cat. No.	Description	P5.8.7.1	P5.8.7.4
474 1896	Erbium Doped Fibre Module, 16 m	1*	
474 5278	Passive mode locker module $1.5 \mu\text{m}$	1*	
474 5279	Optical isolator, $1.5 \mu\text{m}$ , SM Fiber, ST	1*	
474 5254	Collimating Optics, high NA		1
474 1036	Collimating Optics on Carrier		1
474 151	Coupling Optics, XY- Adjustment Holder		1
474 156	Dichroic Beam Combiner 980/1550 nm		1
474 157	Erbium doped fibre 17 m with holder		1
474 321	Si PIN Photodetector		1
474 5304	Diode Laser Head 980 nm		1
474 5305	Diode Laser Head 1550 nm		1
474 211	Adjustment holder, 4 axes, carrier 20 mm		2
474 7120	LIT: Erbium doped Fibre Amplifier (EDFA)		1

\* additionally recommended

In the experiment P5.8.7.1, an Erbium doped fibre is used as active material. Connected to a pump laser, the fluorescence from the erbium fiber is analysed. Before the laser operation at  $1.5 \mu\text{m}$  is studied, the lifetime of the excited state is measured. To form a ring laser a WDM is used to couple the pump light into the fibre and to close the laser ring structure. The ring is opened where a thin glass plate couples a small fraction of the clockwise (cw) and counter clockwise (ccw) laser modes towards a detector.

Experiment P5.8.7.4 realises an optical amplifier. An Erbium doped fiber is pumped below laser level, incoming light from a laser diode triggers stimulated emission and the light intensity is increased.



## P5.8.7

### OPTICAL FIBRES

#### P5.8.7.2

Plastic optical fibre (POF)

#### P5.8.7.3

Glass fibre optics

Glass fibre optics (P5.8.7.3)

Cat. No.	Description	P5.8.7.2	P5.8.7.3
474 5229	Plastic Optical Fibre, 10 m	1	
474 5230	Plastic Optical Fibre, 20 m	1	
474 5231	Plastic Optical Fibre, 30 m	1	
474 5232	F-SMA Connector Mounting Set	1	
474 6425	Coupler F-SMA for POF	2	
474 124	Plastic Fibre Holder with XY-Adjustment	1	
474 304	Dual Channel Receiver	1	
474 125	Dichroic Beam Splitter Unit	1	
474 5464	Oscilloscope, Dual Channel, Digital	1	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	2	1
474 305	Dual Channel LED Transmitter	1	
474 5424	Dual LED - FSMA	1	
474 5442	Profile rail, 500 mm	1	1
474 209	Mounting Plate C25 with Carrier 20 mm	2	
474 251	Transport and Storage Box #01	1	1
474 7118	Manual Plastic Fibre Optics	1	
474 1036	Collimating Optics on Carrier		1
474 151	Coupling Optics, XY- Adjustment Holder		1
474 152	Bare Fibre Holder with Translation Stage		1
474 154	Bare Fibre Holder on Rotation Stage		1
474 5227	Optical Glass Fibre, 1000 m multimode		1
474 6420	Optical Fibre Cleaver and Breaker		1
474 6421	Adjustable Plastic Cover Stripper		1
474 4025	IR converter screen 0.8 - 1.6 $\mu\text{m}$		1
474 306	Photodetector signal conditioning box		1
474 216	SiPIN Photodetector, Mounting Plate C25		1

Cat. No.	Description	P5.8.7.2	P5.8.7.3
501 06	HF-Cable, BNC-BNC, 1.5 m		1
474 302	Controller for Diode Laser		1
474 1022	Diode Laser Head with Adjustment Holder		1
474 7119	Manual Glass Fibre Optics		1
474 5226	Optical Glass Fibre, 1000 m monomode		1*
474 5295	Multimode optical fibre 5000 m, 50/125 $\mu\text{m}$ , on drum		1*

\* additionally recommended

Experiment P5.8.7.2 shows the properties of a plastic optical fiber and shows the basics of wavelength multiplexing and demultiplexing. A red and a blue LED are simultaneously coupled into the fiber and separated at the end.

Experiment P5.8.7.3 introduces to glass fibre optics. Within this experiment the diode laser itself will be characterised with respect to its output power as function of its temperature and injection current. The spatial intensity distribution is measured by means of the provided rotation stage. Cutting and preparing the fibre is part of the practical training. The light of the diode laser is coupled into the fibre by means of adjustable microscope objectives. The coupling efficiency is monitored with the photodetector detecting the light coming out at the end of the fibre. The intensity distribution of the light emerging at the end of the fibre is measured and the numerical aperture determined. By modulating the diode laser by means of the provided microprocessor controlled device the time of flight inside the fibre will be measured. From the results either the length of the fibre or the speed of light is calculated.

P5.8.7

OPTICAL FIBRES

P5.8.7.5  
Optical time domain  
reflectometry (OTDR)

P5.8.7.6  
Signal transmission via glass fibre

P5.8.7.7  
Fibre Optics Workshop



Optical time domain reflectometry (OTDR) (P5.8.7.5)

Cat. No.	Description	P5.8.7.5	P5.8.7.6	P5.8.7.7
474 1036	Collimating Optics on Carrier	1		
474 104	Focussing Optics, f = 60 mm	1		
474 5274	Quarter wave plate, C25	1		
474 230	Beam splitter module	1		
474 151	Coupling Optics, XY- Adjustment Holder	1		
474 150	Bare Fibre Holder with Carrier	1		
474 152	Bare Fibre Holder with Translation Stage	1		
474 5227	Optical Glass Fibre, 1000 m multimode	1		1
474 6420	Optical Fibre Cleaver and Breaker	1		
474 6421	Adjustable Plastic Cover Stripper	1		1
474 4025	IR converter screen 0.8 - 1.6 µm	1		
474 331	Photodetector, Ultrafast with Amplifier	1		
474 5464	Oscilloscope, Dual Channel, Digital	1	1*	
501 06	HF-Cable, BNC-BNC, 1.5 m	1		
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	1		
474 309	Controller for Pulsed Laser Diode	1		
474 5428	Pulsed Diode Laser Module 908 nm	1		
474 5442	Profile rail, 500 mm	1		
474 5444	Profile Rail, 1000 mm	1		
474 209	Mounting Plate C25 with Carrier 20 mm	2		
474 2114	Adjustment holder, 4 axes, rotary insert	1		
474 251	Transport and Storage Box #01	1	1*	
474 7121	Manual OTDR	1		
474 6426	Optical Fibre Transmitter, ST Connectors		1	
474 6427	Optical Fibre Receiver, ST Connectors		1	
474 5228	Optical Multimode Glass Fibre, 5000 m		1	
773 629	Fibre Patch Cable MM, 1.0 m		2	
501 091	BNC T adapter		1	
474 341	BNC load resistor 50 Ohm		1	

Cat. No.	Description	P5.8.7.5	P5.8.7.6	P5.8.7.7
474 5465	CCD Camera, coloured		1	
474 5466	CD player incl. music CD		1	
474 5467	Flat panel TV 19 inch		1	
474 7122	Manual Optical Data Transmission		1	
474 184	Fibre coupling module		1*	
474 5479	Hot Melt Assembly Kit			1
474 5476	HotMelt polishing unit			1
474 5477	Hotmelt ST connector, set of 60			1
474 7123	Manual Fibre Optics Workshop			1
474 6423	Splicing tubes 60 mm, set of 250			1*
474 5473	Fusion Splicer, SM & MM			1*
474 5474	High performance fibre cleaver & breaker			1*

\* additionally recommended

Experiment P5.8.7.5 shows the properties of optical time domain reflectometry (OTDR) in a fiber optic setup. Optical and mechanical imperfections within the fibre and fibre links or mechanical stress all lead to power losses. OTDR is the tool to analyse and locate such imperfections in optical fibres. The basic idea is to feed a light signal into the fibre and monitor the occurrence of light echoes.

Within the experiment P5.8.7.6, a data transmission line will be set-up with an optical fibre with a length of 5 km and the transmission of video as well as audio signals are demonstrated and studied. The set-up comes with a colour CCD video camera and a CD-player as an audio source and a TV screen as a monitor. By means of an optional fibre coupling module the optical signals can be monitored and analyzed and the sensitivity against misalignment studied.

The main goal of the experiment P5.8.7.7 is the connecting of optical glass fibres with ST connectors. Although a variety of other fibre connectors exist the process of connecting however remains the same. Another major technology is the welded connection of bare fibres by means of the fusion splicing technology.

## P5.8.8

### TECHNICAL APPLICATIONS

#### P5.8.8.1

Michelson laser interferometer



Michelson laser interferometer - Basic setup (P5.8.8.1\_b)

Cat. No.	Description	P5.8.8.1 (b)	P5.8.8.1 (c)	P5.8.8.1 (d)
474 5216	Plano-Convex lens f = 40 mm, C25 mount	1	1	1
474 5219	Biconcave Lens f = -5 mm, C25 mount	1	1	1
474 5220	Biconcave Lens f = -10 mm, C25 mount	1	1	1
474 5246	Laser Mirror, flat, M16 Mount	2	2	2
474 113	Laser Mirror Adjustment Holder, left	2	2	2
474 5247	Beam splitter plate 1:1 @ 632 nm, mounted	1	1	1
474 115	Carrier cross-piece, adjustable stage	1	1	1
474 5453	Crossed Hair Target in C25 mount	1	1	1
474 147	EXP10-Modul H, Schirm auf Reiter	1	1	1
474 306	Photodetector signal conditioning box	1	1	1
474 108	SiPIN photodetector	1	1	1
474 5464	Oscilloscope, Dual Channel, Digital	1	1	1
501 06	HF-Cable, BNC-BNC, 1.5 m	1	1	1
474 3034	HeNe Laser High Voltage Supply	1	1	1
474 5421	HeNe Pilot Laser Ø 30 mm	1	1	1
474 5442	Profile rail, 500 mm	4	4	4
474 5454	Laser adjustment holder, soft ring 30 mm, carrier 20 mm	2	2	2
474 209	Mounting Plate C25 with Carrier 20 mm	4	4	4
474 122	Optics cleaning set	1	1	1
671 9700	Ethanol, absolute, 250 ml	1	1	1
474 251	Transport and Storage Box #01	2	2	2
474 7112	LIT: Michelson Interferometer	1	1	1
474 267	Set of spare parts	1*	1*	1*
474 5248	Polarising beam splitter cube		1	1
474 403	Triple Reflector in 1 inch mount		2	2
474 144	Beam Displacer 5 mm		1	1
474 146	Fringe detection unit		1	1

Cat. No.	Description	P5.8.8.1 (b)	P5.8.8.1 (c)	P5.8.8.1 (d)
474 148	Dial Gauge 5 mm / 1 µm, Carrier 20		1	
474 149	Triple Reflector with Translation Stage		1	
474 308	Photodetector Preamplifier		1	1
474 3111	Fringe Counter		1	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m		2	2
474 330	Stepper Motor Controller, 1 Axis, USB			1
474 145	Motorised Translation Stage 50 mm			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1

\* additionally recommended

The setup of the famous Michelson interferometer in experiment P5.8.8.1 is accomplished by the modern modules and components. As light source a two mode Helium Neon laser is used. The contrast function is measured for different path lengths. From these results the coherence length of the probe laser is determined. The interference fringes are either displayed by the translucent white screen or their intensity is measured by means of an oscilloscope. By adding a variety of technical components the interferometer can be upgraded to a technical interferometer to train the calibration of CNC machines.

P5.8.8

TECHNICAL APPLICATIONS

P5.8.8.2  
Laser range finder

P5.8.8.3  
Laser vibrometer



Laser vibrometer (P5.8.8.3)


Cat. No.	Description	P5.8.8.2	P5.8.8.3
474 1031	Module B - Collimating optics on carrier MG-65 with Melles Griot	1	
474 104	Focussing Optics, f = 60 mm	1	1
474 4025	IR converter screen 0.8 - 1.6 $\mu$ m	1	
474 331	Photodetector, Ultrafast with Amplifier	1	
474 5464	Oscilloscope, Dual Channel, Digital	1	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	1	
474 309	Controller for Pulsed Laser Diode	1	
474 5428	Pulsed Diode Laser Module 908 nm	1	
474 5442	Profile rail, 500 mm	2	1
474 209	Mounting Plate C25 with Carrier 20 mm	1	1
474 211	Adjustment holder, 4 axes, carrier 20 mm	1	1
474 251	Transport and Storage Box #01	1	1
474 7113	LIT: Laser Range Finder	1	
474 5220	Biconcave Lens f = -10 mm, C25 mount		1
474 5320	Quarter wave plate, C25		1
474 206	Speaker mounted on carrier 20 mm		1
474 207	Beam Recombiner Assembly		1
474 208	Beam Splitting Assembly		1
474 2071	Rider with tilted mirror		1
474 313	Heterodyne Mixer and AOM driver		1
474 411	Acoustic optic modulator AOM		1
522 621	Function generator S 12		1
474 3312	Fast photodetector for laser use		2
501 06	HF-Cable, BNC-BNC, 1.5 m		3
575 24	Screened cable, BNC/4 mm		1

Cat. No.	Description	P5.8.8.2	P5.8.8.3
474 302	Controller for Diode Laser		1
474 5431	Laser 532 nm, singlemode		1
474 7114	LIT: Laser Vibrometer		1

In experiment P5.8.8.2 a high energy pulsed laser diode is collimated and aimed at the target. The short laser pulse (30 ns) travels with the speed of light. The scattered light is detected by a Si-PIN photodiode after passing the receiver lens. Based on the time of flight and the known speed of light the distance to the target is calculated.

The laser vibrometer in experiment P5.8.8.3 demonstrates the working principle of a contactless measurement of vibrations of a target. Using a heterodyne setup one optical detector is sufficient to realise a quadrature forward / backward counter. The back scattered and Doppler shifted light is coupled back to the interferometer and superimposed with the reference beam. The frequency of the reference beam is shifted by a Bragg cell. The subsequent signal amplifiers can be AC coupled allowing a much higher gain in a simple way than a Michelson interferometer.





## P5.8.8

### TECHNICAL APPLICATIONS

P5.8.8.4

Laser Doppler anemometer

Laser Doppler anemometer (P5.8.8.4)


Cat. No.	Description	P5.8.8.4
474 104	Focussing Optics, $f = 60$ mm	2
474 107	Filter Plate Holder	1
474 187	LDA Beam Splitting Assembly	1
474 1876	LDA Beam Deflection and Focussing	1
474 188	Ultrasonic particle seeder	1
474 315	Ultrasonic Particle Nebuliser	1
474 331	Photodetector, Ultrafast with Amplifier	1
474 5464	Oscilloscope, Dual Channel, Digital	1
501 06	HF-Cable, BNC-BNC, 1.5 m	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	1
474 351	Signal Amplifier	1
474 301	Adaptive Power Supply	1
474 128	Diode laser module	1
474 5442	Profile rail, 500 mm	1
474 251	Transport and Storage Box #01	2
474 7115	LIT: Laser Doppler Anemometer	1

Laser Doppler anemometry is a non-contact optical measurement method to obtain the velocity of a flow (liquid, gas). In the experiment P5.8.8.4 a laser Doppler anemometer is assembled. A laser beam is split into two parts. Focused back to one spot, the laser beams create an interference pattern. Particles in the fluid flow move through the bright and dark zones of the pattern and the scattered light is modulated according to the speed of the particle. Alternatively, the same setup can be explained in terms of Doppler shifted light.

In this experiment, water droplets inside an air stream are used as the scattering particles. The water droplets are too small to stick to surfaces nearby and to wet them. They will just evaporate or bounce off due to surface tension.



Laser Doppler Anemometer Setup




Oscilloscope signal of a water droplet passing through the interference zone

P5.8.8

TECHNICAL APPLICATIONS

P5.8.8.5


HeNe laser gyroscope



HeNe laser gyroscope (P5.8.8.5)

Cat. No.	Description	P5.8.8.5
474 159	Gyroscope base plate	1
474 160	Rotation Unit	1
474 330	Stepper Motor Controller, 1 Axis, USB	1
474 407	Fringe Detection Unit	1
474 308	Photodetector Preamplifier	1
474 311	Fringe Up and Down Counter	1
474 346	Plug-in power supply 12 V - 2.5 A	1
474 5464	Oscilloscope, Dual Channel, Digital	1
501 061	HF-Cable, BNC-Mini BNC, 1.5 m	4
728 950	USB port isolator	1
474 303	HeNe Laser High Voltage supply, adjustable	1
474 5423	Alignment laser 532 nm	1
474 122	Optics cleaning set	1
671 9700	Ethanol, absolute, 250 ml	1
474 7116	LIT: HeNe Laser Gyroscope	1
	additionally required: PC with Windows XP or higher	1

Experiment P5.8.8.5 shows the basics of a laser gyroscope. A three-mirror mono mode HeNe ring laser is set up; using a green pilot laser adjusting the mirrors is made easy. Both CW and CCW light beams are coupled out of the resonator and fed into an interferometer, creating interference patterns. Rotating this laser, the CW and CCW cavity length is no longer equal due to special relativity, and mixing both light rays will generate a beat frequency in the kHz range. This way the absolute rotation of the laser can be measured. This is not a Sagnac style interferometer, having the laser medium inside the moving system creates a much more sensitive setup with beat frequency instead of the phase shift of a rotating Sagnac interferometer.



Optical Setup of the ring laser gyro with interferometer

P5.8.8

TECHNICAL APPLICATIONS

P5.8.8.6

Laser beam analysis



Laser beam analysis (P5.8.8.6)

Cat. No.	Description	P5.8.8.6
474 1036	Collimating Optics on Carrier	1
474 5266	Collimating cylindrical lens $f = 20$ mm	1
474 5267	Collimating cylindrical lens $f = 80$ mm	1
474 5263	Beam expander 6x	1
474 5470	Beammaster BM-7S	1
474 5418	Diode Laser Module, 532 nm	1
474 5420	Dimo diode laser module, 630 nm (red)	1
474 5442	Profile rail, 500 mm	1
474 209	Mounting Plate C25 with Carrier 20 mm	3
474 211	Adjustment holder, 4 axes, carrier 20 mm	1
474 251	Transport and Storage Box #01	1
474 7117	Manual Laser beam analysis	1
	additionally required: PC with Windows XP or higher	1

In experiment P5.8.8.6 two different visible lasers plus collimating and expanding optics are used to demonstrate a variety of beam profile measurements with various beam shaping arrangements to show the possibilities of laser beam forming. The most important property of a laser beam is the intensity distribution in the beam. The beam profiler is no CDD type, but a knife edge type: A precise knife edge is moved through the cross-section of the laser beam. As the blade moves across the beam, it is cut from reaching the photodetector.

Using several knives with different angles on one rotating drum, the beam profile can be calculated from the time responses when the knives pass through the beam.



Optical setup of a semiconductor laser diode

P5.3.4.2  
Determining the wavelength of the  
light of an He-Ne laser using a Mi-  
chelson interferometer



# P6 ATOMIC AND NUCLEAR PHYSICS



P6.1	INTRODUCTORY EXPERIMENTS	231
P6.2	ATOMIC SHELL	239
P6.3	X-RAY PHYSICS	250
P6.4	RADIOACTIVITY	261
P6.5	NUCLEAR PHYSICS	265
P6.6	QUANTUM PHYSICS	272

# P6 ATOMIC AND NUCLEAR PHYSICS



## P6.1 INTRODUCTORY EXPERIMENTS

P6.1.1 Oil-spot experiment	231
P6.1.2 Millikan experiment	232
P6.1.3 Specific electron charge	233
P6.1.4 Planck's constant	234-236
P6.1.5 Dual nature of wave and particle	237
P6.1.6 Paul trap	238

## P6.2 ATOMIC SHELL

P6.2.1 Balmer series of hydrogen	239-240
P6.2.2 Emission and absorption spectra	241-243
P6.2.3 Inelastic collisions of electrons	244
P6.2.4 Franck-Hertz experiment	245-246
P6.2.6 Electron spin resonance	247
P6.2.7 Normal Zeeman effect	248
P6.2.8 Optical pumping (anomalous Zeeman effect)	249

## P6.3 X-RAY PHYSICS

P6.3.1 Detection of X-rays	250-252
P6.3.2 Attenuation of X-rays	253
P6.3.3 Physics of the atomic shell	254
P6.3.5 X-ray energy spectroscopy	255
P6.3.6 Structure of X-ray spectrums	256-257
P6.3.7 Compton effect at X-rays	258
P6.3.8 X-ray tomography	259-260

## P6.4 RADIOACTIVITY

P6.4.1 Detecting radioactivity	261
P6.4.2 Poisson distribution	262
P6.4.3 Radioactive decay and half-life	263
P6.4.4 Attenuation of $\alpha$ -, $\beta$ - and $\gamma$ radiation	264

## P6.5 NUCLEAR PHYSICS

P6.5.1 Demonstrating paths of particles	265
P6.5.2 Rutherford scattering	266
P6.5.3 Nuclear magnetic resonance	267
P6.5.4 $\alpha$ spectroscopy	268
P6.5.5 $\gamma$ spectroscopy	269
P6.5.6 Compton effect	270
P6.5.7 Properties of radiation particles	271

## P6.6 QUANTUM PHYSICS

P6.6.1 Quantum optics	272
P6.6.2 Particles	273



**P6.1.1**  
**OIL-SPOT EXPERIMENT**

P6.1.1.1  
Estimating the size of oil molecules

Estimating the size of oil molecules (P6.1.1.1)


Cat. No.	Description	P6.1.1.1
664 179	Crystallisation dish, 230 mm diam., 3500 ml	1
665 844	Burette, amber glass, 10 ml	1
664 110	Beaker, DURAN, 50 ml, tall	1
665 751	Measuring cylinder, 10 ml, with plastic base	1
665 754	Measuring cylinder, 100 ml, with plastic base	1
300 02	Stand base, V-shaped, small	1
300 43	Stand rod, 75 cm, 12 mm diam.	1
301 01	Leybold multiclamp	1
666 555	Universal clamp, 0...80 mm	1
675 3410	Water, pure, 5 l	1
672 1240	Glycerin trioleate, 100 ml	1
674 2220	Petroleum ether, 40...70 °C, 1 l	1
670 6920	Lycopodium spores, 25 g	1

One important issue in atomic physics is the size of the atom. An investigation of the size of molecules makes it easier to come to a usable order of magnitude by experimental means. This is estimated from the size of an oil spot on the surface of water using simple means.

In the experiment P6.1.1.1, a drop of glycerin trioleate is added to a grease-free water surface dusted with *Lycopodium* spores. Assuming that the resulting oil spot has a thickness of one molecule, we can calculate the size  $d$  of the molecule according to

$$d = \frac{V}{A}$$

from the volume  $V$  of the oil droplet and the area  $A$  of the oil spot. The volume of the oil spot is determined from the number of drops needed to fill a volume of 1 cm<sup>3</sup>. The area of the oil spot is determined using graph paper.



Determining the area  $A$  of the oil spot

### P6.1.2

#### MILLIKAN EXPERIMENT

##### P6.1.2.1

Determining the electric unit charge after Millikan and verifying the charge quantification - Measuring the suspension voltage and the falling speed

##### P6.1.2.2

Determining the electric unit charge after Millikan and verifying the charge quantification - Measuring the rising and falling speed

##### P6.1.2.3

Determining the electric unit charge after Millikan and verifying the charge quantification - Measuring the suspension voltage and the falling speed with CASSY

##### P6.1.2.4

Determining the electric unit charge after Millikan and verifying the charge quantification - Measuring the rising and falling speed with CASSY



Determining the electric unit charge after Millikan and verifying the charge quantification - Measuring the suspension voltage and the falling speed (P6.1.2.1\_b)

Cat. No.	Description	P6.1.2.1 (b)	P6.1.2.2	P6.1.2.3	P6.1.2.4
559 412	Millikan apparatus	1	1	1	1
559 421	Millikan supply unit	1	1	1	1
575 471	Counter S	1			
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	3	4	3	3
313 033	Electronic time clock P		2		
<b>524 013</b>	<b>Sensor-CASSY 2</b>			<b>1</b>	<b>1</b>
524 220	CASSY Lab 2			1	1
524 034	Timer box			1	1
501 461	Connecting leads, 19 A, 100 cm, black, pair			1	1
500 421	Connecting lead 19 A, 50 cm, red				1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1	1

With his famous oil-drop method, *R. A. Millikan* succeeded in demonstrating the quantum nature of minute amounts of electricity in 1910. He caused charged oil droplets to be suspended in the vertical electric field of a plate capacitor and, on the basis of the radius  $r$  and the electric field  $E$ , determined the charge  $q$  of a suspended droplet:

$$q = \frac{4\pi}{3} \cdot r^3 \cdot \frac{\rho \cdot g}{E}$$

$\rho$ : density of oil

$g$ : gravitanational acceleration

He discovered that  $q$  only occurs as a whole multiple of an electron charge  $e$ . His experiments are produced here in two variations.

In the variation P6.1.2.1 and P6.1.2.3, the electric field

$$E = \frac{U}{d}$$

$d$ : plate spacing

is calculated from the voltage  $U$  at the plate capacitor at which the observed oil droplet just begins to hover. The constant falling velocity  $v_1$  of the droplet when the electric field is switched off is subsequently measured to determine the radius. From the equilibrium between the force of gravity and Stokes friction, we derive the equation


$$\frac{4\pi}{3} \cdot r^3 \cdot \rho \cdot g = 6\pi \cdot r \cdot \eta \cdot v_1$$

$\eta$ : viscosity

In the variant P6.1.2.2 and P6.1.2.4, the oil droplets are observed which are not precisely suspended, but which rise with a low velocity  $v_2$ . The following applies for these droplets:

$$q \cdot \frac{U}{d} = \frac{4\pi}{3} \cdot r^3 \cdot \rho \cdot g + 6\pi \cdot r \cdot \eta \cdot v_2$$

Additionally, the falling speed  $v_1$  is measured, as in the variant P6.1.2.1 and P6.1.2.3. The measuring accuracy for the charge  $q$  can be increased by causing the oil droplet under study to rise and fall over a given distance several times in succession and measuring the total rise and fall times.



The histogram reveals the quantum nature of the charge





**P6.1.3**  
**SPECIFIC ELECTRON  
CHARGE**

P6.1.3.1  
Determining the specific charge  
of the electron

Determining the specific charge of the electron (P6.1.3.1)

Cat. No.	Description	P6.1.3.1
555 571	Fine beam tube	1
555 581	Helmholtz coils with holder and measuring device for fine beam tube	1
531 120	Multimeter LDanalog 20	2
521 65	Tube power supply, 0...500 V	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1
311 77	Steel tape measure, 2 m	1
500 614	Safety connecting lead, 25 cm, black	3
500 624	Safety connecting lead, 50 cm, black	3
500 644	Safety connecting lead, 100 cm, black	7
531 835	Universal measuring instrument, Physics	1*
524 0382	Axial B sensor S, ±1000 mT	1*
501 11	Extension cable, 15 pin	1*

\* additionally recommended

The mass  $m_e$  of the electron is extremely difficult to determine in an experiment. It is much easier to determine the specific charge of the electron

$$\epsilon = \frac{e}{m_e}$$

from which we can calculate the mass  $m_e$  for a given electron charge  $e$ . In the experiment P6.1.3.1, a tightly bundled electron beam is diverted into a closed circular path using a homogeneous magnetic field in order to determine the specific electron charge. The magnetic field  $B$  which diverts the electrons into the path with the given radius  $r$  is determined as a function of the acceleration voltage  $U$ . The Lorentz force caused by the magnetic field acts as a centripetal force. It depends on the velocity of the electrons, which in turn is determined by the acceleration voltage. The specific electron charge can thus be determined from the measurement quantities  $U$ ,  $B$  and  $r$  according to the formula

$$\frac{e}{m_e} = 2 \cdot \frac{U}{B^2 \cdot r^2}$$



Circular electron path in fine beam tube

### P6.1.4


#### PLANCK'S CONSTANT

##### P6.1.4.1

Determining Planck's constant -  
Measuring in a compact assembly

##### P6.1.4.5

Determining Planck's constant -  
Recording the current-voltage  
characteristics, measuring in  
a compact assembly



Determining Planck's constant - Measuring in a compact assembly (P6.1.4.1)

Cat. No.	Description	P6.1.4.1	P6.1.4.5 (a)
558 77	Photocell for determining Planck's constant	1	1
558 79	Compact arrangement for determining Planck's constant	1	1
451 15	High pressure mercury lamp	1	1
451 195	Power supply unit for high-pressure mercury lamp	1	1
532 14	Electrometer amplifier	1	
562 791	Plug-in power supply, 12 V AC	1	
578 22	Capacitor, 100 pF, STE 2/19	1	
579 10	Push button (NO), STE 2/19	1	
590 011	Clamping plug	2	
531 120	Multimeter LDanalog 20	1	2
575 24	Screened cable, BNC/4 mm	1	
502 04	Distribution box with earthing socket	1	
500 414	Connecting lead, 19 A, 25 cm, black	1	
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	2
501 461	Connecting leads, 19 A, 100 cm, black, pair	1	1
500 440	Connecting lead, 19 A, 100 cm, yellow/green	1	
532 00	I-measuring amplifier D		1
576 74	Plug-in board, DIN A4, STE		1
576 86	Monocell holder STE 2/50		3
685 48ET5	Batteries 1.5 V (D, mono), set of 5		1
577 93	Potentiometer, 1 kΩ, 10-turn, STE 4/50		1
579 13	Toggle switch, STE 2/19		1
501 48	Bridging plugs, STE 2/19, set of 10		1
501 02	BNC cable, 1 m		1
500 444	Connecting lead 19 A, 100 cm, black		1

When light with the frequency  $\nu$  falls on the cathode of a photocell, electrons are released. Some of the electrons reach the anode where they generate a current in the external circuit, which is compensated to zero by applying a voltage with opposite sign  $U = -U_0$ . The applicable relationship

$$e \cdot U_0 = h \cdot \nu - W \quad W: \text{electronic work function}$$

was first used by *R. A. Millikan* to determine Planck's constant  $h$ .

In the experiment P6.1.4.1, a compact arrangement is used to determine  $h$ , in which the light from a high-pressure mercury vapour lamp is spectrally dispersed in a direct-vision prism. The light of precisely one spectral line at a time falls on the cathode. A capacitor is connected between the cathode and the anode of the photocell which is charged by the anode current, thus generating the opposing voltage  $U$ . As soon as the opposing voltage reaches the value  $-U_0$ , the anode current is zero and the charging of the capacitor is finished.  $U_0$  is measured without applying a current by means of an electrometer amplifier.

In the experiment P6.1.4.5 light from a mercury gas discharge lamp is deflected by a direct view prism, one wavelength selected and focused onto the photocathode. The countervoltage of the anode is varied and the resulting current is measured with high sensitivity. The variation of the characteristic curves under irradiation with different wavelengths leads to the determination of Planck's constant  $h$ .



**P6.1.4**  
**PLANCK'S CONSTANT**

P6.1.4.3  
Determining Planck's constant -  
Selection of wavelengths using  
interference filters on the  
optical bench

Determining Planck's constant - Selection of wavelengths using interference filters on the optical bench (P6.1.4.3\_a)

Cat. No.	Description	P6.1.4.3 (a)
558 77	Photocell for determining Planck's constant	1
558 791	Holder for photocell	1
460 335	Optical bench with standardised profile, 0.5 m	1
460 374	Optics rider, 90/50	2
460 375	Optics rider, 120/50	3
558 792	Filter wheel with diaphragm	1
468 401	Interference filter, 578 nm	1
468 402	Interference filter, 546 nm	1
468 403	Interference filter, 436 nm	1
468 404	Interference filter, 405 nm	1
460 03	Lens in frame, $f=100$ mm	1
460 26	Iris diaphragm	1
451 15	High pressure mercury lamp	1
451 195	Power supply unit for high-pressure mercury lamp	1
532 14	Electrometer amplifier	1
562 791	Plug-in power supply, 12 V AC	1
578 22	Capacitor, 100 pF, STE 2/19	1
579 10	Push button (NO), STE 2/19	1
590 011	Clamping plug	2
531 120	Multimeter LDanalog 20	1
501 10	BNC adapter, straight	1
501 09	BNC/4-mm plug adapter, single-pole	1
340 89ET5	Coupling plugs, 4 mm, set of 5	1
502 04	Distribution box with earthing socket	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1
500 440	Connecting lead, 19 A, 100 cm, yellow/green	2

In determining Planck's constant using the photoelectric effect, it must be ensured that only the light of a single spectral line of the high-pressure mercury vapour lamp falls on the cathode of the photocell at any one time. As an alternative to a prism, it is also possible to use narrow-band interference filters to select the wavelength. This simplifies the subsequent optical arrangement, and it is no longer necessary to darken the experiment room. Also, the intensity of the light incident on the cathode can be easily varied using an iris diaphragm.

In the experiment P6.1.4.3, the capacitor method described previously (see P6.1.4.1) is used to generate the opposing voltage  $U$  between the cathode and the anode of the photocell. The voltage at the capacitor is measured without current using the electrometer amplifier.

*Note:* The opposing voltage  $U$  can alternatively be tapped from a DC voltage source. The I-measuring amplifier  $D$  is recommended for sensitive measurements of the anode current (see P 6.1.4.4).

### P6.1.4

#### PLANCK'S CONSTANT

##### P6.1.4.4

Determining Planck's constant -  
Recording the current-voltage  
characteristics, selection of  
wavelengths using interference  
filters on the optical bench




*Determining Planck's constant - Recording the current-voltage characteristics, selection of wavelengths using interference filters on the optical bench (P6.1.4.4\_b)*

Cat. No.	Description	P6.1.4.4 (a)	P6.1.4.4 (b)
558 77	Photocell for determining Planck's constant	1	1
558 791	Holder for photocell	1	1
460 335	Optical bench with standardised profile, 0.5 m	1	1
460 374	Optics rider, 90/50	2	2
460 375	Optics rider, 120/50	3	3
558 792	Filter wheel with diaphragm	1	1
468 401	Interference filter, 578 nm	1	1
468 402	Interference filter, 546 nm	1	1
468 403	Interference filter, 436 nm	1	1
468 404	Interference filter, 405 nm	1	1
468 406	Interference filter, 365 nm	1	1
460 03	Lens in frame, f=100 mm	1	1
460 26	Iris diaphragm	1	1
451 15	High pressure mercury lamp	1	1
451 195	Power supply unit for high-pressure mercury lamp	1	1
532 00	I-measuring amplifier D	1	1
531 120	Multimeter LDanalog 20	2	
576 74	Plug-in board, DIN A4, STE	1	1
576 86	Monocell holder STE 2/50	3	3
685 48ET5	Batteries 1.5 V (D, mono), set of 5	1	1
577 93	Potentiometer, 1 k $\Omega$ , 10-turn, STE 4/50	1	1
579 13	Toggle switch, STE 2/19	1	1
501 48	Bridging plugs, STE 2/19, set of 10	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	2
500 444	Connecting lead 19 A, 100 cm, black	1	1

Cat. No.	Description	P6.1.4.4 (a)	P6.1.4.4 (b)
524 013	Sensor-CASSY 2		1
524 220	CASSY Lab 2		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

In the experiment P6.1.4.4 one of the emission lines from a mercury gas discharge lamp is selected by interference filters and focused onto the photocathode. The countervoltage of the anode is varied and the resulting current is measured with high sensitivity. The variation of the characteristic curves under irradiation with different wavelengths leads to the determination of Planck's constant  $h$ .



*Current-voltage characteristics of the photocell*

## P6.1.5

### DUAL NATURE OF WAVE AND PARTICLE

#### P6.1.5.1

Diffraction of electrons at a polycrystalline lattice (Debye-Scherrer diffraction)


#### P6.1.5.2

Optical analogy to electron diffraction at a polycrystalline lattice



Diffraction of electrons at a polycrystalline lattice (Debye-Scherrer diffraction) (P6.1.5.1)

Cat. No.	Description	P6.1.5.1	P6.1.5.2
555 626	Electron diffraction tube	1	
555 600	Tube stand	1	
521 70	High-voltage power supply, 10 kV	1	
311 54	Precision vernier callipers	1	
500 611	Safety connecting lead, 25 cm, red	1	
500 621	Safety connecting lead, 50 cm, red	1	
500 641	Safety connecting lead, 100 cm, red	1	
500 642	Safety connecting lead, 100 cm, blue	1	
500 644	Safety connecting lead, 100 cm, black	2	
555 629	Cross grating, 5000/cm, rotatable		1
450 63	Halogen bulb, 12 V/100 W, G6.35		1
450 64	Halogen lamp, 12 V, 50/100 W		1
450 66	Picture slider		1
521 25	Transformer, 2...12 V, 120 W		1
460 03	Lens in frame, f=100 mm		1
460 22	Holder with spring clips		1
441 53	Screen, translucent		1
311 77	Steel tape measure, 2 m		1
460 43	Small optical bench		1
301 01	Leybold multiclamp		5
300 01	Stand base, V-shaped, large		1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair		1



Optical analogon of Debye-Scherrer diffraction (P6.1.5.2)

In 1924, *L. de Broglie* first hypothesized that particles could have wave properties in addition to their familiar particle properties, and that their wavelength depends on the linear momentum  $p$

$$\lambda = \frac{h}{p} \quad h: \text{Planck's constant}$$

His conjecture was confirmed in 1927 by the experiments of *C. Davisson* and *L. Germer* on the diffraction of electrons at crystalline structures.

The experiment P6.1.5.1 demonstrates diffraction of electrons at polycrystalline graphite. As in the Debye-Scherrer method with x-rays, we observe diffraction rings in the direction of radiation which surround a central spot on a screen. These are caused by the diffraction of electrons at the lattice planes of microcrystals which fulfill the Bragg condition

$$2 \cdot d \cdot \sin \vartheta = n \cdot \lambda$$

$\vartheta$ : angular aperture of diffraction ring

$d$ : spacing of lattice planes

As the graphite structure contains two lattice-plane spacings, two diffraction rings in the first order are observed. The electron wavelength

$$\lambda = \frac{h}{\sqrt{2 \cdot m_e \cdot e \cdot U}}$$

$m_e$ : mass of electron,  $e$ : elementary charge

is determined by the acceleration voltage  $U$ , so that for the angular aperture of the diffraction rings we can say

$$\sin \vartheta \propto \frac{1}{\sqrt{U}}$$

The experiment P6.1.5.2 illustrates the Debye-Scherrer method used in the electron diffraction tube by means of visible light. Here, parallel, monochromatic light passes through a rotating cross grating. The diffraction pattern of the cross grating at rest, consisting of spots of light arranged around the central beam in a network-like pattern, is deformed by rotation into rings arranged concentrically around the central spot.

P6.1.6

PAUL TRAP

P6.1.6.1

Observing individual lycopod spores in a Paul trap



Observing individual lycopod spores in a Paul trap (P6.1.6.1\_a)

Cat. No.	Description	P6.1.6.1 (a)
558 80	Paul trap	1
471 830	He-Ne Laser, linearly polarised	1
460 01	Lens in frame, f=5 mm	1
460 335	Optical bench with standardised profile, 0.5 m	1
460 373	Optics rider, 60/50	3
522 27	Power supply, 450 V	1
521 35	Variable extra-low voltage transformer S	1
562 11	U-core with yoke	1
562 121	Clamping device with spring clip	1
562 18	Coil, 50 turns, extra-low voltage	1
562 16	Coil, 10 000 turns	1
531 120	Multimeter LDanalog 20	1
536 211	Measuring resistor, 10 MΩ	1
502 04	Distribution box with earthing socket	1
500 624	Safety connecting lead, 50 cm, black	2
500 641	Safety connecting lead, 100 cm, red	1
500 642	Safety connecting lead, 100 cm, blue	1
500 644	Safety connecting lead, 100 cm, black	1
500 98	Safety adapter sockets, black, set of 6	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2
500 440	Connecting lead, 19 A, 100 cm, yellow/green	1

Spectroscopic measurements of atomic energy levels are normally impaired by the motion of the atoms under study with respect to the radiation source. This motion shifts and broadens the spectral lines due to the Doppler effect, which becomes strongly apparent in high-resolution spectroscopy. The influence of the Doppler effect is reduced when individual atoms are enclosed in a small volume for spectroscopic measurements. For charged particles (ions), this can be achieved using the ion trap developed by *W. Paul* in the 1950's. This consists of two rotationally symmetrical cover electrodes and one ring electrode. The application of an AC voltage generates a time-dependent, parabolic potential with the form

$$U(r, z, t) = U_0 \cdot \cos \omega t \cdot \frac{r^2 - 2z^2}{2 \cdot r_0^2}$$

$z$ : coordinate on the axis of symmetry

$r$ : coordinate perpendicular to axis of symmetry

$r_0$ : inside radius of ring electrode

An ion with the charge  $q$  and the mass  $m$  remains trapped in this potential when the conditions

$$0.4 \cdot \alpha < \frac{q}{m} < 1.2 \alpha \quad \text{where} \quad \alpha = \frac{r_0^2 \cdot \omega^2}{U_0}$$

are fulfilled.

The experiment P6.1.6.1 demonstrates how a Paul trap works using a model which can be operated with no special requirements at standard air pressure and with 50 Hz AC. When a suitable voltage amplitude  $U_0$  is set, it is possible to trap lycopod spores for several hours and observe them under laser light. Tilting of the entire ion trap causes the trapped particles to move radially within the ring electrode. When a voltage is applied between the cover electrodes, it is possible to shift the potential along the  $z$ -axis.





## P6.2.1 BALMER SERIES OF HYDROGEN

P6.2.1.1  
Determining the wavelengths  $H_{\alpha}$ ,  $H_{\beta}$  and  $H_{\gamma}$  from the Balmer series of hydrogen

P6.2.1.2  
Observing the Balmer series of hydrogen using a prism spectrometer

P6.2.1.4  
Observing the Balmer series of hydrogen using a compact spectrometer

Determining the wavelengths  $H_{\alpha}$ ,  $H_{\beta}$  and  $H_{\gamma}$  from the Balmer series of hydrogen (P6.2.1.1)

Cat. No.	Description	P6.2.1.1	P6.2.1.2	P6.2.1.4
451 13	Balmer lamp	1	1	1
451 141	Power supply unit for Balmer lamp	1	1	1
471 23	Ruled grating, 6000/cm (Rowland)	1		
311 77	Steel tape measure, 2 m	1		
460 02	Lens in frame, $f=50$ mm	1		
460 03	Lens in frame, $f=100$ mm	1		
460 14	Adjustable slit	1		
460 22	Holder with spring clips	1		
441 53	Screen, translucent	1		
460 43	Small optical bench	1		
300 01	Stand base, V-shaped, large	1		
301 01	Leybold multiclamp	6		
467 112	School spectroscope		1	
467 251	Compact spectrometer, physics (spectral photometer)			1
460 251	Fibre holder			1
300 11	Saddle base			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1

In the visible range, the emission spectrum of atomic hydrogen has four lines,  $H_{\alpha}$ ,  $H_{\beta}$ ,  $H_{\gamma}$  and  $H_{\delta}$ ; this sequence continues into the ultraviolet range to form a complete series. In 1885, *Balmer* empirically worked out a formula for the frequencies of this series

$$\nu = R_{\infty} \cdot \left( \frac{1}{2^2} - \frac{1}{m^2} \right), \quad m: 3, 4, 5, \dots$$

$$R_{\infty}: 3.2899 \cdot 10^{15} \text{ s}^{-1}: \text{Rydberg constant}$$

which could later be explained using Bohr's model of the atom.

In the experiment P6.2.1.1, the emission spectrum is excited using a Balmer lamp filled with water vapor, in which an electric discharge splits the water molecules into an excited hydrogen atom and a hydroxyl group. The wavelengths of the lines  $H_{\alpha}$ ,  $H_{\beta}$  and  $H_{\gamma}$  are determined using a high-resolution grating. In the first diffraction order of the grating, we can find the following relationship between the wavelength  $\lambda$  and the angle of observation  $\vartheta$ :


$$\lambda = d \cdot \sin \vartheta$$

$$d: \text{grating constant}$$

The measured values are compared with the values calculated according to the Balmer formula.

In the experiment P6.2.1.2 the Balmer series is studied by means of a prism spectroscope (complete device).

In the experiment P6.2.1.4 the Balmer series is studied by means of a computer attached compact spectrometer. The optical adjustment is extremely simple, just placing the optical fiber close to the discharge lamp. As a result you get the complete infrared and visible range of the emission spectrum.



Emission spectrum of atomic hydrogen

### P6.2.1

#### BALMER SERIES OF HYDROGEN

##### P6.2.1.3

Observing the splitting of the Balmer series on deuterated hydrogen (isotope splitting)



Observing the splitting of the Balmer series on deuterated hydrogen (isotope splitting) (P6.2.1.3\_b)

Cat. No.	Description	P6.2.1.3 (b)
451 41	Balmer lamp, deuterated	1
451 141	Power supply unit for Balmer lamp	1
460 02	Lens in frame, f=50 mm	1
460 08	Lens in frame, f=150 mm	1
460 09	Lens in frame, f=300 mm	1
460 14	Adjustable slit	1
471 27	Holographic grating in frame	1
337 47USB	VideoCom	1
460 32	Optical bench with standardised profile, 1 m	1
460 335	Optical bench with standardised profile, 0.5 m	1
460 341	Swivel joint with circular scale	1
460 374	Optics rider, 90/50	6
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The Balmer series of the hydrogen atom is given by the electron transitions to the second energy level (principal quantum number  $n = 2$ ) from higher states ( $m: 3, 4, 5, \dots$ ). The wavelength of the emitted photons is given by

$$\frac{c}{\lambda} = R \left( \frac{1}{n^2} - \frac{1}{m^2} \right) \quad R = \text{Rydberg constant}$$

Here, one assumes that the mass of the nucleus is much bigger than the mass of the electron. For a more exact calculation, the Rydberg constant has to be corrected employing the reduced mass. Therefore, the Rydberg constants  $R_H$  for hydrogen and  $R_D$  for the isotope deuterium whose nucleus consists of a proton and a neutron differ. The spectral lines of the Balmer series of deuterium are shifted to somewhat smaller wavelengths compared to the spectral lines of hydrogen. This effect is called isotopic shift.

In the experiment P6.2.1.3 the Balmer series is studied by means of a high resolution spectrometer setup. A holographic grating with the grating constant  $g$  is used. The wavelength splitting is calculated from the angle  $\beta$  of the 1. order maximum and the angle splitting  $\Delta\beta$ :

$$\Delta\lambda = g \cos\beta \Delta\beta$$





## P6.2.2 EMISSION AND ABSORPTION SPECTRA

P6.2.2.1  
Displaying the spectral lines  
of inert gases and metal vapors

P6.2.2.2  
Qualitative investigation of the  
absorption spectrum of sodium

Displaying the spectral lines of inert gases and metal vapors (P6.2.2.1)

Cat. No.	Description	P6.2.2.1	P6.2.2.2
451 011	Spectral lamp, Ne	1	
451 041	Spectral lamp, Cd	1	
451 062	Spectral lamp, Hg 100	1	
451 111	Spectral lamp, Na	1	1
451 16	Housing for spectral lamps	1	1
451 30	Universal choke, 230 V, 50 Hz	1	1
471 23	Ruled grating, 6000/cm (Rowland)	1	
311 77	Steel tape measure, 2 m	1	
460 02	Lens in frame, f=50 mm	1	
460 03	Lens in frame, f=100 mm	1	
460 14	Adjustable slit	1	
460 22	Holder with spring clips	1	
441 53	Screen, translucent	1	1
460 43	Small optical bench	1	
300 01	Stand base, V-shaped, large	1	
301 01	Leybold multiclamp	6	2
450 60	Lamp housing with cable		1
450 511	Bulbs, 6 V/30 W, E14, set of 2		1
521 210	Transformer, 6/12 V		1
300 02	Stand base, V-shaped, small		2
300 11	Saddle base		1
300 42	Stand rod, 47 cm, 12 mm diam.		2
666 711	Butane gas burner		1
666 712ET3	Butane cartridge, 190 g, set of 3		1
666 962	Double-ended spatula, stainless steel, 150 mm		1
673 0840	Magnesia rods, 25 pieces		1
673 5700	Sodium chloride, 250 g		1

When an electron in the shell of an atom or atomic ion drops from an excited state with the energy  $E_2$  to a state of lower energy  $E_1$ , it can emit a photon with the frequency


$$\nu = \frac{E_2 - E_1}{h}$$

$h$ : Planck's constant

In the opposite case, a photon with the same frequency is absorbed. As the energies  $E_1$  and  $E_2$  can only assume discrete values, the photons are only emitted and absorbed at discrete frequencies. The totality of the frequencies which occur is referred to as the spectrum of the atom. The positions of the spectral lines are characteristic of the corresponding element.

The experiment P6.2.2.1 disperses the emission spectra of metal vapors and inert gases (mercury, sodium, cadmium and neon) using a high-resolution grating and projects these on the screen for comparison purposes.

In the experiment P6.2.2.2, the flame of a Bunsen burner is alternately illuminated with white light and sodium light and observed on a screen. When sodium is burned in the flame, a dark shadow appears on the screen when illuminating with sodium light. From this it is possible to conclude that the light emitted by a sodium lamp is absorbed by the sodium vapor, and that the same atomic components are involved in both absorption and emission.



Emission spectra

## P6.2.2

### EMISSION AND ABSORPTION SPECTRA

#### P6.2.2.3

Investigating the spectrum of a high pressure mercury lamp



Investigating the spectrum of a high pressure mercury lamp (P6.2.2.3\_c)

Cat. No.	Description	P6.2.2.3 (c)
451 15	High pressure mercury lamp	1
451 195	Power supply unit for high-pressure mercury lamp	1
460 02	Lens in frame, f=50 mm	1
460 08	Lens in frame, f=150 mm	1
460 09	Lens in frame, f=300 mm	1
460 14	Adjustable slit	1
471 27	Holographic grating in frame	1
441 531	Screen	1
337 47USB	VideoCom	1
460 335	Optical bench with standardised profile, 0.5 m	1
460 32	Optical bench with standardised profile, 1 m	1
460 341	Swivel joint with circular scale	1
460 373	Optics rider, 60/50	1
460 374	Optics rider, 90/50	4
460 382	Tilting rider, 90/50	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

Spectral lines arise by the transition of electrons from higher to lower energy states in the shell excited atoms. The wavelength of the emitted light depends on this energy difference:

$$\Delta E = h \cdot \nu = \frac{h \cdot c}{\lambda}$$

The multiple energy states in the term scheme of mercury results in a large number of lines with different intensities (transition probabilities). These lines can be observed in the visible range resp. measured in the near UV range.

In the experiment P6.2.2.3 the spectral lines of a high pressure mercury lamp are investigated with a high-resolution spectrometer assembly using a holographic grating. The grating works in reflection, leading to a high intensity of the lines. Different lines are observed and their wavelengths determined, especially the yellow, green, blue, violet and also the ultraviolet line. Some lines are investigated closely, e.g. the yellow double line, and the splitting of the wavelengths is determined.



**P6.2.2**  
**EMISSION AND ABSORPTION SPECTRA**

P6.2.2.4  
Recording the emission spectra of flame colouration

P6.2.2.5  
Recording Fraunhofer lines with a compact spectrometer

P6.2.2.6  
Recording the spectra of gas discharge lamps with a compact spectrometer

Recording the emission spectra of flame colouration (P6.2.2.4)


Cat. No.	Description	P6.2.2.4	P6.2.2.5	P6.2.2.6
467 251	Compact spectrometer, physics (spectral photometer)	1	1	1
460 251	Fibre holder	1	1*	1
300 11	Saddle base	1	1*	1
666 711	Butane gas burner	1		
666 712ET3	Butane cartridge, 190 g, set of 3	1		
666 731	Gas igniter, mechanical	1		
673 0840	Magnesia rods, 25 pieces	1		
604 5681	Powder spatula, steel, 150 mm	1		
667 089	Spotting tile	1		
661 088	Salts for flame tests, set of 9	1		
674 6950	Hydrochloric acid, 0.1 mol/l, 500 ml	1		
467 63	Spectral tube, Hg (with Ar)			1
467 67	Spectral tube, He			1
467 68	Spectral tube, Ar			1
467 69	Spectral tube, Ne			1
467 81	Holder for spectral tubes			1
521 70	High-voltage power supply, 10 kV			1
536 251	Measuring resistor, 100 kΩ			1
300 02	Stand base, V-shaped, small			1
300 40	Stand rod, 10 cm, 12 mm diam.			1
301 01	Leybold multiclamp			1
500 621	Safety connecting lead, 50 cm, red			1
500 622	Safety connecting lead, 50 cm, blue			1
500 611	Safety connecting lead, 25 cm, red			1
500 610	Safety connecting lead, 25 cm, yellow/green			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

\* additionally recommended

In the experiment P6.2.2.4 flame tests with metal salts are performed. A compact spectrometer at the USB port of a PC enables the easy recording of such transient processes and analyses the different emission lines. In contrast to classical observation with the eye, the spectrometer records also lines in the IR region, identifying potassium for example.

In the experiment P6.2.2.5 Fraunhofer absorption lines in the solar spectrum are recorded with a compact spectrometer. The presence of several elements in the solar photosphere is shown.

Experiment P6.2.2.6 records the spectra of gas discharge lamps using a compact and easy to use spectrometer.



Spectra of gas discharge lamps (P6.2.2.6)

P6.2.3

INELASTIC COLLISIONS  
OF ELECTRONS

P6.2.3.1  
Discontinuous energy emission  
of electrons in a gas-filled triode




Discontinuous energy emission of electrons in a gas-filled triode (P6.2.3.1)

Cat. No.	Description	P6.2.3.1
555 614	Gas triode	1
555 600	Tube stand	1
521 65	Tube power supply, 0...500 V	1
531 120	Multimeter LDanalog 20	3
500 622	Safety connecting lead, 50 cm, blue	1
500 641	Safety connecting lead, 100 cm, red	6
500 642	Safety connecting lead, 100 cm, blue	4

In inelastic collision of an electron with an atom, the kinetic energy of the electron is transformed into excitation or ionization energy of the atom. Such collisions are most probable when the kinetic energy is exactly equivalent to the excitation or ionization energy. As the excitation levels of the atoms can only assume discrete values, the energy emission in the event of inelastic electron collision is discontinuous.

The experiment P6.2.3.1 uses a tube triode filled with helium to demonstrate this discontinuous emission of energy. After acceleration in the electric field between the cathode and the grid, the electrons enter an opposing field which exists between the grid and the anode. Only those electrons with sufficient kinetic energy reach the anode and contribute to the current  $I$  flowing between the anode and ground. Once the electrons in front of the grid have reached a certain minimum energy, they can excite the gas atoms through inelastic collision. When the acceleration voltage  $U$  is continuously increased, the inelastic collisions initially occur directly in front of the grid, as the kinetic energy of the electrons reaches its maximum value here. After collision, the electrons can no longer travel against the opposing field. The anode current  $I$  is thus greatly decreased. When the acceleration voltage  $U$  is increased further, the excitation zone moves toward the cathode, the electrons can again accumulate energy on their way to the grid and the current  $I$  again increases. Finally, the electrons can excite gas atoms a second time, and the anode current drops once more.



Anode current  $I$  as a function of the acceleration voltage  $U$  for He



**P6.2.4**  
**FRANCK-HERTZ EXPERIMENT**

P6.2.4.1  
Franck-Hertz experiment with mercury - Recording with the oscilloscope and point by point


P6.2.4.2  
Franck-Hertz experiment with mercury - Recording and evaluation with CASSY

Franck-Hertz experiment with mercury - Recording with the oscilloscope (P6.2.4.1\_b)

Cat. No.	Description	P6.2.4.1 (a)	P6.2.4.1 (b)	P6.2.4.2
555 854	Hg Franck-Hertz tube	1	1	1
555 864	Socket for Hg Franck-Hertz tube, with DIN connector	1	1	1
555 81	Electric oven for tubes, 230 V	1	1	1
555 880	Franck-Hertz supply unit	1	1	1
666 193	Temperature probe, NiCr-Ni, 1.5 mm	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous		1	
575 24	Screened cable, BNC/4 mm		2	
524 013	Sensor-CASSY 2			1
524 220	CASSY Lab 2			1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair			2
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1

In 1914, *J. Franck* and *G. Hertz* reported observing discontinuous energy emission when electrons passed through mercury vapor, and the resulting emission of the ultraviolet spectral line ( $\lambda = 254 \text{ nm}$ ) of mercury. A few months later, *Niels Bohr* recognized that their experiment supported his model of the atom.

This experiment is offered in two variations, experiments P6.2.4.1 and P6.2.4.2, which differ only in the means of recording and evaluating the measurement data. The mercury atoms are enclosed in a tetrode with cathode, grid-type control electrode, acceleration grid and target electrode. The control grid ensures a virtually constant emission current of the cathode. An opposing voltage is applied between the acceleration grid and the target electrode. When the acceleration voltage  $U$  between the cathode and the acceleration grid is increased, the target current  $I$  corresponds closely to the tube characteristic once it rises above the opposing voltage. As soon as the electrons acquire sufficient kinetic energy to excite the mercury atoms through inelastic collisions, the electrons can no longer reach the target, and the target current drops. At this acceleration voltage, the excitation zone is directly in front of the excitation grid. When the acceleration voltage is increased further, the excitation zone moves toward the cathode, the electrons can again accumulate energy on their way to the grid and the target current again increases. Finally, the electrons can excite the mercury atoms once more, the target current drops again, and so forth. The  $I(U)$  characteristic thus demonstrates periodic variations, whereby the distance between the minima  $\Delta U = 4.9 \text{ V}$  corresponds to the excitation energy of the mercury atoms from the ground state  $^1S_0$  to the first  $^3P_1$  state.



Franck-Hertz curve for mercury

### P6.2.4

#### FRANCK-HERTZ EXPERIMENT

##### P6.2.4.3

Franck-Hertz experiment with neon -  
Recording with the oscilloscope and  
point by point

##### P6.2.4.4

Franck-Hertz experiment with neon -  
Recording and evaluation with CASSY




*Franck-Hertz experiment with neon - Recording and evaluation with CASSY (P6.2.4.4)*

Cat. No.	Description	P6.2.4.3 (a)	P6.2.4.3 (b)	P6.2.4.4
555 870	Ne Franck-Hertz tube	1	1	1
555 871	Holder with socket and screen	1	1	1
555 872	Connecting cable for Ne Franck-Hertz tube	1	1	1
555 880	Franck-Hertz supply unit	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous		1	
575 24	Screened cable, BNC/4 mm		2	
524 013	Sensor-CASSY 2			1
524 220	CASSY Lab 2			1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair			2
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)			1

When neon atoms are excited by means of inelastic electron collision at a gas pressure of approx. 10 hPa, excitation is most likely to occur to states which are 18.7 eV above the ground state. The de-excitation of these states can occur indirectly via intermediate states, with the emission of photons. In this process, the photons have a wavelength in the visible range between red and green. The emitted light can thus be observed with the naked eye and e.g. measured using the school spectroscope Kirchhoff/Bunsen (467 112).

The Franck-Hertz experiment with neon is offered in two variations, experiments P6.2.4.3 and P6.2.4.4, which differ only in the means of recording and evaluating the measurement data. In both variations, the neon atoms are enclosed in a glass tube with four electrodes: the cathode  $K$ , the grid-type control electrode  $G_1$ , the acceleration grid  $G_2$  and the target electrode  $A$ . Like the Franck-Hertz experiment with mercury, the acceleration voltage  $U$  is continuously increased and the current  $I$  of the electrons which are able to overcome the opposing voltage between  $G_2$  and  $A$  and reach the target is measured. The target current is always lowest when the kinetic energy directly in front of grid  $G_2$  is just sufficient for collision excitation of the neon atoms, and increases again with the acceleration voltage. We can observe clearly separated luminous red layers between grids  $G_1$  and  $G_2$ ; their number increases with the voltage. These are zones of high excitation density, in which the excited atoms emit light in the visible spectrum.



*Luminous layers between control electrode and acceleration grid*

P6.2.6  
ELECTRON SPIN RESONANCE

P6.2.6.2  
Electron spin resonance at DPPH - determining the magnetic field as a function of the resonance frequency

P6.2.6.3  
Resonance absorption of a passive RF oscillator circuit



Electron spin resonance at DPPH - determining the magnetic field as a function of the resonance frequency (P6.2.6.2)

Cat. No.	Description	P6.2.6.2	P6.2.6.3
514 55	ESR basic unit	1	1
514 571	ESR supply unit	1	1
555 604	Pair of Helmholtz coils	1	
575 214	Oscilloscope 30 MHz, two-channel, analogous	1	1
501 02	BNC cable, 1 m	2	
300 11	Saddle base	3	2
501 20	Connecting lead, 32 A, 25 cm, red	1	
501 25	Connecting lead, 32 A, 50 cm, red	1	
501 26	Connecting lead, 32 A, 50 cm, blue	1	
531 120	Multimeter LDanalog 20		1
575 24	Screened cable, BNC/4 mm		1
501 644	Two-way adapters, black, set of 6		1
590 13	Stand rod with bore holes		1

The magnetic moment of the unpaired electron with the total angular momentum  $j$  in a magnetic field assumes the discrete energy states

$$E_m = -g_j \cdot \mu_B \cdot m \cdot B \quad \text{where } m = -j, -j + 1, \dots, j$$

$$\mu_B = 9.274 \cdot 10^{-24} \frac{\text{J}}{\text{T}} : \text{Bohr's magneton}$$

$g_j$ :  $g$  factor

When a high-frequency magnetic field with the frequency  $\nu$  is applied perpendicularly to the first magnetic field, it excites transitions between the adjacent energy states when these fulfill the resonance condition

$$h \cdot \nu = E_{m+1} - E_m$$

$h$ : Planck's constant

This fact is the basis for electron spin resonance, in which the resonance signal is detected using radio-frequency technology. The electrons can often be regarded as free electrons. The  $g$ -factor then deviates only slightly from that of the free electron ( $g = 2.0023$ ), and the resonance frequency  $\nu$  in a magnetic field of 1 mT is about 27.8 MHz. The actual aim of electron spin resonance is to investigate the internal magnetic fields of the sample substance, which are generated by the magnetic moments of the adjacent electrons and nuclei.

The experiment P6.2.6.2 verifies electron spin resonance in diphenylpicryl-hydrazyl (DPPH). DPPH is a radical, in which a free electron is present in a nitrogen atom. In the experiment, the magnetic field  $B$  which fulfills the resonance condition the resonance frequencies  $\nu$  can be set in a continuous range from 13 to 130 MHz. The aim of the evaluation is to determine the  $g$  factor.

The object of the experiment P6.2.6.3 is to verify resonance absorption using a passive oscillator circuit.




Diagram of resonance condition of free electrons

P6.2.7

NORMAL ZEEMAN EFFECT

P6.2.7.3

Observing the normal Zeeman effect in transverse and longitudinal configuration - spectroscopy using a Fabry-Perot etalon

P6.2.7.4

Measuring the Zeeman splitting of the red cadmium line as a function of the magnetic field - spectroscopy using a Fabry-Perot etalon



Measuring the Zeeman splitting of the red cadmium line as a function of the magnetic field - spectroscopy using a Fabry-Perot etalon (P6.2.7.4\_b)

Cat. No.	Description	P6.2.7.3 (b)	P6.2.7.4 (b)
451 12	Spectral lamp, Cd, with holding plate	1	1
451 30	Universal choke, 230 V, 50 Hz	1	1
562 11	U-core with yoke	1	1
562 131	Coil, 480 turns, 10 A	2	2
560 315	Pole pieces with large bore, pair	1	1
521 55	High current power supply	1	1
471 221	Fabry-Perot etalon mirror, in holder	1	1
460 08	Lens in frame, f=150 mm	2	2
472 601	Quarter-wavelength plate, 140 nm	1	
472 401	Polarisation filter	1	
468 41	Holder for interference filters	1	1
468 400	Interference filter, 644 nm	1	1
460 135	Ocular with scale	1	
460 32	Optical bench with standardised profile, 1 m	1	1
460 381	Rider base with thread	1	1
460 373	Optics rider, 60/50	7	5
501 33	Connecting lead, 32 A, 100 cm, black	3	3
337 47USB	VideoCom		1
524 005	Mobile-CASSY 2		1
524 0381	Combi B sensor S		1
501 11	Extension cable, 15 pin		1
300 02	Stand base, V-shaped, small		1
300 42	Stand rod, 47 cm, 12 mm diam.		1
301 01	Leybold multiclamp		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

The Zeeman effect is the name for the splitting of atomic energy levels in an external magnetic field and, as a consequence, the splitting of the transitions between the levels. The effect was predicted by *H. A. Lorentz* in 1895 and experimentally confirmed by *P. Zeeman* one year later. In the red spectral line of cadmium ( $\lambda = 643.8 \text{ nm}$ ), Zeeman observed a line triplet perpendicular to the magnetic field and a line doublet parallel to the magnetic field, instead of just a single line. Later, even more complicated splittings were discovered for other elements, and were collectively designated the anomalous Zeeman effect. It eventually became apparent that the normal Zeeman effect is the exception, as it only occurs at transitions between atomic levels with the total spin  $S = 0$ .

In the experiment P6.2.7.3, the Zeeman effect is observed at the red cadmium line perpendicular and parallel to the magnetic field, and the polarization state of the individual Zeeman components is determined. The observations are explained on the basis of the radiating characteristic of dipole radiation. The so-called  $\pi$  component corresponds to a Hertzian dipole oscillating parallel to the magnetic field, i.e. it cannot be observed parallel to the magnetic field and radiates linearly polarized light perpendicular to the magnetic field. Each of the two  $\sigma$  components corresponds to two dipoles oscillating perpendicular to each other with a phase differential of  $90^\circ$ . They radiate circularly polarized light in the direction of the magnetic field and linearly polarized light parallel to it.


In the experiment P6.2.7.4, the Zeeman splitting of the red cadmium line is measured as a function of the magnetic field  $B$ . The energy interval of the triplet components

$$\Delta E = \frac{h}{4\pi} \cdot \frac{e}{m_e} \cdot B$$

$m_e$ : mass of electron,  $e$ : electron charge


$h$ : Planck's constant

$B$ : magnetic induction



Zeeman splitting as a function of magnetic field





## P6.2.8 OPTICAL PUMPING (ANOMALOUS ZEEMAN EFFECT)

- P6.2.8.1  
Optical pumping: observing the relaxation signal
- P6.2.8.2  
Optical pumping: measuring and observing the Zeeman transitions in the ground states of Rb-87 with  $\sigma^+$ - and  $\sigma^-$ -pumped light
- P6.2.8.3  
Optical pumping: measuring and observing the Zeeman transitions in the ground states of Rb-85 with  $\sigma^+$ - and  $\sigma^-$ -pumped light
- P6.2.8.4  
Optical pumping: measuring and observing the Zeeman transitions in the ground states of Rb-87 as a function of the magnetic flux density  $B$
- P6.2.8.5  
Optical pumping: measuring and observing the Zeeman transitions in the ground states of Rb-85 as a function of the magnetic flux density  $B$
- P6.2.8.6  
Optical pumping: measuring and observing two-quantum transitions

Optical pumping: measuring and observing the Zeeman transitions in the ground states of Rb-87 with  $\sigma^+$ - and  $\sigma^-$ -pumped light (P6.2.8.2)

Cat. No.	Description	P6.2.8.1	P6.2.8.2-6
558 823	Rubidium high-frequency lamp	1	1
558 826	Helmholtz coils on rider	1	1
558 833	Absorption chamber with Rb cell	1	1
558 835	Silicon photodetector	1	1
558 836	I/U converter for silicon photodetector	1	1
530 88	Plug-in power unit, 230 V/9.2 V DC	1	1
558 814	Supply unit for optical pumping	1	1
521 45	DC power supply 0...±15 V	1	1
501 02	BNC cable, 1 m	2	3
575 230	Digital storage oscilloscope 722	1	1
531 282	Multimeter Metrahit Pro	1	1
504 48	Two-way switch	1	1
468 000	Line filter, 795 nm	1	1
472 410	Polarisation filter for red radiation	1	1
472 611	Quarter-wavelength plate, 200 nm	1	1
460 021	Lens in frame, f=50 mm	1	1
460 031	Lens in frame, f=100 mm	1	1
460 32	Optical bench with standardised profile, 1 m	1	1
460 370	Optics rider, 60/34	6	6
460 374	Optics rider, 90/50	1	1
666 7681	Circulation thermostat SC 100-S5P	1	1
688 115	Silicone tubing 6 x 2 mm, 5.0 m	1	1
501 28	Connecting lead, 32 A, 50 cm, black	4	4
501 38	Connecting lead, 32 A, 200 cm, black	2	2
675 3410	Water, pure, 5 l	2	2
522 551	Function generator, 1 mHz to 12 MHz, internal sweep		1
501 022	BNC cable, 2 m		1

The two hyperfine structures of the ground state of an alkali atom with the total angular momentums

$$F_+ = I + \frac{1}{2}, \quad F_- = I - \frac{1}{2}$$

split in a magnetic field  $B$  into  $2F_{\pm} + 1$  Zeeman levels having an energy which can be described using the Breit-Rabi formula

$$E = \frac{-\Delta E}{2(2I+1)} + \mu_K g_I m_F \pm \frac{\Delta E}{2} \sqrt{1 + \frac{4m_F}{2I+1} \xi + \xi^2}$$

$$\text{where } \xi = \frac{g_J \mu_B - g_I \mu_K}{\Delta E} \cdot B$$

$\Delta E$ : hyperfine structure interval

$I$ : nuclear spin,  $m_F$ : magnetic quantum number

$\mu_B$ : Bohr's magneton,  $\mu_K$ : nuclear magneton

$g_J$ : shell g factor,  $g_I$ : nuclear g factor

Transitions between the Zeeman levels can be observed using a method developed by A. Kastler. When right-handed or left-handed circularly polarized light is directed parallel to the magnetic field, the population of the Zeeman level differs from the thermal equilibrium population, i.e. optical pumping occurs, and RF radiation forces transitions between the Zeeman levels.

The change in the equilibrium population when switching from right-handed to left-handed circular pumped light is verified in the experiment P6.2.8.1.

The experiments P6.2.8.2 and P6.2.8.3 measure the Zeeman transitions in the ground state of the isotopes Rb-87 and Rb-85 and determine the nuclear spin  $I$  from the number of transitions observed. The observed transitions are classified through comparison with the Breit-Rabi formula.

In the experiments P6.2.8.4 and P6.2.8.5, the measured transition frequencies are used for precise determination of the magnetic field  $B$  as a function of the magnet current  $I$ . The nuclear g factors  $g_I$  are derived using the measurement data.

In the experiment P6.2.8.6, two-quantum transitions are induced and observed for a high field strength of the irradiating RF field.

### P6.3.1

#### DETECTION OF X-RAYS

##### P6.3.1.1

Fluorescence of a luminescent screen due to X-rays

##### P6.3.1.2

X-ray photography: Exposure of film stock due to X-rays

##### P6.3.1.5

Investigation of an implant model

##### P6.3.1.6

Influence of a contrast medium on the absorption of X-rays



X-ray photography: Exposure of film stock due to X-rays (P6.3.1.2)

Cat. No.	Description	P6.3.1.1	P6.3.1.2	P6.3.1.5	P6.3.1.6
554 800	X-ray apparatus	1	1	1	1
554 861	X-ray tube, Mo	1	1	1	1
554 838	Film holder, X-ray		1	1	
554 896	X-ray film, Agfa Dentus M2		1		
554 8971	Developer and fixer for X-ray film		1		
554 8931	Changing bag with developer tank		1*		
554 8391	Implant model			1	
554 839	Blood vessel model for contrast medium				1
602 023	Beaker, Boro 3.3, 150 ml, squat				1
602 295	Bottle brown glass wide treath with cap				1
602 783	Glass rod, 200 mm, diam. 6 mm				1
672 6610	Potassium iodide, 100 g				1

\* additionally recommended



Screen of the implant model



Screen of the blood vessel


Soon after the discovery of X-rays by *W. C. Röntgen*, physicians began to exploit the ability of this radiation to pass through matter which is opaque to ordinary light for medical purposes. The technique of causing a luminescent screen to fluoresce with X-ray radiation is still used today for screen examinations, although image amplifiers are used additionally. The exposure of a film due to X-ray radiation is used both for medical diagnosis and materials testing, and is the basis for dosimetry with films.

The experiment P6.3.1.1 demonstrates the transillumination with X-rays using simple objects made of materials with different absorption characteristics. A luminescent screen of zinc-cadmium sulfate is used to detect X-rays; the atoms in this compound are excited by the absorption of X-rays and emit light quanta in the visible light range. This experiment investigates the effect of the emission current  $I$  of the X-ray tube on the brightness and the effect of the high voltage  $U$  on the contrast of the luminescent screen.

The experiment P6.3.1.2 records the transillumination of objects using X-ray film. Measuring the exposure time required to produce a certain degree of exposure permits quantitative conclusions regarding the intensity of the X-rays.

The experiment P6.3.1.5 demonstrates the use of radioscopy to detect hidden objects. A metal rod inside a block of wood is visually invisible, but can be seen by X-ray fluorescence and its dimensions measured.

The experiment P6.3.1.6 demonstrates the use of contrast medium. The radiopaque iodine solution is flowing through channels inside a plate and is clearly visible in the X-ray fluorescence image, but pure water is not.



**P6.3.1**  
**DETECTION OF X-RAYS**

P6.3.1.3  
Detecting X-rays using an ionization chamber

P6.3.1.4  
Determining the ion dose rate of the X-ray tube with molybdenum anode

Detecting X-rays using an ionization chamber (P6.3.1.3)


Cat. No.	Description	P6.3.1.3-4
554 800	X-ray apparatus	1
554 861	X-ray tube, Mo	1
554 840	Plate capacitor, X-ray	1
522 27	Power supply, 450 V	1
532 14	Electrometer amplifier	1
577 02	Resistor, 1 GΩ, STE2/19	1
531 120	Multimeter LDanalog 20	2
575 24	Screened cable, BNC/4 mm	1
501 451	Connecting leads, 19 A, 50 cm, black, pair	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2

As X-rays ionize gases, they can also be measured via the ionization current of an ionization chamber.

The aim of the experiments P6.3.1.3 and P6.3.1.4 is to detect X-rays using an ionization chamber. First, the ionization current is recorded as a function of the voltage at the capacitor plates of the chamber and the saturation range of the characteristic curves is identified. Next, the mean ion dose rate

$$\bar{j} = \frac{I_{\text{ion}}}{m}$$

is calculated from the ionization current  $I_{\text{ion}}$  which the X-radiation generates in the irradiated volume of air  $V$ , and the mass  $m$  of the irradiated air. The measurements are conducted for various emission currents  $I$  and high voltages  $U$  of the X-ray tube.



Mean ion dose rate  $j$  as a function of the tube high voltage  $U$ ,  $I = 1.0 \text{ mA}$

### P6.3.1

#### DETECTION OF X-RAYS

P6.3.1.7  
Digital X-ray photography,  
measurement of the attenuation

P6.3.1.8  
Influence of the high voltage  
and the anode current on the  
contrast (hardening)

P6.3.1.9  
Resolution of the digital X-ray  
photography

P6.3.1.10  
Determining the size of the  
focal spot



Resolution of the digital X-ray photography (P6.3.1.9\_b)

Cat. No.	Description	P6.3.1.7 (a)	P6.3.1.7 (b)	P6.3.1.8 (a)	P6.3.1.8 (b)	P6.3.1.9 (a)	P6.3.1.9 (b)	P6.3.1.10
554 800	X-ray apparatus	1	1	1	1	1	1	1
554 866	X-ray tube, Au	1	1	1	1	1	1	1
554 821	Computed tomography module	1		1		1		
554 838	Film holder, X-ray	1		1		1		
554 8382	Object holder X-ray 1	1		1		1		
554 834	Absorption accessory, X-ray	1	1	1	1			
554 828	X-ray image sensor		1		1		1	1
554 829	Precision rails for X-ray image sensor		1		1		1	1
554 8292	Object holder X-ray 2		1		1		1	1
554 822	Test patterns, 2D						1	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1	1	1	1


In X-ray radiography and image processing it is necessary to know the properties of the radiation source as well as the properties of the irradiated object. The minimum resolution in an image is determined by the size of the focal spot and the visibility of details in an image depends on the different attenuation properties of various materials and the radiation spectrum of the X-ray-source.

The experiment P6.3.1.7 measures the attenuation of X-rays through the quantitative analysis of a digital X-ray image on the computer.


The contrast of a digital X-ray image as a function of the energy and intensity of x-radiation is measured in experiment P6.3.1.8.

In the experiment P6.1.3.9 the resolution digital radiographs is determined using various test images.

The focal spot of the anode of the x-ray tube is measured by a simple geometric projection on the X-ray image sensor in the experiment P6.3.1.10.



Test patterns for determination of the resolution of the X-ray image sensor



**P6.3.2**  
**ATTENUATION OF X-RAYS**

P6.3.2.1  
Investigating the attenuation of X-rays as a function of the absorber material and absorber thickness

P6.3.2.2  
Investigating the wavelength dependency of the attenuation coefficient

P6.3.2.3  
Investigating the relationship between the attenuation coefficient and the atomic number  $Z$

P6.3.2.4  
Investigating the attenuation of X-rays as a function of distance

Investigating the attenuation of X-rays as a function of the absorber material and absorber thickness (P6.3.2.1)

Cat. No.	Description	P6.3.2.1	P6.3.2.2	P6.3.2.3	P6.3.2.4
554 800	X-ray apparatus	1	1	1	1
554 861	X-ray tube, Mo	1	1	1	
554 831	Goniometer	1	1	1	
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1	1	1	
554 834	Absorption accessory, X-ray	1			
554 78	NaCl crystal for Bragg reflection		1	1	
554 832	Set of absorber foils		1	1	
554 866	X-ray tube, Au				1
554 828	X-ray image sensor				1
554 829	Precision rails for X-ray image sensor				1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1		1

The attenuation of X-rays on passing through an absorber with the thickness  $d$  is described by Lambert's law for attenuation:

$$I = I_0 \cdot e^{-\mu d}$$

$I_0$ : intensity of primary beam

$I$ : transmitted intensity

Here, the attenuation is due to both absorption and scattering of the X-rays in the absorber. The linear attenuation coefficient  $\mu$  depends on the material of the absorber and the wavelength  $\lambda$  of the X-rays. An absorption edge, i.e. an abrupt transition from an area of low absorption to one of high absorption, may be observed when the energy  $h \cdot \nu$  of the X-ray quantum just exceeds the energy required to move an electron out of one of the inner electron shells of the absorber atoms.

The object of the experiment P6.3.2.1 is to confirm Lambert's law using aluminium and to determine the attenuation coefficients  $\mu$  for six different absorber materials averaged over the entire spectrum of the X-ray apparatus.


The experiment P6.3.2.2 records the transmission curves

$$T(\lambda) = \frac{I(\lambda)}{I_0(\lambda)}$$

for various absorber materials. The aim of the evaluation is to confirm the  $\lambda^3$  relationship of the attenuation coefficients for wavelengths outside of the absorption edges.

In the experiment P6.3.2.3, the attenuation coefficient  $\mu(\lambda)$  of different absorber materials is determined for a wavelength  $\lambda$  which lies outside of the absorption edge. This experiment reveals that the attenuation coefficient is closely proportional to the fourth power of the atomic number  $Z$  of the absorbers.

In the experiment P6.3.2.4 the intensity of the X-radiation is measured at different distances from the X-ray tube. The inverse square distance law is confirmed.



### P6.3.3

#### PHYSICS OF THE ATOMIC SHELL


P6.3.3.1  
Bragg reflection: diffraction of X-rays at a monocrystal

P6.3.3.2  
Investigating the energy spectrum of an X-ray tube as a function of the high voltage and the emission current

P6.3.3.3  
Duane–Hunt relation and determination of Planck's constant

P6.3.3.5  
Edge absorption: filtering X-rays

P6.3.3.6  
Moseley's law and determination of the Rydberg constant



Investigating the energy spectrum of an X-ray tube as a function of the high voltage and the emission current (P6.3.3.2)

Cat. No.	Description	P6.3.3.1-3	P6.3.3.5	P6.3.3.6
554 801	X-ray apparatus, Mo, complete	1	1	1
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1	1	1
554 832	Set of absorber foils			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

The radiation of an X-ray tube consists of two components: continuous bremsstrahlung radiation is generated when fast electrons are decelerated in the anode. Characteristic radiation consisting of discrete lines is formed by electrons dropping to the inner shells of the atoms of the anode material from which electrons were liberated by collision.

To confirm the wave nature of X-rays, the experiment P6.3.3.1 investigates the diffraction of the characteristic  $K_\alpha$  and  $K_\beta$  lines of the molybdenum anode at an NaCl monocrystal and explains these using Bragg's law of reflection.

The experiment P6.3.3.2 records the energy spectrum of the X-ray apparatus as a function of the high voltage and the emission current using a goniometer in the Bragg configuration. The aim is to investigate the spectral distribution of the continuum of bremsstrahlung radiation and the intensity of the characteristic lines.

The experiment P6.3.3.3 measures how the limit wavelength  $\lambda_{\min}$  of the continuum of bremsstrahlung radiation depends on the high voltage  $U$  of the X-ray tube. When we apply the Duane-Hunt relationship

$$e \cdot U = h \cdot \frac{c}{\lambda_{\min}}$$

$e$ : electron charge

$c$ : velocity of light


to the measurement data, we can derive Planck's constant  $h$ .

The object of the experiment P6.3.3.5 is to filter X-rays using the absorption edge of an absorber, i. e. the abrupt transition from an area of low absorption to one of high absorption.

The experiment P6.3.3.6 determines the wavelengths  $\lambda_K$  of the absorption edges as a function of the atomic number  $Z$ . When we apply Moseley's law

$$\frac{1}{\lambda_K} = R \cdot (Z - \sigma)^2$$

to the measurement data we obtain the Rydberg constant  $R$  and the mean screening  $\sigma$ .



Bragg spectrum of a molybdenum anode

P6.3.5

X-RAY ENERGY SPECTROSCOPY

P6.3.5.1  
Recording and calibrating an X-ray energy spectrum

P6.3.5.2  
Recording the energy spectrum of a molybdenum anode

P6.3.5.3  
Recording the energy spectrum of a copper anode

P6.3.5.4  
Investigation of the characteristic spectra as a function of the element's atomic number: K-lines

P6.3.5.5  
Investigation of the characteristic spectra as a function of the element's atomic number: L-lines

P6.3.5.6  
Energy-resolved Bragg reflection in different orders of diffraction



Recording and calibrating an X-ray energy spectrum (P6.3.5.1)

Cat. No.	Description	P6.3.5.1-2	P6.3.5.3	P6.3.5.4	P6.3.5.5	P6.3.5.6
554 800	X-ray apparatus	1	1	1	1	1
554 861	X-ray tube, Mo	1		1	1	
554 831	Goniometer	1	1	1	1	1
559 938	X-ray energy detector	1	1	1	1	1
524 013	Sensor-CASSY 2	1	1	1	1	1
524 220	CASSY Lab 2	1	1	1	1	1
524 058	MCA box	1	1	1	1	1
501 02	BNC cable, 1 m	1	1	1	1	1
554 862	X-ray tube, Cu		1			1
554 844	Set of targets for K-line fluorescence			1		
554 846	Set of targets for L-line fluorescence				1	
554 78	NaCl crystal for Bragg reflection					1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1	1



Recording and calibrating an X-ray energy spectrum

The X-ray energy detector enables recording of the energy spectrum of X-rays. The detector is a Peltier-cooled photodiode where in the incoming X-rays produce electron-hole pairs. The number of electron-hole pairs and thus the voltage pulse height after amplification is proportional to the X-ray energy. The pulse height analysis is carried out with CASSY used as a multichannel analyzer (MCA-Box), which is connected to a computer (PC).


The object of the experiment P6.3.5.1 is to record the X-ray fluorescence spectrum of a target and to use the known energies for calibration of the energy axis. The target is made of a zincplated steel and emits several fluorescent lines.

The experiments P6.3.5.2 and P6.3.5.3 use the calibrated detector to record emission spectra of either a molybdenum anode or a copper anode. The resulting spectrum shows the characteristic lines of the anode material and the bremsstrahlung continuum.

The experiment P6.3.5.4 demonstrates differences in the characteristic fluorescent K-lines (transitions to K-shell) within the X-ray spectra of different elements. These are used to confirm Moseley's law and show aspects of material analysis.

The experiment P6.3.5.5 shows similar characteristic fluorescent L-lines for heavier elements, demonstrating the X-ray emission from transitions to the L-shell.

In the experiment P6.3.5.6 using the X-ray energy detector in Bragg geometry it is possible to observe different X-ray energies simultaneously, because Bragg condition is fulfilled for different orders.



X-ray fluorescence of different elements (P6.3.5.4/5)

### P6.3.6


#### STRUCTURE OF X-RAY SPECTRUMS

- P6.3.6.1  
Fine structure of the characteristic X-ray radiation of a molybdenum anode
- P6.3.6.2  
Fine structure of the characteristic X-ray radiation of a copper anode
- P6.3.6.3  
Fine structure of the characteristic X-ray radiation of an iron anode
- P6.3.6.4  
Fine structure of the characteristic X-ray radiation of a silver anode
- P6.3.6.5  
Fine structure of the characteristic X-ray radiation of a tungsten anode
- P6.3.6.6  
Determining the binding energy of individual subshells by selective excitation
- P6.3.6.7  
Fine structure of the characteristic X-ray radiation of a gold anode



Fine structure of the characteristic X-ray radiation of a tungsten anode (P6.3.6.5)

Cat. No.	Description	P6.3.6.1	P6.3.6.2	P6.3.6.3	P6.3.6.4	P6.3.6.5-6	P6.3.6.7
554 800	X-ray apparatus	1	1	1	1	1	1
554 861	X-ray tube, Mo	1					
554 831	Goniometer	1	1	1	1	1	1
554 78	NaCl crystal for Bragg reflection	1			1		
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1	1	1	1	1	1
554 862	X-ray tube, Cu		1				
554 791	KBr crystal for Bragg reflection		1				
554 863	X-ray tube, Fe			1			
554 77	LiF crystal for Bragg reflection			1		1	1
554 865	X-ray tube, Ag				1		
554 864	X-ray tube, W					1	
554 866	X-ray tube, Au						1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1	1	1



Bragg spectrum of the tungsten anode

The structure and fine-structure of X-ray spectra gives valuable information on the position of the atomic energy levels. The systematics of X-ray transitions are presented. Starting with molybdenum and completed with other anode materials like copper and iron the K-shell transitions of light and medium elements are investigated.

In contrast to these materials the heavy elements like tungsten show characteristic emission from the L-shell with a lot of details, because the lower level of the transition consists of several sublevels which can also be selectively excited.

The experiment P6.3.6.1 investigates the X-ray spectrum of a molybdenum anode and the fine structure of the  $K_{\alpha}$  line.

The experiments P6.3.6.2 and P6.3.6.3 observe the low-energy characteristic radiation from a copper or iron anode and the fine structure of the  $K_{\alpha}$  line.

The experiment P6.3.6.4 observes the high-energy characteristic radiation of silver and the fine structure splitting due to spin-orbit coupling.

The experiment P6.3.6.5 demonstrates the fine structure of the tungsten L-lines. Due to the splitting of the energy levels there are approximately 11 transitions visible ( $L_{\alpha 1-2}$ ,  $L_{\beta 1-4}$ ,  $L_{\gamma 1-3}$ ,  $L_{\delta}$ ,  $L_{\nu}$ ), which can be used to evaluate the position of the energy levels and to demonstrate allowed and forbidden transitions.

In addition to experiment P6.3.6.5, the experiment P6.3.6.6 measures directly the splitting of the L-shell. At a low acceleration voltage only the L3 level can be excited, with raising voltages transitions to L2 and later L1 become observable. The absolute binding energies of the L-sublevels can be measured directly.

The experiment P6.3.6.7 demonstrates the fine structure of the gold L-lines. Due to the splitting of the energy levels there are approximately 10 transitions visible ( $L_{\alpha 1-2}$ ,  $L_{\beta 1-4}$ ,  $L_{\gamma 1-3}$ ,  $L_{\delta}$ ,  $L_{\nu}$ ), which can be used to evaluate the position of the energy levels and to demonstrate allowed and forbidden transitions.





**P6.3.6**  
**STRUCTURE OF X-RAY SPECTRUMS**

P6.3.6.11  
High-resolution fine structure of the characteristic X-ray radiation of a molybdenum anode

P6.3.6.12  
High-resolution fine structure of the characteristic X-ray radiation of a copper anode

P6.3.6.13  
High-resolution fine structure of the characteristic X-ray radiation of an iron anode

P6.3.6.14  
High-resolution fine structure of the characteristic X-ray radiation of a silver anode

P6.3.6.15  
High-resolution fine structure of the characteristic X-ray radiation of a tungsten anode

P6.3.6.17  
High-resolution fine structure of the characteristic X-ray radiation of a gold anode

High-resolution fine structure of the characteristic X-ray radiation of a molybdenum anode (P6.3.6.11)

Cat. No.	Description	P6.3.6.11	P6.3.6.12	P6.3.6.13	P6.3.6.14	P6.3.6.15	P6.3.6.17
554 800	X-ray apparatus	1	1	1	1	1	1
554 861	X-ray tube, Mo	1					
554 831	Goniometer	1	1	1	1	1	1
554 835	HD Accessory, X-ray	1	1	1	1	1	1
554 78	NaCl crystal for Bragg reflection	1			1		
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1	1	1	1	1	1
554 862	X-ray tube, Cu		1				
554 77	LiF crystal for Bragg reflection		1	1		1	1
554 863	X-ray tube, Fe			1			
554 865	X-ray tube, Ag				1		
554 864	X-ray tube, W					1	
554 866	X-ray tube, Au						1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1	1	1


The experiment P6.3.6.11 investigates the high-resolution X-ray spectrum of a molybdenum anode and the fine structure of the  $K_{\alpha}$ -line.

The experiments P6.3.6.12 and P6.3.6.13 observe the low-energy characteristic radiation from a copper or iron anode and the high-resolution fine structure of the  $K_{\alpha}$ -line.

The experiment P6.3.6.14 observes the high-energy characteristic radiation of silver and the high-resolution fine structure splitting due to spin-orbit coupling.

The experiment P6.3.6.15 demonstrates the high-resolution fine structure of the tungsten L-lines. Due to the splitting of the energy levels there are approximately 10 transitions visible ( $L_{\alpha}1-2$ ,  $L_{\beta}1-5$ ,  $L_{\gamma}1-3$ ), which can be used to evaluate the position of the energy levels and to demonstrate allowed and forbidden transitions.

The experiment P6.3.6.17 demonstrates the fine structure of the gold L-lines. Due to the splitting of the energy levels there are approximately 10 transitions visible ( $L_{\alpha}1-2$ ,  $L_{\beta}1-4$ ,  $L_{\gamma}1-3$ ,  $L_{\delta}1$ ,  $L_{\delta}2$ ), which can be used to evaluate the position of the energy levels and to demonstrate allowed and forbidden transitions.



Separation of the  $K_{\alpha 1}$  and  $K_{\alpha 2}$  lines in the Bragg spectrum

P6.3.7

COMPTON EFFECT AT X-RAYS

P6.3.7.1

Compton effect: verifying the energy loss of the scattered X-ray quantum

P6.3.7.2

Compton effect: Measurement the energy of the scattered photons as a function of the scattering angle



Compton effect: Measurement the energy of the scattered photons as a function of the scattering angle (P6.3.7.2)

Cat. No.	Description	P6.3.7.1	P6.3.7.2
554 800	X-ray apparatus	1	1
554 861	X-ray tube, Mo	1	1
554 831	Goniometer	1	1
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1	
554 836	Compton accessory, X-ray	1	
554 8371	Compton accessory, X-ray II		1
559 938	X-ray energy detector		1
524 013	Sensor-CASSY 2		1
524 058	MCA box		1
524 220	CASSY Lab 2		1
501 02	BNC cable, 1 m		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1


At a time (early 1920's) when the particle nature of light (photons) suggested by the photoelectric effect was still being debated, the Compton experiment, the scattering of X-rays on weakly bound electrons, in 1923 gave another evidence of particle-like behaviour of X-rays in this process. Compton investigated the scattering of X-rays passing through matter. According to classical physics the frequency of the radiation should not be changed by the scattering process. However, *A. H. Compton* observed a frequency change for scattered X-rays. He interpreted this in the particle model as a collision of the X-ray photon and an electron of the scattering material. Assuming total energy and momentum to be conserved, energy is transferred from the photon to the electron, so the energy of the scattered photon depends on the scattering angle  $\vartheta$ .

The experiment P6.3.7.1 verifies the Compton shift using the endwindow counter. The change of frequency or wavelength due to the scattering process is apparent as a change of the attenuation of an absorber, which is placed either in front of or behind the scattering body.

The object of the experiment P6.3.7.2 is to record directly the energy spectra of the scattered X-rays with the X-ray energy detector as a function of the scattering angle  $\vartheta$ . The energy  $E(\vartheta)$  of the scattered photons at different angles is determined and compared with the calculated energy obtained from conservation of energy and momentum by using the relativistic expression for the energy:

$$E(\vartheta) = \frac{E_0}{1 + \frac{E_0}{m \cdot c^2} \cdot (1 - \cos \vartheta)}$$

$E_0$ : energy of the photon before the collision  
 $m$ : mass of electron at rest  
 $c$ : velocity of light



Energy shift of the scattered X-rays at different angles (P6.3.7.2)

## P6.3.8

### X-RAY TOMOGRAPHY

#### P6.3.8.1

Measurement and presentation of a computed tomogram with the computed tomography module

#### P6.3.8.2

Computed tomography of simple geometrical objects with the computed tomography module

#### P6.3.8.3

Medical basics of computed tomography with the computed tomography module

#### P6.3.8.4

Determining absorption coefficients and Hounsfield units with computed tomography with the computed tomography module

#### P6.3.8.5

Computed tomography of biological samples with the computed tomography module



Measurement and presentation of a computed tomogram with the computed tomography module (P6.3.8.1)

Cat. No.	Description	P6.3.8.1	P6.3.8.2-4	P6.3.8.5
554 800	X-ray apparatus	1	1	1
554 831	Goniometer	1	1	1
554 866	X-ray tube, Au	1	1	1
554 821	Computed tomography module	1	1	1
554 823	Phantom, 3D		1	
554 825	LEGO® adapter		1	
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

In 1972 the first computed tomographic scanner was built by Godfrey Hounsfield who, together with Allan Cormack, was awarded the Nobel Prize in Physiology or Medicine in 1979. The basic idea of computed tomography (CT) is the illumination of an object by X-rays from numerous different angles. Our educational X-ray apparatus allows the illumination of objects by X-rays. The resulting 2D-projections are visualised at the fluorescence screen or with the X-ray image sensor. By turning an object using the built-in goniometer of the X-ray apparatus, and recording the 2D-projections from each angular step, the computer can reconstruct the object illuminated by X-rays. Our e-learning software visualises the back projection (necessary for reconstructing the computed tomography) concurrently with the scanning process. The 3D-model is then displayed on the PC screen.


Experiment P6.3.8.1 discusses the basics of computed tomography. The computed tomographies of simple geometrical objects are recorded and displayed.

Experiment P6.3.8.2 shows the CT of simple geometrical objects to demonstrate the basic properties and the resolution power of tomography.

Experiment P6.3.8.3 shows the CT of medical related objects to demonstrate the basic properties of tomography and to investigate possible problems with artefacts and image errors.

Experiment P6.3.8.4 analyses the absorption coefficient of water inside a plastic body to demonstrate the capabilities of CT in distinguishing different kinds of tissues and to calibrate Hounsfield units and it discusses hardening effects of the X-rays.

Experiment P6.3.8.5 analyses the CT of real biological specimens and applies to the results of the previous experiments.



Computed tomography of a Lego figure (P6.3.8.2)

**P6.3.8**

**X-RAY TOMOGRAPHY**

P6.3.8.11

Measurement and presentation of a computed tomogram with the X-ray image sensor

P6.3.8.12

Computed tomography of simple geometrical objects with the X-ray image sensor

P6.3.8.13

Medical basics of computed tomography with the X-ray image sensor

P6.3.8.14

Determining absorption coefficients and Hounsfield units with computed tomography with the X-ray image sensor

P6.3.8.15

Computed tomography of biological samples with the X-ray image sensor



Measurement and presentation of a computed tomogram with the X-ray image sensor (P6.3.8.11)

Cat. No.	Description	P6.3.8.11	P6.3.8.12-4	P6.3.8.15
554 800	X-ray apparatus	1	1	1
554 831	Goniometer	1	1	1
554 866	X-ray tube, Au	1	1	1
554 820P	Computed Tomography Pro package	1	1	1
554 823	Phantom, 3D		1	
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1


Experiment P6.3.8.11 discusses the basics of computed tomography. The computed tomographies of simple geometrical objects are recorded and displayed.

Experiment P6.3.8.12 shows the CT of simple geometrical objects to demonstrate the basic properties and the resolution power of tomography.


Experiment P6.3.8.13 shows the CT of medical related objects to demonstrate the basic properties of tomography and to investigate possible problems with artefacts and image errors.

Experiment P6.3.8.14 analyses the absorption coefficient of water inside a plastic body to demonstrate the capabilities of CT in distinguishing different kinds of tissues and to calibrate Hounsfield units and it discusses hardening effects of the X-rays.

Experiment P6.3.8.15 analyses the CT of real biological specimens and applies to the results of the previous experiments.



High-resolution computed tomograms of a snail shell



## P6.4.1

### DETECTING RADIOACTIVITY

#### P6.4.1.1

Ionization of air through radioactivity

#### P6.4.1.3

Demonstrating radioactive radiation with a Geiger counter

#### P6.4.1.4

Recording the characteristic of a Geiger-Müller (end-window) counter tube

Ionization of air through radioactivity (P6.4.1.1)

Cat. No.	Description	P6.4.1.1	P6.4.1.3	P6.4.1.4
559 821	Am-241 preparation, 74 kBq	1		
546 311	Zinc and grid electrodes	1		
532 14	Electrometer amplifier	1		
532 16	Connecting rod	1		1
577 03	Resistor, 10 GΩ, STE 2/19	1		
531 120	Multimeter LDanalog 20	1		
522 27	Power supply, 450 V	1		
500 412	Connecting lead, 19 A, 25 cm, blue	1		
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	1	
501 451	Connecting leads, 19 A, 50 cm, black, pair	1		
546 282	Geiger counter with adapter		1	
559 435	Ra 226 preparation, 5 kBq		1	1
521 70	High-voltage power supply, 10 kV		1	
575 302	Oszilloscope 30 MHz, digital, PT1265		1	
575 24	Screened cable, BNC/4 mm		1	
666 555	Universal clamp, 0...80 mm		1	
301 01	Leybold multiclamp		1	
300 41	Stand rod, 25 cm, 12 mm Ø		1	
300 11	Saddle base		1	2
500 610	Safety connecting lead, 25 cm, yellow/green		1	
559 01	End-window counter with cable for α, β, γ and X-rays			1
575 48	Digital counter			1
590 13	Stand rod with bore holes			1
591 21	Clip plug, large			1

In 1895, *H. Becquerel* discovered radioactivity while investigating uranium salts. He found that these emitted a radiation which was capable of fogging light-sensitive photographic plates even through black paper. He also discovered that this radiation ionizes air and that it can be identified by this ionizing effect.

In the experiment P6.4.1.1, a voltage is applied between two electrodes, and the air between the two electrodes is ionized by radioactivity. The ions created in this way cause a charge transport which can be detected using an electrometer as a highly sensitive ammeter.

The experiment P6.4.1.3 uses a Geiger counter to detect radioactivity. A potential is applied between a cover with hole which serves as the cathode and a fine needle as the anode; this potential is just below the threshold of the disruptive field strength of the air. As a result, each ionizing particle which travels within this field initiates a discharge collision.

The experiment P6.4.1.4 records the current-voltage characteristic of a Geiger-Müller counter tube. Here too, the current increases proportionally to the voltage for low voltage values, before reaching a saturation value which depends on the intensity or distance of the preparation.

P6.4.2

POISSON DISTRIBUTION

P6.4.2.1

Statistical variations in determining counting rates



Statistical variations in determining counting rates (P6.4.2.1)

Cat. No.	Description	P6.4.2.1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
524 0331	Geiger-Müller counter tube S	1
559 835	Radioactive preparations, set of 3	1
591 21	Clip plug, large	1
590 02ET2	Clip plugs, small, set of 2	1
532 16	Connecting rod	2
300 11	Saddle base	1*
587 07	Tweeter	1*
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

\* additionally recommended

For each individual particle in a radioactive preparation, it is a matter of coincidence whether it will decay over a given time period  $\Delta t$ . The probability that any particular particle will decay in this time period is extremely low. The number of particles  $n$  which will decay over time  $\Delta t$  thus shows a Poisson distribution around the mean value  $\mu$ . In other words, the probability that  $n$  decays will occur over a given time period  $\Delta t$  is


$$W_{\mu}(n) = \frac{\mu^n}{n!} e^{-\mu}$$

$\mu$  is proportional to the size of the preparation and the time  $\Delta t$ , and inversely proportional to the half-life  $T_{1/2}$  of the radioactive decay.

Using a computer-assisted measuring system, the experiment P6.4.2.1 determines multiple pulse counts  $n$  triggered in a Geiger-Müller counter tube by radioactive radiation over a selectable gate time  $\Delta t$ . After a total of  $N$  counting runs, the frequencies  $h(n)$  are determined at which precisely  $n$  pulses were counted, and displayed as histograms. For comparison, the evaluation program calculates the mean value  $\mu$  and the standard deviation

$$\sigma = \sqrt{\mu}$$

of the measured intensity distribution  $h(n)$  as well as the Poisson distribution  $w_{\mu}(M)$ .



Measured and calculated Poisson distribution Histogram:  $h(n)$ , curve:  $N \cdot w_{\mu}(n)$



**P6.4.3**  
**RADIOACTIVE DECAY AND HALF-LIFE**

P6.4.3.3  
Determining the half-life of Cs-137 - Point-by-point recording of a decay curve

P6.4.3.4  
Determining the half-life of Cs-137 - Recording and evaluating the decay curve with CASSY

Determining the half-life of Cs-137 - Recording and evaluating the decay curve with CASSY (P6.4.3.4)

Cat. No.	Description	P6.4.3.3 (b)	P6.4.3.4
559 8150Z	Cs/Ba 137 m isotope generator	1	1
524 0331	Geiger-Müller counter tube S	1	1
524 005	Mobile-CASSY 2	1	
300 02	Stand base, V-shaped, small	1	1
300 42	Stand rod, 47 cm, 12 mm diam.	1	1
301 01	Leybold multiclamp	2	2
666 555	Universal clamp, 0...80 mm	2	2
664 043	Test tubes, Fiolax, 16 x 160 mm, set of 10	1	1
664 103	Beaker, DURAN, 250 ml, squat	1	1
524 013	Sensor-CASSY 2		1
524 220	CASSY Lab 2		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1

For the activity of a radioactive sample, we can say:

$$A(t) = \left| \frac{dN}{dt} \right|$$

Here,  $N$  is the number of radioactive nuclei at time  $t$ . It is not possible to predict when an individual atomic nucleus will decay. However, from the fact that all nuclei decay with the same probability, it follows that over the time interval  $dt$ , the number of radioactive nuclei will decrease by

$$dN = -\lambda \cdot N \cdot dt$$

$\lambda$ : decay constant

Thus, for the number  $N$ , the law of radioactive decay applies:

$$N(t) = N_0 \cdot e^{-\lambda \cdot t}$$


$N_0$ : number of radioactive nuclei at time  $t = 0$

Among other things, this law states that after the half-life

$$t_{1/2} = \frac{\ln 2}{\lambda}$$

the number of radioactive nuclei will be reduced by half.

To determine the half-life of Ba-137m in the experiments P6.4.3.3 and P6.4.3.4, a plastic bottle with Cs-137 stored at salt is used. The metastable isotope Ba-137m arising from the  $\beta$ -decay is released by an elution solution. The half-time amounts to 2.6 minutes approx.



Decaying radioactivity after elution of Ba-137m

### P6.4.4

#### ATTENUATION OF $\alpha$ -, $\beta$ - AND $\gamma$ RADIATION

##### P6.4.4.2

Attenuation of  $\beta$  radiation  
when passing through matter

##### P6.4.4.3

Confirming the inverse-square  
law of distance for  $\beta$  radiation

##### P6.4.4.4

Absorption of  $\gamma$  radiation  
through matter



Attenuation of  $\beta$  radiation when passing through matter (P6.4.4.2)

Cat. No.	Description	P6.4.4.2	P6.4.4.3	P6.4.4.4
559 835	Radioactive preparations, set of 3	1	1	1
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1	1	
575 471	Counter S	1	1	
559 18	Holder with absorber foils	1		
590 02ET2	Clip plugs, small, set of 2	1	1	1
591 21	Clip plug, large	1	1	
532 16	Connecting rod	2	2	1
300 11	Saddle base	2	2	
460 97	Metal rule, 0.5 m		1	
667 9183	Geiger counter with ticker			1
559 94	Absorbers and targets, set			1
666 555	Universal clamp, 0...80 mm			1
666 572	Stand ring with stem, 70 mm diam.			1
300 02	Stand base, V-shaped, small			1
300 42	Stand rod, 47 cm, 12 mm diam.			1
301 01	Leybold multiclamp			3
559 855	Co-60 preparation, 74 kBq			1*

\* additionally recommended

High-energy  $\alpha$  and  $\beta$  particles release only a part of their energy when they collide with an absorber atom. For this reason, many collisions are required to brake a particle completely. Its range  $R$

$$R \propto \frac{E_0^2}{n \cdot Z}$$

depends on the initial energy  $E_0$ , the number density  $n$  and the atomic number  $Z$  of the absorber atoms. Low-energy  $\alpha$  and  $\beta$  particles or  $\gamma$  radiation are braked to a certain fraction when passing through a specific absorber density  $dx$ , or are absorbed or scattered and thus disappear from the beam. As a result, the radiation intensity  $I$  decreases exponentially with the absorption distance  $x$

$$I = I_0 \cdot e^{-\mu \cdot x} \quad \mu: \text{attenuation coefficient}$$

The experiment P6.4.4.2 examines the attenuation of  $\beta$  radiation from Sr-90 in aluminum as a function of the absorber thickness  $d$ . This experiment shows an exponential decrease in the intensity.

As a comparison, the absorber is removed in the experiment P6.4.4.3 and the distance between the  $\beta$  preparation and the counter tube is varied. As one might expect for a point-shaped radiation source, the following is a good approximation for the intensity:

$$I(d) \propto \frac{1}{d^2}$$

The experiment P6.4.4.4 examines the attenuation of  $\gamma$  radiation in matter. Here too, the decrease in intensity is a close approximation of an exponential function. The attenuation coefficient  $\mu$  depends on the absorber material and the  $\gamma$  energy.





**P6.5.1**  
**DEMONSTRATING PATHS OF PARTICLES**


P6.5.1.1  
 Demonstrating the tracks of  $\alpha$  particles in a Wilson cloud chamber

*Demonstrating the tracks of  $\alpha$  particles in a Wilson cloud chamber (P6.5.1.1)*

Cat. No.	Description	P6.5.1.1
559 57	Wilson cloud chamber	1
559 595	Ra-226 radium source for Wilson chamber, 5 kBq	1
450 60	Lamp housing with cable	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
460 20	Condenser with diaphragm holder	1
522 27	Power supply, 450 V	1
521 210	Transformer, 6/12 V	1
301 06	Bench clamp	1
300 11	Saddle base	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1
671 9720	Ethanol, denaturated, 1 l	1

In a Wilson cloud chamber, a saturated mixture of air, water and alcohol vapor is briefly caused to assume a supersaturated state due to adiabatic expansion. The supersaturated vapor condenses rapidly around condensation seeds to form tiny mist droplets. Ions, which are formed e.g. through collisions of  $\alpha$  particles and gas molecules in the cloud chamber, make particularly efficient condensation seeds.

In the experiment P6.5.1.1, the tracks of  $\alpha$  particles are observed in a Wilson cloud chamber. Each time the pump is vigorously pressed, these tracks are visible as traces of droplets in oblique light for one to two seconds. An electric field in the chamber clears the space of residual ions.



*Droplet traces in the Wilson cloud chamber*

P6.5.2

RUTHERFORD SCATTERING

P6.5.2.1

Rutherford scattering: measuring the scattering rate as a function of the scattering angle and the atomic number




Rutherford scattering: measuring the scattering rate as a function of the scattering angle and the atomic number (P6.5.2.1)

Cat. No.	Description	P6.5.2.1
559 820Z	Am-241 preparation, 330 kBq	1
559 56	Rutherford scattering chamber	1
559 52	Aluminium foil in frame	1
559 931	Discriminator preamplifier	1
562 791	Plug-in power supply, 12 V AC	1
575 471	Counter S	1
378 73	Rotary-vane vacuum pump S 1.5	1
378 005	T-piece, DN 16 KF	1
378 040ET2	Centring rings (adapters), DN 10/16 KF, set of 2	1
378 045ET2	Centring rings, DN 16 KF, set of 2	1
378 050	Clamping ring, DN 10/16 KF	2
378 771	Air inlet valve, DN 10 KF	1
378 031	Hose nozzle, DN 16 KF	1
667 186	Vacuum rubber tubing, 8 mm diam.	1
501 01	BNC cable, 0.25 m	1
575 24	Screened cable, BNC/4 mm	1

The fact that an atom is "mostly empty space" was confirmed by *Rutherford, Geiger and Marsden* in one of the most significant experiments in the history of physics. They caused a parallel beam of  $\alpha$  particles to fall on an extremely thin sheet of gold leaf. They discovered that most of the  $\alpha$  particles passed through the gold leaf virtually without deflection, and that only a few were deflected to a greater degree. From this, they concluded that atoms consist of a virtually massless extended shell, and a practically point-shaped massive nucleus.

The experiment P6.5.2.1 reproduces these observations using an Am-241 preparation in a vacuum chamber. The scattering rate  $N(\vartheta)$  is measured as a function of the scattering angle  $\vartheta$  using a Geiger-Müller counter tube. As scattering materials, a sheet of gold leaf ( $Z = 80$ ) and aluminum foil ( $Z = 13$ ) are provided. The scattering rate confirms the relationship

$$N(\vartheta) \propto \frac{1}{\sin^4 \frac{\vartheta}{2}} \text{ and } N(\vartheta) \propto Z^2$$



Scattering rate  $N$  as a function of the scattering angle  $\vartheta$



**P6.5.3**  
**NUCLEAR MAGNETIC RESONANCE**

P6.5.3.1  
Nuclear magnetic resonance in polystyrene, glycerin and Teflon

Nuclear magnetic resonance in polystyrene, glycerin and Teflon (P6.5.3.1\_a)

Cat. No.	Description	P6.5.3.1 (a)
514 602	NMR supply unit	1
514 606	NMR probe	1
562 11	U-core with yoke	1
562 131	Coil, 480 turns, 10 A	2
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1
501 02	BNC cable, 1 m	2
500 622	Safety connecting lead, 50 cm, blue	1
500 641	Safety connecting lead, 100 cm, red	1
500 642	Safety connecting lead, 100 cm, blue	1
531 835	Universal measuring instrument, Physics	1*
524 0381	Combi B sensor S	1*
501 11	Extension cable, 15 pin	1*

\* additionally recommended




Diagram of resonance condition of hydrogen

The magnetic moment of the nucleus entailed by the nuclear spin  $I$  assumes the energy states

$$E_m = -g_l \cdot \mu_K \cdot m \cdot B \quad \text{with } m = -I, -I+1, \dots, I$$

$$\mu_K = 5.051 \cdot 10^{-27} \frac{\text{J}}{\text{T}} : \text{ nuclear magneton}$$

$g_l$ : g factor of nucleus

in a magnetic field  $B$ . When a high-frequency magnetic field with the frequency  $\nu$  is applied perpendicularly to the first magnetic field, it excites transitions between the adjacent energy states when these fulfill the resonance condition

$$h \cdot \nu = E_{m+1} - E_m$$

$h$ : Planck's constant

This fact is the basis for nuclear magnetic resonance, in which the resonance signal is detected using radio-frequency technology. For example, in a hydrogen nucleus the resonance frequency in a magnetic field of 1 T is about 42.5 MHz. The precise value depends on the chemical environment of the hydrogen atom, as in addition to the external magnetic field  $B$  the local internal field generated by atoms and nuclei in the near vicinity also acts on the hydrogen nucleus. The width of the resonance signal also depends on the structure of the substance under study.

The experiment P6.5.3.1 verifies nuclear magnetic resonance in polystyrene, glycerine and Teflon. The evaluation focuses on the position, width and intensity of the resonance lines.

Additionally, the relaxation time of the spin system can be observed by a beat frequency measurement.

P6.5.4

α SPECTROSCOPY


P6.5.4.1  
α spectroscopy of radioactive samples

P6.5.4.2  
Determining the energy loss of  
α radiation in air

P6.5.4.3  
Determining the energy loss of  
α radiation in aluminum and in gold

P6.5.4.4  
Determining age using a  
Ra-226 sample

P6.5.4.5  
Fine structure of Am-241 α radiation



α spectroscopy of radioactive samples (P6.5.4.1)

Cat. No.	Description	P6.5.4.1	P6.5.4.2	P6.5.4.3	P6.5.4.4	P6.5.4.5
559 565	Alpha spectroscopy chamber	1	1	1	1	1
559 921	Semiconductor detector	1	1	1	1	1
559 825	Am-241 preparation, open, 3.7 kBq	1	1	1		1
559 435	Ra 226 preparation, 5 kBq	1		1	1	
524 013	Sensor-CASSY 2	1	1	1	1	1
524 058	MCA box	1	1	1	1	1
524 220	CASSY Lab 2	1	1	1	1	1
559 931	Discriminator preamplifier	1	1	1	1	1
501 16	Multi-core cable, 6-pole, 1.5 m	1	1	1	1	1
501 02	BNC cable, 1 m	1	1	1	1	1
501 01	BNC cable, 0.25 m	1	1	1	1	1
378 73	Rotary-vane vacuum pump S 1.5	1	1	1	1	1
378 005	T-piece, DN 16 KF	1		1	1	1
378 040ET2	Centring rings (adapters), DN 10/16 KF, set of 2	1		1	1	1
378 771	Air inlet valve, DN 10 KF	1		1	1	1
378 045ET2	Centring rings, DN 16 KF, set of 2	1	2	1	1	1
378 050	Clamping ring, DN 10/16 KF	2	3	2	2	2
378 031	Hose nozzle, DN 16 KF	1	1	1	1	1
667 186	Vacuum rubber tubing, 8 mm diam.	1	1	1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1*				1*
378 015	Cross piece, DN 16 KF		1			
378 776	Variable leak valve, DN 16 KF		1			
378 510	Pointer manometer, DN 16 KF		1			
311 77	Steel tape measure, 2 m		1			
559 521	Gold and aluminium foil in holder			1		
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1	1

\* additionally recommended

Up until about 1930, the energy of α rays was characterized in terms of their range in air. For example, a particle of 5.3 MeV (Po-210) has a range of 3.84 cm. Today, α energy spectra can be studied more precisely using semiconductor detectors. These detect discrete lines which correspond to the discrete excitation levels of the emitting nuclei.

The aim of the experiment P6.5.4.1 is to record and compare the α energy spectra of the two standard preparations Am-241 and Ra-226. To improve the measuring accuracy, the measurement is conducted in a vacuum chamber.

In the experiment P6.5.4.2, the energy  $E$  of α particles is measured as a function of the air pressure  $p$  in the vacuum chamber. The measurement data is used to determine the energy per unit of distance  $dE/dx$  which the α particles lose in the air. Here,

$$x = \frac{p}{p_0} \cdot x_0$$

$x_0$ : actual distance

$p_0$ : standard pressure

is the apparent distance between the preparation and the detector.

The experiment P6.5.4.3 determines the amount of energy of α particles lost per unit of distance in gold and aluminum as the quotient of the change in the energy  $\Delta E$  and the thickness  $\Delta x$  of the metal foils.

In the experiment P6.5.4.4, the individual values of the decay chain of Ra-226 leading to the α energy spectrum are analyzed to determine the age of the Ra-226 preparation used here. The activities  $A_1$  and  $A_2$  of the decay chain "preceding" and "following" the longer-life isotope Pb-210 are used to determine the age of the sample from the relationship

$$A_2 = A_1 \cdot \left(1 - e^{-\frac{t}{\tau}}\right)$$

$\tau = 32.2$  a: lifetime of Pb-210

The aim of the experiment P6.5.4.5 is to record the fine structure in the α spectrum of Am-241. An alpha decay can end in several excited states of the daughter nucleus, revealing nuclear excitation levels of the nucleus. Experimentally, this can be recorded using an unsealed radioactive source.



**P6.5.5**  
**γ SPECTROSCOPY**

P6.5.5.1  
Detecting γ radiation with a scintillation counter

P6.5.5.2  
Recording and calibrating a γ spectrum

P6.5.5.3  
Absorption of γ radiation

P6.5.5.4  
Identifying and determining the activity of radioactive samples

P6.5.5.5  
Recording a β spectrum with a scintillation counter

P6.5.5.6  
Coincidence and γ-γ angular correlation in positron decay

P6.5.5.7  
Coincidence at γ decay of cobalt

Absorption of γ radiation (P6.5.5.3)

Cat. No.	Description	P6.5.5.1	P6.5.5.2	P6.5.5.3	P6.5.5.4	P6.5.5.5	P6.5.5.6	P6.5.5.7
559 845	Mixed nuclide preparation, α, β, γ	1					1	1
559 901	Scintillation counter	1	1	1	1	1	2	2
559 891	Socket for scintillator shielding	1	1	1	1	1	1	1
559 912	Detector output stage	1	1	1	1	1	2	2
521 68	High-voltage power supply, 1.5 kV	1	1	1	1	1	2	2
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>
524 058	MCA box	1	1	1	1	1	2	2
524 220	CASSY Lab 2	1	1	1	1	1	1	1
300 42	Stand rod, 47 cm, 12 mm diam.	1	1	1		1	1	1
301 01	Leybold multiclamp	1	1	1		1	1	1
666 555	Universal clamp, 0...80 mm	1	1	1		1	1	1
575 214	Oscilloscope 30 MHz, two-channel, analogous	1*						
501 02	BNC cable, 1 m	1*						
559 835	Radioactive preparations, set of 3		1	1		1		
559 855	Co-60 preparation, 74 kBq		1*	1*				1
559 94	Absorbers and targets, set			1		1		
559 89	Scintillator shielding				1	1		
559 88	Marinelli beaker				2			
559 885	Calibrating preparation CS-137, 5 KBq				1			
672 5210	Potassium chloride, 250 g				4			
559 865	Na-22 preparation, 74 kBq						1	
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1	1	1	1	1

\* additionally recommended

γ-spectra recorded with the scintillation counter allow to identify different nuclei and give insight into fundamental aspects of nuclear physics and the interaction of radiation with matter, like Compton scattering or photoeffect.

In the experiment P6.5.5.1, the output pulses of the scintillation counter are investigated using the oscilloscope and the MCA-Box with CASSY. The total absorption peak and the Compton distribution are identified in the pulse-amplitude distribution generated with monoenergetic γ radiation.

The aim of the experiment P6.5.5.2 is to record and compare the γ energy spectra of standard preparations. The total absorption peaks are used to calibrate the energy of the scintillation counter and to identify the preparations.

The experiment P6.5.5.3 examines the attenuation of γ radiation in various absorbers. The aim here is to show how the attenuation coefficient μ depends on the absorber material and the γ energy.

A Marinelli beaker is used in the experiment P6.5.5.4 for quantitative measurements of weakly radioactive samples. This apparatus encloses the scintillator crystal virtually completely, ensuring a defined measurement geometry. Lead shielding considerably reduces the interfering background from the laboratory environment.

The experiment P6.5.5.5 records the continuous spectrum of a pure β radiator (Sr-90/Y-90) using the scintillation counter. To determine the energy loss  $dE/dx$  of the β particles in aluminum, aluminum absorbers of various thicknesses  $x$  are placed in the beam path between the preparation and the detector.

In the experiment P6.5.5.6, the spatial correlation of the two γ quanta in electron-positron pair annihilation is demonstrated. The conservation of momentum requires emission of the two quanta at an angle of 180°. Selective measurement of a coincidence spectrum leads to the suppression of non-correlated lines.


The experiment P6.5.5.7 shows the decay of Cobalt-60 in detail and proves the existence of a decay chain by coincidence measurements.

P6.5.6

COMPTON EFFECT


P6.5.6.1

Quantitative observation of the Compton effect



Quantitative observation of the Compton effect (P6.5.6.1)

Cat. No.	Description	P6.5.6.1
559 800	Equipment set for Compton scattering	1
559 809	Cs-137 preparation, 3.7 MBq	1
559 845	Mixed nuclide preparation, $\alpha$ , $\beta$ , $\gamma$	1
559 901	Scintillation counter	1
559 912	Detector output stage	1
521 68	High-voltage power supply, 1.5 kV	1
524 013	Sensor-CASSY 2	1
524 058	MCA box	1
524 220	CASSY Lab 2	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1



Measuring arrangement

In the Compton effect, a photon transfers a part of its energy  $E_0$  and its linear momentum

$$p_0 = \frac{E_0}{c}$$

$c$ : speed of light in a vacuum

to a free electron by means of elastic collision. Here, the laws of conservation of energy and momentum apply just as for the collision of two bodies in mechanics. The energy

$$E(\vartheta) = \frac{E_0}{1 + \frac{E_0}{m \cdot c^2} \cdot (1 - \cos \vartheta)}$$

$m$ : mass of electron at rest

and the linear momentum

$$p = \frac{E}{c}$$

of the scattered photon depend on the scattering angle  $\vartheta$ . The effective cross-section depends on the scattering angle and is described by the *Klein-Nishina* formula:

$$\frac{d\sigma}{d\Omega} = \frac{1}{2} \cdot r_0^2 \cdot \frac{p^2}{p_0^2} \cdot \left( \frac{p_0}{p} + \frac{p}{p_0} - \sin^2 \vartheta \right)$$

$r_0$ :  $2.5 \cdot 10^{-15}$  m: classic electron radius

In the experiment P6.5.6.1, the Compton scattering of  $\gamma$  quanta with the energy  $E_0 = 667$  keV at the quasi-free electrons of an aluminium scattering body is investigated. For each scattering angle  $\vartheta$ , a calibrated scintillation counter records one  $\gamma$  spectrum with and one without aluminum scatterer as a function of the respective scattering angle. The further evaluation utilizes the total absorption peak of the differential spectrum. The position of this peak gives us the energy  $E(\vartheta)$ . Its integral counting rate  $N(\vartheta)$  is compared with the calculated effective cross-section.

P6.5.7  
PROPERTIES OF  
RADIATION PARTICLES

P6.5.7.1  
Deflection of beta radiation  
in a magnetic field



Deflection of beta radiation in a magnetic field (P6.5.7.1)

Cat. No.	Description	P6.5.7.1
559 835	Radioactive preparations, set of 3	1
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1
575 471	Counter S	1
LDS 00001	Stopwatch, digital	1
562 11	U-core with yoke	1
562 13	Coil, 250 turns	2
560 31	Bored pole pieces, pair	1
559 23	Swivelling clamp	1
559 18	Holder with absorber foils	1
531 120	Multimeter LDanalog 20	1
521 35	Variable extra-low voltage transformer S	1
300 11	Saddle base	1
300 41	Stand rod, 25 cm, 12 mm $\emptyset$	1
301 01	Leybold multiclamp	1
501 25	Connecting lead, 32 A, 50 cm, red	2
501 26	Connecting lead, 32 A, 50 cm, blue	2
524 005	Mobile-CASSY 2	1*
524 0381	Combi B sensor S	1*
501 11	Extension cable, 15 pin	1*

\* additionally recommended

Historically, it was easy to see that there are different kinds of radioactive radiation. But to identify which particles were involved took some time. The behaviour of those particles in an magnetic field was and still is the most valuable tool to characterize charged particles.

In the experiment P6.5.7.1  $\beta$  radiation passes through a magnetic field and the resulting angular distribution is recorded. From these data, the energy of the decay can be estimated.

P6.6.1

QUANTUM OPTICS

P6.6.1.1  
Quantum eraser



Quantum eraser (P6.6.1.1)

Cat. No.	Description	P6.6.1.1
473 40	Laser optics base plate	1
471 830	He-Ne Laser, linearly polarised	1
473 411	Laser mount	1
473 421	Optics base	9
473 431	Holder for beam divider	2
473 432	Beam divider, 50 %	2
473 461	Planar mirror with fine adjustment	2
473 471	Spherical lens, $f = 2.7$ mm	2
473 49	Polarising filter for laser optics base plate	3
441 53	Screen, translucent	2
300 11	Saddle base	2
311 02	Metal rule, 1 m	1

Quantum optics is a field of research in physics, dealing with the application of quantum mechanics to phenomena involving light and its interactions with matter.

A basic principle of quantum mechanics is complementarity: each quantum-mechanical object has both wave-like and particle-like properties. In the experiment P6.6.1.1 an analogue experiment to a quantum eraser is built up. It shows the complementarity of which-way information and interference.



P6.6.2  
PARTICLES

P6.6.2.1  
Detection of Muons




Detection of Muons (P6.6.2.1)

Cat. No.	Description	P6.6.2.1
524 013	Sensor-CASSY 2	1
524 220	CASSY Lab 2	1
524 033	GM box	2
559 012	Pancake GM counter tube	2
300 11	Saddle base	1*
587 07	Tweeter	1*
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1*
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

\* additionally recommended

The origin of high-energy particles is not only from radioactive decay of unstable nuclei on earth, they can also be found in the cosmic radiation. Such particles are an interesting topic on their own.

Experiment P6.6.2.1 uses the coincidence detection of two GM detectors to identify natural muons passing through this muon telescope. The specific properties of the muons can be shown this way and recording the muon flux over the day shows their origin is related to the sun.




Count rate distribution for the muon detector oriented in different directions

P6.1.1.1  
Estimating the size  
of oil molecules



# P7 SOLID-STATE PHYSICS



P7.1	PROPERTIES OF CRYSTALS	277
P7.2	CONDUCTION PHENOMENA	281
P7.3	MAGNETISM	288
P7.4	SCANNING PROBE MICROSCOPY	290
P7.5	APPLIED SOLID-STATE PHYSICS	292

# P7 SOLID-STATE PHYSICS



## P7.1 PROPERTIES OF CRYSTALS

P7.1.1 Crystal structure	277
P7.1.2 X-ray scattering	278-279
P7.1.4 Elastic and plastic deformation	280

## P7.2 CONDUCTION PHENOMENA

P7.2.1 Hall effect	281-282
P7.2.2 Electrical conductivity in solids	283
P7.2.3 Photoconductivity	284
P7.2.4 Luminescence	285
P7.2.5 Thermoelectricity	286
P7.2.6 Superconductivity	287

## P7.3 MAGNETISM

P7.3.1 Dia-, para- and ferromagnetism	288
P7.3.2 Ferromagnetic hysteresis	289

## P7.4 SCANNING PROBE MICROSCOPY

P7.4.1 Scanning tunneling microscope	290-291
--------------------------------------	---------

## P7.5 APPLIED SOLID-STATE PHYSICS

P7.5.1 X-ray fluorescence analysis	292
------------------------------------	-----

P7.1.1

CRYSTAL STRUCTURE

P7.1.1.1

Structure of a body-centered cubic and face-centered cubic lattice



Structure of a body-centered cubic and face-centered cubic lattice (P7.1.1.1)

Cat. No.	Description	P7.1.1.1
554 60	Field emission microscope	1
554 605	Connection panel for FEM	1
301 339	Stand bases, pair	1
521 70	High-voltage power supply, 10 kV	1
521 39	Variable extra-low voltage transformer	1
531 130	Multimeter LDanalog 30	1
500 614	Safety connecting lead, 25 cm, black	2
500 624	Safety connecting lead, 50 cm, black	2
500 641	Safety connecting lead, 100 cm, red	1
500 642	Safety connecting lead, 100 cm, blue	1
500 644	Safety connecting lead, 100 cm, black	2

In the field emission microscope, the extremely fine tip of a tungsten monocrystal is arranged in the center of a spherical luminescent screen. In the vicinity of the tip, the electric field between the crystal and the luminescent screen reaches such a high field strength that the conducting electrons can “tunnel” out of the crystal and travel radially to the luminescent screen. Here, an image of the emission distribution of the crystal tip is created, magnified by a factor of

$$V = \frac{R}{r}$$

$R = 5 \text{ cm}$ : radius of luminescent screen

$r = 0.1 - 0.2 \text{ }\mu\text{m}$ : radius of tip

In the first part of the experiment P7.1.1.1, the tungsten tip is purified by heating it to a white glow. The structure which appears on the luminescent screen after the electric field is applied corresponds to the body-centered cubic lattice of tungsten, which is observed in the (110) direction, i.e. the direction of one of the diagonals of a cube face. Finally, a minute quantity of barium is vaporized in the tube, so that individual barium atoms can precipitate on the tungsten tip to produce bright spots on the luminescent screen. When the tungsten tip is heated carefully, it is even possible to observe the thermal motion of the barium atoms.




Image of tungsten tip: hot electrode

P7.1.2

X-RAY SCATTERING

P7.1.2.1

Bragg reflection: determining the lattice constants of monocrystals

P7.1.2.2

Laue diagrams: investigating the lattice structure of monocrystals

P7.1.2.3

Debye-Scherrer photography: determining the lattice plane spacings of polycrystalline powder samples

P7.1.2.4

Debye-Scherrer Scan: determining the lattice plane spacings of poly-crystalline powder samples



Laue diagrams: investigating the lattice structure of monocrystals (P7.1.2.2)

Cat. No.	Description	P7.1.2.1	P7.1.2.2	P7.1.2.3	P7.1.2.4
554 800	X-ray apparatus	1	1	1	1
554 861	X-ray tube, Mo	1	1	1	
554 831	Goniometer	1			1
559 01	End-window counter with cable for $\alpha$ , $\beta$ , $\gamma$ and X-rays	1			1
554 77	LiF crystal for Bragg reflection	1			
554 78	NaCl crystal for Bragg reflection	1			
554 838	Film holder, X-ray		1	1	
554 896	X-ray film, Agfa Dentus M2		1	1	
554 87	LiF crystal for Laue diagrams		1		
554 88	NaCl crystal for Laue diagrams		1		
554 8971	Developer and fixer for X-ray film		1	1	
554 8931	Changing bag with developer tank		1*	1*	
673 5700	Sodium chloride, 250 g			1	1
673 0520	Lithium fluoride, analytically pure, 10 g			1	1
667 091	Pestle, 96 mm long			1	1
667 092	Mortar, porcelain, 70 mm $\emptyset$			1	1
666 960	Powder spatula, stainless steel, 150 mm			1	1
311 54	Precision vernier callipers			1	
554 862	X-ray tube, Cu				1
554 842	Crystal powder holder				1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1			1

\* additionally recommended


X-rays are an essential tool to determine the structure of crystals. The lattice planes inside a crystal are identified by their Miller indices  $h, k, l$  and reflect the X-rays only if the Laue or Bragg conditions are fulfilled. The distribution of reflexes allows to calculate the lattice constant and crystal structure of the investigated crystal.

In the experiment P7.1.2.1, the Bragg reflection of Mo- $K_{\alpha}$  radiation ( $\lambda = 71.080 \text{ pm}$ ) at NaCl and LiF monocrystals is used to determine the lattice constant. The  $K_{\beta}$  component of the X-ray radiation can be suppressed using a zirconium filter.

Laue diagrams from NaCl and LiF monocrystals are created in the experiment P7.1.2.2 using the „bremsstrahlung“ of the X-ray apparatus as „white“ X-radiation. The positions of the „multicolored“ reflections on an X-ray film behind the crystal and their intensities can be used to determine the crystal structure and the lengths of the crystal axes through application of the Laue condition.

In the experiment P7.1.2.3, Debye-Scherrer photographs are produced by irradiating samples of a fine crystal powder with Mo- $K_{\alpha}$  radiation. Some of the randomly oriented crystallites conform to the Bragg condition and diffract X-rays into cones for which the aperture angles  $\vartheta$  can be derived from a photograph. This experiment determines the lattice spacing corresponding to  $\vartheta$  as well as its Laue indices  $h, k, l$ , and thus the lattice structure of the crystallite.

The experiment P7.1.2.4 records the Debye-Scherrer pattern with an end window counter instead of X-ray film. The diffracted reflections of a fine powder sample are recorded as a function of diffraction angle. The intensity peaks in the diffraction spectrum allow the calculation of the separations of adjacent lattice planes.




Laue diagram of NaCl and Debye-scherrer photograph of NaCl

P7.1.2

X-RAY SCATTERING

P7.1.2.5


Digital Laue diagrams: investigating the lattice structure of monocrystals



Digital Laue diagrams: investigating the lattice structure of monocrystals (P7.1.2.5)

Cat. No.	Description	P7.1.2.5
554 800	X-ray apparatus	1
554 866	X-ray tube, Au	1
554 8381	Pinhole collimator with Laue crystals	1
554 828	X-ray image sensor	1
554 829	Precision rails for X-ray image sensor	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

To create *Laue* diagrams of NaCl and LiF monocrystals, the bremsstrahlung of the X-ray apparatus is used in the experiment P7.1.2.5 as a „white“ spectrum of X-rays. The positions of the „multicolored“ diffraction spots on an X-ray image sensor behind the crystal and their intensities can be used to determine the crystal structure and the dimensions of the unit cell through application of the *Laue* condition. The X-ray image sensor allows capturing the *Laue* pattern in one minute (much faster than conventional film) and the digital evaluation eases the determination of the diffraction angles on the computer.



Laue diagram of NaCl and Debye-scherrer photograph of NaCl



The X-ray image sensor is sensitive enough to expose the discrete reflexes behind a crystal. The zero order is shielded by a metal plate.

P7.1.4

ELASTIC AND PLASTIC  
DEFORMATION

P7.1.4.1

Investigating the elastic and plastic extension of metal wires

P7.1.4.2

Investigating the elastic and plastic extension of metal wires - Recording and evaluating with CASSY



Investigating the elastic and plastic extension of metal wires (P7.1.4.1)

Cat. No.	Description	P7.1.4.1	P7.1.4.2 (a)	P7.1.4.2 (b)
550 35	Copper resistance wire, 0.2 mm diam., 100 m	1	1	1
550 51	Iron resistance wire, 0.2 mm diameter, 100 m	1	1	1
342 61	Weights, 50 g, set of 12	2		
340 911ET2	Pulley 50 mm Ø, plug-in, pair, set of 2	1		
381 331	Pointer for linear expansion	1		
340 82	Double scale	1		
314 04ET5	Support clip, for plugging in, set of 5	1		
301 07	Simple bench clamp	2	2	2
301 01	Leybold multiclamp	4	3	3
301 25	Support block	3		
301 26	Stand rod, 25 cm, 10 mm diam.	3	2	2
301 27	Stand rod, 50 cm, 10 mm diam.	1		
300 44	Stand rod, 100 cm, 12 mm diam.	1	1	1
524 013	Sensor-CASSY 2		1	
524 220	CASSY Lab 2		1	
524 042	Force sensor S, ±50 N		1	1
524 082	Rotary motion sensor S		1	1
311 77	Steel tape measure, 2 m		1	1
524 005	Mobile-CASSY 2			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)		1	

The shape of a crystalline solid is altered when a force is applied. We speak of elastic behavior when the solid resumes its original form once the force ceases to act on it. When the force exceeds the elastic limit, the body is permanently deformed. This plastic behavior is caused by the migration of discontinuities in the crystal structure.

In the experiments P7.1.4.1 and P7.1.4.2, the extension of iron and copper wires is investigated by hanging weights from them. A precision pointer indicator or the rotary motion sensor S attached to a CASSY measures the change in length  $\Delta s$ , i. e. the extension

$$\epsilon = \frac{\Delta s}{s}$$

s: length of wire

After each new tensile load

$$\sigma = \frac{F}{A}$$


F: weight of load pieces  
A: wire cross-section

the students observe whether the pointer or the rotary motion sensor returns to the zero position when the strain is relieved, i.e. whether the strain is below the elasticity limit  $\sigma_e$ . Graphing the measured values in a tension-extension diagram confirms the validity of Hooke's law

$$\sigma = E \cdot \epsilon$$

E: modulus of elasticity

up to a proportionality limit  $\sigma_p$ .



Load-extension diagram for a typical metal wire





## P7.2.1

### HALL EFFECT

#### P7.2.1.1

Investigating the Hall effect in silver

#### P7.2.1.2

Investigating the anomalous Hall effect in tungsten

Investigating the Hall effect in silver (P7.2.1.1\_b)

Cat. No.	Description	P7.2.1.1 (b)	P7.2.1.2 (b)
586 81	Hall effect apparatus (silver)	1	
<b>524 005</b>	<b>Mobile-CASSY 2</b>	<b>1</b>	<b>1</b>
524 0381	Combi B sensor S	1	1
501 11	Extension cable, 15 pin	1	1
532 13	Microvoltmeter	1	1
531 130	Multimeter LDanalog 30	1	1
521 55	High current power supply	1	1
521 39	Variable extra-low voltage transformer	1	1
562 11	U-core with yoke	1	1
560 31	Bored pole pieces, pair	1	1
562 13	Coil, 250 turns	2	2
300 41	Stand rod, 25 cm, 12 mm Ø	1	1
301 01	Leybold multiclamp	1	1
300 02	Stand base, V-shaped, small	1	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	4	4
501 33	Connecting lead, 32 A, 100 cm, black	2	2
586 84	Hall effect apparatus (tungsten)		1

In the case of electrical conductors or semiconductors within a magnetic field  $B$ , through which a current  $I$  is flowing perpendicular to the magnetic field, the Hall effect results in an electric potential difference

$$U_H = R_H \cdot B \cdot I \cdot \frac{1}{d} \quad d: \text{thickness of sample}$$

The Hall coefficient

$$R_H = \frac{1}{e} \cdot \frac{p \cdot \mu_p^2 - n \cdot \mu_n^2}{(p \cdot \mu_p + n \cdot \mu_n)^2} \quad e: \text{elementary charge}$$

depends on the concentrations  $n$  and  $p$  of the electrons and holes as well as their mobilities  $\mu_n$  and  $\mu_p$ , and is thus a quantity which depends on the material and the temperature

The experiments P7.2.1.1 and P7.2.1.2 determine the Hall coefficient  $R_H$  of two electrical conductors by measuring the Hall voltage  $U_H$  for various currents  $I$  as a function of the magnetic field  $B$ . A negative value is obtained for the Hall coefficient of silver, which indicates that the charge is being transported by electrons. A positive value is found as the Hall coefficient of tungsten. Consequently, the holes are mainly responsible for conduction in this metal.

## P7.2.1

### HALL EFFECT

#### P7.2.1.3

Determining the density and mobility of charge carriers in n-Germanium

#### P7.2.1.4

Determining the density and mobility of charge carriers in p-Germanium

#### P7.2.1.5

Determining the band gap of germanium



Determining the density and mobility of charge carriers in p-Germanium (P7.2.1.4)


Cat. No.	Description	P7.2.1.3	P7.2.1.4	P7.2.1.5
586 850	Base unit for Hall effect	1	1	1
586 853	N-doped germanium on plug-in board	1		
562 11	U-core with yoke	1	1	
562 13	Coil, 250 turns	2	2	
560 31	Bored pole pieces, pair	1	1	
521 536	DC Power Supply 2 x 0...16 V/2 x 0...5 A	1	1	
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	<b>1</b>	<b>1</b>
524 220	CASSY Lab 2	1	1	1
524 0381	Combi B sensor S	1	1	
501 11	Extension cable, 15 pin	1	1	
300 02	Stand base, V-shaped, small	1	1	1
300 41	Stand rod, 25 cm, 12 mm Ø	1	1	
301 01	Leybold multiclamp	1	1	
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	7	7	4
586 852	P-doped germanium on plug-in board		1	
586 851	Undoped germanium on plug-in board			1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A			1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1	1

The experiments P7.2.1.3 and P7.2.1.4 explore the temperature-dependency of the Hall voltage and the electrical conductivity


$$\sigma = e \cdot (p \cdot \mu_p + n \cdot \mu_n)$$

using doped germanium samples. The concentrations of the charge carriers and their mobilities are determined under the assumption that, depending on the doping, one of the concentrations  $n$  or  $p$  can be ignored.

In the experiment P7.2.1.5, the electrical conductivity of undoped germanium is measured as a function of the temperature to provide a comparison. The measurement data permits determination of the band gap between the valence band and the conduction band in germanium.



Hall voltage when heating up the p-Ge sample



## P7.2.2

### ELECTRICAL CONDUCTIVITY IN SOLIDS

#### P7.2.2.1

Measuring the temperature-dependency of a noble-metal resistor

#### P7.2.2.2

Measuring the temperature-dependency of a semiconductor resistor

Measuring the temperature-dependency of a noble-metal resistor (P7.2.2.1)

Cat. No.	Description	P7.2.2.1	P7.2.2.2
586 80	Noble metal resistor	1	
555 81	Electric oven for tubes, 230 V	1	1
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
524 0673	NiCr-Ni adapter S, type K	1	1
529 676	Temperature probe, NiCr-Ni, 1.5 mm, type K	1	1
524 031	Current source box	1	1
502 061	Safety connecting box, with earth	1	1
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	1	1
586 821	Semiconductor resistor, 5 kΩ		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

The temperature-dependency of the specific resistance  $\rho$  is a simple test for models of electric conductivity of conductors and semiconductors. In electrical conductors,  $\rho$  increases with the temperature, as the collisions of the quasi-free electrons from the conduction band with the atoms of the conductor play an increasingly important role. In semiconductors, on the other hand, the specific resistance decreases as the temperature increases, as more and more electrons move from the valence band to the conduction band, thus contributing to the conductivity.

The experiments P7.2.2.1 and P7.2.2.2 measure the resistance values as a function of temperature using a Wheatstone bridge. The computer-assisted CASSY measured-value recording system is ideal for recording and evaluating the measurements. For the noble metal resistor, the relationship

$$R = R_0 \cdot \frac{T}{\Theta}$$


$$\Theta = 240 \text{ K: Debye temperature of platinum}$$

is verified with sufficient accuracy in the temperature range under study. For the semiconductor resistor, the evaluation reveals a dependency with the form

$$R \propto e^{\frac{\Delta E}{2kT}}$$

$$k = 1.38 \cdot 10^{-23} \frac{\text{J}}{\text{K}} : \text{Boltzmann constant}$$

with the band spacing  $E = 0.5 \text{ eV}$ .




Resistance against temperature

### P7.2.3

#### PHOTOCONDUCTIVITY

##### P7.2.3.1

Recording the current-voltage characteristics of a CdS photoresistor




Recording the current-voltage characteristics of a CdS photoresistor (P7.2.3.1)

Cat. No.	Description	P7.2.3.1
578 02	Photoresistor LDR 05, STE 2/19	1
450 511	Bulbs, 6 V/30 W, E14, set of 2	1
450 60	Lamp housing with cable	1
460 20	Condenser with diaphragm holder	1
460 14	Adjustable slit	1
472 401	Polarisation filter	2
460 08	Lens in frame, f=150 mm	1
460 32	Optical bench with standardised profile, 1 m	1
460 374	Optics rider, 90/50	6
460 21	Holder for plug-in elements	1
521 546	DC Power Supply 0 ... 16 V, 0 ... 5 A	1
521 210	Transformer, 6/12 V	1
531 282	Multimeter Metrahit Pro	1
531 303	Multimeter Metrahit X-tra	1
500 422	Connecting lead 19 A, 50 cm, blue	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	2

Photoconductivity is the phenomenon in which the electrical conductivity  $\sigma$  of a solid is increased through the absorption of light. In CdS, for example, the absorbed energy enables the transition of activator electrons to the conduction band and the reversal of the charges of traps, with the formation of electron holes in the valence band. When a voltage  $U$  is applied, a photocurrent  $I_{ph}$  flows.

The object of the experiment P7.2.3.1 is to determine the relationship between the photocurrent  $I_{ph}$  and the voltage  $U$  at a constant radiant flux  $\Phi_e$  as well as between the photocurrent  $I_{ph}$  and the radiant flux  $\Phi_e$  at a constant voltage  $U$  in the CdS photoresistor.



## P7.2.4

### LUMINESCENCE

#### P7.2.4.1

Exciting luminescence through irradiation with ultraviolet light and electrons

Exciting luminescence through irradiation with ultraviolet light and electrons (P7.2.4.1)

Cat. No.	Description	P7.2.4.1
555 618	Luminescence tube	1
555 600	Tube stand	1
521 70	High-voltage power supply, 10 kV	1
451 15	High pressure mercury lamp	1
451 195	Power supply unit for high-pressure mercury lamp	1
469 79	Filter, ultra-violet	1
500 611	Safety connecting lead, 25 cm, red	1
500 621	Safety connecting lead, 50 cm, red	1
500 641	Safety connecting lead, 100 cm, red	1
500 642	Safety connecting lead, 100 cm, blue	1
500 644	Safety connecting lead, 100 cm, black	2

Luminescence is the emission of light following the absorption of energy. This energy can be transmitted in the form of e.g. high-energy electrons or photons which have an energy greater than that of the emitted photons. Depending on the type of decay, we distinguish between fluorescence and phosphorescence. In fluorescence, the emission of photons fades exponentially very rapidly when excitation is switched off (i.e. about  $10^{-8}$  s). Phosphorescence, on the other hand, can persist for several hours.

In the experiment P7.2.4.1, the luminescence of various solids following irradiation with ultraviolet light or electrons is demonstrated. These samples include yttrium vanadate doped with europium (red fluorescent), zinc silicate doped with manganese (green fluorescent) and barium magnesium aluminate doped with europium (blue fluorescent).

*Note:* It is possible to recognize individual emission lines within the band spectrum using a pocket spectroscope.



## P7.2.5

### THERMOELECTRICITY

#### P7.2.5.1

Seebeck effect: Determining the thermoelectric voltage as a function of the temperature differential



Seebeck effect: Determining the thermoelectric voltage as a function of the temperature differential (P7.2.5.1\_a)

Cat. No.	Description	P7.2.5.1 (a)	P7.2.5.1 (b)
557 01	Thermocouples, set	1	1
590 011	Clamping plug	2	2
532 13	Microvoltmeter	1	
382 34	Thermometer, -10...+110 °C/0.2 K	1	1
666 767	Hotplate, 1500 W, 180 mm diam.	1	1
664 104	Beaker, DURAN, 400 ml, squat	1	1
524 005	Mobile-CASSY 2		1
524 040	µV box		1

When two metal wires with different Fermi energies  $E_F$  touch, electrons move from one to the other. The metal with the lower electronic work function  $W_A$  emits electrons and becomes positive. The transfer does not stop until the contact voltage

$$U = \frac{W_{A,1} - W_{A,2}}{e}$$

$e$ : elementary charge

is reached. If the wires are brought together in such a way that they touch at both ends, and if the two contact points have a temperature differential  $T = T_1 - T_2$ , an electrical potential, the thermoelectric voltage

$$U_T = U(T_1) - U(T_2)$$

is generated. Here, the differential thermoelectric voltage


$$\alpha = \frac{dU_T}{dT}$$

depends on the combination of the two metals.

In the experiment P7.2.5.1, the thermoelectric voltage  $U_T$  is measured as a function of the temperature differential  $T$  between the two contact points for thermocouples with the combinations iron/constantan, copper/constantan and chrome-nickel/constantan. One contact point is continuously maintained at room temperature, while the other is heated in a water bath. The differential thermoelectric voltage is determined by applying a best-fit straight line

$$U_T = \alpha \cdot T$$

to the measured values.



Thermoelectric voltage as a function of the temperature. Top: chrome-nickel/constantan, Middle: iron/constantan, Bottom: copper/constantan



## P7.2.6 SUPERCONDUCTIVITY

P7.2.6.1  
Determining the transition temperature of a high-temperature superconductor

P7.2.6.2  
Meissner-Ochsenfeld effect for a high-temperature superconductor


Determining the transition temperature of a high-temperature superconductor (P7.2.6.1)

Cat. No.	Description	P7.2.6.1	P7.2.6.2
667 552	Transition temperature and electrical resistance, experiment kit	1	
<b>524 013</b>	<b>Sensor-CASSY 2</b>	<b>1</b>	
524 220	CASSY Lab 2	1	
501 45	Connecting lead, 19 A, 50 cm, red/blue, pair	2	
667 551	Meissner-Ochsenfeld effect, experiment kit		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	

In 1986, *K. A. Müller* and *J. G. Bednorz* succeeded in demonstrating that the compound  $\text{YBa}_2\text{Cu}_3\text{O}_7$  becomes superconducting at temperatures far greater than any known up to that time. Since then, many high-temperature superconductors have been found which can be cooled to their transition temperature using liquid nitrogen. Like all superconductors, high-temperature superconductors have no electrical resistance and demonstrate the phenomenon known as the Meissner-Ochsenfeld effect, in which magnetic fields are displaced out of the superconducting body.

The experiment P7.2.6.1 determines the transition temperature of the high-temperature superconductor  $\text{YBa}_2\text{Cu}_3\text{O}_{7-x}$ . For this purpose, the substance is cooled to below its critical temperature of  $T_c = 92 \text{ K}$  using liquid nitrogen. In a four-point measurement setup, the voltage drop across the sample is measured as a function of the sample temperature using the computer-assisted measured value recording system CASSY.

In the experiment P7.2.6.2, the superconductivity of the  $\text{YBa}_2\text{Cu}_3\text{O}_{7-x}$  body is verified with the aid of the Meissner-Ochsenfeld effect. A low-weight, high field-strength magnet placed on top of the sample begins to hover when the sample is cooled to below its critical temperature so that it becomes superconducting and displaces the magnetic field of the permanent magnet.



Resistance against temperature

### P7.3.1

#### DIA-, PARA- AND FERROMAGNETISM

##### P7.3.1.1

Dia-, para- and ferromagnetic materials in an inhomogeneous magnetic field




*Dia-, para- and ferromagnetic materials in an inhomogeneous magnetic field (P7.3.1.1)*

Cat. No.	Description	P7.3.1.1
560 41	Set of rods for para- and diamagnetism	1
562 11	U-core with yoke	1
562 13	Coil, 250 turns	2
560 31	Bored pole pieces, pair	1
521 39	Variable extra-low voltage transformer	1
300 02	Stand base, V-shaped, small	1
300 41	Stand rod, 25 cm, 12 mm Ø	2
301 01	Leybold multiclamp	1
500 422	Connecting lead 19 A, 50 cm, blue	1
501 46	Connecting leads, 19 A, 100 cm, red/blue, pair	1

Diamagnetism is the phenomenon in which an external magnetic field causes magnetization in a substance which is opposed to the applied magnetic field in accordance with Lenz's law. Thus, in an inhomogeneous magnetic field, a force acts on diamagnetic substances in the direction of decreasing magnetic field strength. Paramagnetic materials have permanent magnetic moments which are aligned by an external magnetic field. Magnetization occurs in the direction of the external field, so that these substances are attracted in the direction of increasing magnetic field strength. Ferromagnetic substances in magnetic fields assume a very high magnetization which is orders of magnitude greater than that of paramagnetic substances.

In the experiment P7.3.1.1, three 9 mm long rods with different magnetic behaviors are suspended in a strongly inhomogeneous magnetic field so that they can easily rotate, allowing them to be attracted or repelled by the magnetic field depending on their respective magnetic property.



*Placement of a sample in the magnetic field*





**P7.3.2**  
**FERROMAGNETIC HYSTERESIS**

P7.3.2.1  
Recording the initial magnetization curve and the hysteresis curve of a ferromagnet

Recording the initial magnetization curve and the hysteresis curve of a ferromagnet (P7.3.2.1\_a)

Cat. No.	Description	P7.3.2.1	
		(a)	(b)
562 11	U-core with yoke	1	1
562 121	Clamping device with spring clip	1	1
562 14	Coil, 500 turns	2	2
522 621	Function generator S 12	1	
524 013	Sensor-CASSY 2	1	1
524 220	CASSY Lab 2	1	1
577 20	Resistor, 10 Ω, STE 2/19	1	
576 71	Plug-in board section, STE	1	
500 424	Connecting lead 19 A, 50 cm, black	1	
500 444	Connecting lead 19 A, 100 cm, black	7	4
524 011USB	Power-CASSY USB		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

In a ferromagnet, the magnetic induction

$$B = \mu_r \cdot \mu_0 \cdot H$$

$$\mu_0 = 4\pi \cdot 10^{-7} \frac{\text{Vs}}{\text{Am}} : \text{magnetic field constant}$$

reaches a saturation value  $B_s$  as the magnetic field  $H$  increases. The relative permeability  $\mu_r$  of the ferromagnet depends on the magnetic field strength  $H$ , and also on the previous magnetic treatment of the ferromagnet. Thus, it is common to represent the magnetic induction  $B$  in the form of a hysteresis curve as a function of the rising and falling field strength  $H$ . The hysteresis curve differs from the magnetization curve, which begins at the origin of the coordinate system and can only be measured for completely demagnetized material.

In the experiment P7.3.2.1, a current  $I_1$  in the primary coil of a transformer which increases (or decreases) linearly over time generates the magnetic field strength

$$H = \frac{N_1}{L} \cdot I_1$$

$L$ : effective length of iron core

$N_1$ : number of windings of primary coil


The corresponding magnetic induction value  $B$  is obtained through integration of the voltage  $U_2$  induced in the secondary coil of a transformer:

$$B = \frac{1}{N_2 \cdot A} \cdot \int U_2 \cdot dt$$

$A$ : cross-section of iron core

$N_2$ : number of windings of secondary coil

The computer-assisted measurement system CASSY is used to control the current and to record and evaluate the measured values. The aim of the experiment is to determine the relative permeability  $\mu_r$  in the magnetization curve and the hysteresis curve as a function of the magnetic field strength  $H$ .



Magnetic flux in the iron core against the coil current

P7.4.1

SCANNING TUNNELING  
MICROSCOPE

P7.4.1.1  
Investigating a graphite surface  
using a scanning tunneling microscope

P7.4.1.2  
Investigating a gold surface  
using a scanning tunneling microscope

P7.4.1.3  
Investigating a MoS<sub>2</sub> sample  
using a scanning tunneling microscope



Investigating a graphite surface using a scanning tunneling microscope (P7.4.1.1)

Cat. No.	Description	P7.4.1.1-2	P7.4.1.3
554 581	Scanning tunnel microscope	1	1
554 584	Molybdenum disulphide (MoS <sub>2</sub> ) sample		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

The scanning tunneling microscope was developed in the 1980's by G. Binnig and H. Rohrer. It uses a fine metal tip as a local probe; the probe is brought so close to an electrically conductive sample that the electrons "tunnel" from the tip to the sample due to quantum-mechanical effects. When an electric field is applied between the tip and the sample, an electric current, the tunnel current, can flow. As the tunnel current varies exponentially with the distance, even an extremely minute change in distance of 0.01 nm results in a measurable change in the tunnel current. The tip is mounted on a platform which can be moved in all three spatial dimensions by means of piezoelectric control elements. The tip is scanned across the sample to measure its topography. A control circuit maintains the distance between tip and sample extremely precisely at a constant distance by maintaining a constant tunnel current value. The controlled motions performed during the scanning process are recorded and imaged using a computer. The image generated in this manner is a composite in which the sample topography and the electrical conductivity of the sample surface are superimposed.

The experiments P7.4.1.1, P7.4.1.2 and P7.4.1.3 use a scanning tunneling microscope specially developed for practical experiments, which operates at standard air pressure. At the beginning of the experiment, a measuring tip is made from platinum wire. The graphite sample is prepared by tearing off a strip of tape. When the gold sample is handled carefully, it requires no cleaning; the same is valid for the MoS<sub>2</sub> probe. The investigation of the samples begins with an overview scan. In the subsequent procedure, the step width of the measuring tip is reduced until the positions of the individual atoms of the sample with respect to each other are clearly visible in the image.





### P7.4.1

#### SCANNING TUNNELING MICROSCOPE


##### P7.4.1.4

Investigating of surfaces  
using an atomic force  
microscope (AFM)

*Investigating of surfaces using an atomic force microscope (AFM) (P7.4.1.4)*

Cat. No.	Description	P7.4.1.4
554 5863	LD AFM extended	1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1

The experiment P7.4.1.4 uses an AFM to investigate several different surface structures from a CD stamper to blood cells. In contrast to a STM non-conductive samples can be used. A microscopic „needle“ is touching the surface and scans it. In static mode, that needle just scratches along the surface and records the height variations. Using extended modes, not only topography can be recorded, but also local properties like elasticity or magnetism. For example the needle is set into vibration and slightly touches the surface. Hard and soft areas of the sample will dampen the oscillation of the needle differently and can be imaged.



*AFM image of an uncoated CD with data structures in 2D and 3D*

## P7.5.1

### X-RAY FLUORESCENCE ANALYSIS

#### P7.5.1.1

Application of X-ray fluorescence for the non-destructive analysis of the chemical composition

#### P7.5.1.2

Determination of the chemical composition of a brass sample by X-ray fluorescence analysis




Application of X-ray fluorescence for the non-destructive analysis of the chemical composition (P7.5.1.1)

Cat. No.	Description	P7.5.1.1	P7.5.1.2
554 800	X-ray apparatus	1	1
554 861	X-ray tube, Mo	1	1
554 831	Goniometer	1	1
559 938	X-ray energy detector	1	1
554 848	Set of target alloys	1	1
524 013	Sensor-CASSY 2	1	1
524 058	MCA box	1	1
524 220	CASSY Lab 2	1	1
501 02	BNC cable, 1 m	1	1
554 844	Set of targets for K-line fluorescence		1
554 846	Set of targets for L-line fluorescence		1
	additionally required: PC with Windows XP/Vista/7/8/10 (x86 or x64)	1	1

X-ray fluorescence is a very useful tool for a non-destructive analysis of the chemical composition of a target alloy. When irradiating a sample with X-rays, all the different elements it contains emit characteristic X-rays due to fluorescence, which are fingerprints of every single element.

In the experiment P7.5.1.1, X-ray fluorescence is used to do qualitative analysis by identifying the substances in four alloy samples, made from chrome-nickel steel, two different kinds of brass and rare earth magnet.

In the experiment P7.5.1.2, the composition of one brass alloy is analysed quantitatively. The weight percentage of each component in the alloy is calculated from the strength of different fluorescence lines.



X-Ray fluorescence spectrum of a brass sample

# P8 REGISTER



$\alpha, \beta, \gamma$ 

$\alpha$ radiation .....	264
$\alpha$ spectrum .....	268
$\beta$ radiation .....	264
$\beta$ spectrum .....	269
$\gamma$ radiation .....	264
$\gamma$ spectrum .....	269

## 3

3D .....	259
----------	-----

## A

Abbe refractometer .....	209
Aberration,chromatic .....	169
Aberration, lens .....	169
Absorber, saturable .....	218
Absorption .....	206, 212, 214
Absorption edge .....	253
Absorption of $\beta$ radiation .....	269
Absorption of $\gamma$ radiation .....	264
Absorption of light .....	173
Absorption of microwaves .....	137
Absorption of X-rays .....	253
Absorption spectra .....	174
Absorption spectrum .....	241
AC power meter .....	118
Acceleration .....	15
AC-DC generator .....	123
Acousto-optic modulator .....	54
Action = reaction .....	19
Active power .....	132
Activity determination .....	269
AD converter .....	164
Adder .....	160
Additive colour mixing .....	171
Adiabatic exponent .....	79
Aerodynamics .....	60-62
AFM .....	291

Air resistance .....	62
Airfoil .....	62
Alloy composition .....	292
Amontons' law .....	78
Ampere, definition of .....	113
Amplifier .....	157
Amplifier, optical .....	220
Amplitude hologram .....	187
Amplitude modulation (AM) .....	135
AND-Gate .....	162
Angle of inclination .....	121
Angled projection .....	23
Angular acceleration .....	26
Angular velocity .....	26
Anharmonic oscillation .....	36
Annihilation radiation .....	269
Anomalous Hall effect .....	281
Anomalous Zeeman effect .....	249
Anomaly of water .....	67
Antenna .....	139, 170
Apparent power .....	132
Archimedes' principle .....	56
Astigmatism .....	169
Atom, size of .....	231
Atomic force microscope .....	291
Atomic size .....	278
Attenuation .....	252
Attenuation X-rays .....	253
Attenuation of $\alpha, \beta$ and $\gamma$ radiation .....	264
Autocollimation .....	167

## B

Babinet's theorem .....	176
Balmer series .....	239
Balmer-Serie .....	239
Band gap .....	281
Barrel distortions .....	169
Beam analysis, Laser .....	227

Beam deflection .....	211	Capacitive impedance.....	126, 128-129
Beam profiler.....	205	Capacitor .....	101-103, 126
Beam, Gaussian .....	205	Cathode rays .....	147
Beamsplitter .....	211	Cavendish hemispheres .....	99
Beats.....	44-45, 51	CCD Camera .....	212
Bell .....	133	Center of gravity .....	24
Bending.....	6	Central force .....	24
Bending radius.....	3	Centrifugal and centripetal force.....	27
Bernoulli equation .....	62	Chaotic oscillation .....	36
Bessel method.....	167	Characteristic radiation .....	254
Biot-Savart's law.....	114	Characteristic(s) of a diode.....	153
Bipolar transistors.....	156	Characteristic(s) of a field-effect transistor.....	156
Biprism .....	180	Characteristic(s) of a glow lamp.....	153
Birefringence.....	190, 192-193	Characteristic(s) of a light-emitting diode.....	153
Birefringent tuner .....	216	Characteristic(s) of a photoresistor .....	284
Black body.....	196	Characteristic(s) of a phototransistor .....	158
Block and tackle .....	9	Characteristic(s) of a solar battery.....	152
Bohr's magneton.....	248	Characteristic(s) of a transistor .....	156
Bohr's model of the atom.....	244-246	Characteristic(s) of a tube diode .....	140
Boyle-Mariotte's law .....	78	Characteristic(s) of a tube triode .....	141
Bragg reflection .....	254	Characteristic(s) of a varistor.....	153
Braun tube .....	143	Characteristic(s) of a Z-diode.....	153
Break-away method.....	59	Charge carrier concentration .....	281
Breit-Rabi formula.....	249	Charge distribution.....	99
Bremsstrahlung .....	254	Charge transport.....	104
Brewster angle .....	189	Charge, electric.....	89-93, 140-144
Bridge rectifier.....	156	Chromatic aberration.....	169
Brightness control .....	160	Circular motion .....	24, 26-27
Brownian motion of molecules.....	77	Circular polarization.....	41
Building materials.....	68	Circular waves.....	42
Buoyancy.....	57	Cloud chamber.....	265
<b>C</b>		Coercive force .....	289
Calcite .....	190	Coherence .....	179, 184, 208
Calliper gauge .....	3	Coherence length.....	184
Canal rays.....	147	Coherence time.....	184
Capacitance of a plate capacitor .....	101	Coil .....	127
Capacitance of a sphere .....	100	Coincidence .....	269
		Collision .....	18-19, 25

# REGISTER

Color mixing .....	212	Curve form factor .....	130
Colour filter .....	174	Cushion distortions.....	169
Colour mixing.....	171	CW value .....	61
Coma.....	169	Czerny-Turner monochromator.....	211
Comparator.....	160		
Complementary colours.....	171	<b>D</b>	
Composition of forces .....	7	DA converter.....	164
Compton effect .....	254	Damped oscillation.....	35
Compton scattering .....	258	Daniell element.....	110
Condensation heat.....	74	De Broglie wavelength.....	237
Conductivity.....	281-283	Debye temperature.....	283
Conductor, electric .....	99, 105-107, 283-284	Debye-Scherrer diffraction of electrons.....	237
Conoscopic image.....	206	Debye-Scherrer photograph .....	278
Conoscopic ray path .....	193	Debye-Sears effect.....	54
Conservation of angular momentum .....	27	Decimeter waves .....	135
Conservation of energy.....	18, 25, 27, 31	Decomposition of forces.....	7
Conservation of linear momentum.....	18-19, 25	Decomposition of white light.....	171
Constant-current source.....	151	Deflection of electrons in a magnetic field .....	142-144
Constant-voltage source .....	151	Deflection of electrons in an electric field.....	143
Control, closed-loop.....	160	Density balance .....	4
Control, open-loop .....	160	Density maximum of water.....	67
Cork-powder method.....	47	Density measuring .....	4
Coulomb's law.....	91-93	Density of air.....	4
Counter .....	163	Density of liquids .....	4
Counter tube .....	261	Density of solids.....	4
Counting rates, determination of.....	262	Detection of radioactivity.....	261
Coupled pendulums.....	37	Detection of X-rays.....	250
Coupling of oscillations.....	37	Deuterium spectrum .....	240
Cp, CV.....	79	D-Flip-Flop.....	163
Crest factor.....	132	Diamagnetism .....	288
Critical point.....	76	Dielectric constant .....	101
Cross grating .....	176	Dielectric constant of water .....	135
Crystal lattice .....	277	Differentiator .....	160
Crystal structure .....	279	Diffraction.....	208
CT.....	259	Diffraction at a crossed grating .....	176
Current source .....	151	Diffraction at a double slit.....	43, 51, 137, 176-178
Current transformation of a transformer .....	119	Diffraction at a grating.....	43, 51, 176
Curvature of image.....	169	Diffraction at a half-plane.....	178



Diffraction at a multiple grating.....	43, 51, 176-177	Eddy currents .....	118
Diffraction at a pinhole diaphragm.....	176	EDFA.....	220
Diffraction at a post.....	176	Edge absorption.....	254
Diffraction at a single slit.....	43, 51, 137, 176-178	Edge, diffraction at.....	178
Diffraction at a standing wave.....	54	Edison effect.....	140
Diffraction of electrons.....	237	Effective voltage.....	132
Diffraction of light.....	176-178	Efficiency of a heat pump.....	85
Diffraction of microwaves.....	137	Efficiency of a hot air engine.....	83
Diffraction of ultrasonic waves.....	51	Efficiency of a solar collector.....	69
Diffraction of water waves.....	43	Efficiency of a transformer.....	119
Diffraction of X-rays.....	254	Einstein coefficients.....	214
Diode.....	140, 154-156	Elastic collision.....	18-19, 25
Diode characteristic.....	140, 153-154	Elastic deformation.....	280
Diode laser.....	211	Elastic rotational collision.....	27
Directional characteristic.....	135	Elastic strain constant.....	6
Directional characteristic of antennas.....	139	Electric charge.....	89-93, 99, 140-144
Dispersion.....	206	Electric conductor.....	99, 105-107, 283-284
Dispersion of glasses.....	171	Electric current as charge transport.....	104
Dispersion of liquids.....	171	Electric energy.....	73, 131-132
Distortions.....	169	Electric field.....	94-98, 103
Doping.....	281	Electric generator.....	123
Doppler effect.....	42, 52, 213	Electric motor.....	124
Doppler effect, optical.....	225	Electric oscillator circuit.....	53
Dosimetry.....	251	Electric potential.....	96
Double mirror.....	180	Electric power.....	131
Double pendulum.....	37	Electric work.....	131
Double slit, diffraction at.....	43, 51, 137, 176-178	Electrical machines.....	122-124
DPSS Laser.....	218	Electrochemistry.....	110
Dualism of wave and particle.....	237	Electrolysis.....	109
Duane and Hunt's law.....	254	Electromagnet.....	111
Dynamic pressure.....	60	Electromagnetic oscillations.....	53
<b>E</b>		Electromechanical devices.....	133
E, determination of.....	232	Electrometer.....	90
E/m, determination of.....	144	Electron charge.....	232
Earth inductor.....	121	Electron diffraction.....	237
Earth's magnetic field.....	121	Electron holes.....	281-282, 284
Echo sounder.....	50	Electron spin.....	247-249
		Electron spin resonance.....	247

Electrostatic induction .....	89-90, 99	Fieldmill.....	96
Electrostatics.....	90	Fine beam tube .....	233
Elliptical polarization .....	190	Fine structure.....	256
Emission.....	206	Fine structure, X-ray.....	256
Emission spectra .....	243	Fixed pulley.....	9
Emission spectrum.....	241	Flame colouration.....	243
Energy loss of $\alpha$ radiation.....	268	Flame probe.....	96
Energy spectrum of X-rays.....	254	Fletcher's trolley.....	12
Energy, electrical.....	73, 131-132	Flip-Flop.....	163
Energy, heat.....	73	Flow velocity .....	213
Energy, mechanical.....	9-10, 16-19, 24, 27, 72	Fluorescence.....	174
Energy, conservation of.....	18-19, 25, 27, 31	Fluorescent screen.....	250
Energy, mechanical.....	31	Focal point, focal length.....	167
Energy-band interval .....	281	Force.....	6-9, 11, 13
Equilibrium.....	8	Force along the plane.....	10
Equilibrium of angular momentum.....	8	Force in an electric field.....	98
Equipotential surface.....	96	Force normal to the plane.....	10
ESR.....	247	Force, measuring on current-carrying conductors.....	113
Evaporation heat.....	74	Forced oscillation.....	35
Excitation of atoms.....	244-246	Foucault-Michelson method.....	198
Expansion coefficient.....	65	Fourier transformation.....	53
<b>F</b>		Franck-Hertz experiment.....	245
Fabry Perot.....	215	Fraunhofer lines.....	243
Fabry-Perot Interferometer .....	186	Free fall.....	21-23
Falling-ball viscosimeter.....	57	Frequency.....	32, 35-47, 50-53, 134-135, 137
Faraday constant .....	109	Frequency doubling.....	218
Faraday effect.....	194	Frequency modulation (FM).....	135
Faraday's cup .....	99	Frequency response.....	130
Feedback.....	134	Fresnel biprism.....	180
Ferromagnetism.....	288	Fresnel mirror.....	208
FET.....	156	Fresnel zone plate.....	208
Fiber cleaving.....	221	Fresnel's laws .....	189
Fiber laser.....	220	Fresnel's mirror .....	180
Fibre.....	174	Friction.....	10
Field effect transistor .....	156	Friction coefficient.....	11
Field emission microscope.....	277	Fringe counter.....	223
		Full-wave rectifier .....	156

<b>G</b>	
Gain profile.....	216
Galilean telescope.....	170
Galvanic element.....	110
Gas discharge.....	145
Gas discharge spectra.....	243
Gas elastic resonance apparatus.....	79
Gas laser.....	205
Gas laws.....	78
Gas thermometer.....	78
Gaussian beam.....	205
Gay-Lussac's law.....	78
Geiger counter.....	261
Geiger-Müller counter tube.....	261
Generator circuits.....	157
Generator, electric.....	123
Geometrical optics.....	167-170
Glassfiber.....	221
Glowing layer.....	146
Golden rule of mechanics.....	9
Graetz circuit.....	154
Grating spectrometer.....	203-205
Grating, diffraction at.....	43, 51, 176
Gratings.....	211
Gravitation torsion balance after Cavendish.....	5
Gravitational acceleration.....	21-22, 32-33
Gravitational constant.....	5
Gyroscope.....	28
Gyroscope, laser.....	226
<b>H</b>	
H determination.....	269
$H_{\alpha}$ -line.....	239
H, determination of.....	234-236, 234, 254
Half-life.....	126-127, 263
Half-plane, diffraction at.....	178
Half-shadow polarimeter.....	191
Half-wave rectifier.....	156
Hall effect.....	281
Hammer interrupter.....	133
Harmonic oscillation.....	33-36
Heat capacity.....	71
Heat conduction.....	68
Heat energy.....	73
Heat engine.....	81, 83-84
Heat equivalent, electric.....	73
Heat equivalent, mechanical.....	72
Heat insulation.....	68
Heat pump.....	82-83, 85
Helical spring.....	6
Helical spring after Wilberforce.....	38
Helical spring waves.....	39
Helium Neon Laser.....	216
Helmholtz coils.....	114
He-Ne laser.....	205
Heterodyne mixing.....	224
High voltage.....	120
High-temperature superconductor.....	287
Hologram.....	187
Holographic grating.....	205, 240, 242
Holography.....	210
Homogeneous electric field.....	97
Hooke's law.....	6
Hot-air engine.....	81-84
Huygens principle.....	208
Huygens' principle.....	42
Hydrogen spectrum.....	240
Hydrostatic pressure.....	55
Hyperfine structure.....	249
Hysteresis.....	289
<b>I</b>	
Ideal gas.....	78
Image charge.....	98
Image distortion.....	169
Impedance.....	126-128

# REGISTER

Implant model.....	250	Keplerian telescope .....	170
Inclined plane.....	10	Kerr effect.....	192
Independence principle .....	23	Kinetic energy.....	16
Induction .....	115-117, 122	Kinetic theory of gases.....	77-79
Inductive impedance.....	127-129	Kirchhoff's law of radiation.....	196
Inelastic collision.....	18-19, 25	Kirchhoff's laws.....	106
Inelastic electron collision.....	244-246	Kirchhoff's voltage balance .....	98
Inelastic rotational collision .....	27	Klein-Nishina formula .....	270
Integrator .....	160	Kundt's tube .....	47
Interference .....	179, 208, 215, 272		
Interference of light.....	180	<b>L</b>	
Interference of microwaves .....	137	Lambert-Beer law .....	206
Interference of ultrasonic waves.....	51	Lambert's law of radiation .....	195
Interference of water waves.....	43	Laser.....	205, 213-214, 216-220, 226
Interferometer.....	183-185	Laser beam analysis .....	227
Interferometer.....	209	Laser doppler Anemometer.....	225
Interferometer, Fabry-Perot.....	186	Laser doppler anemometry.....	213
Interferometer, heterodyne.....	224	Laser gyroscope.....	226
Interferometer, Mach-Zehnder .....	209	Laser resonator .....	216
Interferometer, Michelson.....	209	Laser safety .....	214
Internal resistance.....	108, 151-152	Laser vibrometer.....	224
Intrinsic conduction .....	281	Laser, fiber.....	220
Inverse-square law of distance.....	264	Laser, pulse .....	224
Inverting operational amplifier.....	160	Laser, pulsed .....	218
Ion dose rate .....	251	Laser, ring.....	220
Ion trap .....	238	Laserdiode.....	217
Ionization chamber.....	251	Latent heat.....	74
IR position detector.....	5	Laue.....	279
Isoelectric lines.....	95	Laue diagram.....	278
Isotope splitting .....	240	Laws of images .....	167
		Laws of radiation .....	196
<b>J</b>		LDA.....	225
JK-Flip-Flop .....	163	Leaf spring.....	6
Jones Matrix.....	206	Lecher line.....	136
		LED .....	153-155
<b>K</b>		Length measurement .....	3
$K_{\alpha}$ -line.....	254	Lens aberration.....	169
K-edge.....	253	Lenses.....	206

Leslie's cube.....	196	Magnetic moment .....	112
Lever .....	8	Magnetization curve.....	289
Lever with unequal sides .....	8	Magnetostriktion .....	183
Light emitting diode.....	153-155	Magnets.....	111
Light guide .....	174	Magnifier .....	170
Light waveguide .....	158	Maltese-cross tube .....	142
Light, velocity of.....	198-199, 201	Malus' law .....	189
Light, velocity of.....	200	Mathematical pendulum.....	32
Line spectrum.....	202-203, 239, 241	Maxwell measuring bridge.....	129
Linear air track.....	15-18	Maxwell's wheel.....	31
Linear expansion .....	65	Measuring bridge, Maxwell.....	129
Linear motion .....	11-13, 15-17	Measuring bridge, Wheatstone.....	106
Lines of force.....	94	Measuring bridge, Wien.....	129
Lines of magnetic force.....	111	Measuring range, expanding.....	108
Littrow condition .....	205	Mechanical energy .....	9, 16-19, 24-25, 27, 72
Littrow prism.....	216	Meissner-Ochsenfeld effect .....	287
Lloyd's experiment.....	43	Melde's law .....	41
Logic gate.....	162-164	Melting heat.....	74
Longitudinal laser modes.....	215	Mercury spectrum .....	242
Longitudinal waves.....	39	Metallic conductor .....	283
Loose pulley.....	9	Michelson interferometer.....	183-184, 209, 223, 272
Luminescence.....	285	Micrometer screw.....	3
Luminous intensity .....	195	Microscope .....	170
Luminous zone .....	146	Microwaves.....	137
<b>M</b>			
Machine .....	123-125	Millikan experiment .....	232
Machine(s).....	122	Mirror, Fresnel.....	208
Machine(s), electrical.....	125	Mixing temperature .....	70
Machine(s), simple .....	9	Mobility of charge carriers.....	281
Mach-Zehnder Interferometer.....	209	Mode locking.....	220
Mach-Zehnder-Interferometer .....	185	Modulation of light .....	193
Magnetic field of a coil.....	114	Modulus of elasticity .....	6
Magnetic field of Helmholtz coils.....	114	Molecular motion.....	77
Magnetic field the Earth.....	121	Molecule, size of.....	231
Magnetic field, earth .....	121	Mollier diagram .....	85
Magnetic focusing.....	142	Moment of inertia.....	28
Magnetic force .....	113	Monochromator.....	211
		Moseley's law.....	254
		Motions with reversal of direction .....	15-17

# REGISTER

Motions, one-dimensional.....	11-13, 15-17	Optical amplifier.....	220
Motions, two-dimensional.....	24	Optical analogon.....	237
Motions, uniform.....	12-13, 15-17, 26-27	Optical cavity.....	205
Motions, uniformly accelerated.....	11-13, 15-17, 26-27	Optical pumping.....	214
Motor, electric.....	124	Optical resonator.....	215
Multimeter.....	130	Optical time domain reflectometry.....	222
Multiple slit, diffraction at.....	43, 51, 176-178	Optical transmission line.....	158
Multiplexer.....	163	Optoelectronics.....	158
Muon.....	273	Orbital spin.....	248
<b>N</b>			
Nd:YAG Laser.....	218	OR-Gate.....	162
n-doped germanium.....	281	Oscillation of a string.....	46
Newton rings.....	181	Oscillation period.....	32, 34-38, 79, 134
Newton, definition of.....	13	Oscillations.....	32-38, 46, 53
Newton's experiments with white light.....	171	Oscillator.....	157
Newton's law.....	25	Oscillator circuit.....	53
Newtons rings.....	208	OTDR.....	222
Night vision.....	212	<b>P</b>	
NMR.....	267	Parallel connection of capacitors.....	101
Non-inverting operational amplifier.....	160	Parallel connection of resistors.....	106
Non-self-maintained gas discharge.....	145	Parallelogram of forces.....	7
Normal Hall effect.....	281	Paramagnetism.....	288
Normal Zeeman effect.....	248	Paths of particles.....	265
NTC resistor.....	153	Path-time diagram.....	11-13, 15-17, 26
Nuclear magnetic resonance.....	267	Paul trap.....	238
Nuclear magneton.....	249	p-doped germanium.....	281
Nuclear spin.....	249	Peak voltage.....	132
Nutation.....	28	Pendulum, amplitude.....	33
<b>O</b>			
Ohmic resistance.....	105-108	Pendulums, coupled.....	37
Ohm's law.....	105	Pendulums, mathematical and physical.....	32
Oil spot experiment.....	231	Performance number.....	85
One-sided lever.....	8	Permanent magnets.....	111
Operational amplifier.....	159	Perrin tube.....	143
Opposing force.....	25	Phase hologram.....	187
Optical activity.....	191	Phase transition.....	74-76
		Phase velocity.....	39, 41-42
		Phosphorescence.....	206
		Photoconductivity.....	284

Photodiode.....	158	pV diagram.....	81–82, 84
Photoelasticity.....	190	Pyknometer.....	4
Photoelectric effect.....	234–236	<b>Q</b>	
Photoresistor.....	153	Q-switch.....	218
Phototransistor.....	158	Quantum eraser.....	272
Physical pendulum.....	32	Quantum nature.....	187, 189–190, 196–197
PID controller.....	160	Quantum nature of charges.....	232
Pinhole diaphragm, diffraction at.....	176	Quartz, right-handed and left-handed polarization.....	191
Planck's constant.....	234–236, 254	<b>R</b>	
Plastic deformation.....	280	Radiant flux density.....	195
Plastic fibers.....	221	Radioactive dating.....	268
Plate capacitor.....	101–103	Radioactive decay.....	263
PMMA fibre.....	174	Radioactivity.....	261
Pockels effect.....	193	Reactance.....	126–128
Poisson distribution.....	262	Reactive power.....	132
Polarimeter.....	191	Real gas.....	76
Polarisation.....	207	Recoil.....	18
Polarity of electrons.....	143	Rectification.....	140, 154–155
Polarization of decimeter waves.....	135	Redox pairs.....	110
Polarization of light.....	189–194	Reflection.....	207
Polarization of microwaves.....	137	Reflection hologram.....	210
Post, diffraction at.....	176	Reflection of light.....	167
Potentiometer.....	106	Reflection of microwaves.....	137
Power plant generator.....	123	Reflection of ultrasonic waves.....	50
Power transmission of a transformer.....	120	Reflection of water waves.....	42
Precession.....	28	Reflection spectra.....	175
Pressure.....	55	Reflection, law of.....	42, 50, 167
Primary colours.....	171	Refraction.....	206
Prism spectrometer.....	202	Refraction of light.....	167
Prism, Littrow.....	216	Refraction of microwaves.....	137
Prisms.....	206	Refraction of water waves.....	42
Projection parabola.....	23	Refraction, law of.....	42
Propagation of electrons.....	142	Refractive index.....	42, 171, 185, 189, 200–201
Propagation of water waves.....	42	Refractometer.....	209
Propagation velocity of voltage pulses.....	199	Refrigerating machine.....	82
Propagation velocity of waves.....	40–42	Relay.....	133
PTC resistor.....	153		
Pulsed Laser.....	219		

# REGISTER

Remanence.....	289	SHG.....	218
Resistance.....	105	Shift register.....	163
Resistors, special.....	153	Simple machines.....	9
Resonance.....	35	Single slit, diffraction at.....	43, 51, 137, 177-178
Resonance absorption.....	247	Single slit, diffraction at a.....	176
Resonator, Laser.....	216	Slide gauge.....	3
Reversing pendulum.....	32	Sliding friction.....	11
Revolving-field generator.....	123	Slit, diffraction at.....	43, 51, 137, 176-178
Rigid body.....	24	Snellius' law.....	42, 167, 206
Ring laser.....	220	Sodium D-lines.....	203
Rocket principle.....	18	Solar battery.....	152
Rolling friction.....	11	Solar collector.....	69
Rotating the plane of polarization.....	191	Solid state Laser.....	214
Rotating-mirror method.....	198	Sound.....	53
Rotational motion.....	24-26	Sound waves.....	44-45, 47-49
Rotational oscillation.....	35	Sound, velocity of in air.....	48
Rotor.....	122-125	Sound, velocity of in gases.....	48
Rüchardt.....	79	Sound, velocity of in solids.....	49
Rutherford scattering.....	266	Spatial coherence.....	179
Rydberg constant.....	254	Special resistors.....	153
<b>S</b>		Specific conductivity.....	283
Saccharimeter.....	191	Specific electron charge.....	144, 233, 248
Safety, Laser.....	214	Specific heat.....	71
Sagnac Interferometer.....	226	Specific resistance.....	105
Saturable Absorber.....	218	Spectra, absorption.....	174
Scanning tunnelling microscope.....	290	Spectra, reflection.....	175
Scattering of $\gamma$ quanta.....	270	Spectrograph.....	211
Scintillation counter.....	269	Spectrometer.....	174-175, 202, 239, 243
Second Harmonic generation.....	218	Spectrum.....	202
Secondary colours.....	171	Speech analysis.....	53
Seebeck effect.....	286	Speed of light.....	221
Self-excited generator.....	123	Spherical aberration.....	206
Self-maintained gas discharge.....	145	Spherical aberration.....	169
Semiconductor detector.....	268	Spherometer.....	3
Semiconductors.....	283	Spin.....	247-249, 267
Series connection of capacitors.....	101	Spring.....	6
Series connection of resistors.....	106	Spring pendulum.....	34
		Standard potentials.....	110



Standing wave .....	39, 43, 47, 136-138	Tomography .....	259
Standing-wave ratio .....	138	Torsion balance .....	91
Static friction .....	10	Torsion collision .....	27
Static pressure .....	60	Torsion pendulum, Pohl .....	35
Stationary-field generator .....	123	Total pressure .....	60
Stator .....	122-125	Total reflection .....	206
Stefan-Boltzmann's law .....	196	Total reflection of microwaves .....	137
Steiner's law .....	30	Traffic-light control system .....	160
Stirling process .....	81-84	Transformer .....	119
Straight waves .....	42	Transformer under load .....	119
Subtractive colour mixing .....	171	Transistor .....	156
Subtractor .....	160	Transit time measurement .....	199
Sugar solution, concentration of .....	191	Transition temperature .....	287
Superconductivity .....	287	Translational motion .....	24
Superpositioning principle .....	23	Transmission hologram .....	188
Surface tension .....	59	Transmission of filters .....	204
Synchronous motor .....	124	Transmitter .....	135
<b>T</b>		Transversal waves .....	39
Telescope .....	170	Transverse modes .....	205
TEM modes .....	205	Triode .....	141
Temperature .....	70	Tube diode .....	140
Temperature variations .....	68	Tube triode .....	141
Terrestrial telescope .....	170	Tuning fork .....	44
Thermal emission in a vacuum .....	143	Two-beam interference .....	44
Thermal expansion of liquids .....	66	Two-dimensional motion .....	24
Thermal expansion of solid bodies .....	65	Two-point regulator .....	160
Thermal expansion of water .....	67	Two-pole rotor .....	124
Thermodynamic cycle .....	81-85	Two-pronged lightning rod .....	120
Thermoelectric voltage .....	286	Two-quantum transitions .....	249
Thermoelectricity .....	286	Two-sided lever .....	8
Thomson tube .....	144	Tyndall effect .....	189
Thread waves .....	39	<b>U</b>	
Three-phase generator .....	125	Ultrasonic waves .....	50-52
Three-phase machine .....	125	Ultrasound .....	54
Three-pole rotor .....	124	Ultrasound in liquids .....	54
Time constant L/R .....	127	Uniform acceleration .....	11-13, 15-17, 24, 26
Time constant RC .....	126	Uniform motion .....	11-13, 15-17, 24, 26-27

Universal motor.....	124	Wheatstone measuring bridge.....	106
<b>V</b>		Wheel and axle.....	8
Vapour pressure.....	75	White light.....	171
Velocity.....	11, 13, 15-17	White light reflection hologram.....	187
Velocity filter for electrons.....	144	Wien measuring bridge.....	129
Venturi tube.....	60	Wilberforce, helical spring.....	38
Verdet's constant.....	194	Wilson cloud chamber.....	265
Vernier.....	3	Wind speed.....	60
Vibrometer.....	224	Wind tunnel.....	62
VideoCom.....	13, 17, 19, 22, 37, 178, 204, 240, 242	Work, electrical.....	73, 81-84, 131-132
Viscosity.....	57	Work, mechanical.....	9-10, 16-17, 72, 81-84
Voltage amplification with a tube triode.....	141	<b>X</b>	
Voltage balance.....	98	XOR-Gate.....	162
Voltage control.....	160	X-ray contrast medium.....	250
Voltage divider.....	106	X-ray fine structure.....	256
Voltage pulse.....	115	X-ray fluorescence.....	255
Voltage series.....	110	X-ray hardening.....	252
Voltage source.....	151	X-ray photography.....	250
Voltage transformation in a transformer.....	119	X-ray scattering.....	258
Voltage-current characteristic.....	154	X-ray spectra.....	255
Volume flow.....	60	X-ray structural analysis.....	278
Volume measurement.....	4	X-ray tomography.....	259
Volumetric expansion.....	66	X-rays.....	250-255, 278-279, 292
Volumetric expansion coefficient.....	66	<b>Y</b>	
Vowel analysis.....	53	Young's experiment.....	43, 51, 137, 176-179
<b>W</b>		<b>Z</b>	
Wagner interrupter.....	133	Z-diode.....	153-155
Waltenhofen's pendulum.....	118	Zeeman effect.....	248
Water.....	67	Zone plate.....	208
Water waves.....	42		
Wave machine.....	40		
Waveguide.....	138		
Wavelength.....	39-42, 46-47, 183-184		
Wavelength multiplex.....	221		
Waveplates.....	207		
Waves.....	39-53, 135-139, 176-178, 180-185, 187-188		

P2.5.2.1  
Pressure-dependency of  
the volume of a gas at  
a constant temperature  
(Boyle-Mariotte's law)







# CONTACT

## GERMANY:

LD DIDACTIC GmbH  
Leyboldstr. 1  
50354 Huerth  
Germany  
Tel.: +49 2233 604 0  
Fax: +49 2233 604 222  
E-Mail: [info@ld-didactic.de](mailto:info@ld-didactic.de)  
[www.ld-didactic.com](http://www.ld-didactic.com)



[WWW.LD-DIDACTIC.COM](http://WWW.LD-DIDACTIC.COM)

BRANDS OF THE LD DIDACTIC GROUP

LEYBOLD® Feedback ELWE® TECHNIK